


PushTracker E2 — Battery Life Tips & Tricks



Note: Whilst Permobil has done internal testing on the battery life of the PushTracker E2 with SmartDrive, its performance can of course be affected by other functions that the device is performing (e.g., apps with sensors, updates being performed in the background, use of location services and GPS, etc). During the set up and daily use of the PushTracker E2, a few considerations to obtain maximal battery life are listed below.

During initial setup

1. Fully charge the battery the first time.
2. Check to make sure the PushTracker E2 has the most up-to-date operating system. **Perform all updates on the charger** before use.
 - Settings App > System > About > System Updates
3. Consider turning off the “Location” setting.
 - Settings App > Connectivity > Location
4. Consider choosing a Mobvoi watch face with static features (e.g., no second hand).



Other tips for battery life performance

1. Turn off the “Always-On” display.
 - Settings App > Display > Always-On Screen
2. Turn off “Tilt-To-Wake” & “Wrist Gestures”.
 - Settings App > Gestures
3. Turn off Wi-Fi when not performing updates (otherwise the PushTracker E2 will continually search for networks).
 - Settings App > Connectivity > Wi-Fi
4. Maintain a screen brightness of level 3 or less when possible.
 - Settings App > Display > Adjust Brightness
5. Turn “Ok Google” off.
 - Settings App > Personalisation > Ok Google Detection
6. Turn on “Airplane Mode” (the PushTracker E2 will prompt you to reconnect Bluetooth when going to use your SmartDrive).
 - Settings App > Connectivity > Airplane Mode
7. Monitor battery activity through the Wear OS mobile app on your phone.
 - Open Wear OS app > Scroll down to Advanced Settings > Select Watch Battery
8. Utilise a power bank.