
permobil | SEATING + POSITIONING

ROHO® DRY FLOATATION® Cushions Operation Manual

Single Compartment, Sensor Ready, ENHANCER, SELECT

Contact Information

ROHO, Inc.
100 N. Florida Avenue
Belleville, IL 62221-5429 U.S.A.

U.S.A.:
800-851-3449
Fax 888-551-3449

Customer Support:
orders.roho@permobil.com

Outside the U.S.A.:
618-277-9150
Fax 618-277-6518
intlorders.roho@permobil.com

Warning

Indicates that not following the specified procedure could lead to potentially hazardous conditions resulting in serious injury.

Caution

Indicates that not following the specified procedure could lead to potentially hazardous conditions resulting in minor to moderate injury or damage to the equipment or other property.

Incident Reporting

In the event of a serious incident or injury, contact Customer Support.

Contents

Intended Use	2
Important Safety Information	
Parts Detail	3
Product Specifications	
Cushion Setup	3-5
Troubleshooting	6
Cleaning and Disinfecting	6-7
Storage, Transport, Disposal and Recycling	7
Limited Warranty	
Español	8
Français	15
Deutsch	26
Italiano	33
Nederlands	40
Dansk	47
Svenska	54
Norsk	61
Suomi	68
Português	75
hrvatski	82
Čeština	89
Eesti	96
Ελληνικά	103
magyar	110
Íslenska	117
Lietuviškai	124
Polski	131
Română	138
Slovenčina	145
Slovenščina	152
Türkçe	159

Symbols Glossary

	Manufacturer		Fits		Size (width x depth x height)		Front of product
	Date of manufacture		Model number		Weight limit		Back of product
	Consult instructions for use		Serial number		Medical Device		Bottom of product
	DO NOT use the ISO FLO Memory Control as a handle; the product may tear.		Quantity		Authorized representative in the European Community		

Please Note

Supplier: This manual must be given to the user of this product.
Operator (Individual or Caregiver): Before using this product, read the instructions and save for future reference.

Products may be covered by one or more U.S. and foreign patents and trademarks, including CONTOUR SELECT®, DRY FLOATATION®, ENHANCER®, HIGH PROFILE®, ISO FLO®, LOW PROFILE®, MID PROFILE®, QUADTRO SELECT®, ROHO®, Sensor Ready®, shape fitting technology®, Smart Check® and Yellow Rope®. Zytel® is a registered trademark of DuPont.

Product information is changed as needed; current product information is available at permobilus.com.

ROHO, Inc. is part of Permobil Seating and Positioning.

© 2007, 2020 Permobil

Rev. 2020-07-24

The following products are featured in this operation manual:

LOW PROFILE Single Compartment Cushion

HIGH PROFILE Single Compartment Cushion

HIGH PROFILE Sensor Ready Cushion

ENHANCER Cushion

CONTOUR SELECT Cushion

QUADTRO SELECT LOW PROFILE Cushion

QUADTRO SELECT MID PROFILE Cushion

QUADTRO SELECT HIGH PROFILE Cushion

Intended Use

Each ROHO DRY FLOATATION product in this manual is an adjustable, air-filled, cellular-design wheelchair support surface. There is no weight limit, yet the cushion must be properly sized to the individual.

Each of the following products utilizes DRY FLOATATION Technology and is intended to conform to an individual's seated shape to provide skin/soft tissue protection, positioning, and an environment to facilitate wound healing.

SELECT and ENHANCER Cushions:

ROHO QUADTRO SELECT LOW PROFILE Cushion (QUADTRO SELECT LOW PROFILE); ROHO QUADTRO SELECT MID PROFILE Cushion (QUADTRO SELECT MID PROFILE); ROHO QUADTRO SELECT HIGH PROFILE Cushion (QUADTRO SELECT HIGH PROFILE); ROHO CONTOUR SELECT Cushion (CONTOUR SELECT); ROHO ENHANCER Cushion (ENHANCER)

Each of the following products utilizes DRY FLOATATION Technology and is intended to conform to an individual's seated shape to provide skin/soft tissue protection and an environment to facilitate wound healing.

Single Compartment Cushions:

ROHO LOW PROFILE Single Compartment Cushion (LOW PROFILE Single Compartment Cushion); ROHO HIGH PROFILE Single Compartment Cushion (HIGH PROFILE Single Compartment Cushion); ROHO HIGH PROFILE Single Compartment Cushion with Sensor Ready Technology (HIGH PROFILE Sensor Ready Cushion)
The Sensor Ready Cushion is compatible with Smart Check Cushion Monitoring Device.

Contraindication: Based on clinical, scientific or engineering evidence, these products may not be suitable for individuals who require accommodation of pelvic asymmetry greater than ____ in. (____ cm), as specified per product:

1 in. (2.5cm) - QUADTRO SELECT LOW PROFILE, ENHANCER, CONTOUR SELECT, LOW PROFILE Single Compartment Cushion

2 in. (5cm) - QUADTRO SELECT MID PROFILE

3 in. (7.5cm) - QUADTRO SELECT HIGH PROFILE, HIGH PROFILE Sensor Ready Cushion, HIGH PROFILE Single Compartment Cushion

A clinician who is experienced in seating and positioning should determine whether the product is appropriate for the individual's particular seating needs.

The product's intended use statement alone should not be relied upon to make that determination.

The medical products featured in this document are intended to be only one part of an overall care regimen that includes all seating and mobility equipment and therapeutic interventions. A clinician who is experienced in seating and positioning should identify this care regimen after evaluating the individual's physical needs and overall medical condition.

A clinician should also evaluate for visual, reading, and cognitive impairment to determine the need for caregiver assistance or the need for other assistive technology, such as large-print instructions, to ensure proper product use.

Important Safety Information

⚠ Warnings:

- Skin/soft tissue breakdown can occur due to a number of factors, which vary by individual. Check skin frequently, at least once a day. Redness, bruising, or darker areas (when compared to normal skin) may indicate superficial or deep tissue injury and should be addressed. If there is any discoloration to skin/soft tissue, STOP USE immediately. If the discoloration does not disappear within 30 minutes after disuse, immediately consult a healthcare professional.
- DO NOT use the product on top of, or in conjunction with, any other products or materials, except otherwise indicated in this manual. Doing so may cause the individual to become unstable and vulnerable to falling. If your cushion feels unstable, see "Troubleshooting" or consult with a clinician for proper use.
- DO NOT place obstructions between the individual and the cushion. Except for the compatible covers and accessories listed in "Product Specifications" in this manual, any items placed between the individual and the cushion, whether inside or outside the cover, will make the cushion less effective.
- The cushion and the cover MUST be compatible sizes and MUST be used as directed in this manual. IF THEY ARE NOT: 1) the benefits provided by the cushion may be reduced or eliminated, increasing the risk to skin or to other soft tissue, and 2) the individual may become unstable and vulnerable to falling.
- If you are unable to perform any tasks described in this manual, seek assistance by contacting your clinician, equipment provider, distributor, or Customer Support.
- DO NOT expose your cushion to high heat, open flames or hot ashes. Testing or certification claims, including for flammability, may no longer apply to this device when it is combined with other products or materials. Check testing and certification claims for all products used in combination with this device.

⚠ Warnings:

- The cushion surface will acclimate to the surrounding temperature. Take appropriate precautions, especially when the cushion will contact unprotected skin.
- Use caution when handling the cushion: Metal grommets on the corners may be a finger trap hazard.
- DO NOT use the cushion as a water flotation device (e.g., a Life Preserver). It will NOT support you in water.

⚠ Cautions:

- Keep the cushion away from sharp objects.
- Changes in altitude may require adjustment to the cushion.
- If the cushion has been in temperatures lower than 32°F (0°C) and exhibits unusual stiffness, allow the cushion to warm to room temperature.
- DO NOT use a pump, cover, or repair kit other than those compatible with your product.
- DO NOT modify the cushion or any components. Doing so may result in damage to the product and will void the warranty.
- DO NOT allow your cushion to come into contact with oil-based lotions or lanolin. They may degrade the material.
- Prolonged exposure to ozone may degrade materials used in the cushion, affect the performance of the cushion, and void the product warranty.
- Periodically, check the cover and components for damage and replace if necessary. See "Troubleshooting".
- DO NOT use any of the following as a handle for carrying or pulling the cushion; the product may tear: inflation valve(s), hose(s), quick disconnect (Sensor Ready Cushion), or the ISOFLO Memory Control (SELECT Cushion). Carry the product by the cushion base or by the carrying handle on the cover.


Parts Detail

Package contents: cushion, cover*; repair kit, hand pump, operation manual, limited warranty supplement, product registration card.

* Instructions are provided with the cover.

To order replacement parts, contact your product provider or Customer Support.

The serial number label is located on the bottom of the cushion.

Product Specifications

Cushion Dimensions and Weights

	SINGLE COMPARTMENT Cushions		QUADTRO SELECT Cushions			Contoured Cushions	
	LOW PROFILE	HIGH PROFILE**	LOW PROFILE	MID PROFILE	HIGH PROFILE	CONTOUR SELECT	ENHANCER
Height	2¼ in. (5.5cm)	4¼ in. (10.5cm)	2¼ in. (5.5cm)	3¼ in. (8.5cm)	4¼ in. (10.5cm)	2¼ to 4¼ in. (5.5 to 10.5cm)	2¼ to 4¼ in. (5.5 to 10.5cm)
Average Weight*	2.5 lbs. (1.1 kg)	3.25 lb. (1.5 kg)	2.5 lbs. (1.1 kg)	3.0 lbs. (1.4 kg)	3.25 lb. (1.5 kg)	3.5 lbs. (1.6 kg)	3.5 lbs. (1.6 kg)
Min. Width	6 cells; 11½ in. (29cm)	6 cells; 11½ in. (29cm)	5 cells; 11 in. (28cm)	8 cells; 16 in. (41cm)	6 cells; 12¾ in. (32.5cm)	7 cells; 14½ in. (36.5cm)	8 cells; 14¾ in. (37.5cm)
Min. Depth	6 cells; 11½ in. (29cm)	6 cells; 11½ in. (29cm)	6 cells; 11¾ in. (30cm)	8 cells; 15¼ in. (38.5cm)	6 cells; 11¼ in. (30cm)	8 cells; 15¼ in. (38.5cm)	8 cells; 14¾ in. (37.5cm)
Max. Width	13 cells; 23½ in. (59.5cm)	15 cells; 26¾ in. (68cm)	13 cells; 24½ in. (62.5cm)	11 cells; 21¼ in. (54cm)	14 cells; 26¼ in. (67cm)	13 cells; 24½ in. (62.5cm)	11 cells; 20 in. (50.5cm)
Max. Depth	11 cells; 20 in. (50.5cm)	14 cells; 23½ in. (59.5cm)	12 cells; 22 in. (56cm)	11 cells; 20¼ in. (51.5cm)	14 cells; 23¾ in. (60cm)	12 cells; 22 in. (56cm)	11 cells; 20 in. (50.5cm)

* based on 9 x 9 cell cushions: 16½ in. x 16½ in. (42cm x 42cm)

1 in. = 2.54 cm = 25.4 mm

** applies to the Sensor Ready HIGH PROFILE Cushion, except Max. Depth is 12 cells; 21¼ in (55cm).

Note: Contact Customer Support for specific cushion sizes.

Covers and Accessories Compatible with ROHO DRY FLOATATION Cushions:

Not available in all markets. Contact Customer Support.

- ROHO Heavy Duty Cushion Cover
- ROHO Standard Cushion Cover
- ROHO SELECT Cushion Cover
- ROHO Planar Solid Seat Insert
- ROHO Contour Base

Cushion materials: cushion: neoprene rubber, made without latex; inflation valve: Zytel; ISO FLO Memory Control: polypropylene, polycarbonate, aluminum; Quick Disconnect: acetal

Cover materials - Refer to the care label on the cover.

Cushion expected lifetime: 5 years

IF YOU OWN A SMART CHECK: The Sensor Ready Cushion is equipped with Sensor Ready Technology and can be used with Smart Check, which is a cushion feedback system. Smart Check can be used to find the individual's inflation level and help check their cushion. For more information, contact Customer Support.

Note: If you own Smart Check, STOP HERE and refer to the cushion setup instructions provided with Smart Check.

Cushion Setup

⚠ Warnings:

- DO NOT use a product that is underinflated or overinflated, because 1) the product benefits will be reduced or eliminated, resulting in an increased risk to skin and other soft tissue, and 2) the individual may become unstable and vulnerable to falling. Carefully follow the instructions for inflation, placement and hand check. If the product does not appear to be holding air, or if you are not able to inflate or deflate the product, see "Troubleshooting". Immediately contact your equipment provider, distributor, or Customer Support if the problem persists.
- When possible, use hook and loop fasteners to further secure the cushion to the seating surface. Hook fasteners are available upon request.
- After setting up the cushion the first time, perform a hand check frequently, at least once a day. When the cushion is used by a different individual or in a different wheelchair, repeat cushion setup. Follow the hand check instructions in this manual.



Cushion Setup is continued on the next page.

Cushion Setup, continued

Inflate, cover, and transfer.



IMPORTANT - For SELECT Cushions:

Before inflating, you must unlock the ISOFLLO Memory Control so that all of the air cells will inflate. Push the green knob to unlock the ISOFLLO Memory Control.

- 1 Inflate:** Slide the hand pump nozzle over the inflation valve. Open the inflation valve. Inflate the cushion until all of the air cells feel firm. Close the inflation valve until it comes to a stop. Remove the hand pump. Make sure the inflation valve is closed. Repeat for all inflation valves.

Note: To open, turn the inflation valve at least one full rotation. The inflation valve will spin freely when open.

- 2 Cover:** Place the cushion inside the cover. Refer to the instructions provided with the cover.

Note: Keeping a cover on your cushion will help to protect the cushion components.

3 Transfer:

- a. Place the cushion on the wheelchair with the air cells facing up and the non-skid surface of the cover facing down. The symbol showing the seated figure on the cover may be used to help orient cushion. (Fig.1)

Note: Consult your prescriber about alternate positions of inflation valve(s) and ISOFLLO Memory Control (for SELECT Cushions).

- b. The individual should transfer to the wheelchair and sit in their usual sitting position. Make sure that the cushion is centered underneath the individual. The inflation valve should be in the front, on the left side of the seated individual. (Fig. 2)


Fig. 1

Fig. 2

Perform a hand check.

- 4** To perform a hand check, use the width of your fingers to check the amount of air in the cushion. When the cushion has the correct amount of air*, the individual is both immersed in the cushion and supported by air:
- a. Slide your hand between the cushion and the individual, and feel for the lowest bony prominences. *Suggestion: To help locate the bony prominences while performing a hand check, lift and lower the leg.*
- b. Refer to the Hand Check Guide below to adjust the cushion.


* The recommended distance between the lowest bony prominence and the seat surface is
½ - 1 in. (1.5 - 2.5 cm).

IMPORTANT! Perform cushion setup with all valves. SELECT and ENHANCER Cushions may require additional positioning and adjustment. See additional instructions on the next page.

HAND CHECK GUIDE

UNDERINFLATED!

Not able to move fingers.
Not enough air in the cushion.
Risk: Not supported by air.


If the cushion is UNDERINFLATED:

Add some air. Close the inflation valve. Perform a hand check again. Repeat until the cushion has the correct amount of air.

OVERINFLATED!

Can move fingers too freely.
Too much air in the cushion.
Risk: Not immersed in the cushion.


If the cushion is OVERINFLATED:

Remove some air. Close the inflation valve. Perform a hand check again. Repeat until the cushion has the correct amount of air.

CORRECT AMOUNT OF AIR!

Able to slightly move fingers.
Immersed in the cushion and supported by air.

If your cushion has the CORRECT AMOUNT OF AIR, it is ready to use.

Cushion Setup, continued


ENHANCER

- Follow the instructions on the previous page for Cushion Setup. Explanation of inflation valves:

"A" - inflation valve for outer cells/abductor (air cells shown in white)

"B" - inflation valve for posterior contour section (air cells shown in gray)

- Make sure that both inflation valves are closed when setup is complete.


SELECT Cushion Adjustment for a Fixed or Flexible Pelvic Obliquity or Tilt

Fixed Pelvic Obliquity or Tilt:

To adjust the QUADTRO SELECT Cushion or the CONTOUR SELECT Cushion to accommodate a fixed pelvic obliquity or a fixed posterior pelvic tilt:

- Follow the Cushion Setup instructions on the previous page to make sure the cushion has the correct amount of air.
- The individual should sit in their normal position while the air transfers through the ISOFLO Memory Control. When the air transfer is complete, push the red knob toward the locked position.
- Perform a hand check for each compartment of the cushion again, as instructed on the previous page.

GREEN KNOB: Unlock

Flexible Pelvic Obliquity or Tilt:

To adjust the QUADTRO SELECT Cushion to correct a flexible pelvic obliquity or a flexible posterior pelvic tilt:

- Follow the Cushion Setup instructions on the previous page to make sure the cushion has the correct amount of air.
- Make sure the ISOFLO Memory Control is in the unlocked position
- Follow the instructions below.

RED KNOB: Lock

Flexible pelvic obliquity:

- Manually correct to the desired pelvic position. (For example, if the left side of the pelvis is lower than the right side, either have the individual lean toward the right, or place downward pressure on the right hip until the desired position is achieved.)
- Push the red knob toward the locked position.
- Perform a hand check as instructed on the previous page.

Flexible posterior pelvic tilt:

- Check to make sure there is approximately 1" (2.5 cm) of space between the individual's bony prominences and the seating surface. Add more air if needed.
- Lift both of the individual's knees to put more downward pressure on the pelvis and to force air to the front of the cushion. This will create a "well" (a lower, concave area) for the pelvis and will help prevent the pelvis from sliding forward.
- Push the red knob toward the locked position.
- Perform a hand check as instructed on the previous page.

Troubleshooting

For additional assistance, contact your equipment provider, distributor, or Customer Support.

The cushion is not holding air.	Inflate the cushion. Inspect the quick disconnect (Sensor Ready Cushion), the inflation valve(s), and hose(s) for damage. Confirm that the inflation valve(s) are completely closed. Look for holes in the cushion. If very small holes or no holes are visible, follow the instructions in the Repair Kit provided with your product. For damage to the inflation valve, quick disconnect, or hose, or for large holes or leaks in the cushion, see the Limited Warranty supplement, or contact Customer Support.
Not able to inflate or deflate the cushion.	Confirm that the inflation valve(s) are open. Inspect the inflation valve(s) and hose(s) for damage. When inflating SELECT Cushions: Confirm that the green knob of the ISOFL0 Memory Control is pushed toward the center of the ISOFL0 in the UNLOCKED position so that all chambers will inflate. For the Sensor Ready Cushion: also check the quick disconnect for damage.
The cushion is uncomfortable or unstable.	- Make sure that the cushion is not underinflated or overinflated (See Cushion Setup). - The cushion cells and the top of the cover should face up. - The cushion should be centered beneath the individual. - The cushion should be properly sized to the individual and the wheelchair. Allow the individual at least one hour to adjust to the cushion.
The cushion slides on the wheelchair.	Make sure that the non-skid bottom of the cushion cover faces away from the individual. Use directional labels on the cushion and on the cover to orient the cushion. When possible, use hook and loop fasteners to further secure the cushion to the wheelchair. Hook fasteners are available upon request. Refer to the instructions provided with the cover.
A component is damaged.	Contact Customer Support for a replacement.

Cleaning and Disinfecting

Note: For cover care instructions, refer to the manual provided with the cover.

⚠ Warnings:

- DO NOT allow water or cleaning solution to enter the cushion. Make sure the inflation valve(s) are closed.
- Sensor Ready Cushion: Make sure the Smart Check is detached from the cushion. Water in a cushion with Sensor Ready (SR) technology can result in damage to the Smart Check and cause a false reading when using Smart Check.
- Cleaning and disinfecting are separate processes. Cleaning must precede disinfection. If the product becomes soiled, or between use by different individuals: clean, disinfect, and check the product for proper functioning.
- Clean the product regularly and when it becomes soiled.
- Follow all manufacturer bleach container and germicidal disinfectant safety guidelines.
- Note about sterilization: High temperatures accelerate aging and will damage the product assembly. The products featured in this manual are not packaged in a sterile condition, nor are they intended or required to be sterilized prior to use. If institution protocol demands sterilization: First, follow the cleaning and disinfecting instructions; then, open the inflation valve(s) and use the lowest possible sterilizing temperature, but not to exceed 70°C (158°F), for the shortest time possible. DO NOT autoclave.

⚠ Cautions:

- Use of the following may cause damage to the cushion: abrasives (steel wool, scouring pads); caustic, automatic dishwasher detergents; cleaning products that contain petroleum or organic solvents, including acetone, toluene, Methyl Ethyl Ketone (MEK), naphtha, dry cleaning fluid, adhesive remover; oil-based lotions, lanolin, ozone gas; ultraviolet light cleaning methods.
- DO NOT expose to direct sunlight.

Recommendation: When washing by hand, use room-temperature water.


Note: Deflated air cells may be easier to thoroughly clean. To deflate the air cells: with the inflation valve(s) open, roll up and/or compress the cushion, and then completely close the inflation valve(s).

Hand Wash and Hand Disinfect the Cushion (Preferred Method)

Hand Wash the Cushion: Remove the cover, close the inflation valve(s), and place in a large sink. Hand wash, using mild liquid hand soap, hand dishwashing liquid, laundry detergent, or a multipurpose detergent (follow product label instructions). Gently scrub all cushion surfaces, using a soft plastic bristle brush, a sponge or a cloth. Rinse with clean water. Air dry.

Hand Disinfect the Cushion: Follow the hand wash instructions, using 1 part household liquid bleach per 9 parts water. Keep the product wet with the bleach solution for 10 minutes. Rinse thoroughly with clean water. Air dry.

⚠ Caution: Rinse thoroughly. Detergent residue may cause air cells to stick together.

Cleaning and Disinfecting, continued

Machine Wash and Disinfect the Cushion

⚠ IMPORTANT: Machine washing can cause wear and damage, including cuts and tears to the cushion and the detachment of required product information. Washing and disinfecting the cushion by hand and air drying is strongly recommended.

⚠ Warnings:

- DO NOT allow water or cleaning solution to enter the cushion. Make sure the inflation valve(s) are closed.
- Sensor Ready Cushion: Make sure the Smart Check is detached from the cushion. Water in a cushion with Sensor Ready (SR) technology can result in damage to the Smart Check and cause a false reading when using Smart Check.
- High water temperature may cause cushion discoloration and label detachment.
- DO NOT expose cushions to temperatures above 70°C (158°F), which may damage cushion components.
- For machine washing or disinfecting, make sure the washing machine capacity is sufficient. A washing machine that isn't large enough may not remove all soil or contaminants from a cushion.
- Using a damaged cushion will reduce or eliminate the benefits of the cushion and increase risk to skin and other soft tissue.

Prepare the Cushion for Machine Cleaning:

- 1) Remove the cover and close the inflation valve(s). Examine the cushion for damage, and patch any holes (refer to Repair Kit instructions provided with the product).
- 2) Remove stains or odors as needed: Place the cushion in warm water and a disinfectant cleaning agent. Gently scrub, using a soft plastic bristle brush, and rinse with clean water.
- 3) Deflate the cushion: Open the inflation valve and roll up the cushion, starting at the corner opposite the valve. After you have removed as much air as possible, close the inflation valve.
- 4) Inspect for any sharp edges in the washing machine that may damage the cushion.

Machine Wash the Cushion:

Place the cushion in a laundry bag and then into the washing machine. When folding the cushion, make sure the air cells face outward. Use warm water (85° - 105°F / 30° - 40°C) and laundry detergent or machine- and cushion-safe cleaner/biocide, placed in the detergent dispenser. Gentle cycle: Use a machine setting that does not include a high-speed spin cycle.

Machine Disinfect the Cushion: Follow the machine washing instructions and use hot water 60°C (140°F). Use biocide (refer to the manufacturer instructions). Bleach with care: Only use a dilution of 1 part household liquid bleach per 9 parts water. Rinse with clean water.

⚠ Caution: Rinse thoroughly. Detergent residue may cause air cells to stick together.

Drying the Cushion: Remove the cushion from the laundry bag and air dry. For faster drying, use a hand-held hair dryer, hold no closer than 6 in. (15 cm) from the cushion, and use a sweeping motion.

⚠ Warning: Overheating the neoprene rubber can cause permanent product damage.

Care Symbols

Hand wash.	Line dry in the shade.	Bleach (1 part bleach: 9 parts water).	Machine wash, warm, gentle.
Machine wash hot, normal, at the temperature shown.	Use a handheld dryer.		

Storage, Transport, Disposal and Recycling

Storage and Transport: Clean and disinfect the product before storing. Open the inflation valve(s). Deflate the product. Store the product in a container that protects the product from moisture, contaminants and damage. After storing, follow instructions to prepare the product for use. After transport, inspect the product to ensure it was not damaged. (see Troubleshooting).

Disposal: The components of the products in this manual are associated with no known environmental hazards when properly used and when disposed of in accordance with all local/regional regulations. At the end of its useful life, treat as healthcare waste and dispose of according to local healthcare-waste guidelines. Incinerate with other healthcare waste. Incineration must be performed by a qualified, licensed waste management facility.

Recycling: Contact a local recycling agency to determine recycling options for the product.

Limited Warranty

Limited Warranty Term from the date the product was originally purchased: SELECT Cushions: 36 months; all other DRY FLOATATION products: 24 months. The warranty does not apply to punctures, tears, burns, or misuse. See also the Limited Warranty supplement provided with your product, or contact Customer Support.

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LOW PROFILE Single Compartment-dynan

HIGH PROFILE Single Compartment-dynan

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2,5 cm (1 tum) – QUADTRO SELECT LOW PROFILE, ENHANCER, CONTOUR SELECT, LOW PROFILE Single Compartment-dyna

5 cm (2 tum) – QUADTRO SELECT MID PROFILE

7,5 cm (3 tum) – QUADTRO SELECT HIGH PROFILE, HIGH PROFILE Sensor Ready-dyna, HIGH PROFILE Single Compartment-dyna

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Single Compartment-dynor:

ROHO LOW PROFILE Single Compartment-dynan (LOW PROFILE Single Compartment-dynan); ROHO HIGH PROFILE Single Compartment-dynan (HIGH PROFILE Single Compartment-dynan);

ROHO HIGH PROFILE Single Compartment-dynan med Sensor Ready-teknologi (HIGH PROFILE Sensor Ready-dynan)

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