

New! Stander Size 1



permobil

Distributed in Canada by Permobil www.permobil.com
Vancouver Toronto Montreal Tel: 800-265-9830 Fax: 877-636-8944

NEW



Stander

It's finally here – the Stander that has it all. We're thrilled to introduce the first, size of a total redesign of our Supine and Prone Stander families. This new Stander can be configured as supine, prone or multi-position.

All models allow hip abduction up to 30 degrees to promote hip joint development. Independent adjustments for leg length accommodate contractures and discrepancies. A broad range of positioning components allow users to remain comfortably upright and active for longer. The frame can be raised and lowered to position a child at eye level with peers, or in easy reach of activities.

All this function is contained within an elegant and approachable design that centers the child while opening worlds of possibility.



Find letters of medical necessity and studies on the therapeutic benefits of standing at:
www.rifton.com/newstander

"The ability to fully customize joint positioning through the lower extremities for maximal comfort and alignment is a game changer! I have a student who only tolerated five to ten minutes of static standing in a less adjustable stander. When we tried the new Rifton Stander the student was happy and comfortable for nearly an hour, and we only got him out because we ran out of time. *"*

*Linda K. Miller, PT, DPT
MOVE International Trainer
Prospect Center, NY*

● **Supine Stander configuration** (S410)

Stander features



Supine Stander



The Stander has abduction up to 30°.



The Stander has 0°–90° gas-spring-assisted tilt adjustment.



Prone Stander

Multi-Position Stander



The Stander can be configured in prone or supine. If the multi-position configuration is ordered, the conversion from one to the other is quick and easy.

Quick and easy conversion



Prone

- 1 **Add/Remove** the upper trunk support.
- 2 **Rotate** the tray.
- 3 **Swap** the supine hip strap and prone pelvic harness.
- 4 **Add/Remove** the supine knee straps.
- 5 **Turn** the sandals 180°.



Supine



Hip abduction

The Stander supports abduction up to 30 degrees. Research has shown that weight bearing in an abducted position can improve hip joint health and development during childhood, particularly for children at risk of hip displacement.



When the frame is lowered, users can interact at eye level with their peers.



With the frame raised, users can access activities and environments.



Sit-to-stand transfers are easy in the prone configuration.



Pads and supports swing away for easy lateral transfers in the supine configuration.



Components



Knee cuffs (Prone only)

The knee cuffs provide lateral stabilization and rotational control and prevent knee hyperextension.



Head laterals (Supine only)

The head laterals mount on the headrest to provide lateral head support. They adjust independently in height, and width.



Butterfly harness (Supine only)

The butterfly harness can be used instead of the trunk strap for greater upper trunk stabilization.



Standard



Rotation control

Hip strap (Supine only)

This strap supports the user at the level of the hips. The rotation control version provides rotational stabilization.



Standard



Rotation control

Trunk strap

The trunk strap supports the user's chest. The rotation control version wraps around the user and provides rotational stabilization.





Tray

The tray can be used in prone and supine configuration. Its position and angle adjust easily, and it swings away for transfer.



Additional trunk laterals

The laterals attach at chest or hip level and swing away for transfer. The pads rotate for individual height adjustment.




Medial/lateral thigh supports

The medial/lateral thigh supports provide abduction or adduction. The mounting bracket is off-center, so the pad has two possible height positions depending which side it is installed on.



Stander dimensions

User dimensions (inches)		S410 Size 1
Height		28–42
<div></div> <div>Key user dimension: height Select the appropriate stander by the user's overall height. Choose the model that allows for growth. Important: User's weight must not exceed the maximum working load.</div>		
Item dimensions (inches)		S410 Size 1
Base length x width		32½ x 21½
Height of horizontal body support above floor		22–28
Sandal to knee height		6¾–11
Knee to hip height		6¼–10½
Sandal to top of headrest		31½–43
Width between trunk laterals		5½–10
Width between medial/lateral thigh supports		2½–4
Basic item weight PRONE (lb)		47
Basic item weight SUPINE (lb)		52
Max. working load (lb)		75



Range of adjustment



Rifton Equipment
PO Box 260
Rifton NY 12471-0260



permobil

Distributed in Canada by Permobil www.permobil.com
Vancouver Toronto Montreal Tel: 800-265-9830 Fax: 877-636-8944