

Permobil

Heavy duty products



permobil



Permobil M300 Corpus® HD

The M300 Corpus HD takes all the things you love about Permobil’s mid-wheel drive and amplifies them. Outfitted with stronger actuators, more powerful motors, and a reinforced Corpus seating system, it’s built to withstand the rigors of an active heavy duty lifestyle. Designed specifically for users up to 450 lb, the suspension has been strengthened to provide you with more stability and completed with elegant features like fully adjustable legrests and up to four power seat functions.




Specifications


Max user weight	450 lb
Max speed	5 mph
Range	13 mi (Grp 34) 20 mi (Grp 24)
Base width	26"
Base length	38"
Seat-to-floor height	18", 19", 20" (elevator) 18", 19", 20" (fixed)
Base turning radius	22"
Ground clearance	3"
Recommended battery type	Group 34 or 24
Weight including batteries	425 lb (Group 34) 447 lb (Group 24)
Drive electronics	R-Net 120A


Optional seat elevator	8"
Tilt options	Posterior: 0° – 45°
Backrest angle - power recline	85° – 150°
Legrest angle - power	85° – 175°
Seat depth	19" – 25" (by 1" incr.)
Seat width	19", 21", 23", 25"
Backrest height	21", 24" – 29" (by 1" incr.)


** Actual driving range will vary based on driving and battery conditions.
All specifications may vary based upon configuration.
Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.*


Color options



Volcano red



Pop Star pink



Sublime green


Cobalt blue


Purple rain


Platinum silver


Phantom grey


Onyx black

Permobil M3 Corpus

Meticulously engineered to help you maintain better positioning, feel more secure and reduce driving fatigue so you can conquer daily activities with ease. Its compact base offers maneuverability in tight spaces without sacrificing performance. Featuring the innovative FlexLink suspension geometry with independent caster swing-arms helping to keep you stable and in control as you traverse the obstacles in your world at up to 6 mph.




Specifications

Max user weight	330 lb
Max speed	6 mph
Range	16 mi (Grp 34) 20 mi (Grp 24)
Base width	24"
Base length	35.5"
Seat-to-floor height	17.5", 18.5", 19.5"
Base turning radius	20.5"
Ground clearance	3"
Recommended battery type	Group 34 or 24
Weight including batteries	385 lb (Group 34) 405 lb (Group 24)
Drive electronics	R-Net 120A
Optional seat elevator	12"


Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10° & 20° (up to 300 lbs)
Backrest angle - power recline	85° – 180°
Legrest angle - power	85° – 180°
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)
Backrest height	20", 23" – 28" (1" incr.)

** Actual driving range will vary based on driving and battery conditions.
All specifications may vary based upon configuration.
Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.*

Color options



Premium options



Permobil M5 Corpus

Take control with high-performance motors enabling a top speed of 7.5 mph and exceptional low-end torque that lets you pivot like a pro. Experience unsurpassed stability and ride comfort with FlexLink Pro independent drive wheel suspension. The M5 Corpus will take you to new heights with increased stability and full suspension in seat elevation.

Up to 14" of ActiveHeight™ & 20° of ActiveReach™

Base width: 25"
Base length: 37.5"



Weight capacity— 330 lb

Specifications

Max user weight	330 lb
Max speed	7.5 mph
Range	20 mi
Base width	25"
Base length	37.5"
Seat-to-floor height	17.5", 18.5", 19.5"
Base turning radius	21.5"
Ground clearance	3"
Recommended battery type	Group 24
Weight including batteries	423 lb
Drive electronics	R-Net 120A
Optional seat elevator	14"

Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10° & 20° (up to 300 lbs)
Backrest angle - power recline	85° – 180°
Legrest angle - power	85° – 180°
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)
Backrest height	20", 23" – 28" (1" incr.)

* Actual driving range will vary based on driving and battery conditions.
All specifications may vary based upon configuration.
Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.

Color options

Chili

Sunset

Canary

Spring

Aqua

Sky

Violet

Flamingo

Arctic

Midnight

Gold

Platinum

Bronze

Premium options

Permobil F3 Corpus

The F3 Corpus with a 6 mph top-end speed features optimized battery weight distribution, lower center of mass, and fully adjustable ComfortRide suspension that improve stability while traversing inclines or declines. Experience the enhanced functionality that combines outstanding chassis performance with Permobil's industry-leading Corpus seating system.



Specifications

Max user weight	330 lb	Drive electronics	R-Net 120A
Max speed	6 mph 5 mph <i>without support wheels</i>	Optional seat elevator	12"
Range	16 mi (Grp 34) 20 mi (Grp 24)	Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20° & 30°
Base width	24"	Backrest angle - power recline	85° – 180°
Base length	40" with support wheels 36" without support wheels	Legrest angle - power	85° – 180°
Seat-to-floor height	17.5", 18.5", 19.5"	Seat depth	14" – 22" (1" incr.)
Base turning radius	26.5"	Seat width	17" – 23" (2" incr.)
Ground clearance	3"	Backrest height	20", 23" – 28" (1" incr.)
Recommended battery type	Group 34 or 24	<i>* Actual driving range will vary based on driving and battery conditions. All specifications may vary based upon configuration. Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.</i>	
Weight including batteries	386 lb (Group 34) 407 lb (Group 24)		

Color options



Premium options

Permobil F5 Corpus

Power meets performance on the F5 Corpus with a 7.5 mph top-end speed and upgraded ComfortRide Pro suspension. All-wheel independent, fully adjustable, oil dampened shocks provide the ultimate in ride comfort and reduced vibration while driving in all seat positions, including an unprecedented 14” ActiveHeight seat elevation. Climb and reach higher with increased obstacle negotiation, exceptional motor torque, and enhanced stability allowing up to 45° of ActiveReach forward tilt functionality.




Specifications

Max user weight	330 lb
Max speed	7.5 mph 6 mph without support wheels
Range	20 mi
Base width	25.5"
Base length	42" with support wheels 38" without support wheels
Seat-to-floor height	17.5", 18.5", 19.5"
Base turning radius	30"
Ground clearance	3"
Recommended battery type	Group 24
Weight including batteries	412 lb
Drive electronics	R-Net 120A


Optional seat elevator	14"
Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20°, & 45°
Backrest angle - power recline	85° – 180°
Legrest angle - power	85° – 180°
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)
Backrest height	20", 23" – 28" (1" incr.)

** Actual driving range will vary based on driving and battery conditions.
All specifications may vary based upon configuration.
Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.*

Color options



Premium options



TiLite Aero T

Our most advanced aluminum chair is uniquely modeled and handmade to each user's specifications. The dual-tube frame of the Aero T provides superior rigidity, while minimizing weight and maximizing performance. Thanks to the TiFit process, you get an optimized fit plus adjustability that meets your needs.



Note: The ROHO® AGILITY® Active Carbon back supports up to 265 lb.



















Specifications

Max user weight	265 lb standard 350 lb heavy duty	Camber	0°, 2°, 4°, or 6°
Frame features	Titanium dual-tube rigid, TiFit, full adjustability	Center of gravity adjustment	5"
Total chair width*	23"	Front seat width (seat taper)	0", 1", 2", 3"
Seat width	12" – 20"	Precision lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat depth	12" – 20"	Ergonomic seat	6", 7", 8"
Front seat height	15" – 21"	Transit options	Q-Straint lockdown approved
Rear seat height	12.5" – 21"	Transport weight*	12 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components








* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

Frame colors

High gloss paint			Candy paint			Metallic paint		
								
TiLite red	Neon pink	High gloss black	Ruby red	Azure blue	Electric plum	Tangerine metallic	Acid green metallic	Ocean blue metallic
Matte textured paint			Matte paint		Pearlescent paint			
								
Flame orange	OD green	Hyper blue	Titanium grey	Black lava	Midnight blue pearl	White pearl		

Component colors

Anodized

		
Red	Orange	Gold
		
Green	Blue	Purple
		
Black		



TiLite TRA

With the power and style of its titanium dual-tube frame, plus TiFit precision for an exact fit, the TRA provides performance and adjustability made just for you.





















Weight capacity— 350 lb

Specifications

Max user weight	265 lb standard 350 lb heavy duty	Camber	0°, 2°, 4°, or 6°
Frame features	Titanium dual-tube rigid, TiFit, full adjustability	Center of gravity adjustment	5"
Total chair width*	23"	Front seat width (seat taper)	0", 1", 2", 3"
Seat width	12" – 20"	Precision lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat depth	12" – 20"	Ergonomic seat	6", 7", 8"
Front seat height	15" – 21"	Transit options	Q-Strait lockdown approved
Rear seat height	12.5" – 21"	Transport weight*	11.1 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components








* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

Frame colors

High gloss paint			Candy paint		Metallic paint			
								
TiLite red	Neon pink	High gloss black	Ruby red	Azure blue	Electric plum	Tangerine metallic	Acid green metallic	Ocean blue metallic
Matte textured paint			Matte paint		Pearlescent paint		Natural titanium	
								
Flame orange	OD green	Hyper blue	Titanium grey	Black lava	Midnight blue pearl	White pearl	Satin	Polished

Component colors

Anodized

		
Red	Orange	Gold
		
Green	Blue	Purple
		
Black		

TiLite ZRA

With full adjustability and a minimalist mono-tube frame, the ZRA dishes out style without boundaries. With its TiFit titanium frame, you have the complete package: strength, style, and the ability to set up and configure your chair just for you.



Weight capacity— 300 lb




Specifications

Max user weight	265 lb standard 300 lb heavy duty	Camber	0°, 2°, 4°, or 6°
Frame features	Titanium mono-tube rigid, TiFit, full adjustability	Center of gravity adjustment	6"
Total chair width*	23"	Front seat width (seat taper)	0", 1.5", 2.5"
Seat width	10" – 20"	Precision Lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat depth	10" – 20"	Ergonomic seat	6", 7", 8"
Front seat height	16" – 21"	Transit options	Transit tie-down
Rear seat height	12.5" – 20"	Transport weight*	10.7 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components


*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

Frame colors


High gloss paint




Matte textured paint




Candy paint



Matte paint



Metallic paint



Pearlescent paint




Natural titanium



Component colors

Anodized



TiLite TR

The sweeping curves of the titanium dual-tube frame are the perfect fusion of striking design and driven performance. Lighter and more efficient, but with strength and a vibration-dampening ride too, you can see why the TR is our flagship rigid chair.











Specifications








Max user weight	265 lb standard	Center of gravity adjustment	2.25"
Frame features	Titanium dual-tube rigid, TiFit, limited adjustability	Front seat width (seat taper)	0", 1.5", 2.5"
Total chair width*	23"	Precision lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat width	10" – 20"	Fixed seat back	Fixed height: 8" – 20" Fixed angle: 80° – 101°
Seat depth	10" – 20"	Ergonomic seat	4" – 10"
Front seat height	16" – 21"	Transit options	Transport bracket
Rear seat height	12.5" – 20"	Transport weight*	9 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components
Camber	0°, 2°, 4°, 6°, 8°, or 12°		

*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

Frame colors

High gloss paint	Candy paint	Metallic paint	
			
TiLite red	Neon pink	High gloss black	Azure blue
			Electric plum
			Tangerine metallic
			Acid green metallic
			Ocean blue metallic
Matte textured paint	Matte paint	Pearlescent paint	Natural titanium
			
Flame orange	OD green	Hyper blue	Black lava
			Midnight blue pearl
			White pearl
			Satin
			Polished

Component colors

Anodized		
		
Red	Orange	Gold
		
Green	Blue	Purple
		
Black		

TiLite Aero X

The Aero X features an innovative modular design with clean modern lines and configurability. Featuring multiple configurations that include heavy duty, hemi-height, and amputee options. Add superior rolling dynamics and full adjustability, and you get a chair that excels in both performance and price.



















Specifications

Max user weight	265 lb standard 350 lb heavy duty	Front frame angle	Fixed front: 70°, 80°, or 85° Swing-away: 70°, 80°, or 85°
Frame features	Aluminum box folding frame, modular design, full adjustability	Camber	0°, 2°, 4°, or 6°
Total chair width folded*	10.5"	Center of gravity adjustment	6"
Total chair width*	22.5"	Fixed seat back	Fixed height: 11" – 20"
Seat width	12" – 22"	Transit options	Transit tie-down, transport bracket
Seat depth	12" – 20"	Transport weight*	Fixed front: 15.6 lb Swing-away: 18 lb
Front seat height	13.5" – 21.5"	Warranty	Five years on frame One year on TiLite components
Rear seat height	13" – 20"		

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

Frame colors

High gloss paint			Candy paint		Metallic paint			
								
TiLite red	Neon pink	High gloss black	Ruby red	Azure blue	Electric plum	Tangerine metallic	Acid green metallic	Ocean blue metallic
Matte textured paint			Matte paint		Pearlescent paint			
								
Flame orange	OD green	Hyper blue	Titanium grey	Black lava	Midnight blue pearl	White pearl		



TiLite 2GX

The TiLite 2GX combines the advantages of a modular folding frame with the high strength and vibration-dampening characteristics of titanium. Add in the Tru-Fit full adjustability that meets your needs and you are ready to roll in a chair designed to empower you.



Specifications

Max user weight	265 lb standard 350 lb heavy duty	Front frame angle	Fixed front: 70°, 80°, or 85° Swing-away: 70°, 80°, or 85°
Frame features	Titanium box folding frame, modular design, full adjustability	Camber	0°, 2°, 4°, or 6°
Total chair width folded*	10.75"	Center of gravity adjustment	6"
Total chair width*	23"	Fixed seat back	Fixed height: 11" – 20"
Seat width	18" – 22"	Transit options	Transport bracket
Seat depth	16" – 20"	Transport weight*	Fixed front: 13.1 lb Swing-away: 15.5 lb
Front seat height	13.5" – 21"	Warranty	Five years on frame One year on TiLite components
Rear seat height	13" – 20"		

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

Frame colors

High gloss paint			Candy paint			Metallic paint		
TiLite red	Neon pink	High gloss black	Ruby red	Azure blue	Electric plum	Tangerine metallic	Acid green metallic	Ocean blue metallic
Matte textured paint			Matte paint		Pearlescent paint	Natural titanium		
Flame orange	OD green	Hyper blue	Titanium grey	Black lava	Midnight blue pearl	White pearl	Satin	Polished



SmartDrive MX2+

SmartDrive’s intuitive controls add a power assist for a variety of manual wheelchair scenarios like cruising up hills or over thick carpets. Its lightweight design attaches quickly and its unique OmniWheel design moves with you as you pivot, turn, wheelie, and navigate curbs in style. SmartDrive gives you freedom to do more.



Weight capacity— 331 lb

Specifications

Max user weight	30 – 331 lb
Product dimensions	Length: 15.3” Width: 5.7” Height: 9.5”
Product weight	13.2 lb drive unit
Driving range*	12.3 mi level/no resistance
Motor power	250 W brushless DC
Maximum speed	5.5 mph
Operating temp	-13°F – 122°F
Chair type compatibility	Rigid, folding, tilt-in space, one-arm drive, and standing

Drive wheel diameter	22”, 24”, 25”, 26”
Warranty	2-Year limited warranty
Bluetooth® LE	FCC ID QQQBT113 IC ID 5123A-BGTBLE113

**Performance measurements will vary based on user and chair characteristics, driving and battery conditions.*

The SmartDrive has been tested and conforms to all applicable requirements of ANSI/RESNA Standards for Wheelchairs - Volume 1 and Volume 2 and EN 12184.

A wheelchair’s specifications provided by its manufacturer may be slightly affected by the addition of the SmartDrive.

Cushions & back supports



ROHO Contour Select

No weight limit when properly sized to cushion and individual.



ROHO® Hybrid Elite® SR® Single Compartment cushion*



ROHO Hybrid Elite Dual Compartment cushion*



ROHO CONTOUR SELECT® cushion

Ideal for individuals who

May require accommodation of pelvic asymmetry, up to approx. 2"
May require pelvic and lower extremity positioning and a more stable environment for transfers

May require correction of pelvic asymmetry (medial lateral), up to approx. 2"
May require a more stable environment for transfers

May require accommodation of pelvic asymmetry, up to approx. 1"
May need correction of pelvic asymmetry (medial lateral AND anterior posterior)
May require different positioning of the left and right femur

Dimensions	Width: 13.75" – 23.75"	Width: 13.75" – 23.75"	Width: 14.5" – 24.25"
	Depth: 14.75" – 22.75"	Depth: 14.75" – 22.75"	Depth: 15.25" – 22"
	Height: 4.25" without load	Height: 4.25" without load	Height: 2.25" – 4.25"
Weight limit	Available up to 500 lb	Available up to 500 lb	No weight limit on properly sized cushions
Medicare codes	E2622, E2623	E2622, E2623	E2624, E2625

*Hybrid Elite images do not show the cushion cover, which is required for use.



ROHO QUADTRO SELECT® HIGH PROFILE® cushion



ROHO QUADTRO SELECT® MID PROFILE® cushion



ROHO ENHANCER® cushion

Ideal for individuals who

May be at risk for or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
May require accommodation or correction of medial lateral and/or anterior posterior pelvic asymmetries, up to approx. 3"
May require different positioning of the left and right femur

May be at risk for or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
May require accommodation or correction of medial lateral and/or anterior posterior pelvic asymmetries, up to approx. 2"
May require different positioning of the left and right femur

May require accommodation of pelvic asymmetry, up to approx. 1"
Require positioning of the pelvis and lower extremities
May need increased lateral support
May need to partially offload the ischial tuberosities

Dimensions	Width: 11" – 26.25"	Width: 11" – 26.25"	Width: 14.75" – 20"
	Depth: 11.75" – 23.75"	Depth: 11.75" – 23.75"	Depth: 14.75" – 20"
	Height: 4.25"	Height: 3.25"	Height: 2.25" – 4.25"
Weight limit	No weight limit on properly sized cushions	No weight limit on properly sized cushions	No weight limit on properly sized cushions
Medicare codes	E2624, E2625	E2624, E2625	E2624, E2625



ROHO HIGH PROFILE Dual
Compartment cushion



ROHO HIGH PROFILE
Single Compartment
cushion



ROHO MID PROFILE
Single Compartment
cushion

Ideal for individuals who	May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury	May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury	May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
	May require accommodation of pelvic asymmetry, up to approx. 3"	May require accommodation of pelvic asymmetry, up to approx. 3"	May require accommodation of pelvic asymmetry, up to approx. 2"
	May require correction of pelvic asymmetry (medial lateral and/or anterior posterior)	May require maximum amount of immersion and envelopment	May require immersion and envelopment
Dimensions	Width: 14.75" – 20" Depth: 14.75" – 21.75" Height: 4.25"	Width: 11.5" – 26.75" Depth: 11.5" – 21.75" Height: 4.25"	Width: 11.5" – 26.75" Depth: 11.5" – 21.75" Height: 3.25"
Weight limit	No weight limit on properly sized cushions	No weight limit on properly sized cushions	No weight limit on properly sized cushions
Medicare codes	E2624, E2625	E2622, E2623	E2622, E2623



Comfort Support Pro 7 Series
cushion



Comfort Saddle 7 Series
cushion

Ideal for individuals who	May need relief in the coccyx area	May need postural alignment
	Requires/wants cushion to strap to chair	May need relief in the coccyx area
Dimensions	Width: 22" – 30"	May require minimal to moderate skin protection
	Depth: 18"– 20"	Requires/wants cushion to strap to chairs
Weight limit	750 lb	
Medicare codes	E2602	E2608



Comfort Embrace cushion




Comfort Embrace cushion
with ATI



Comfort Ascent cushion

Ideal for individuals who	May need postural alignment	May need postural alignment	May need support for the greater trochanters
	May need support to the greater trochanters	May need support to the greater trochanters	May require mild-moderate skin protection
Dimensions	Width: 18" – 24" Depth: 14" – 24" Height: 4" for Zero elevation, 3.75" for Anti-Thrust	Width: 18" – 24" Depth: 14" – 24" Height: 4" for Zero elevation, 3.75" for Anti-Thrust	Width: 16" – 24" Depth: 16" – 20" Height: 4"
	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb
Medicare codes	E2607, E2608	E2624, 2625	E2603, E2604



Comfort M2 cushion



Comfort M2 cushion with
ATI



Vicair® Adjuster 02
cushion (standard)

Ideal for individuals who	May need coccyx relief	May need coccyx relief	May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
	May need postural alignment	May need postural alignment	Have normal or impaired sensation
Dimensions	May need support to the greater trochanters	May need support to the greater trochanters	Experience discomfort from prolonged sitting
		May need dynamic positioning	May require different positioning of the left and right femur
Weight limit	Width: 18" – 24" Depth: 14" – 24" Height: 4.25" for Zero elevation, 3.75" for Anti-Thrust	Width: 18" – 24" Depth: 14" – 24" Height: 4.25" for Zero elevation, 3.75" for Anti-Thrust	Width: 14" – 24" Depth: 14" – 22" Height: 4"
	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	550 lb
Medicare codes	E2607, E2608	E2624, E2625	E2622, E2623



Vicair Adjuster O2 cushion (low)



Vicair Active O2 cushion (low)



Vicair Vector O2 cushion (standard)

Ideal for individuals who	May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury	May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury	May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
	Have normal or impaired sensation	Have normal or impaired sensation	Have normal or impaired sensation
	Experience discomfort from prolonged sitting	May need support of the greater trochanters	May need support of the greater trochanters
	May require different positioning of the left and right femur	May require alignment/correction of the pelvis	May require alignment/correction of the pelvis
Dimensions	Width: 16" – 20"	Width: 14" – 24"	Width: 14" – 24"
	Depth: 16" – 20"	Depth: 14" – 22"	Depth: 14" – 22"
	Height: 2.4" (Kwik Strap only)	Height: 2.4"	Height: 4"
Weight limit	550 lb	440 lb	550 lb
Medicare codes	E2622, E2623	E2624, E2625	E2624, E2625



Comfort ACTA-BACK back support



Comfort ACTA-BACK Deep back support



Comfort ACTA-BACK Contour back support

Ideal for individuals who	May require a basic amount of postural support	May require a significant amount of postural support	
	Those who are more independent with transfers	Those who are more dependent with transfers	May require lower thoracic trunk support
	Those who have good trunk control	May need the maximum amount of trunk support	
Dimensions	Width: 20" – 24"	Width: 20" – 24"	Width: 20" – 24"
	Length: 14" – 22"	Length: 14" – 22"	Length: 16" – 22"
	20" wide = 300 lb	20" wide = 300 lb	20" wide = 300 lb
Weight limit	22" wide = 350 lb	22" wide = 350 lb	22" wide = 350 lb
	24" wide = 400 lb	24" wide = 400 lb	24" wide = 400 lb
Medicare codes	E2613, E2614	E2615, E2616	E2620, E2621



Comfort ACTA-BACK LTS back support



Comfort ACTA-EMBRACE back support



Comfort ACTA-RELIEF back support

Ideal for individuals who	May require a basic amount of postural support		May require a significant amount of postural support
	May require additional lateral trunk support	Complex user that truly needs a custom fit	Those who have moderate to severe kyphosis
	May require laterals that are width and height adjustable		May need the maximum amount of trunk support
Dimensions	Width: 20" – 24"	Width: 20" – 24"	Width: 20"
	Length: 14" – 22"	Length: 14" – 22"	Length: 10" – 22"
	20" wide = 300 lb	20" wide = 300 lb	
Weight limit	22" wide = 350 lb	22" wide = 350 lb	20" wide = 300 lb
	24" wide = 400 lb	24" wide = 400 lb	
Medicare codes	E2620, E2621	E2613, E2614	E2615



Not seeing your size?

The Comfort Inception cushion (E2609) and Custom back (E2617) can be built up to 28" wide. The Inception cushion can also be made thicker as needed.

Contact your Permobil rep for more details.

©2022 Permobil. The following are trademarks and registered trademarks of Permobil: ROHO®, Hybrid Elite®, CONTOUR SELECT®, QUADTRO SELECT®, HIGH PROFILE®, MID PROFILE®, ENHANCER®, ACTA-BACK®, ACTA-EMBRACE®. Vicair® is a registered trademark of Vicair.

Permobil is the distributor of Vicair O2 cushions.

Permobil has a policy of continual product improvement and reserves the right to amend this document. The current version of this document is available at permobil.com.

Medicare Coding: Decisions about how to code a product in the U.S. are the responsibility of the provider/supplier. Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process.



Terry Giles

*Custom Fabrication Technician—Power Products
U.S. Navy—Machinery Repairman First Class
(Surface Warfare Specialist/Seabee Combat Warfare Specialist)*

Proudly built in America

Behind every Permobil product is the hard work of hundreds of dedicated Americans. From the world's most advanced power and manual wheelchairs, to proven seating and positioning solutions, each Permobil product in this brochure is assembled with passion and care by employees

in one of our U.S. production locations. Permobil is proud to employ many of our nation's finest veterans and active military personnel who each day continue our mission that EVERYONE has the right to freedom, independence, and mobility.