

Permobil

TiLite TRA



Titanium



Rigid



Adjustable



Minimized
weight

permobil

Power and style combined with full adjustability

TiLite TRA is a light, efficient active wheelchair made of a titanium dual-tube frame capturing both power and style. The center of gravity can be moved forward and backward. All backrest options are adjustable. With the precision lock components, adjustments are made fast and easy.



TiLite TRA

Seat width	30,5 to 51 cm
Seat depth	30,5 to 51 cm
Seat height	Front from 38 to 53,5 cm Rear from 33 to 53,5 cm (adjustable)
Backrest height	Folding 21,5 to 53,5 cm (adjustable)
Front frame angle	70° - 80° - 85° - 90°
Backrest angle	Folding 80° - 101° (adjustable)
Camber	0° - 2° - 4° - 6° - 8° - 12°
Weight without rear wheels	Starting from 5,0 kg approx.
Wheelchair weight	Starting from 8,3 kg approx.
Maximum load	120,2 kg

Personalise your ride*

NEON
PINK

TiLite
RED

ACID GREEN
METALLIC

ELECTRIC
PLUM

BLACK
LAVA

Permobil reserves the right to change specifications and design without prior notice. *More color options available.

Actual product colors may vary.

permobil

permobil.com | customer_support_uk@permobil.com | T. +44 (0)1484 722 888



TILITETRADATASHEETEN210909