

Permobil

# Heavy duty products



**permobil**





Up to 8" of seat elevation

Base width: 26"  
Base length: 38"

## Permobil M300 Corpus® HD

The M300 Corpus HD takes all the things you love about Permobil's mid-wheel drive and amplifies them. Outfitted with stronger actuators, more powerful motors, and a reinforced Corpus seating system, it's built to withstand the rigors of an active heavy duty lifestyle. Designed specifically for users up to 450 lb, the suspension has been strengthened to provide you with more stability and completed with elegant features like fully adjustable legrests and up to four power seat functions.



Weight capacity— 450 lb


### Specifications

<b>Max user weight</b>	450 lb
<b>Max speed</b>	5 mph
<b>Range</b>	13 mi (Grp 34)   20 mi (Grp 24)
<b>Base width</b>	26"
<b>Base length</b>	38"
<b>Seat-to-floor height</b>	18", 19", 20" (elevator) 18", 19", 20" (fixed)
<b>Base turning radius</b>	22"
<b>Ground clearance</b>	3"
<b>Recommended battery type</b>	Group 34 or 24
<b>Weight including batteries</b>	425 lb (Group 34) 447 lb (Group 24)
<b>Drive electronics</b>	R-Net 120A

<b>Optional seat elevator</b>	8"
<b>Tilt options</b>	Posterior: 0° – 45°
<b>Backrest angle - power recline</b>	85° – 150°
<b>Legrest angle - power</b>	85° – 175°
<b>Seat depth</b>	19" – 25" (by 1" incr.)
<b>Seat width</b>	19", 21", 23", 25"
<b>Backrest height</b>	21", 24" – 29" (by 1" incr.)

\* Actual driving range will vary based on driving and battery conditions.  
All specifications may vary based upon configuration.  
Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.

### Color options



Volcano red   Pop Star pink   Sublime green   Cobalt blue   Purple rain   Platinum silver   Phantom grey   Onyx black



# Permobil M3 Corpus

Meticulously engineered to help you maintain better positioning, feel more secure and reduce driving fatigue so you can conquer daily activities with ease. Its compact base offers maneuverability in tight spaces without sacrificing performance. Featuring the innovative FlexLink suspension geometry with independent caster swing-arms helping to keep you stable and in control as you traverse the obstacles in your world at up to 6 mph.



## Specifications

<b>Max user weight</b>	330 lb
<b>Max speed</b>	6 mph
<b>Range</b>	16 mi (Grp 34)   20 mi (Grp 24)
<b>Base width</b>	24"
<b>Base length</b>	35.5"
<b>Seat-to-floor height</b>	17.5", 18.5", 19.5"
<b>Base turning radius</b>	20.5"
<b>Ground clearance</b>	3"
<b>Recommended battery type</b>	Group 34 or 24
<b>Weight including batteries</b>	385 lb (Group 34) 405 lb (Group 24)
<b>Drive electronics</b>	R-Net 120A
<b>Optional seat elevator</b>	12"


<b>Tilt options</b>	Posterior: 0° – 50° ActiveReach: 5°, 10° & 20° (up to 300 lb)
<b>Backrest angle - power recline</b>	85° – 180°
<b>Legrest angle - power</b>	85° – 180°
<b>Seat depth</b>	14" – 22" (1" incr.)
<b>Seat width</b>	17" – 23" (2" incr.)
<b>Backrest height</b>	20", 23" – 28" (1" incr.)

*\* Actual driving range will vary based on driving and battery conditions. All specifications may vary based upon configuration. Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.*

## Color options



## Premium options





# Permobil M5 Corpus

Take control with high-performance motors enabling a top speed of 7.5 mph and exceptional low-end torque that lets you pivot like a pro. Experience unsurpassed stability and ride comfort with FlexLink Pro independent drive wheel suspension. The M5 Corpus will take you to new heights with increased stability and full suspension in seat elevation.

Up to 14" of ActiveHeight™ & 20° of ActiveReach™

Base width: 25"  
Base length: 37.5"



Weight capacity— 330 lb

## Specifications

Max user weight	330 lb
Max speed	7.5 mph
Range	20 mi
Base width	25"
Base length	37.5"
Seat-to-floor height	17.5", 18.5", 19.5"
Base turning radius	21.5"
Ground clearance	3"
Recommended battery type	Group 24
Weight including batteries	423 lb
Drive electronics	R-Net 120A
Optional seat elevator	14"


Tilt options	Posterior: 0° – 50°
	ActiveReach: 5°, 10° & 20° (up to 300 lb)
Backrest angle - power recline	85° – 180°
Legrest angle - power	85° – 180°
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)
Backrest height	20", 23" – 28" (1" incr.)

\* Actual driving range will vary based on driving and battery conditions. All specifications may vary based upon configuration. Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.

## Color options



## Premium options





# Permobil F3 Corpus

The F3 Corpus with a 6 mph top-end speed features optimized battery weight distribution, lower center of mass, and fully adjustable ComfortRide suspension that improve stability while traversing inclines or declines. Experience the enhanced functionality that combines outstanding chassis performance with Permobil's industry-leading Corpus seating system.



## Specifications

<b>Max user weight</b>	330 lb	<b>Drive electronics</b>	R-Net 120A
<b>Max speed</b>	6 mph 5 mph <i>without support wheels</i>	<b>Optional seat elevator</b>	12"
<b>Range</b>	16 mi (Grp 34)   20 mi (Grp 24)	<b>Tilt options</b>	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20° & 30°
<b>Base width</b>	24"	<b>Backrest angle - power recline</b>	85° – 180°
<b>Base length</b>	40" with support wheels 36" without support wheels	<b>Legrest angle - power</b>	85° – 180°
<b>Seat-to-floor height</b>	17.5", 18.5", 19.5"	<b>Seat depth</b>	14" – 22" (1" incr.)
<b>Base turning radius</b>	26.5"	<b>Seat width</b>	17" – 23" (2" incr.)
<b>Ground clearance</b>	3"	<b>Backrest height</b>	20", 23" – 28" (1" incr.)
<b>Recommended battery type</b>	Group 34 or 24	<i>* Actual driving range will vary based on driving and battery conditions. All specifications may vary based upon configuration. Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.</i>	
<b>Weight including batteries</b>	386 lb (Group 34) 407 lb (Group 24)		

## Color options




## Premium options



# Permobil F5 Corpus

Power meets performance on the F5 Corpus with a 7.5 mph top-end speed and upgraded ComfortRide Pro suspension. All-wheel independent, fully adjustable, oil dampened shocks provide the ultimate in ride comfort and reduced vibration while driving in all seat positions, including an unprecedented 14" ActiveHeight seat elevation. Climb and reach higher with increased obstacle negotiation, exceptional motor torque, and enhanced stability allowing up to 45° of ActiveReach forward tilt functionality.




## Specifications

<b>Max user weight</b>	330 lb
<b>Max speed</b>	7.5 mph 6 mph without support wheels
<b>Range</b>	20 mi
<b>Base width</b>	25.5"
<b>Base length</b>	42" with support wheels 38" without support wheels
<b>Seat-to-floor height</b>	17.5", 18.5", 19.5"
<b>Base turning radius</b>	30"
<b>Ground clearance</b>	3"
<b>Recommended battery type</b>	Group 24
<b>Weight including batteries</b>	412 lb
<b>Drive electronics</b>	R-Net 120A

<b>Optional seat elevator</b>	14"
<b>Tilt options</b>	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20°, & 45°
<b>Backrest angle - power recline</b>	85° – 180°
<b>Legrest angle - power</b>	85° – 180°
<b>Seat depth</b>	14" – 22" (1" incr.)
<b>Seat width</b>	17" – 23" (2" incr.)
<b>Backrest height</b>	20", 23" – 28" (1" incr.)


*\* Actual driving range will vary based on driving and battery conditions. All specifications may vary based upon configuration. Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.*

## Color options



Chili Sunset Canary Spring Aqua Sky Violet Flamingo Arctic Midnight

## Premium options



Gold Platinum Bronze



# TiLite Aero T

Our most advanced aluminum chair is uniquely modeled and handmade to each user's specifications. The dual-tube frame of the Aero T provides superior rigidity, while minimizing weight and maximizing performance. Thanks to the TiFit process, you get an optimized fit plus adjustability that meets your needs.



Note: The ROHO® AGILITY™ Active Carbon back supports up to 265 lb.




## Specifications

<b>Max user weight</b>	265 lb standard 350 lb heavy duty	<b>Camber</b>	0°, 2°, 4°, or 6°
<b>Frame features</b>	Titanium dual-tube rigid, TiFit, full adjustability	<b>Center of gravity adjustment</b>	5"
<b>Total chair width*</b>	23"	<b>Front seat width (seat taper)</b>	0", 1", 2", 3"
<b>Seat width</b>	12" – 20"	<b>Precision lock folding seat back</b>	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
<b>Seat depth</b>	12" – 20"	<b>Ergonomic seat</b>	6", 7", 8"
<b>Front seat height</b>	15" – 21"	<b>Transit options</b>	Q-Straint lockdown approved
<b>Rear seat height</b>	12.5" – 21"	<b>Transport weight*</b>	12 lb
<b>Front frame angle</b>	70°, 75°, 80°, 85°, or 90°	<b>Warranty</b>	Five years on frame One year on TiLite components

\*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


## Frame colors

### High gloss paint




TiLite red Neon pink High gloss black

### Candy paint



Ruby red Azure blue Electric plum

### Metallic paint



Tangerine metallic Acid green metallic Ocean blue metallic

### Matte textured paint



Flame orange OD green Hyper blue

### Matte paint



Titanium grey Black lava


### Pearlescent paint



Midnight blue pearl White pearl

## Component colors

### Anodized



Red Orange Gold



Green Blue Purple



Black





# TiLite TRA

With the power and style of its titanium dual-tube frame, plus TiFit precision for an exact fit, the TRA provides performance and adjustability made just for you.



Weight capacity— 350 lb


## Specifications

<b>Max user weight</b>	265 lb standard 350 lb heavy duty	<b>Camber</b>	0°, 2°, 4°, or 6°
<b>Frame features</b>	Titanium dual-tube rigid, TiFit, full adjustability	<b>Center of gravity adjustment</b>	5"
<b>Total chair width*</b>	23"	<b>Front seat width (seat taper)</b>	0", 1", 2", 3"
<b>Seat width</b>	12" – 20"	<b>Precision lock folding seat back</b>	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
<b>Seat depth</b>	12" – 20"	<b>Ergonomic seat</b>	6", 7", 8"
<b>Front seat height</b>	15" – 21"	<b>Transit options</b>	Q-Strait lockdown approved
<b>Rear seat height</b>	12.5" – 21"	<b>Transport weight*</b>	11.1 lb
<b>Front frame angle</b>	70°, 75°, 80°, 85°, or 90°	<b>Warranty</b>	Five years on frame One year on TiLite components

\* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


## Frame colors

### High gloss paint




TiLite red Neon pink High gloss black

### Candy paint




Ruby red Azure blue Electric plum

### Metallic paint




Tangerine metallic Acid green metallic Ocean blue metallic

### Matte textured paint




Flame orange OD green Hyper blue

### Matte paint




Titanium grey Black lava Midnight blue pearl

### Pearlescent paint



White pearl Satin Polished


### Natural titanium



Green Blue Purple

## Component colors

### Anodized



Red Orange Gold  
Green Blue Purple  
Black



# TiLite ZRA

With full adjustability and a minimalist mono-tube frame, the ZRA dishes out style without boundaries. With its TiFit titanium frame, you have the complete package: strength, style, and the ability to set up and configure your chair just for you.



Weight capacity— 300 lb




## Specifications

<b>Max user weight</b>	265 lb standard 300 lb heavy duty	<b>Camber</b>	0°, 2°, 4°, or 6°
<b>Frame features</b>	Titanium mono-tube rigid, TiFit, full adjustability	<b>Center of gravity adjustment</b>	6"
<b>Total chair width*</b>	23"	<b>Front seat width (seat taper)</b>	0", 1.5", 2.5"
<b>Seat width</b>	10" – 20"	<b>Precision Lock folding seat back</b>	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
<b>Seat depth</b>	10" – 20"	<b>Ergonomic seat</b>	6", 7", 8"
<b>Front seat height</b>	16" – 21"	<b>Transit options</b>	Transit tie-down
<b>Rear seat height</b>	12.5" – 20"	<b>Transport weight*</b>	10.7 lb
<b>Front frame angle</b>	70°, 75°, 80°, 85°, or 90°	<b>Warranty</b>	Five years on frame One year on TiLite components

\*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


## Frame colors

### High gloss paint




TiLite red Neon pink High gloss black

### Candy paint




Ruby red Azure blue Electric plum

### Metallic paint




Tangerine metallic Acid green metallic Ocean blue metallic

### Matte textured paint




Flame orange OD green Hyper blue

### Matte paint




Titanium grey Black lava Midnight blue pearl

### Pearlescent paint



White pearl Satin Polished


### Natural titanium



Titanium grey Black lava Midnight blue pearl

## Component colors

### Anodized



Red Orange Gold  
Green Blue Purple  
Black





## TiLite ZR

With the ZR's minimalist mono-tube frame, less truly is more. The clean aesthetic of the ZR turns heads while the titanium frame makes sure you have the strength and vibration-dampening ride to go with your beauty. Style and strength to take on whatever life throws at you.



Weight capacity— 300 lb


## Specifications

<b>Max user weight</b>	265 lb standard 300 lb heavy duty	<b>Center of gravity adjustment</b>	2.25"
<b>Frame features</b>	Titanium dual-tube rigid, TiFit, limited adjustability	<b>Front seat width (seat taper)</b>	0", 1.5", 2.5"
<b>Total chair width*</b>	23"	<b>Precision lock folding seat back</b>	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
<b>Seat width</b>	10" – 20"	<b>Fixed seat back</b>	Fixed height: 8" – 20" Fixed angle: 80° – 101°
<b>Seat depth</b>	10" – 20"	<b>Ergonomic seat</b>	4" – 10"
<b>Front seat height</b>	16" – 21"	<b>Transit options</b>	Transport bracket
<b>Rear seat height</b>	12.5" – 20"	<b>Transport weight*</b>	9 lb
<b>Front frame angle</b>	70°, 75°, 80°, 85°, or 90°	<b>Warranty</b>	Five years on frame One year on TiLite components
<b>Camber</b>	0°, 2°, 4°, 6°, 8°, or 12°		

\* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


## Frame colors

### High gloss paint




TiLite red Neon pink High gloss black

### Candy paint




Ruby red Azure blue Electric plum

### Metallic paint




Tangerine metallic Acid green metallic Ocean blue metallic

### Matte textured paint




Flame orange OD green Hyper blue

### Matte paint




Titanium grey Black lava Midnight blue pearl

### Pearlescent paint



White pearl Satin Polished


### Natural titanium



Black

## Component colors

### Anodized



Red Orange Gold  
Green Blue Purple  
Black



# TiLite Aero X

The Aero X features an innovative modular design with clean modern lines and configurability. Featuring multiple configurations that include heavy duty, hemi-height, and amputee options. Add superior rolling dynamics and full adjustability, and you get a chair that excels in both performance and price.




## Specifications

<b>Max user weight</b>	265 lb standard 350 lb heavy duty	<b>Front frame angle</b>	Fixed front: 70°, 80°, or 85° Swing-away: 70°, 80°, or 85°
<b>Frame features</b>	Aluminum box folding frame, modular design, full adjustability	<b>Camber</b>	0°, 2°, 4°, or 6°
<b>Total chair width folded*</b>	10.5"	<b>Center of gravity adjustment</b>	6"
<b>Total chair width*</b>	22.5"	<b>Fixed seat back</b>	Fixed height: 11" – 20"
<b>Seat width</b>	12" – 22"	<b>Transit options</b>	Transit tie-down, transport bracket
<b>Seat depth</b>	12" – 20"	<b>Transport weight*</b>	Fixed front: 15.6 lb Swing-away: 18 lb
<b>Front seat height</b>	13.5" – 21.5"	<b>Warranty</b>	Five years on frame One year on TiLite components
<b>Rear seat height</b>	13" – 20"		

\* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.


## Frame colors

### High gloss paint




TiLite red Neon pink High gloss black

### Candy paint




Ruby red Azure blue Electric plum

### Metallic paint




Tangerine metallic Acid green metallic Ocean blue metallic

### Matte textured paint



Flame orange OD green Hyper blue

### Matte paint



Titanium grey Black lava Midnight blue pearl

### Pearlescent paint



White pearl





# TiLite 2GX

The TiLite 2GX combines the advantages of a modular folding frame with the high strength and vibration-dampening characteristics of titanium. Add in the Tru-Fit full adjustability that meets your needs and you are ready to roll in a chair designed to empower you.



Weight capacity— 350 lb

## Specifications

<b>Max user weight</b>	265 lb standard 350 lb heavy duty	<b>Front frame angle</b>	Fixed front: 70°, 80°, or 85° Swing-away: 70°, 80°, or 85°
<b>Frame features</b>	Titanium box folding frame, modular design, full adjustability	<b>Camber</b>	0°, 2°, 4°, or 6°
<b>Total chair width folded*</b>	10.75"	<b>Center of gravity adjustment</b>	6"
<b>Total chair width*</b>	23"	<b>Fixed seat back</b>	Fixed height: 11" – 20"
<b>Seat width</b>	18" – 22"	<b>Transit options</b>	Transport bracket
<b>Seat depth</b>	16" – 20"	<b>Transport weight*</b>	Fixed front: 13.1 lb Swing-away: 15.5 lb
<b>Front seat height</b>	13.5" – 21"	<b>Warranty</b>	Five years on frame One year on TiLite components
<b>Rear seat height</b>	13" – 20"		

\* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

## Frame colors

<b>High gloss paint</b>			<b>Candy paint</b>			<b>Metallic paint</b>					
TiLite red	Neon pink	High gloss black	Ruby red	Azure blue	Electric plum	Tangerine metallic	Acid green metallic	Ocean blue metallic			
<b>Matte textured paint</b>			<b>Matte paint</b>			<b>Pearlescent paint</b>			<b>Natural titanium</b>		
Flame orange	OD green	Hyper blue	Titanium grey	Black lava	Midnight blue pearl	White pearl	Satin	Polished			





## SmartDrive MX2+

SmartDrive's intuitive controls add a power assist for a variety of manual wheelchair scenarios like cruising up hills or over thick carpets. Its lightweight design attaches quickly and its unique OmniWheel design moves with you as you pivot, turn, wheelie, and navigate curbs in style. SmartDrive gives you freedom to do more.



Weight capacity— 331 lb

### Specifications

<b>Max user weight</b>	30 – 331 lb
<b>Product dimensions</b>	Length: 15.3" Width: 5.7" Height: 9.5"
<b>Product weight</b>	13.2 lb drive unit
<b>Driving range*</b>	12.3 mi level/no resistance
<b>Motor power</b>	250 W brushless DC
<b>Maximum speed</b>	5.5 mph
<b>Operating temp</b>	-13°F – 122°F
<b>Chair type compatibility</b>	Rigid, folding, tilt-in-space, one-arm drive, and standing

<b>Drive wheel diameter</b>	22", 24", 25", 26"
<b>Warranty</b>	2-Year limited warranty
<b>Bluetooth® LE</b>	FCC ID QOQBT113 IC ID 5123A-BGTBLE113

*\*Performance measurements will vary based on user and chair characteristics, driving and battery conditions.*

*The SmartDrive has been tested and conforms to all applicable requirements of ANSI/RESNA Standards for Wheelchairs - Volume 1 and Volume 2 and EN 12184.*

*\*\*A wheelchair's specifications provided by its manufacturer may be slightly affected by the addition of the SmartDrive.\*\**



# Cushions & back supports



## ROHO Contour Select

No weight limit when properly sized to cushion and individual.



## ROHO® Hybrid Elite® SR® Single Compartment cushion\*

**Ideal for individuals who** May require accommodation of pelvic asymmetry, up to approx. 2"  
May require pelvic and lower extremity positioning and a more stable environment for transfers

<b>Dimensions</b>	Width: 13.75" – 23.75"
	Depth: 14.75" – 22.75"
	Height: 4.25" without load
<b>Weight limit</b>	Available up to 500 lb
<b>Medicare codes</b>	E2622, E2623



## ROHO Hybrid Elite Dual Compartment cushion\*

May require correction of pelvic asymmetry (medial lateral), up to approx. 2"  
May require a more stable environment for transfers

<b>Dimensions</b>	Width: 13.75" – 23.75"
	Depth: 14.75" – 22.75"
	Height: 4.25" without load
<b>Weight limit</b>	Available up to 500 lb
<b>Medicare codes</b>	E2622, E2623



## ROHO CONTOUR SELECT® cushion

May require accommodation of pelvic asymmetry, up to approx. 1"  
May need correction of pelvic asymmetry (medial lateral AND anterior posterior)  
May require different positioning of the left and right femur

<b>Dimensions</b>	Width: 14.5" – 24.25"
	Depth: 15.25" – 22"
	Height: 2.25" – 4.25"
<b>Weight limit</b>	No weight limit on properly sized cushions
<b>Medicare codes</b>	E2624, E2625

\*Hybrid Elite images do not show the cushion cover, which is required for use.



## ROHO QUADTRO SELECT® HIGH PROFILE® cushion

**Ideal for individuals who** May be at risk for or currently have skin/soft tissue breakdown and/or any category/stage pressure injury  
May require accommodation or correction of medial lateral and/or anterior posterior pelvic asymmetries, up to approx. 3"  
May require different positioning of the left and right femur

<b>Dimensions</b>	Width: 11" – 26.25"
	Depth: 11.75" – 23.75"
	Height: 4.25"
<b>Weight limit</b>	No weight limit on properly sized cushions
<b>Medicare codes</b>	E2624, E2625



## ROHO QUADTRO SELECT® MID PROFILE® cushion

May be at risk for or currently have skin/soft tissue breakdown and/or any category/stage pressure injury  
May require accommodation or correction of medial lateral and/or anterior posterior pelvic asymmetries, up to approx. 2"  
May require different positioning of the left and right femur

<b>Dimensions</b>	Width: 11" – 26.25"
	Depth: 11.75" – 23.75"
	Height: 3.25"
<b>Weight limit</b>	No weight limit on properly sized cushions
<b>Medicare codes</b>	E2624, E2625



## ROHO ENHANCER® cushion

May require accommodation of pelvic asymmetry, up to approx. 1"  
Require positioning of the pelvis and lower extremities  
May need increased lateral support  
May need to partially offload the ischial tuberosities

<b>Dimensions</b>	Width: 14.75" – 20"
	Depth: 14.75" – 20"
	Height: 2.25" – 4.25"
<b>Weight limit</b>	No weight limit on properly sized cushions
<b>Medicare codes</b>	E2624, E2625





**ROHO HIGH PROFILE Dual Compartment cushion**



**ROHO HIGH PROFILE Single Compartment cushion**



**ROHO MID PROFILE Single Compartment cushion**

**Ideal for individuals who**

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury  
 May require accommodation of pelvic asymmetry, up to approx. 3"  
 May require correction of pelvic asymmetry (medial lateral and/or anterior posterior)

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury  
 May require accommodation of pelvic asymmetry, up to approx. 3"  
 May require maximum amount of immersion and envelopment

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury  
 May require accommodation of pelvic asymmetry, up to approx. 2"  
 May require immersion and envelopment

<b>Dimensions</b>	Width: 14.75" – 20" Depth: 14.75" – 21.75" Height: 4.25"	Width: 11.5" – 26.75" Depth: 11.5" – 21.75" Height: 4.25"	Width: 11.5" – 26.75" Depth: 11.5" – 21.75" Height: 3.25"
<b>Weight limit</b>	No weight limit on properly sized cushions	No weight limit on properly sized cushions	No weight limit on properly sized cushions
<b>Medicare codes</b>	E2624, E2625	E2622, E2623	E2622, E2623



**Comfort Support Pro 7 Series cushion**



**Comfort Saddle 7 Series cushion**

**Ideal for individuals who**

May need relief in the coccyx area  
 Requires/wants cushion to strap to chair

May need postural alignment  
 May need relief in the coccyx area  
 May require minimal to moderate skin protection  
 Requires/wants cushion to strap to chairs

<b>Dimensions</b>	Width: 22" – 30" Depth: 18" – 20"	Width: 22" – 30" Depth: 18" – 20"
<b>Weight limit</b>	750 lb	750 lb
<b>Medicare codes</b>	E2602	E2608



**Comfort Embrace cushion**



**Comfort Embrace cushion with ATI**



**Comfort Ascent cushion**

**Ideal for individuals who**

May need postural alignment  
 May need support to the greater trochanters

May need postural alignment  
 May need support to the greater trochanters  
 May need dynamic positioning

May need support for the greater trochanters  
 May require mild-moderate skin protection

<b>Dimensions</b>	Width: 18" – 24" Depth: 14" – 24" Height: 4" for Zero elevation, 3.75" for Anti-Thrust	Width: 18" – 24" Depth: 14" – 24" Height: 4" for Zero elevation, 3.75" for Anti-Thrust	Width: 16" – 24" Depth: 16" – 20" Height: 4"
<b>Weight limit</b>	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb
<b>Medicare codes</b>	E2607, E2608	E2624, 2625	E2603, E2604



**Comfort M2 cushion**



**Comfort M2 cushion with ATI**



**Vicair® Adjuster O2 cushion (standard)**

**Ideal for individuals who**

May need coccyx relief  
 May need postural alignment  
 May need support to the greater trochanters

May need coccyx relief  
 May need postural alignment  
 May need support to the greater trochanters  
 May need dynamic positioning

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury  
 Have normal or impaired sensation  
 Experience discomfort from prolonged sitting  
 May require different positioning of the left and right femur

<b>Dimensions</b>	Width: 18" – 24" Depth: 14" – 24" Height: 4.25" for Zero elevation, 3.75" for Anti-Thrust	Width: 18" – 24" Depth: 14" – 24" Height: 4.25" for Zero elevation, 3.75" for Anti-Thrust	Width: 14" – 24" Depth: 14" – 22" Height: 4"
<b>Weight limit</b>	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	18" wide = 300 lb 20" wide = 350 lb 22" wide or > = 400 lb	550 lb
<b>Medicare codes</b>	E2607, E2608	E2624, E2625	E2622, E2623





**Vicair Adjuster O2 cushion (low)**

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury  
 Have normal or impaired sensation  
 Experience discomfort from prolonged sitting  
 May require different positioning of the left and right femur

**Ideal for individuals who**

<b>Dimensions</b>	Width: 16" – 20" Depth: 16" – 20" Height: 2.4" (Kwik Strap only)
<b>Weight limit</b>	550 lb
<b>Medicare codes</b>	E2622, E2623



**Vicair Active O2 cushion (low)**

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury  
 Have normal or impaired sensation  
 May need postural alignment  
 May need support of the greater trochanters  
 May require alignment/correction of the pelvis  
 May need dynamic positioning  
 Those who want easy transfers

<b>Dimensions</b>	Width: 14" – 24" Depth: 14" – 22" Height: 2.4"
<b>Weight limit</b>	440 lb
<b>Medicare codes</b>	E2624, E2625



**Vicair Vector O2 cushion (standard)**

May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury  
 Have normal or impaired sensation  
 May need postural alignment  
 May need support of the greater trochanters  
 May require alignment/correction of the pelvis  
 May need dynamic positioning

<b>Dimensions</b>	Width: 14" – 24" Depth: 14" – 22" Height: 4"
<b>Weight limit</b>	550 lb
<b>Medicare codes</b>	E2624, E2625



**Comfort ACTA-BACK LTS back support**

May require a basic amount of postural support  
 May require additional lateral trunk support  
 May require laterals that are width and height adjustable

**Ideal for individuals who**

<b>Dimensions</b>	Width: 20" – 24" Length: 14" – 22"	Width: 20" – 24" Length: 14" – 22"	Width: 20" Length: 10" – 22"
<b>Weight limit</b>	20" wide = 300 lb 22" wide = 350 lb 24" wide = 400 lb	20" wide = 300 lb 22" wide = 350 lb 24" wide = 400 lb	20" wide = 300 lb
<b>Medicare codes</b>	E2620, E2621	E2613, E2614	E2615



**Comfort ACTA-EMBRACE® back support**

Complex user that truly needs a custom fit



**Comfort ACTA-RELIEF® back support**

May require a significant amount of postural support  
 Those who have moderate to severe kyphosis  
 May need the maximum amount of trunk support



**Comfort ACTA-BACK® back support**

May require a basic amount of postural support  
 Those who are more independent with transfers  
 Those who have good trunk control

**Ideal for individuals who**

<b>Dimensions</b>	Width: 20" – 24" Length: 14" – 22"	Width: 20" – 24" Length: 14" – 22"	Width: 20" – 24" Length: 16" – 22"
<b>Weight limit</b>	20" wide = 300 lb 22" wide = 350 lb 24" wide = 400 lb	20" wide = 300 lb 22" wide = 350 lb 24" wide = 400 lb	20" wide = 300 lb 22" wide = 350 lb 24" wide = 400 lb
<b>Medicare codes</b>	E2613, E2614	E2615, E2616	E2620, E2621




**Comfort ACTA-BACK Deep back support**

May require a significant amount of postural support  
 Those who are more dependent with transfers  
 May need the maximum amount of trunk support



**Comfort ACTA-BACK Contour back support**

May require lower thoracic trunk support



**Not seeing your size?**

The Comfort Inception cushion (E2609) and Custom back (E2617) can be built up to 28" wide. The Inception cushion can also be made thicker as needed.

Contact your Permobil rep for more details.

©2022 Permobil. The following are trademarks and registered trademarks of Permobil: ROHO®, Hybrid Elite®, CONTOUR SELECT®, QUADTRO SELECT®, HIGH PROFILE®, MID PROFILE®, ENHANCER®, ACTA-BACK®, ACTA-EMBRACE®. Vicair® is a registered trademark of Vicair.

Permobil is the distributor of Vicair O2 cushions.

Permobil has a policy of continual product improvement and reserves the right to amend this document. The current version of this document is available at permobil.com.

Medicare Coding: Decisions about how to code a product in the U.S. are the responsibility of the provider/supplier. Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process.





**Terry Giles**

*Custom Fabrication Technician—Power Products  
U.S. Navy—Machinery Repairman First Class  
(Surface Warfare Specialist/Seabee Combat Warfare Specialist)*

## Proudly built in America

Behind every Permobil product is the hard work of hundreds of dedicated Americans. From the world's most advanced power and manual wheelchairs, to proven seating and positioning solutions, each Permobil product in this brochure is assembled with passion and care by employees

in one of our U.S. production locations. Permobil is proud to employ many of our nation's finest veterans and active military personnel who each day continue our mission that EVERYONE has the right to freedom, independence, and mobility.