

TiLite CR1

Effective Aug 2025 (New Items in Red)

MONO-TUBE RIGID CARBON FIBER FRAME

Page 1 of 5

PROVIDER INFORMATION		ATP INFORMATION	
Provider:	Account #:	ATP:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:		Name:	Ship Via:
TiLite Sales Rep:		Address:	
TiLite Information:	Promo Code:	City:	State: Zip:
Ship Date:	Order #:	Phone:	Fax:

Please Email to the Permobil Manual Customer Experience Team at orders@permobil.com or Fax to 509.586.2413

FRAME STYLE	
CFS1	☐ TiLite CR1 Carbon Fiber Frame (Minimal Adjustability. 231 lb. (105 kg) Weight Limit. For Rear Seat Width A and Seat Depth B of 15" or greater, use K0005. For Rear Seat Width A and Seat Depth B of less than 15", use E1235.....\$8,320

SMARTDRIVE POWER ASSIST SYSTEM	
CMAXSD1	☐ Permobil SmartDrive MX2+ (requires CCB2 SmartDrive Ready Carbon Fiber Camber Tube).....\$7,529 E0986
CMAXSD18	☐ SwitchControl Buttons.....NCO
CMAXSD8	☐ SwitchControl Button with Mono Jack Option.....NCO
CMAXSD19	☐ NEW Buddy Button.....NCO
CMAXSD4	☐ Carrying Bag.....\$84
CMAXSD10	☐ Universal Padded Handle.....\$22

*SmartDrive apps available on Watch OS (Apple Watch Series 5 or newer and Samsung Galaxy Watch 4 or newer).

COLOR ANODIZE PACKAGE	
(Color Anodize Package includes: Rear Strut Band-Clamps, Footrest Clamps, Caster Center Caps, Side Guard Caps. If no anodize color is selected, the default will be Charcoal Grey.)	
☐ CANO16 Charcoal Grey.....STD	
☐ CANO12 Gold ☐ CANO15 Green ☐ CANO13 Purple ☐ CANO14 Orange	
☐ CANO10 Blue ☐ CANO11 Red.....NCO	
TRIM RINGS (Includes: Frame and Backrest)	
☐ CANO16 Charcoal Grey ☐ CANO17 Polished Aluminum (with clear anodize).....NCO	

TiFit: TAILORED INDIVIDUALIZED FRAME	
SEAT WIDTH "A" Measure from outside of seat tube to the outside of the opposite seat tube.	
CSW1	☐ 13" (requires CBUP3 Padded Tension Adjustable by strap back upholstery).....STD
CSW1	☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18".....STD


SEAT DEPTH "B" Measure from front of back post (at the intersection with the seat tube) to the front edge of seat sling. Frame bend will start 1" in front of seat sling unless additional **Custom Frame Depth** is specified or the combination of specifications affect the rolling dynamics or functionality of the chair. Please note that 1" is the equivalent to .5" in front of the seat sling on other TiLite frame styles.

Seat Depth	Rear Seat-to-Floor Height "D"	Front Seat-to-Floor Height "C"	Frame Type
☐ 13" Deep	1.6" less than FSTFH	☐ 17.5" ☐ 17.75" (requires 3" casters and 24" rear wheels)	Frame 1
		☐ 18" – 20.75" (in .25" increments) Set at _____"	
	2.3" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 2
	2.9" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 3


2" additional custom frame depth available in .5" increments

Seat Depth	Rear Seat-to-Floor Height "D"	Front Seat-to-Floor Height "C"	Frame Type
☐ 14" Deep	1.7" less than FSTFH	☐ 17.5" ☐ 17.75" (requires 3" casters and 24" rear wheels)	Frame 1
		☐ 18" – 20.75" (in .25" increments) Set at _____"	
	2.4" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 2
	3.1" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 3

2" additional custom frame depth available in .5" increments



Frame 1
(seat angle "J" is 7 degrees from horizontal)



TiLite CR1

Effective Aug 2025 (New Items in Red)

MONO-TUBE RIGID CARBON FIBER FRAME

Page 2 of 5

Seat Depth	Rear Seat-to-Floor Height "D"	Front Seat-to-Floor Height "C"	Frame Type
☐ 15" Deep	1.8" less than FSTFH	☐ 17.5" ☐ 17.75" (requires 3" casters and 24" rear wheels)	Frame 1
		☐ 18" – 20.75" (in .25" increments) Set at _____"	
	2.6" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 2
	3.4" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 3

2" additional custom frame depth available in .5" increments

Seat Depth	Rear Seat-to-Floor Height "D"	Front Seat-to-Floor Height "C"	Frame Type
☐ 16" Deep	1.9" less than FSTFH	☐ 17.5" ☐ 17.75" (requires 3" casters and 24" rear wheels)	Frame 1
		☐ 18" – 20.75" (in .25" increments) Set at _____"	
	2.8" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 2
	3.6" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 3

2" additional custom frame depth available in .5" increments

Seat Depth	Rear Seat-to-Floor Height "D"	Front Seat-to-Floor Height "C"	Frame Type
☐ 17" Deep	2.1" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 1
	3.0" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 2
	3.8" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 3

2" additional custom frame depth available in .5" increments

Seat Depth	Rear Seat-to-Floor Height "D"	Front Seat-to-Floor Height "C"	Frame Type
☐ 18" Deep	2.2" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 1
	3.1" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 2
	4.0" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 3

1" additional custom frame depth available in .5" increments

Seat Depth	Rear Seat-to-Floor Height "D"	Front Seat-to-Floor Height "C"	Frame Type
☐ 19" Deep	2.3" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 1
	3.3" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 2
	4.3" less than FSTFH	☐ 18" – 20.75" (in .25" increments) Set at _____"	Frame 3

No custom frame depth available

Custom Frame Depth ☐ + _____" Frame Depth Longer Than Seat Depth. Measure from front of seat sling to beginning of frame bend and can be given as specific as .5" increments. Cannot be more than 2" additional to the 1" standard (total 3" max) unless noted under the seat depth box. _____NCO

CUSTOM CONFIGURATION OPTIONS

STANDARD FRONT END "I" Measure from inside of front frame tube to inside of opposite front frame tube.

CFE1 2.1" Taper (footrest width is 4.2" shorter than seat width "A") _____STD

FRONT ANGLE

CFA1 85° _____STD

SEAT TO FOOTREST "E" Measure from front edge of seat sling to top middle of footrest. Measurement must be 1" to 3" shorter than Front Seat Height "C" ☐ Set at: _____ (.25" increments. *Nylon Covered Hard Foam Pad (CMD #25) is available for individuals wanting a shorter footrest height.*)

CFTR13 ☐ Carbon Fiber Open Loop Footrest (with center support tube) _____STD

CFTR14 ☐ Carbon Fiber Footrest w/Flat Carbon Fiber Cover _____\$200

CFCT1 ☐ Flat Carbon Fiber Cover Grip Tape _____\$15


2" Sided Aluminum (Custom Sized) Footrest Cover (CMD #49) is available for individuals with smaller feet and wanting to keep their feet centered on the footrest.

FOOTREST TRIM RINGS


CFAR3 ☐ Stainless Steel _____STD

CFAR4 ☐ Titanium _____\$120

Frame 2
(seat angle "J" is 10 degrees from horizontal)



Frame 3
(seat angle "J" is 13 degrees from horizontal)



BACK SUPPORTS**ROHO® AGILITY® ACTIVE ALUMINUM Contour Back Support**Shell Height ☐ 8" ☐ 10"

BRS14	☐ AGILITY Active Aluminum Min Contour with Foam.....	\$738	E2620
BRS15	☐ AGILITY Active Aluminum Mid Contour with Foam.....	\$738	E2620

ROHO® AGILITY® ACTIVE CARBON Contour Back Support (2x2 Twill Weave Carbon Fiber)Shell Height ☐ 10" ☐ 13" ☐ 16"

	☐ AGILITY Carbon Min Contour Back with Foam w/ Air Insert		
BRS7	Foam w/ Air Insert.....	\$738	E2620
ACB1	AGILITY Active Carbon Upgrade.....	\$464	
	Fixed Hardware.....	NCO	
	☐ AGILITY Carbon Mid Contour Back with Foam w/ Air Insert		
BRS9	Foam w/ Air Insert.....	\$738	E2620
ACB2	AGILITY Active Carbon Upgrade.....	\$464	
	Fixed Hardware.....	NCO	

BACK UPHOLSTERY

CBUP3	☐ Padded Nylon, Tension Adjustable by Strap (required when seat width "A" is 13").....	STD
CBUP0	☐ Omit Back Upholstery.....	NCO

TiFit: TAILORED INDIVIDUALIZED FRAME**SEAT BACK TYPE FOLDING, FIXED-HEIGHT "F" CARBON SEAT BACK (CSB15)**

CSBA1	☐ Backrest tubes set at: _____ (10" – 17", can be given as specific as .25" increments.) available range with AGILITY Active solid seatback: (8" – 11", 14" – 17").....	STD
-------	---	-----

SEAT BACK ANGLE "G" Measure angle from front of back post to floor.

CSBA1	☐ Set at: _____ (available range 80° - 101°, in 1° increments).....	STD
-------	--	-----

SIDE GUARDS

CSG1	☐ Downward Pivoting/Fixed Lockout Carbon Fiber Side Guards (Profiled .125" above Rear Wheels or as close as possible depending upon chair configuration).....	STD
CSG2	☐ Downward Pivoting/Fixed Lockout Fendered Carbon Fiber Side Guards (Custom Cut to Profile of Rear Wheels).....	\$460

PUSH HANDLES

CSB10	☐ Integrated Carbon Fiber Push Handles (short).....	\$300
-------	--	-------

CENTER OF GRAVITY "H" Measure from front of back post at the seat tube to the center of the rear axle. Please note that due to its design and materials, the CR1 has a lighter than normal front end compared to other TiLite models which will affect the overall feel when selecting this measurement. There is a maximum 1" of adjustability in either direction from the specified location.

CCOG1	☐ 2.25" ☐ 2.5" ☐ 2.75" ☐ 3" ☐ 3.25" ☐ 3.5" ☐ 3.75" ☐ 4" ☐ 4.25" ☐ 4.5" ☐ 4.75" ☐ 5".....	STD
-------	--	-----

CAMBER TUBES

CCBR1	Aero Full Carbon Fiber Camber Tube with integrated Zero Play Carbon Fiber Camber Plugs (Not available on 13" wide chairs or certain configuration combinations) ☐ 0° ☐ 2° ☐ 4°.....	STD
CCBR2	SmartDrive Ready (round) Carbon Fiber Camber Tube with integrated Zero Play Carbon Fiber Camber Plugs ☐ 0° ☐ 2° ☐ 4°.....	NCO

FRONT WHEELS

CFW24	☐ 3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire (required when FSTFH "C" is 17.5") ☐ Silver Hub OR ☐ Black Hub.....	NCO
CFW25	☐ 4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub.....	STD
CFFW26	☐ 4"x1.25" (10cmx4cm) LiteSpeed-6 Carbon Fiber Wheel w/Soft Roll Tire.....	\$190

FRONT FORKS

CFK2	☐ TiLite Slipstream-CR1 Single-Sided Fork (Shorter length and width, stainless steel caster axle).....	STD
CTFS1	☐ Titanium Fork Stem.....	STD
CTFA1	☐ Titanium Caster Axle.....	\$130

PERFORMANCE REAR WHEELS

CRW14	Spinergy® Spox ☐ Black..... ☐ 24"(540) ☐ 25"(559)	STD
CRW16	Spinergy® Carbon Blade (Not available with Natural-Fit or Natural-Fit LT with Standard Grip or Super Grip.)..... ☐ 24"(540) ☐ 25"(559)	\$1,429

TiLite CR1

Effective Aug 2025 (New Items in Red)

MONO-TUBE RIGID CARBON FIBER FRAME

Page 4 of 5

AXLES

CAXL1	☐ Stainless Quick Release.....	STD
CAXL2	☐ Titanium Quick Release.....	\$277
CAXL3	☐ Quad Release.....	\$169

PNEUMATIC REAR TIRES

CRTR31	☐ MBL® SpeedLite.....	STD
CRTR10	☐ MBL® Gazelle (Kevlar bead, lightest weight option)	\$99
CRTR8	☐ Schwalbe® Marathon Plus Evolution.....	\$152

HANDRIMS

CRIM21	☐ Aluminum – Black Anodized.....	STD
CRIM23	☐ Spinerge® Aluminum – Anodized to match chair color anodize package.....	\$240
CRIM12	☐ Titanium – Satin Finish.....	\$509
CRIM20	☐ Q-Grip (Exclusive High-Friction, Anti-Slip Coating)	\$358

The Natural-Fit® LT (Reduced Weight, Smaller Profile)

CRIM17	☐ Standard Grip (Not available with Short Tabs or Spinerge Carbon Blade Rear Wheels.)	\$440
CRIM18	☐ Super Grip (Not available with Short Tabs or Spinerge Carbon Blade Rear Wheels.)	\$514
CRIM19	☐ Without Thumb Grip.....	\$423

Tab Length

CRIM7	☐ Long Tabs.....	STD
CRIM8	☐ Short Tabs (Trimmed, not available with The Natural-Fit® Handrims with Thumb Grip.)	\$30

WHEEL LOCKS

CWLK9	☐ MBL® Composite Scissor Lock.....	STD
CWLK10	☐ Omit Scissor Lock (Note that selecting this option prevents adding wheel locks at a later date.)	NCO

SEAT UPHOLSTERY

CSUP2	☐ Tension Adjustable by Straps (Black upholstery only.)	STD
-------	--	-----

SEAT CUSHION

CSSH2	☐ Comfort Curve Cushion (Stretch-Air cover)	\$122	E2601
CSSH23	☐ Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion	\$428	E2607
CSSH11	☐ ROHO Single Compartment (Fits seat widths 14" to 19") Profile: ☐ Low ☐ Mid ☐ High.....	\$534	E2622
CSSH14	☐ ROHO Quadtro Select (Fits seat widths 14" to 19") Profile: ☐ Low ☐ Mid ☐ High.....	\$534	E2624

POSITIONING BELT

CBLT1	Velcro Adjustable, Black Nylon Webbing.....	\$61	E0978
-------	---	------	-------

IMPACT GUARDS

CMPCT2	☐ Black Neoprene Impact Guards with Ultrasuede Stripe.....	\$120
--------	---	-------

CALF STRAP

CCLF3	Bodypoint® Padded Velcro Adjustable.....	\$109	K0038
	☐ Medium (Recommended for 13" – 16" Chair Widths.)		
	☐ Large (Recommended for 17" – 18" Chair Widths.)		

BACKPACK

CBCK3	☐ Permobil Medical Necessities Backpack.....	\$103
-------	---	-------

SEAT POUCH

CPCH3	☐ Permobil Medical Necessities Seat Pouch.....	\$60
CPCH4	☐ Bodypoint® Mobility Bag.....	\$109

ANTI-TIPS

CTIP1	☐ Rear Aluminum – Black.....	\$174	E0971
-------	---	-------	-------

Military Patch

CAFP1	4" Logo Patch ☐ Army ☐ Navy ☐ Air Force ☐ Marines.....	\$115
-------	--	-------

GENERAL

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at permobil.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of

TiLite CR1

Effective Aug 2025 (New Items in Red)

MONO-TUBE RIGID CARBON FIBER FRAME

Page 5 of 5

Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. ROHO®, AGILITY®, Low Profile® and Quadtro Select® are registered trademarks of ROHO Inc. Spinergy® is a registered trademark of Spinergy Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® is a registered trademark of Three Rivers Holdings, LLC. The CURVE® is a registered trademark of THE COMFORT COMPANY® and are licensed to Permobil .

NOTES

QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR