



PDG Fuze T50 - AADL W653



Dealer Name*:

Account #*:

City, Province, Post Code:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference*:

Client Weight:

*Required field

Special considerations that need to be addressed here (e.g., diagnosis):

October 1, 2026. Prices are in Canadian Dollars (CAD)

Wheelchair & Frame

☒ **Fuze T50 50° Manual Tilt Wheelchair** Also available: the narrow Fuze T50n with an overall chair width of only seat width + 4". See separate order form. W653 \$ 2,955.00

Frame Type

- ☐ **Fixed Width Frame, 300lbs Weight Capacity** Select a fixed seat width below. No charge
- ☐ **Fixed Width Frame, 400lbs Weight Capacity** Not compatible with transit tie down system. Select a fixed seat width below. No charge
- ☐ **Adjustable Width Frame, 250lbs Weight Capacity** Select an adjustable seat width below. \$ 200.00

Frame Options

- ☐ **Take-Apart Frame** The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. This feature reduces the weight of the heaviest piece of the wheelchair to 29lbs or less, perfect for transporting the chair in the trunk of a car or SUV. Not compatible with power tilt, drum brakes, IV pole, O2 holder, ventilator tray, ventilator battery tray, dynamic rocker back or 13" seat to floor height. \$ 650.00
- ☐ **Transit Tie Down System (TTDS)** Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details. Tested to maximum weight capacity of 300lbs. No charge

Seat Measurements

Seat Width | Overall chair width is +7.5" with 12/16" rear wheels and +8.5" with 20/22/24" rear wheels.

- Fixed:**
- ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" No charge
- ☐ 20" ☐ 21" ☐ 22" W111 \$ 200.00
- ☐ 23" ☐ 24" W111 \$ 500.00
- Adjustable 12"-16":** ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" No charge
- Adjustable 15"-20":** ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" No charge
- ☐ 20" W111 \$ 200.00

Seat Depth

- Select Depth:**
- ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" No charge
- ☐ 20" \$ 100.00
- ☐ 21" ☐ 22" \$ 700.00

21"-22" seat depths are not adjustable and require straight stroller handle back. Not compatible with reclining or fold-down back.

Frame Finishes

Select Color: ☐ High Gloss Black ☐ Hyper Blue

Armrests

Armrest Type

- ☐ **No Armrests** No charge
- ☐ **Single Post, Height Adjustable Armrests** No charge
- ☐ 10-14" Range ☐ 12"-16" Range

☐ **Fixed Height Armrests** 12" height, full length armpads only. No charge

Armrest Position: ☐ Mid (default) ☐ Front

If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.

☐ **Locking Flip-Up Armrests** Not recommended for users with excessive strength or tone. Not compatible with fold-down back. Only with regular armpads. \$ 545.00

☐ Set at 10" ☐ Set at 11" ☐ Set at 12" ☐ Set at 13" ☐ Set at 14"

Armpads

☐ **Regular Armpads** No charge

☐ Desk Length 2" x 10" ☐ Full Length 2" x 14"

☐ **Waterfall Armpads** \$ 275.00

☐ Desk Length 2" x 10" ☐ Full Length 2" x 14"

☐ **Comfort Armpads with Interlocking Hardware** May interfere with wheels in some tilt configurations. \$ 496.00

☐ Desk Length 3.75" x 10" ☐ Full Length 3.75" x 14"

☐ **Gel Armpads** May interfere with wheels in some tilt configurations. \$ 305.00

☐ Narrow 2" x 12" ☐ Wide 3.5" x 12"

☐ **Comfort Standard Arm Trough with Interlocking Hardware, Right** Not compatible with lap tray or locking flip-up armrests. May interfere with wheels in some tilt configurations. \$ 562.00

☐ 6" x 17" ☐ 6" x 22"

☐ **Comfort Standard Arm Trough with Interlocking Hardware, Left** Not compatible with lap tray or locking flip-up armrests. May interfere with wheels in some tilt configurations. \$ 562.00

☐ 6" x 17" ☐ 6" x 22"

Armrest Wideners | Moves armrests outwards by up to 1" each side. Seat size does not change.

☐ **Armrest Wideners, 1" Total Additional Width** \$ 245.00

☐ **Armrest Wideners, 2" Total Additional Width** \$ 395.00

Seat Height

Seat to floor height is adjustable within the range of the chosen axle plate (Low/Standard/High). Chairs with swing-away front rigging are set with 3° of pre-tilt to avoid interference. This increases the actual seat to floor height at the front of the seat by up to 1" compared to the selected value. Easy-Out and standard plug-in front rigging are not affected.

		5" Casters					6" Casters & 6" Forks					6" Casters & 8" Forks					8" Casters & 8" Forks					6" Frog Legs Casters & Suspension Forks		
RW	24"			☐ Low	☐ Low	☐ Low		☐ Low	☐ Low	☐ Low	☐ Std	☐ Low	☐ Low	☐ Low	☐ Std	☐ Std	☐ Low	☐ Low	☐ Std	☐ Std	☐ Std	☐ Low	☐ Low	☐ Std
	22"		☐ Low	☐ Low	☐ Low	☐ Std	☐ Low	☐ Low	☐ Low	☐ Std	☐ Std	☐ Low	☐ Low	☐ Std	☐ Std	☐ Std	☐ Low	☐ Std	☐ Std	☐ Std		☐ Low	☐ Std	☐ Std
	20"	☐ Low	☐ Low	☐ Low	☐ Std	☐ Std	☐ Low	☐ Low	☐ Std	☐ Std	☐ Std	☐ Low	☐ Std	☐ Std	☐ Std		☐ Std	☐ Std	☐ Std		☐ High	☐ Std	☐ Std	☐ Std
	16"		☐ Std	☐ Std	☐ Std		☐ Std	☐ Std	☐ Std		☐ High	☐ Std	☐ Std		☐ High	☐ High	☐ Std		☐ High	☐ High	☐ High	☐ Std		☐ High
	12"	☐ Std	☐ Std		☐ High	☐ High	☐ Std		☐ High	☐ High	☐ High		☐ High	☐ High	☐ High		☐ High	☐ High	☐ High	☐ High	☐ High	☐ High	☐ High	☐ High
		13"	14"	15"	16"	17"	14"	15"	16"	17"	18"	15"	16"	17"	18"	19"	16"	17"	18"	19"	20"	16"	17"	18"
Seat to Floor Height																								

Seat Pan

Seat Pan

☐ **No Seat Pan** No charge

☐ **Solid Aluminium Seat Pan** No charge

Seat Cushions

Cushions over 20" wide subject to increased pricing.

☐ **Comfort M2 Cushion Zero Elevation** \$ 580.00

☐ Comfort-Tek Cover ☐ Stretch Air Cover

☐ **Comfort M2 Cushion Anti-Thrust** \$ 580.00

☐ Comfort-Tek Cover ☐ Stretch Air Cover

☐ **Comfort Embrace Cushion Zero Elevation** \$ 508.00

☐ Comfort-Tek Cover ☐ Stretch Air Cover

☐ **Comfort Embrace Cushion Anti-Thrust** \$ 508.00

☐ Comfort-Tek Cover ☐ Stretch Air Cover

☐ Comfort Saddle 7 Series Foam Cushion Zero Elevation	<i>Solace cover.</i>	\$ 818.00
☐ Comfort Saddle 7 Series Gel Cushion Zero Elevation	<i>Solace cover.</i>	\$ 1,133.00
☐ Comfort Support Pro 7 Series Foam Cushion Zero Elevation	<i>Solace cover.</i>	\$ 445.00
☐ Comfort Support Pro 7 Series Gel Cushion Zero Elevation	<i>Solace cover.</i>	\$ 693.00
☐ ROHO Quadtro Select High Profile Cushion		\$ 777.00
☐ ROHO Sensor Ready High Profile Cushion		\$ 637.00

Front Casters

Urethane Casters

☐ 5" x 1.25" Caster	No charge
☐ 6" x 1.25" Caster	No charge
☐ 8" x 1.25" Caster	No charge
☐ 6" x 1.50" Caster	\$ 230.00
☐ 8" x 2.00" Caster	\$ 305.00

Pneumatic Casters

☐ 6" x 1.25" Pneumatic Caster	\$ 265.00
☐ 8" x 1.25" Pneumatic Caster	\$ 280.00
☐ 8" x 2.00" Pneumatic Caster	\$ 325.00

Soft-Roll Casters

☐ 6" x 1.5" Soft-Roll Caster	<i>Urethane tire and silver anodized aluminum hub.</i>	\$ 515.00
---	--	-----------

Frog Legs Suspension Forks

☐ FrogLegs Suspension Forks	<i>Not compatible with 12" rear wheels. Requires 6" soft-roll caster.</i>	\$ 995.00
--	---	-----------

Rear Wheels

☐ Composite Rear Wheel	No charge
☐ 12" ☐ 16" ☐ 20" ☐ 22" ☐ 24"	
☐ Spoke Rear Wheel	\$ 215.00
☐ 20" ☐ 22" ☐ 24"	
☐ Heavy Duty Spoke Wheel	\$ 625.00
☐ 24"	

Rear Wheel Position: ☐ Rear (More Stable) ☐ Mid (Default) ☐ Front (Less Stable)

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

Axles

☐ Quick Release Rear Axles	No charge
☐ Fixed Rear Axles	No charge

Rear Tires

☐ Urethane Tires	No charge
☐ Pneumatic Tires <i>May result in a finished seat height 1/2" greater than stated.</i>	\$ 190.00
☐ Pneumatic Tires with Flat-Free Inserts <i>May result in a finished seat height 1/2" greater than stated. Not available with spoke wheels.</i>	\$ 245.00

Handrims | *Handrims are only available with 20", 22" and 24" rear wheels.*

☐ No Handrims	No charge
☐ Standard Aluminium Handrims	No charge
☐ Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	No charge
☐ Projection Handrims <i>Composite wheels only. 10 projections, 30°. 1-1.8" length.</i>	\$ 845.00

Wheel Locks

☐ User Wheel Locks	No charge
☐ Push to Lock ☐ Pull to Lock	
☐ Attendant Wheel Locks <i>Not compatible with drum brakes or extensions.</i>	\$ 230.00
☐ Push to Lock ☐ Pull to Lock	

- ☐ **User and Attendant Wheel Locks** *Not compatible with drum brakes.* \$ 290.00
☐ Push to Lock ☐ Pull to Lock
- ☐ **Foot Operated Wheel Lock** *Rear-mounted unilateral single-activation lock. Not available with power tilt in combination with 12" or 16" rear wheels. Not compatible with O2 holder, IV pole or ventilator tray.* \$ 995.00
- ☐ **Drum Brakes** *Available with 12" composite or 20" spoke rear wheels only. Urethane tires only. Not with plastic coated handrims. 12" wheels add 1" to overall chair width. 12" compatible with 13", 14", 17.5", 18.5", 19.5" STF heights only. 20" wheels compatible with 16", 17", 18" STF heights only.* \$ 795.00
- ☐ **User Activated Unilateral Wheel Lock** *Pull to lock, Available up to 20" seat widths only.* \$ 1,295.00
☐ Left Activation ☐ Right Activation

Wheel Lock Options

- ☐ **Wheel Lock Extensions** *Not compatible with attendant wheel lock.* No charge

Anti-Tippers

- ☒ **Rear Anti-Tippers** *Required safety option.* No charge

Backrest

- Back Post Style:** ☐ Stroller Handle Back ☐ Cane Handle Back
Cane handle back not compatible with power tilt with manual override.
- Back Post Type:** ☐ Straight ☐ 8° Bend
- Back Post Height:** ☐ Short (17") ☐ Tall (23")
Measured from seat pan to handle.
- Folding Style:** ☐ Fold-Down ☐ Non-Folding
Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.
- Seat to Back Angle:** ☐ Set at 90° ☐ Set at 98° ☐ Set at 106° ☐ Set at 114° ☐ Set at 122°

Backrest Options

- ☐ **Angle/Height Adjustable Stroller Handle:** *Not compatible with cane handle back.* No charge
☐ Button on Outside ☐ Button on Inside
- ☐ **Endomorph Back** *Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not compatible with flip-up armrests, dynamic rocker back or cane handle back.* \$ 695.00
☐ 1" Offset Per Side ☐ 1.5" Offset Per Side ☐ 2" Offset Per Side
- ☐ **Seating Dynamics Dynamic Rocker Back Interface** *Not compatible with Take-Apart frame, fold-down back, endomorph back, 30° recline or cane handle back. Not available with 21" and 22" seat depths. Only available with stroller handle - 8° bend. Finished back post height is 7" taller than selected back post height. Maximum compatible seat depth is 20". 90° and 98° back angle only.* \$ 2,995.00

Back Supports

Comfort Back Supports | Backs over 20" wide or tall subject to increased pricing.

- ☐ **Comfort Unite 1 Back** *Includes 2-point quick release Compass hardware. Select shell height and cover.* \$ 961.00
☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 1 Back - 5" Wing** *Includes 2-point quick release Compass hardware. Select shell height and cover.* \$ 1,245.00
☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 1 Back w/ROHO Insert - 5" Wing** *Includes 2-point quick release Compass hardware. Select shell height and cover.* \$ 1,499.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 1 Back - 7" Wing** *Includes 2-point quick release Compass hardware. Select shell height and cover.* \$ 1,245.00
☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 1 Back w/ROHO Insert - 7" Wing** *Includes 2-point quick release Compass hardware. Select shell height and cover.* \$ 1,499.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 2 Back** *Includes 2-point quick release Compass hardware. Select shell height and cover.* \$ 1,106.00
☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 2 Back w/ROHO Insert** *Includes 2-point quick release Compass hardware. Select shell height and cover.* \$ 1,254.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover

- ☐ **Comfort Unite 3 Back** Includes 2-point quick release Compass hardware. Select shell height and cover. \$ 1,287.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 3 Back w/ROHO Insert** Includes 2-point quick release Compass hardware. Select shell height and cover. \$ 1,499.00
☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Embrace Back** Includes 2-point quick release Compass hardware. Select shell height and cover. \$ 1,254.00
☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Relief Back** Includes Comfort-Tek cover and 2-point quick release Compass hardware. Select shell height. \$ 1,329.00
☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"

Positioning Accessories

Lateral Supports

- ☐ **BodiLink Lateral Supports** Select size, cover and side.
☐ 3" x 5.5" ☐ 4" x 6.5" ☐ 6" x 6.5"
☐ Comfort-Tek Cover ☐ Stretch Air Cover
☐ Left ☐ Right \$ 220.00
☐ Pair \$ 440.00
- ☐ **BodiLink Lateral Support Mounting Hardware**
☐ Individual \$ 368.00
☐ Pair \$ 736.00

Headrests

- ☐ **BodiLink Headrest Pad** Select size and cover. \$ 260.00
☐ 10" Wide ☐ 12" Wide ☐ 14" Wide ☐ 18" Wide
☐ Comfort-Tek Cover ☐ Stretch Air Cover
- ☐ **BodiLink Headrest Hardware, Removable**
☐ Mid 6 ☐ Mid 8 \$ 373.00
☐ Lat 6 ☐ Lat 8 \$ 489.00
☐ Max 8 \$ 545.00
- ☐ **BodiLink Headrest Hardware, Flip-Down**
☐ Mid 6 ☐ Mid 8 \$ 419.00
☐ Lat 6 ☐ Lat 8 \$ 530.00
☐ Max 8 \$ 591.00

Recline

- ☐ **Reclining Backrest (Manual)** 30° maximum recline. Not compatible with cane handle back or fold-down back. Not available with 21" and 22" seat depths. Finished back post height is 3" taller than selected back post height. \$ 1,445.00

Tilt Activation

Attendant Tilt Only

- ☐ **Attendant Tilt Only, Right Back Post** No charge
- ☐ **Attendant Tilt Only, Left Back Post** No charge

Power Tilt | 40° tilt range. Includes 2 batteries for off-chair charging. Requires stroller handle back.

- ☐ **Standard Power Tilt, 250lb Weight Capacity** \$ 1,795.00
- ☐ **Standard Power Tilt, 350lb Weight Capacity** Requires fixed frame. \$ 1,995.00
- ☐ **Power Tilt with Manual Tilt Override, 250lb Weight Capacity** Allows for instand manual release of tilt in case of power tilt failure. \$ 2,095.00
- ☐ **Power Tilt with Manual Tilt Override, 350lb Weight Capacity** Allows for instand manual release of tilt in case of power tilt failure. Requires fixed frame. \$ 2,295.00
- ☐ **Toggle Switch Control** No charge
☐ Installed on Right ☐ Installed on Left
- ☐ **Button Switch Control** No charge
☐ Installed in Right ☐ Installed in Left

- ☐ **Stealth Egg Switch Control** *Not installed.* No charge
- ☐ **Swing-Away Mounting Hardware** *Only for Stealth egg switch. Not compatible with all headrest makes and models.* \$ 1,040.00

Tilt Limit | *Optional. If no tilt limit is selected the full tilt range will be available.*

Tilt Limit: ☐ Set at 20° ☐ Set at 30° ☐ Set at 40°

Front Rigging

- ☐ **No Front Rigging** No charge

Easy-Out Front Rigging | *Effortless insertion and removal. No hanger brackets remain on the frame.*

- ☐ **Fixed 70°** \$ 95.00
- ☐ **Universal Adjustable Angle (0°-100°)** *Select angle. Angle adjustment by loosening one bolt, in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings.* \$ 545.00
- ☐ 0° ☐ 10° ☐ 20° ☐ 30° ☐ 40° ☐ 50° ☐ 60° ☐ 70° ☐ 80° ☐ 90° ☐ 100°
- ☐ **High Mount Option for Easy-Out Adjustable Angle** *For shorter lower leg lengths of 8"-12.5".* \$ 545.00

Swing Away Front Rigging | *Legacy design. Not compatible with Easy-Out.*

- ☐ **Swing Away 60°** No charge
- ☐ **Swing Away 70°** No charge
- ☐ **Swing Away 90°** *Includes 11"-16" extension tubes. Not compatible with 8" casters.* \$ 75.00
- ☐ **Swing Away Elevating Legrests (ELR)** *Includes calf pads.* W106 \$ 375.00
- ☐ **1" Offset Hangers for Swing Away Front Rigging** \$ 95.00

Plug In Front Rigging | *Legacy design. Not compatible with Easy-Out.*

- ☐ **Plug In 70°** No charge

Footplates

- ☐ **Composite Footplates** No charge
- ☐ **Aluminium Angle Adjustable Footplates** W107 \$ 165.00
- ☐ **One Piece Fixed Footplate** *Available on plug in front rigging only. Requires selection of composite or aluminum angle adjustable footplates.* \$ 625.00
- ☐ **One Piece Flip-Up Footplate** *Requires selection of composite or aluminum angle adjustable footplates.* \$ 695.00
- ☐ **Interlocking Angle Adjustable Footplates** *Available for seat widths up to 21"* \$ 395.00

Extension Tubes | *Consider the height of your seat cushion.*

- ☐ **Standard Extension Tubes, Short** No charge
- ☐ **Standard Extension Tubes, Medium** *Only with 90° swing-away front rigging.* No charge
- ☐ **Standard Extension Tubes, Long** \$ 75.00

Footplate Height: ☐ Specify:

Ext Tube	Short	8" - 13"	11" - 16"	12" - 16"	4" - 9"	13" - 16"
Medium		N/A	N/A	N/A	10" - 12"	N/A
Long		N/A	14" - 22"	14" - 19"	12" - 16"	14" - 17"

High Mount Easy-Out

Standard Easy-Out

60°, 70°, ELR Swing-Away

90° Swing-Away

Standard Plug-In

Legrest

Front Rigging Accessories

- ☐ **Calf Pads, Easy-Out Front Rigging** *Not compatible with gel padding. Includes reinforced premium calf pad mount.* \$ 400.00
- ☐ **Calf Pads, Swing Away Front Rigging** *Not compatible with gel padding.* \$ 305.00
- ☐ **Gel Padding** *Not compatible with calf pads or Easy-Out front rigging.* \$ 250.00
- ☐ **Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging** \$ 400.00
- ☐ **Gel Knee Pads, Swing Away ELR Front Rigging** \$ 250.00
- ☐ **Length Adjustable Heel Loops** No charge
- ☐ **Standard Calf Strap** \$ 50.00
- ☐ **Bodypoint Padded Calf Strap** *Sizes greater than 22" are not padded.* \$ 117.00
- ☐ **Comfort Foot Box, Double** *Without separator. Select seat width*
- ☐ 16 ☐ 18" \$ 454.00
- ☐ 20" ☐ 22" \$ 644.00
- ☐ 24" \$ 776.00

☐ Comfort Foot Box, Single	
☐ Left ☐ Right	\$ 285.00
☐ Pair	\$ 570.00

Positioning Belts

☐ Bodypoint Non-Padded Center Pull Push Button Belt, Large <i>Not compatible with gel padding. Includes reinforced premium calf pad mount.</i>	\$ 125.00
☐ Bodypoint 2 Point Padded Center Pull Push Button Belt <i>Not compatible with gel padding.</i>	\$ 215.00
☐ Bodypoint 4 Point Padded Center Pull Push Button Belt <i>Not compatible with calf pads or Easy-Out front rigging.</i>	\$ 175.00
☐ Bodypoint Non-Padded Chest Belt 46"	\$ 165.00
☐ Q'Straint Lap Belt <i>Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.</i>	\$ 470.00

Miscellaneous Options

☐ One-Arm Drive <i>Available with 22" and 24" composite rear wheels. Minimum 17" and 18" seat to floor height respectively. Large hand setup includes 24" rear wheels with 22" driving handrim. Not compatible with drum brakes, quick release axles or pneumatic tires.</i>	\$ 2,800.00
Installed On: ☐ Right ☐ Left	
Handrim Setup: ☐ Standard ☐ Large Hand	
<i>Large hand setup only available with 24" rear wheels.</i>	
☐ Spoke Guards <i>Available for 20"-24" composite rear wheels. Made from durable and lightweight aluminium.</i>	\$ 395.00
☐ Black ☐ Frame color	
☐ Hip Guide, Bracket Only <i>Installed with receivers. Accessory mount.</i>	\$ 200.00
☐ Gel Hip Guides <i>Reduces available seat width by 2"</i>	\$ 845.00
☐ Hip Belt Mount <i>Provides improved mounting options for Bodypoint hip belt or other third-party belts.</i>	\$ 75.00
☐ IV Pole	\$ 545.00
☐ Installed on Right ☐ Installed on Left	
☐ O2 Holder <i>May add to overall width of chair.</i>	\$ 695.00
☐ Installed on Right ☐ Installed on Left	
☐ Ventilator Tray <i>Reduced ground clearance. Chair may be unstable when unoccupied. Minimum seat width is 16". Not compatible with O2 holder, power tilt or foot operated wheel locks.</i>	\$ 545.00
☐ Clear Lap Tray (Slide On) <i>Available for seat widths 17"-22" only. Not compatible with contour arm trough.</i>	\$ 835.00
☐ Medical Necessities Backpack <i>Black with Permobil logo.</i>	\$ 245.00
☐ Adhesive Back Hook Strip <i>For cushion mounting, 2"x30".</i>	\$ 15.00

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes