

UK (English)
 Date
 User
 End user weight (kg)
 Language user manual

Company stamp
 Responsible

Folding frame wheelchair

£ 2750,00

Retail price - VAT not included

Wheelchair measures - The tolerance on all the dimensions indicated and/or required is ± 10 mm

Seat width (LS) 36 39 42 45 48
 The seat width is measured from the outer sides of the frame
 For seat enlargement see section enlargement of mudguard and armrest with integrated mudguard


Seat depth (PS) 37,5 40 42,5 45 47,5
 The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one. - calculated with angle backrest 90° to the seat


Footplate distance (DP) (cm) ±0.5cm - adjustment steps 1.5cm
(Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)

Setting Standard Prudential Active
 camber 2° 0° 4°
 Backrest Height (TS) Std.Value 33,5(min 24.5-max 47) adjustment steps 1.5cm
 Backrest angle (to the seat) 90° -6° 6°

Front height (STD 49) - Rear height (STD 43)

Front frame

Aluminium standard frame Angle 100°

 Short frame - Ideal with Seat depth (PS) 37,5 - 40
 Long frame - Ideal with Seat depth (PS) 42,5 - 45 - 47,5

Aluminium frame Angle 90° - (max 5" caster)

 Short frame - Ideal with Seat depth (PS) 37,5 - 40
 Long frame - Ideal with Seat depth (PS) 42,5 - 45 - 47,5

Frame colour

<input checked="" type="checkbox"/> Opaque black	<input type="checkbox"/> Tekno anthracite	<input type="checkbox"/> Matt grey
<input type="checkbox"/> POLAR white	<input type="checkbox"/> Matt Midnight blue	<input type="checkbox"/> Matt red
<input type="checkbox"/> Matt orange		

Backrest

Adjustable velcro straps and padded transpiring upholstery

Handle tube

Tubes with standard handle
 Folding handles

Not available with: Push handles for companion, Brake handles for companion (Drum Brake)



02958 **270,00**

Casters

<input checked="" type="checkbox"/> 4" Soft		
<input type="checkbox"/> 3" Sport	<i>Lightest configuration</i>	
<input type="checkbox"/> 5" Soft		
<input type="checkbox"/> 6" Soft	(Available only with 100° angle)	
<input type="checkbox"/> 6" Pneumatic	(Available only with 100° angle)	
<input type="checkbox"/> 3" Active Soft (wide)	(Aluminium alloy rim)	06824 95,00
<input type="checkbox"/> 4" Active Sport (narrow)	(Aluminium alloy rim)	01398 95,00
<input type="checkbox"/> 4" Active Soft (wide)	(Aluminium alloy rim)	03943 95,00
<input type="checkbox"/> 5" Active Soft (wide)	(Aluminium alloy rim)	07518 95,00
<input type="checkbox"/> 5" Pneumatic		10844 68,00

Support front fork

Standard support front fork
 Front fork support V Design 2.0
 Exclude the installation of QUICK RELEASE FRONT FORKS



06850 **158,00**

customer_support_uk@permobil.com - https://www.permobil.com/en-gb

Date Order User Address

Colour set: Includes rear wheel hub and rim, rear wheel plates and front forks
Only standard Progeo 22"/24" wheels have coloured hubs and rims

☆ Black	22" 24" 25" 26"		
○ Bright silver	22" 24"	00833	135,00
○ Bright lime	22" 24"	00833	135,00
○ Bright blue	22" 24"	00833	135,00
○ Bright red	22" 24"	00833	135,00
○ Bright fuchsia	22" 24"	00833	135,00
○ POLAR white	22" 24"	00833	135,00
○ Tekno anthracite	22" 24"	00833	135,00
○ Matt Midnight blue	22" 24"	00833	135,00
○ Matt grey	22" 24"	00833	135,00
○ Matt red	22" 24"	00833	135,00

Rear wheels

☆ 24" Ultra-light Progeo - Standard!			
□ 22" Ultra-light Progeo - Standard!			
□ 25" Ultra-light Progeo - Standard!	(colour: black)		
□ 26" Progeo - Standard!	(colour: black)		
○ 24" ULTRA-LIGHT PROGEO KEY SPORT	<i>Lightest configuration</i>	04383	695,00
○ 24" Ultra-light "Spinergy LX"		01685	695,00
○ 24" Ultra-light "Spinergy LXL"	□ Black hub and white spokes □ Hub and spokes black	07846	695,00
○ 24" Ultra-light "Spinergy Flex rim"	Integrated aluminium/rubber handrim	02270	1.150,00
○ 25" Ultra-light "Spinergy LX"		02101	695,00
○ 25" Ultra-light "Spinergy LXL"	□ Black hub and white spokes	07847	695,00
○ 26" Ultra-light "Spinergy LX"		02102	695,00
○ 24" Off-road (aluminium handrim black - Wheel distance 4cm)		01392	270,00
○ 24" with brake handles for companion (drum brake) (aluminium handrim - Wheel distance 3cm)		03384	386,00
Available only with: camber 0° and Setting Standard			
○ One arm drive - (24" rear wheels Aluminium handrim High pressure Tyre) - camber 0°	○ Right ○ Left	02103	430,00

Rear wheel tyre

☆ High pressure tyre	(22"-24"-25"-26")	<i>Lightest configuration</i>	
□ High pressure profile tyre	(24")		
□ Marathon Plus	(24"-25")	with anti-puncture protection	01669 95,00
□ Anti-puncture - Colour Black	(24")		11083 95,00
Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"			
□ Anti-puncture - Colour Grey	(22")		04181 105,00

Handrim	Wheel distance ☆ 3 cm	□ 4 cm		
☆ Aluminium	(22"-24"-25"-26")	<i>Lightest configuration</i>		
○ Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170	160,00
○ Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684	160,00
○ Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016	290,00
○ Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017	290,00
○ TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023	290,00
Oval countour with a rubber insert				
○ NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11025	290,00
Ergonomic black profile				
○ NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11026	290,00
Ergonomic profile with rubber coating (ideal for quadriplegics)				
○ Silicon cover for standard aluminium handrim		Available for measures (22"- 24" - 26")	04050	172,00


Wheel locks

☆ Bent lever push to lock			
□ Straight lever push to lock			
□ Straight lever pull to lock			
○ Sport brake		06428	135,00
○ Brake with extension lever push to lock		01949	68,00
○ Brake with extension lever pull to lock		03372	68,00


Date	Order	User	Address
------	-------	------	---------

Footplate - See reference tables N°3


- Automatic Closure footplate rake and depth adjustable
 Minimum footplate distance 35cm. Internal 2/3 Internal External
 Minimum footplate distance 30cm. Available only with Multi-function and rotative system - 03945 203,00


11526



- One piece footplate - Flip-up! rake and depth adjustable
 Minimum footplate distance 35cm. Internal 2/3 Internal External
 Minimum footplate distance 30cm. Available only with Multi-function and rotative system - 03945 203,00


11527


- One piece aluminium footplate - Flip-up! rake and depth adjustable
 Minimum footplate distance 35cm. Internal 2/3 Internal 2/3 External External


11528


- Plastic footrest rake adjustable
 Minimum footplate distance 35cm. 2/3 Internal 2/3 External
 Need a stretch calf band.
NOT AVAILABLE WITH FRONT FORK SUPPORT V DESIGN 2.0



06861 93,00



- One piece Carbon fibre footplate rake adjustable, swing away, flip-up.
 Minimum footplate distance 35cm. Internal External


03291 294,00

- Footplate COMPACT, automatic Closure, rake and depth adjustable
 Minimum footplate distance 35cm. Internal External



09957 270,00

Mudguard

- Mudguard**
- ABS mudguard
 - ABS Straight sideguard
 - Carbon fibre mudguard
 - Carbon fibre Straight sideguard
- Mudguard enlargement**
- + 0 cm
 - + 1 cm + 2 cm Not available with Rear wheels with brake handles for companion (Drum Brake)
- 
03946 337,00
08918 337,00
- Removable version**
The distance between the two side guards is 1 cm wider than the seat width
- Removable ABS mudguard
 - Removable ABS fibre straight sideguard
 - Removable carbon fibre mudguard
 - Removable Carbon fibre straight sideguard
- 
08921 135,00
08921 135,00
08919 466,00
08920 466,00


Armrests

- Tubular armrests - Steel support with safety locking system
 (Swing away, removable and height adjustable armrest)
 Not available with anti-tippers if the rear seat height is less than 40cm


00694 270,00

- Removable armrest with integrated mudguard**

 Desk armrests PROGEO (Removable and height adjustable)



09956 240,00


- Enlargement of Desk armrests PROGEO**

 + 0 cm
 + 2 cm Available only with Ultra-light Progeo standard wheels and camber 0° or 2°

12648

- Removable armrest with ABS straight sideguard
The distance between the two side guards is 1 cm wider than the seat width






08881 240,00
- Removable armrest with carbon-fiber straight sideguard
The distance between the two side guards is 1 cm wider than the seat width




08882 486,00

Accessories

- Pad for front frame protectors **07058 56,00**
- Backrest-seat connection cloth **10241 25,00**
- Removable storage bag (with quick coupling) **09811 135,00**

Date	Order	User	Address
------------	-------------	------------	---------------

Optional 1		<u>Only one optional devices for side can be chosen</u>		
<input type="radio"/> Transit wheels	Exclude the installation of : Anti-tippers - Tipping aid - Crutches support		00836	158,00
<input type="radio"/> Anti-tipper with rapid hookup	(inside rotation to reduce the encumbrance) <input type="radio"/> Right <input type="radio"/> Left		00023 00024	145,00 145,00
<input type="radio"/> Tipping aid	<input type="radio"/> Right <input type="radio"/> Left		00025 01046	79,00 79,00
<input type="radio"/> Crutches support	<input type="radio"/> Right <input type="radio"/> Left		00838 00837	79,00 79,00

Optional 2			
<input type="radio"/> Tetraclip	(Quick release axle for quadriplegic)		03947 189,00
<input type="radio"/> Quick release front forks	To be matched only with: Standard support front fork		00030 158,00
<input type="radio"/> Push handles for companion	(Removable and height adjustable)		01048 203,00
<input type="radio"/> Safety pelvic belt			03972 68,00
<input type="radio"/> Harness			04103 145,00
<input type="radio"/> Lumbar support	(to insert in the backrest upholstery)		00013 43,00
<input type="radio"/> Swing away lateral support Not available with Posture backrest PHYSIO	<input type="radio"/> Right <input type="radio"/> Left		02097 270,00 02098 270,00
<input type="radio"/> Spoke guards Black	(22"-24")		01043 158,00
<input type="radio"/> Spoke guards Transparent	(22"-24")		00835 158,00
<input type="radio"/> Spoke guards Progeo design	(22"-24")		00019 189,00








Additional Kit			
<input type="radio"/> 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)		03378 534,00
<input type="radio"/> 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)		09357 895,00

Date Order User Address

PHYSIO - Posture backrest (Std. Backrest upholstery not provided)





Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

- AIR** Carbon shell
not available for seat width 48

BACKREST HEIGHT	28	34	40			
Shell's shape 	Lateral retaining wing A=high - R=regular					
	A	R	A	R	A	R
						
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

1.121,00

- POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape 	Lateral retaining wing R=regular		
	R	R	R
			
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>


1.085,00

Configuration DEEP with side support PAD (height adjustable)

Alluminum lateral with padded, breathable upholstery.



Small (h 13cm)



- Right 10272

49,00

- Left 10272

49,00

Large (h 20cm)



- Right 10273


49,00

- Left 10273

49,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Alluminum lateral with padded, breathable upholstery.



- Large (18x12)

05647 196,00

- Small (14x9)

05643 196,00

- Lateral support extension

Quantity (number pcs).....
10271 74,00

- Lateral support extension

Quantity (number pcs).....

- Small (14x9)

05642 196,00

- Large (18x12)

05646 196,00

- Headrest adjustable - Double ball mounting system

Available only with backrest Physio Posture



- Large (35x12,5)
- Small (25x12,5)

05648 385,00

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888



customer_support_uk@permobil.com - https://www.permobil.com/en-gb

Date	Order	User	Address
------------	-------------	------------	---------------

Note

Wheelchair : Progeo TEKNA ADVANCE - Rev. 2026.0	Prices - VAT excluded <small>to be added to the price list</small>	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th colspan="2" style="text-align: left; padding: 2px;">Send this page only if any optional have been chosen</th> </tr> <tr> <td style="width: 15%; text-align: center; padding: 2px;"><input type="checkbox"/></td> <td style="padding: 2px;"> <input checked="" type="checkbox"/> standard <small>(supplied if no requests are indicated)</small> </td> </tr> <tr> <td style="text-align: center; padding: 2px;"><input type="checkbox"/></td> <td style="padding: 2px;"> <input type="checkbox"/> on request <small>(without extra charge)</small> </td> </tr> <tr> <td style="text-align: center; padding: 2px;"><input type="checkbox"/></td> <td style="padding: 2px;"> <input type="checkbox"/> optional <small>(cost to be added at the basic price)</small> </td> </tr> </table>	Send this page only if any optional have been chosen		<input type="checkbox"/>	<input checked="" type="checkbox"/> standard <small>(supplied if no requests are indicated)</small>	<input type="checkbox"/>	<input type="checkbox"/> on request <small>(without extra charge)</small>	<input type="checkbox"/>	<input type="checkbox"/> optional <small>(cost to be added at the basic price)</small>
Send this page only if any optional have been chosen										
<input type="checkbox"/>	<input checked="" type="checkbox"/> standard <small>(supplied if no requests are indicated)</small>									
<input type="checkbox"/>	<input type="checkbox"/> on request <small>(without extra charge)</small>									
<input type="checkbox"/>	<input type="checkbox"/> optional <small>(cost to be added at the basic price)</small>									

REFERENCE TABLES

1		Front height			
With closed frame (95°) max 5" caster					
With closed frame (90°) max 5" caster					
With closed frame (85°) max 4" caster					
FORK		Casters			
		3"	4"	5"	6"
SMALL		max 49	max 50	max 51	
MEDIUM		max 53	max 54	max 55	max 57
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.					

2		Rear height			
		Rear wheels			
		22"	24"	25"	26"
H min		min. 36	min. 38	min. 39	min. 40
H max		max 44	max 48	max 49	max 51

3	Footplate option	Footplate distance
	Automatic Closure footplate	min. 36
	Automatic Closure footplate , with multifunction rotative system	min. 31
	One piece footplate - Flip-up!	min. 36
	One piece footplate - Flip-up!, with multifunction rotative system	min. 31
	One piece aluminium footplate - Flip-up!	min. 36
	Plastic footrest	min. 36
	One piece Carbon fibre footplate	min. 36
Calculated with the best configuration and with the footplates positioned on inside.		
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all couloers)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)

Date	Order	User	Address
------	-------	------	---------

12" wheel - Motor 250W (750W maximum power)
 14" wheel - Motor 1000W (1500W maximum power)

- Color display with switcher
- Maximum speed limit settable from display
- 5 speed levels selectable
- 10 display-selectable acceleration levels
- Adaptive Cruise Control
- Disc and electronic brake control lever
- Motor brake (adjustable in 10 levels from the display)
- Electronic brake with control bottom (adjustable in 10 levels from the display)
- Power regain when braking
- Error messages management from display (with error history)
- Rear gear beeper
- 48V - 5.8 Ah battery Autonomy up to 25km
- 48V - 11.6 Ah Power plus battery Autonomy up to 50km
- Battery charger 48V
- Rear gear
- USB port
- Easle
- Mudguard
- Front spoiler with predisposition for FULL LED light

(Certified software)

Mototronik

10979 N/A



Handlebar

☆ STANDARD controls

☆ Thumb accelerator

☆ Right
12592

Left
12593

Grip throttle


Right
10913

Left
10914

Disc and electronic brake control lever


☆ Right
12594

Left
12595




QUADRIPLEGIC controls
 Handlebar PUSH&PULL

10905 428,00




☆ **MOTOTRONIK 12" wheel**

Motor 250W (750W maximum power)
 Maximum speed according to law



MOTOTRONIK 14" wheel

Motor 1000W (1500W maximum power)
 Maximum speed according to law



10903 428,00

Battery







☆ 48V - 5.8 Ah battery *Autonomy up to 25 km*

48V - 11.6 Ah Power plus battery *Autonomy up to 50 km*

10904 289,00

Date Order User Address

Accessories

- | | | | | |
|---|---|---|--------|--------|
| <input type="radio"/> FULL LED light (high intensity)
Integrated in the spoiler |  | 10908 | 107,00 | |
| <input type="radio"/> Handlebar extension
Adjustable from 15cm up to 23cm |  | 10907 | 182,00 | |
| <input type="radio"/> Motor unit facilitated unlock lever
Recommended with handlebar for quadriplegic |  | 10915 | 82,00 | |
| <input type="radio"/> Additional FLY 48V - 2.9 Ah battery
Autonomy up to 15km
Approved for air transport | | 10906 | 697,00 | |
| <input type="radio"/> Rearview mirror | <input type="radio"/> Right
<input type="radio"/> Left |  | 10909 | 32,00 |
| | | | 10909 | 32,00 |
| <input type="radio"/> Basket (removable) - (with KLICKfix® attachment system) | |  | 10912 | 212,00 |
| <input type="radio"/> Removable front weight - (with KLICKfix® attachment system)
(it increases the weight over the power wheel to improve traction) | |  | 10910 | 311,00 |
| <input type="radio"/> Travel bag | | | 10916 | 155,00 |

Note

Send this page only if any optional have been chosen

Wheelchair : Progeo TEKNA ADVANCE - Rev. 2026.0

Prices - VAT excluded
to be added to the price list

- standard (supplied if no requests are indicated)
- on request (without extra charge)
- optional (cost to be added at the basic price)