



# PDG Bentley LT - AADL W341



Dealer Name\*: \_\_\_\_\_ Account #\*: \_\_\_\_\_

City, Province, Post Code: \_\_\_\_\_

Phone: \_\_\_\_\_ Email: \_\_\_\_\_ ATP: \_\_\_\_\_

Rep Name: \_\_\_\_\_ P.O. Number: \_\_\_\_\_ Quote Number: \_\_\_\_\_

Client Reference\*: \_\_\_\_\_ Client Weight: \_\_\_\_\_

Special considerations that need to be addressed here (e.g., diagnosis): \_\_\_\_\_ \*Required field

## October 2024, Revised January 2026. Prices are in Canadian Dollars (CAD)

### Wheelchair & Frame

Bentley 20° Manual Tilt Wheelchair Base Model	W341	\$ 2,685.00
Bentley LT Lightweight Aluminum Frame Upgrade <i>Aluminum frame not recommended for agitated users. Consider standard Bentley with steel frame option for agitated or heavier users. See separate order form.</i>		\$ 300.00

#### Frame Type | All frames offer 300lbs weight capacity.

Fixed Width Frame <i>Select a fixed seat width below.</i>	No charge
Adjustable Width Frame <i>Select an adjustable seat width below.</i>	\$ 200.00

#### User Weight | Required for setup of tilt uplift force.

<b>User Weight Range:</b>	Up to 75lbs	75lbs - 100lbs	100lbs - 115lbs	115lbs - 150lbs	150lbs - 200lbs	200lbs - 250lbs
	250lbs - 300lbs					

#### Frame Options

Transit Tie Down System (TTDS) <i>Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See owners manual for complete details.</i>	No charge
---	-----------

### Seat Measurements

#### Seat Width | Overall chair width is +9.5" with 12"/16" rear wheels, +10.5" with 20"/22"/24" rear wheels.

<b>Fixed:</b>	15"    16"    17"    18"    19"	No charge
	20"    21"    22"	W111    \$ 250.00
	23"    24"	W111    \$ 500.00
<b>Adjustable 15"-20":</b>	15"    16"    17"    18"    19"	No charge
	20"	W111    \$ 250.00

#### Seat Depth

<b>Select Depth:</b>	16"    17"    18"    19"	No charge
	20"    21"    22"	W111    \$ 250.00

16" - 20" seat depths are adjustable.

### Frame Finishes

<b>Select Color:</b>	High Gloss Black	Hyper Blue
----------------------	------------------	------------

## Armrests

### Armrest Type

No Armrests No charge

Single Post, Height Adjustable Armrests No charge

8"-11" Range    10-14" Range    12"-16" Range

Armrest Position:    Rear (Default)    Front

*If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.*

### Armpads

Regular Armpads No charge

Desk Length 2" x 10"    Full Length 2" x 14"

### Clothing Guards

Clothing Guards, Single Post Armrests *Not compatible with 8"-11" height range.* \$ 125.00

### Armrest Extensions | *Not compatible with flip-up armrests or desk length armpads.*

Armrest Extensions - 3" Forward \$ 75.00

Armrest Extensions - 3" Rearward \$ 75.00

## Seat Height

### Seat Height - 6" Casters

Low						High						
12"												
16"												
20"												
22"												
24"												
	13"	14"	15"	16"	17"	18"	15"	16"	17"	18"	19"	20"
	Seat to Floor Height											

### Seat Height - 8" Casters

Low				High				
16"								
24"								
	15"	16"	17"	18"	17"	18"	19"	20"
	Seat to Floor Height							

## Seat Pan

No Seat Pan No charge

Solid Aluminium Seat Pan No charge

## Front Casters

### Urethane Casters

6" x 1.25" Caster No charge

8" x 1.25" Caster *Only compatible with 16" and 24" rear wheels.* No charge

## Rear Wheels

Composite Rear Wheel No charge

12"    16"    20"    22"    24"

Rear Wheel Position:    Rear    Front (Default)

*If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.*

### Axles

Quick Release Rear Axles No charge

Fixed Rear Axles No charge

**Rear Tires**

Urethane Tires

No charge

## Handrims

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	\$ 175.00

## Wheel Locks

User Wheel Locks	No charge
Push to Lock   Pull to Lock	

## Wheel Lock Options

Wheel Lock Extensions <i>Not compatible with attendant wheel lock.</i>	No charge
--	-----------

## Anti-Tippers

Rear Anti-Tippers <i>Required safety option.</i>	No charge
--	-----------

## Backrest

### Stroller Handle Back

Back Post Height:	Short (22")	Tall (25")			
	<i>Measured from seat pan to handle.</i>				
Folding Style:	Fold-Down				
	<i>Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.</i>				
Seat to Back Angle:	Set at 90°	Set at 98°	Set at 106°	Set at 114°	Set at 122°
	<i>Back angles &gt;98° in combination with seat depths &gt;19" may limit range of tilt.</i>				

## Backrest Options

Angle/Height Adjustable Stroller Handle:	No charge
Button on Outside   Button on Inside	
Spring Back Suspension <i>Not compatible with dynamic rocker back or reclining backrest.</i>	\$ 105.00

## Tilt Activation

### Attendant Tilt Only

Attendant Tilt Only, Right Back Post	No charge
Attendant Tilt Only, Left Back Post	No charge
Two Trigger Setup <i>A trigger mounted on each side of the stroller handle. Reduces grip force required to activate tilt. Not compatible with reclining backrest or drum brakes.</i>	No charge

## Front Rigging

No Front Rigging	No charge
Easy-Out Front Rigging   <i>Effortless insertion and removal. No hanger brackets remain on the frame.</i>	
Fixed 70°	\$ 95.00
Swing Away Front Rigging   <i>Legacy design. Not compatible with Easy-Out.</i>	
Swing Away 70°	No charge
Swing Away Elevating Legrests (ELR) <i>Includes calf pads.</i>	W106 \$ 375.00
1" Offset Hangers for Swing Away Front Rigging	\$ 95.00
Plug In Front Rigging   <i>Legacy design. Not compatible with Easy-Out.</i>	
Plug In 70°	No charge

## Footplates

Composite Footplates		\$ 0.00
Aluminium Angle Adjustable Footplates	W107	\$ 165.00

### Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short		No charge
Standard Extension Tubes, Long		\$ 75.00

Footplate Height: Specify:

Ext	Short	11" - 16"	12" - 16"	13" - 16"
	Long	14" - 22"	14" - 19"	14" - 17"
		Standard Easy-Out	60°, 70°, ELR Swing-Away	Standard Plug-In
<b>Legrest</b>				

## Front Rigging Accessories

Length Adjustable Heel Loops		No charge
Standard Calf Strap	Basic Seating	\$ 50.00

## General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

## Notes