

Room to grow

Small chair with mini kit at its smallest dimensions

- 1 Backrest height 12½"
- 2 Trunk support width 5½"
- 3 Seat depth 7"
- 4 Seat width 7"
- 5 Armrest height above seat 5"
- 6 Seat height above footboard 7"



Small chair at its largest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 11½"
- 3 Seat depth 12"
- 4 Seat width 12"
- 5 Armrest height above seat 7½"
- 6 Seat height above footboard 12"



Small standard base
R820 Rifton Activity Chair

User dimensions (inches)	
Height	32–48
with mini kit	28–38
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.	
Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.	
Item dimensions (inches)	R820 small Standard base
Frame width	short legs: 21 long legs: 23
	short legs w/ casters: 22½ long legs w/ casters: 23½
Seat height above floor	short legs: 9½–12½ long legs: 18½–21½
	short legs w/ casters: 13½–16½ long legs w/ casters: 18½–21½
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7–12
Seat width with hip guides	7–9
without hip guides	12
Seat depth	8–12
with mini kit	7–11
Armrest height above seat	5–7½
Trunk support width	5½–11½
Backrest height	12½–15½
Headrest height above seat	14½–21
Max. working load (lb)	75



Small hi/lo base
R830 Rifton Activity Chair

User dimensions (inches)	
Height	32–48
with mini kit	28–38
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.	
Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.	
Item dimensions (inches)	R830 small Hi/lo base
Frame width	26
Seat height above floor	10–25
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7–12
Seat width with hip guides	7–9
without hip guides	12
Seat depth	8–12
with mini kit	7–11
Armrest height above seat	5–7½
Trunk support width	5½–11½
Backrest height	12½–15½
Headrest height above seat	14½–21
Max. working load (lb)	75

Retail order form - January 2025

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.



* Required		Standard R820 small	Hi/lo R830 small
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R821 \$1,050 ☐ R831 \$2,029 ☐
		Adjustable w/ back spring	R822 \$1,203 ☐ R832 \$2,251 ☐
* Base	CHOOSE	Without spring	R823 \$610 ☐ R833 \$3,667 ☐
		With spring	R824 \$719 ☐ N/A
* Arm supports (pair)	CHOOSE	Armrests	R827 \$429 ☐ R827 \$429 ☐
		Forearm prompts	R809 \$1,050 ☐ R809 \$1,050 ☐
* Adjustable legs	CHOOSE	Short	R803 \$86 ☐ N/A
		Long	R894 \$148 ☐ N/A
		Short with casters	R805 \$522 ☐ N/A
		Long with casters	R896 \$578 ☐ N/A
* Pads	CHOOSE	Tan	R815 \$429 ☐ R815 \$429 ☐
		Blue	R816 \$429 ☐ R816 \$429 ☐
		Red	R817 \$429 ☐ R817 \$429 ☐
		Pink	R971 \$429 ☐ R971 \$429 ☐
		Purple	R972 \$429 ☐ R972 \$429 ☐
		Green	R973 \$429 ☐ R973 \$429 ☐
Headrest	CHOOSE	Contoured	R897 \$429 ☐ R897 \$429 ☐
		Flat	R898 \$396 ☐ R898 \$396 ☐
		Adjustable winged	R881 \$531 ☐ R881 \$531 ☐
Laterals (pair)	CHOOSE	Small	R828 \$346 ☐ R828 \$346 ☐
		Small with chest strap	R819 \$429 ☐ R819 \$429 ☐
Push handles (recommended for all chairs with wheels)		R878 \$282 ☐ Included	
Tray (requires armrests)		R667 \$429 ☐ R667 \$429 ☐	
Handhold		R649 \$44 ☐ R649 \$44 ☐	
Chest strap (wide)		R802 \$132 ☐ R802 \$132 ☐	
Butterfly harness		R834 \$202 ☐ R834 \$202 ☐	
Thigh belt		R811 \$132 ☐ R811 \$132 ☐	
Pelvic harness		R845 \$266 ☐ R845 \$266 ☐	
Hip guides (pair)		R814 \$357 ☐ R814 \$357 ☐	
Abductor		R884 \$218 ☐ R884 \$218 ☐	
Adductors (pair)		R892 \$346 ☐ R892 \$346 ☐	
Leg prompt		R854 \$382 ☐ R854 \$382 ☐	
Footboard (recommended for chairs with long legs)		R847 \$762 ☐ Included	
Mini kit (backrest insert and trunk support)		R812 \$542 ☐ R812 \$542 ☐	
Ankle straps (pair)		R899 \$148 ☐ R899 \$148 ☐	
Sandals (pair)		R944 \$156 ☐ R944 \$156 ☐	
Wedges (pair, require sandals)		R697 \$124 ☐ R697 \$124 ☐	
Lumbar and seat support kit		R886 \$124 ☐ R886 \$14 ☐	
Backrest filler pad		R849 \$218 ☐ R849 \$218 ☐	
Whitmyer headrest adapter		R804 \$35 ☐ R804 \$35 ☐	
Total		R820 \$ _____ R830 \$ _____	

Room to grow

Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 7¾"




Medium chair at its largest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 15¼"



Medium standard base
R840 Rifton Activity Chair

User dimensions (inches)	
Height	40 – 56
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.	
Item dimensions (inches)	R840 medium Standard base
Frame width	short legs: 23¼ long legs: 24½
	short legs w/ casters: 24½ long legs w/ casters: 25½
Seat height above floor	short legs: 12½ – 16½ long legs: 19 – 23
	short legs w/ casters: 16 – 20 long legs w/ casters: 19 – 23
Seat angle (tilt-in-space) – degrees	15° forward, 15° back
Backrest angle – degrees	5° forward, 20° back
Footboard knee angle – degrees	45° – 110°
Seat height above footboard	7¾ – 15¼
Seat width with hip guides without hip guides	8½ – 11½ 14
Seat depth	11 – 16
Armrest height above seat	7 – 10½
Distance between laterals	6½ – 12
Backrest height	15½ – 19
Headrest height above seat	17½ – 24½
Max. working load (lb)	150



Medium hi/lo base
R850 Rifton Activity Chair

User dimensions (inches)	
Height	40 – 56
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.	
Item dimensions (inches)	R850 medium Hi/lo base
Frame width	27½
Seat height above floor	12 – 23
Seat angle (tilt-in-space) – degrees	15° forward, 25° back
Backrest angle – degrees	5° forward, 20° back
Footboard knee angle – degrees	45° – 110°
Seat height above footboard	7¾ – 15¼
Seat width with hip guides without hip guides	8½ – 11½ 14
Seat depth	11 – 16
Armrest height above seat	7 – 10½
Distance between laterals	6½ – 12
Backrest height	15½ – 19
Headrest height above seat	17½ – 24½
Max. working load (lb)	150

Retail order form - January 2025

Use dimension chart to select appropriate size. All Rifton products are made to order and are non-returnable and non-exchangeable.



* Required		Standard R840 medium		Hi/lo R850 medium		
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R841	\$1,074 ☐	R851	\$2,104 ☐
		Adjustable w/ back spring	R842	\$1,236 ☐	R852	\$2,299 ☐
* Base	CHOOSE	Without spring	R843	\$626 ☐	R853	\$4,522 ☐
		With spring	R844	\$751 ☐	N/A	
* Arm supports (pair)	CHOOSE	Armrests	R868	\$440 ☐	R868	\$440 ☐
		Forearm prompts	R869	\$1,098 ☐	R869	\$1,098 ☐
* Adjustable legs	CHOOSE	Short	R893	\$98 ☐	N/A	
		Long	R894	\$148 ☐	N/A	
		Short with casters	R895	\$531 ☐	N/A	
		Long with casters	R896	\$578 ☐	N/A	
* Pads	CHOOSE	Tan	R835	\$440 ☐	R835	\$440 ☐
		Blue	R836	\$440 ☐	R836	\$440 ☐
		Red	R837	\$440 ☐	R837	\$440 ☐
		Pink	R974	\$440 ☐	R974	\$440 ☐
		Purple	R975	\$440 ☐	R975	\$440 ☐
		Green	R976	\$440 ☐	R976	\$440 ☐
Headrest	CHOOSE	Contoured	R897	\$429 ☐	R897	\$429 ☐
		Flat	R898	\$396 ☐	R898	\$396 ☐
		Adjustable winged	R881	\$531 ☐	R881	\$531 ☐
Laterals (pair)	CHOOSE	Small	R828	\$346 ☐	R828	\$346 ☐
		Small with chest strap	R829	\$440 ☐	R829	\$440 ☐
		Large	R838	\$382 ☐	R838	\$382 ☐
Push handles (recommended for all chairs with wheels)		R878	\$282 ☐	Included		
Tray (requires armrests)		R668	\$487 ☐	R668	\$487 ☐	
Handhold		R649	\$44 ☐	R649	\$44 ☐	
Chest strap (wide)		R863	\$148 ☐	R863	\$148 ☐	
Butterfly harness	CHOOSE	Regular	R864	\$218 ☐	R864	\$218 ☐
		Slim cut	R882	\$218 ☐	R882	\$218 ☐
Thigh belt		R865	\$148 ☐	R865	\$148 ☐	
Pelvic harness		R866	\$282 ☐	R866	\$282 ☐	
Hip guides (pair)		R867	\$372 ☐	R867	\$372 ☐	
Abductor	CHOOSE	Small	R884	\$218 ☐	R884	\$218 ☐
		Large	R891	\$256 ☐	R891	\$256 ☐
Adductors (pair)		R892	\$346 ☐	R892	\$346 ☐	
Leg prompt		R858	\$396 ☐	R858	\$396 ☐	
Footboard (recommended for chairs with long legs)		R848	\$785 ☐	Included		
Ankle straps (pair)		R899	\$148 ☐	R899	\$148 ☐	
Sandals (pair)		R954	\$170 ☐	R954	\$170 ☐	
Wedges (pair, require sandals)		R698	\$124 ☐	R698	\$124 ☐	
Lumbar and seat support kit		R887	\$148 ☐	R887	\$148 ☐	
Backrest filler pad		R849	\$218 ☐	R849	\$218 ☐	
Whitmyer headrest adapter		R804	\$35 ☐	R804	\$35 ☐	
Total		R840	\$	R850	\$	

Room to grow

Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 11¼"



Large chair at its largest dimensions


- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 21¼"



Large standard base
R860 Rifton Activity Chair

User dimensions (inches)

Height 50–74



Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

Item dimensions (inches)

R860 large Standard base

Frame width	short legs: 26¾ long legs: 26¾
	short legs w/ casters: 27¾ long legs w/ casters: 27¾
Seat height above floor	short legs: 16½–20½ long legs: 20–24
	short legs w/ casters: 16½–20½ long legs w/ casters: 20–24
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250



Large hi/lo base
R870 Rifton Activity Chair

User dimensions (inches)

Height 50–74



Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

Item dimensions (inches)

R870 large Hi/lo base

Frame width	29¾
Seat height above floor	16–24
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250

Retail order form - January 2025

Use dimension chart to select appropriate size. All Rifton products are made to order and are non-returnable and non-exchangeable.



* Required		Standard R860 large		Hi/lo R870 large	
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R861 \$1,222 ☐	R871 \$2,283 ☐	
		Adjustable w/ back spring	R862 \$1,410 ☐	R872 \$2,448 ☐	
* Base	CHOOSE	Without spring	R873 \$751 ☐	R883 \$5,075 ☐	
		With spring	R874 \$877 ☐	N/A	
* Arm supports (pair)	CHOOSE	Armrests	R868 \$440 ☐	R868 \$440 ☐	
		Forearm prompts	R879 \$1,146 ☐	R879 \$1,146 ☐	
* Adjustable legs	CHOOSE	Short	R808 \$124 ☐	N/A	
		Long	R894 \$148 ☐	N/A	
		Short with casters	R895 \$531 ☐	N/A	
		Long with casters	R896 \$578 ☐	N/A	
* Pads	CHOOSE	Tan	R875 \$464 ☐	R875 \$464 ☐	
		Blue	R876 \$464 ☐	R876 \$464 ☐	
		Red	R877 \$464 ☐	R877 \$464 ☐	
		Pink	R977 \$464 ☐	R977 \$464 ☐	
		Purple	R978 \$464 ☐	R978 \$464 ☐	
		Green	R979 \$464 ☐	R979 \$464 ☐	
Headrest	CHOOSE	Contoured	R897 \$429 ☐	R897 \$429 ☐	
		Flat	R898 \$396 ☐	R898 \$396 ☐	
		Adjustable winged	R881 \$531 ☐	R881 \$531 ☐	
Laterals (pair)	CHOOSE	Small	R828 \$346 ☐	R828 \$346 ☐	
		Large	R838 \$382 ☐	R838 \$382 ☐	
		Large with chest strap	R839 \$464 ☐	R839 \$464 ☐	
Push handles (recommended for all chairs with wheels)		R878 \$282 ☐	Included		
Tray (requires armrests)		R669 \$510 ☐	R669 \$510 ☐		
Handhold		R649 \$44 ☐	R649 \$44 ☐		
Chest strap (wide)		R807 \$170 ☐	R807 \$170 ☐		
Butterfly harness	CHOOSE	Regular	R889 \$256 ☐	R889 \$256 ☐	
		Slim cut	R826 \$256 ☐	R826 \$256 ☐	
Thigh belt		R813 \$170 ☐	R813 \$170 ☐		
Pelvic harness		R846 \$314 ☐	R846 \$314 ☐		
Hip guides (pair)		R867 \$372 ☐	R867 \$372 ☐		
Abductor	CHOOSE	Small	R884 \$218 ☐	R884 \$218 ☐	
		Large	R891 \$256 ☐	R891 \$256 ☐	
Adductors (pair)		R892 \$346 ☐	R892 \$346 ☐		
Leg prompt		R859 \$414 ☐	R859 \$414 ☐		
Footboard (recommended for chairs with long legs)		R857 \$843 ☐	Included		
Ankle straps (pair)		R855 \$170 ☐	R855 \$170 ☐		
Sandals (pair)		R964 \$188 ☐	R964 \$188 ☐		
Wedges (pair, require sandals)		R699 \$132 ☐	R699 \$132 ☐		
Lumbar and seat support kit		R888 \$156 ☐	R888 \$156 ☐		
Backrest filler pad		R856 \$218 ☐	R856 \$218 ☐		
Whitmyer headrest adapter		R804 \$35 ☐	R804 \$35 ☐		
Total		R860 \$	R870 \$		