

Manual wheelchairs

Tailored performance



permobil



	Titanium performance		Titanium adjustable		Aluminum adjustable		Folding wheelchairs	
	TiLite TR	TiLite ZR	TiLite TRA	TiLite ZRA	TiLite Aero T	TiLite Aero Z	TiLite Aero X	TiLite 2GX
Frame features	1" Titanium Dual-Tube Ultimate TiFit	1.25" Titanium Mono-Tube Ultimate TiFit	1" Titanium Dual-Tube Adjustable TiFit	1.25" Titanium Mono-Tube Adjustable TiFit	1" Aluminum Dual-Tube Adjustable TiFit	1.25" Aluminum Mono-Tube Adjustable TiFit	1" Aluminum Box Frame Modular Design	1" Titanium Box Frame Modular Design
Seat width	12" – 20"	10" – 20"	12" – 20"	10" – 20"	12" – 20"	10" – 20"	12" – 22"	12" – 22"
Seat depth	12" – 20"	10" – 20"	12" – 20"	10" – 20"	12" – 20"	10" – 20"	12" – 20"	12" – 20"
Front seat height	15" – 21"	16" – 21"	15" – 21"	16" – 21"	15" – 21"	16" – 21"	13.5" – 21.5"	13.5" – 21"
Rear seat height	12.5" – 21"	12.5" – 20"	13" – 21"	12.5" – 20"	13" – 21"	12.5" – 20"	13" – 20"	13" – 20"
Front frame angle	60°, 65°, 70°, 75°, 80°, 85°, or 90°	70°, 75°, 80°, 85°, or 90°	70°, 75°, 80°, 85°, or 90°	70°, 75°, 80°, 85°, or 90°	70°, 75°, 80°, 85°, or 90°	70°, 75°, 80°, 85°, or 90°	Fixed Front — 70°, 80°, or 85° Swing Away — 70°, 80°, 85°, or 90°	Fixed Front — 70°, 80°, or 85° Swing Away — 70°, 80°, 85°, or 90°
Seat back type	Folding Seat Back <i>Adj. Height 9.5" – 22"</i> <i>Adj. Angle 80° – 101°</i> Fixed Seat Back <i>Fixed Height 5" – 20"</i> <i>Fixed Angle 80° – 101°</i>	Folding Seat Back <i>Adj. Height 9.5" – 22"</i> <i>Adj. Angle 80° – 101°</i> Fixed Seat Back <i>Fixed Height 8" – 20"</i> <i>Fixed Angle 80° – 101°</i>	Folding Seat Back <i>Adj. Height 9.5" – 22"</i> <i>Adj. Angle 80° – 101°</i>	Folding Seat Back <i>Adj. Height 9.5" – 22"</i> <i>Adj. Angle 80° – 101°</i>	Folding Seat Back <i>Adj. Height 9.5" – 22"</i> <i>Adj. Angle 80° – 101°</i>	Folding Seat Back <i>Adj. Height 9.5" – 22"</i> <i>Adj. Angle 80° – 101°</i>	Folding Seat Back <i>Adj. Height 11" – 20.5"</i> <i>Adj. Angle 80° – 101°</i> Fixed Seat Back <i>Adj. Height 11" – 20"</i>	Folding Seat Back <i>Adj. Height 11" – 20.5"</i> <i>Adj. Angle 80° – 101°</i> Fixed Seat Back <i>Adj. Height 11" – 20"</i>
Center of gravity adjustment	5"	2.25"	5"	6"	5"	6"	6"	6"
Transit options*	Q-Straint Lock Down Approved	Transport Bracket	Q-Straint Lock Down Approved	Transit Tie-Down	Q-Straint Lock Down Approved	Transit Tie-Down	Transit Tie-Down Transport Bracket	Transport Bracket
Transport weight**	9.3 lb.	9.3 lb.	11.1 lb.	10.7 lb.	12 lb.	12.4 lb.	Fixed Front — 15.6 lb. Swing Away — 18 lb.	Fixed Front — 13.1 lb. Swing Away — 15.5 lb.
Weight limit	265 lb.	Standard — 265 lb. Heavy Duty — 300 lb.	Standard — 265 lb. Heavy Duty — 350 lb.	Standard — 265 lb. Heavy Duty — 300 lb.	Standard — 265 lb. Heavy Duty — 350 lb.	265 lb.	Standard — 265 lb. Heavy Duty — 350 lb.	Standard — 265 lb. Heavy Duty — 350 lb.

*See Transit Option Flyers and Owners Manuals at www.permobilus.com for specifications, limitations and precautions. **Weights based on 16" x 16" wheelchairs (12" x 12" for youth chairs) with all standard features except rear wheels. For reference, a pair of standard 24" [540] rear wheels weighs 7.2 lb. For youth chairs, a pair of standard 22" [501] rear wheels weighs 7.0 lb. A pair of standard wheel locks weighs 0.9 lb.



Pediatric

TiLite Twist

TiLite Pilot

Frame features	1" Aluminum Dual-Tube Adjustable TiFit	1" Carbon/Aluminum Mono-Tube Growable
Seat width	8" – 15"	8" – 15"
Seat depth	8" – 16"	9" – 15"
Front seat height	13" – 20"	12" – 16.5"
Rear seat height	12" – 19.5"	12" – 16.5"
Front frame angle	70°, 80°, or 90°	–
Seat back type	Folding Seat Back Adj. Height 8.5" – 21" Adj. Angle 80° – 101°	Non-Folding Adjustable-Height Rigid Back Adj. Height 2" Adj. Angle 80° – 101°
Center of gravity adjustment	6"	3"
Transit options*	Transit Tie-Down	Transit Tie-Down
Transport weight**	12 lb.	11 lb.
Weight limit	165 lb.	120 lb.

*See Transit Option Flyers and Owners Manuals at www.permobilus.com for specifications, limitations and precautions. **Weights based on 16" x 16" wheelchairs (12" x 12" for youth chairs) with all standard features except rear wheels. For reference, a pair of standard 24" (540) rear wheels weighs 7.2 lb. For youth chairs, a pair of standard 22" (501) rear wheels weighs 7.0 lb. A pair of standard wheel locks weighs 0.9 lb.