



PDG Eclipse - AADL W320



Dealer Name*: _____ Account #*: _____

City, Province, Post Code: _____

Phone: _____ Email: _____ ATP: _____

Rep Name: _____ P.O. Number: _____ Quote Number: _____

Client Reference*: _____ Client Weight: _____

Special considerations that need to be addressed here (e.g., diagnosis): _____ *Required field

October 2024, Revised January 2026. Prices are in Canadian Dollars (CAD)

Wheelchair & Frame

Eclipse Bariatric Manual Wheelchair W320 \$ 3,100.00

Frame Type

600lbs Weight Capacity No charge

1000lbs Weight Capacity *Not compatible with Transit Tie Down System* \$ 300.00

Frame Options

Transit Tie Down System (TTDS) *TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. weight capacity of 600lbs. Not compatible with O2 holder or IV pole.* \$ 395.00

Seat Measurements

Seat Width | Overall chair width is +7.5"

| | | | | | | |
|---------------|----------------|----------------|----------------|----------------|------|-------------|
| Select Width: | 20" Seat Width | 22" Seat Width | 24" Seat Width | 26" Seat Width | W111 | \$ 675.00 |
| | 28" Seat Width | 30" Seat Width | | | W111 | \$ 875.00 |
| | 32" Seat Width | 36" Seat Width | | | W111 | \$ 1,100.00 |
| | 40" Seat Width | 44" Seat Width | 48" Seat Width | | W111 | \$ 1,300.00 |

Seat Depth (Adjustable)

| | | | | | | |
|---------------|----------------|----------------|----------------|----------------|----------------|-----------|
| Select Depth: | 16" Seat Depth | 17" Seat Depth | 18" Seat Depth | 19" Seat Depth | 20" Seat Depth | No charge |
| | 21" Seat Depth | 22" Seat depth | | | | |
| | 23" Seat Depth | 24" Seat Depth | | | W111 | \$ 400.00 |

24" seat depth not available with TTDS.

Frame Finishes

Select Color: High Gloss Black Hyper Blue

Armrests

Armrest Type

| | |
|--|-----------|
| No Armrests | No charge |
| Dual Post, Height Adjustable Armrests | \$ 295.00 |
| 10-14" Range 12"-16" Range | |
| Dual Post, Fixed Height Armrests (12") | No charge |

Armpads

| | |
|--|--------------|
| Regular Armpads | No charge |
| Desk Length 2" x 10" Full Length 2" x 14" | |
| Comfort Armpads with Interlocking Hardware <i>May interfere with wheels in some tilt configurations.</i> | A1 \$ 200.00 |
| Desk Length 3.75" x 10" Full Length 3.75" x 14" | |
| Gel Armpads <i>May interfere with wheels in some tilt configurations.</i> | \$ 135.00 |
| Narrow 2" x 12" Wide 3.5" x 12" | A1 |

Armrest Options

| | |
|---|-----------|
| Forward Guards | \$ 200.00 |
| Armrest Extensions - 3' Forward <i>Not compatible with desk length armpads.</i> | \$ 75.00 |

Seat Height

Seat Height | *Front and rear height can be maximum 1" apart.*

| | | | | | | |
|--------------------|---------------|---------------|---------------|---------------|---------------|-----------|
| Front Seat Height: | 14" STF Front | 15" STF Front | 16" STF Front | 17" STF Front | 18" STF Front | |
| | 19" STF Front | 20" STF Front | | | | \$ 200.00 |
| Rear Seat Height: | 14" STF Rear | 15" STF Rear | 16" STF Rear | 17" STF Rear | 18" STF Rear | |
| | 19" STF Rear | 20" STF Rear | | | | \$ 200.00 |

Seat Pan

| | |
|--------------------------|-----------|
| Solid Aluminium Seat Pan | \$ 165.00 |
|--------------------------|-----------|

Front Casters

| | |
|---------------------------------------|-----------|
| 5" Casters with Heavy Duty Steel Fork | No charge |
|---------------------------------------|-----------|

Rear Wheels

| | |
|-----------------------------|-----------|
| 24" Heavy Duty Spoke Wheels | No charge |
|-----------------------------|-----------|

| | | | | | |
|----------------------|---|------------|------------|--------------------|--|
| Rear Wheel Position: | Position 1 (Rear) | Position 2 | Position 3 | Position 4 (Front) | |
| | <i>Positions 3 and 4 not compatible with TTDS. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.</i> | | | | |

Axles

| | |
|--------------------------|-----------|
| Quick Release Rear Axles | No charge |
|--------------------------|-----------|

Rear Tires

| | |
|----------------|-----------|
| Urethane Tires | No charge |
|----------------|-----------|

Handrims

| | |
|-----------------------------|-----------|
| No Handrims | No charge |
| Standard Aluminium Handrims | No charge |
| Plastic Coated Handrims | \$ 175.00 |

Wheel Locks

| | |
|------------------------------|-----------|
| User Wheel Locks | No charge |
| Push to Lock Pull to Lock | |
| Wheel Lock Extensions | No charge |

Anti-Tippers

Rear Anti-Tippers *Required safety option.*

\$ 75.00

Backrest

Back Post Style: Stroller Handle Back Cane Handle Back
Cane handle not compatible with headrest mounting plate.

Back Post Height: Short (22") Tall (25")
Measured from seat pan to handle.

Seat to Back Angle: Set at 85° Set at 90° Set at 95° Set at 100°

Backrest Options

Endomorph Back *Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not available with reclining back.* \$ 395.00

1" Offset Per Side 2" Offset Per Side \$ 395.00

Removable Backrest \$ 0.00

Back Supports

Back Upholstery

No Back Upholstery No charge

Tension Adjustable Back Upholstery \$ 375.00

Position of Back Upholstery

| Uph Type | 22" Back Post | | | | | | 25" Back Post | | | | | | | | |
|----------|-------------------|-----|-----|-----|-----|-----|---------------|-----|-----|-----|-----|-----|-----|-----|-----|
| | Regular | | | | | | | | | | | | | | |
| Tall | | | | | | | | | | | | | | | |
| | 17" | 18" | 19" | 20" | 21" | 22" | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" |
| | Upholstery Height | | | | | | | | | | | | | | |

Recline

Reclining Backrest (Manual) *30° maximum recline. Not compatible with cane handle back or endomorph back. Reclining backrest is 1" narrower than seat width. Back frame is not removable. Finished back post height is 3" taller than selected back post height. Not recommended for user weight above 450lbs. Maximum seat depth is 23".* \$ 895.00

Front Rigging

No Front Rigging No charge

Center Mount Front Rigging Bracket \$ 300.00

Easy-Out Front Rigging | Effortless insertion and removal. No hanger brackets remain on the frame.

Fixed 60° \$ 95.00

Universal Adjustable Angle (0°-60°) *Select Angle. Angle adjustment by loosening one bolt. Tool-free leg length adjustment.* \$ 295.00

0° 10° 20° 30° 40° 50° 60°

AA High Mount Option *For shorter lower leg lengths of 8"-12.5".* \$ 295.00

Swing Away Front Rigging | Legacy design. Not compatible with Easy-Out.

Swing Away 60° No charge

Swing Away Elevating Legrests (ELR) *Includes calf pads.* W106 \$ 375.00

1" Offset Hangers for Swing Away Front Rigging \$ 75.00

Footplates

Aluminium Angle Adjustable Footplates W107 \$ 165.00

Extension Tubes | No selection required for Easy-Out AA High-Mount option.

Standard Extension Tubes, Short No charge

Standard Extension Tubes, Long \$ 75.00

Footplate Height: Specify:

| Ext | Short | 11" - 16" | 12" - 16" |
|-----|---------------------|-------------------|--------------------------|
| | Long | N/A | 14" - 22" |
| | High Mount Easy-Out | Standard Easy-Out | 60°, 70°, ELR Swing-Away |
| | Legrest | | |

Front Rigging Accessories

| | | |
|---|---------------|-----------|
| Oversized 9"x10" Calf Pad <i>Not compatible with seat widths less than 25"</i> | Basic Seating | \$ 695.00 |
| Calf Pads <i>Not compatible with Easy-Out front rigging.</i> | Basic Seating | \$ 135.00 |
| Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging | Basic Seating | \$ 200.00 |
| Gel Knee Pads, Swing Away ELR Front Rigging | Basic Seating | \$ 100.00 |
| Length Adjustable Heel Loops | | No charge |
| Standard Calf Strap | Basic Seating | \$ 50.00 |
| Bodypoint Padded Calf Strap <i>Available up to a maximum size of 32"</i> | Basic Seating | \$ 0.00 |

Miscellaneous Options

| | | |
|---|---------------|-----------|
| Q'Straint Lap Belt <i>Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.</i> | Basic Seating | \$ 245.00 |
| IV Pole <i>Not compatible with TTDS.</i> | A1 | \$ 295.00 |
| O2 Holder <i>Not compatible with TTDS.</i> | W110 | \$ 395.00 |
| Medical Necessities Backpack <i>Black with Permobil logo.</i> | A1 | \$ 95.00 |
| Adhesive Back Hook Strip <i>For cushion mounting, 2"x30".</i> | A1 | \$ 10.00 |

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes