



# TiLite Twist - Order Form



Dealer Name\*: \_\_\_\_\_ Account #\*: \_\_\_\_\_

City, Province, Post Code: \_\_\_\_\_

Phone: \_\_\_\_\_ Email: \_\_\_\_\_ ATP: \_\_\_\_\_

Rep Name: \_\_\_\_\_ P.O. Number: \_\_\_\_\_ Quote Number: \_\_\_\_\_

Client Reference\*: \_\_\_\_\_ Client Weight: \_\_\_\_\_

Special considerations that need to be addressed here (e.g., diagnosis): \_\_\_\_\_ \*Required field

## Frame Style

TWFS2	<b>TiLite Twist</b> Full adjustability. 165lb weight limit. Seat width 'A' grows 2" and seat depth 'B' grows 3".	\$ 3,976.00
TWFSW7	<b>Abduction Frame, 2" Outward Flare</b> Requires a No Taper Footrest.	\$ 181.00

## Transit Options

TWTTDO	<b>Transit Tie-Downs</b> Tested to WC-19. Consists of four tie-down brackets. Transit option is not available with user weights less than 51lb (23kg). Please refer to the Twist Transit Tie-Down Option Brochure for important information regarding transit tie-down options.	\$ 427.00
--------	---	-----------

## Frame Finishes

<b>Frame Color:</b>	TWFF8 High Gloss Black	TWFF10 Red Pepper	TWFF49 Neon Pink
	TWFF14 Ocean Blue Metallic	TWFF15 Tangerine Metallic	TWFF30 Acid Green Metallic
	TWFF17 Midnight Blue Pearl	TWFF18 White Pearl	TWFF19 Black Lava
	TWFF27 Titanium Grey	TWFF59 Pomegranate Red	TWFF62 Blue Mist
	TWFF57 Indigo Purple	TWFF40 Hyper Blue	TWFF41 Flame Orange
	TWFF42 OD Green	TWFF60 Silver Streak	TWFF63 Lavendar Dusk
	TWFF58 Black Galaxy	TWFF61 Seafoam Green	TWFF64 Daffodil Yellow
	<b>Caster Arm Color:</b>	TWFF8 High Gloss Black	TWFF10 Red Pepper
TWFF14 Ocean Blue Metallic		TWFF15 Tangerine Metallic	TWFF30 Acid Green Metallic
TWFF17 Midnight Blue Pearl		TWFF18 White Pearl	TWFF19 Black Lava
TWFF27 Titanium Grey		TWFF45 Pomegranate Red	TWFF48 Blue Mist
TWFF43 Indigo Purple		TWFF40 Hyper Blue	TWFF41 Flame Orange
TWFF42 OD Green		TWFF46 Silver Streak	TWFF49 Lavendar Dusk
TWFF44 Black Galaxy		TWFF47 Seafoam Green	TWFF50 Daffodil Yellow

## Color Anodize Package

<b>Standard:</b>	ZAN01 Black						Standard
<b>Anodized Package:</b>	ZAN12 Gold	ZAN15 Green	ZAN13 Purple	ZAN14 Orange	ZAN10 Blue	ZAN11 Red	\$ 99.00

*Includes caster mount, caster fork and caster caps. If no anodize colour is selected, default will be black.*

## Frame Measurements

Measurements shown below are our most popular selections. Not seeing what you want? Contact customer service or view our Creative Mobility Division (CMD) catalog [here](#) to view additional options. Please refer to Reference Diagrams at end of form for more information.

**Seat Width 'A'** | Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post. Seat width grows 2". All upholstery and components (exception cushion) needed to grow the seat width are included. TIShaft back release bar available at no charge upon request.

**TWSW5:** 8"-10" 9"-11" 10"-12" 11"-13" 12"-14" 13"-15"

**TWSW5:** Set at

Default will be smallest seat width for selected range. Please specify for other settings (1" increments).

**Seat Depth 'B'** | Measure from front of back post to the front edge of seat sling. Seat depth grows 3". Front edge of seat sling will be 1.5" from the beginning of the bend on no-taper and tapered frames, and 2" on abduction frames. All upholstery and components needed to grow the seat depth are included (except cushion).

**TWSD4:** 8"-11" 9"-12" 10"-13" 11"-14" 12"-15" 13"-16"

**TWSD4:** Set at

8"-11" range requires 18" or 20" rear wheels. Default will be smallest seat depth for selected range. Please specify for other settings (1" increments).

**Front Seat Height 'C'** | Measure from floor to top of seat tube at front edge of seat sling. Twist adjusts up and down in 1/2" increments. 2" taper footrest width 'I' available with 12"-15" or 13"-16" seat depths only. For 2" taper footrest width 'I' add 1.5" to minimum front seat height.

**TWFSH1:** 13" 13.5" 14" 14.5" 15" 15.5" 16" 16.5" 17" 17.5" 18" 18.5"  
19" 19.5" 20"

<b>Front Wheel</b>	3"	13"-14.5"	13.5"-15"	13"-16"	14"-17"	14.5"-17.5"	15.5"-18.5"
	3" Glide	14.5"-15.5"	15"-16.5"	14"-17"	15"-18"	16"-19"	16.5"-19.5"
	4"	13.5"-15"	14"-15.5"	13.5"-16.5"	14"-17.5"	15"-18"	16"-18.5"
	4" Glide	14.5"-16"	15.5"-17"	14.5"-17.5"	15.5"-18.5"	16"-19"	17"-20"
	5"			14"-17"	14.5"-17.5"	15.5"-18.5"	16.5"-19.5"
		8"-11" Depth	9"-12" Depth	10"-13" Depth	11"-14" Depth	12"-15" Depth	13"-16" Depth
<b>Front Seat Height Range</b>							

**Rear Seat Height 'D'** | Measure from floor to top of seat tube at back post. Twist adjusts up and down in 1/2" increments.

**TWRSH1:** 12" 12.5" 13" 13.5" 14" 14.5" 15" 15.5" 16" 16.5" 17" 17.5"  
18" 18.5" 19" 19.5"

<b>Rear Wheel</b>	18"	12"-17"
	20"	13"-17.5"
	22"	14"-18.5"
	24"	15"-19.5"

**Rear Seat Height Range**

**Front Angle 'L'** | Measure from frame front to floor.

**TWFA1:** 70° 80° 90°

Standard

## Footrest

### Standard Footrests

**Seat to Footrest 'E':** 9" 9.5" 10" 10.5" 11" 11.5" 12" 12.5" 13" 13.5" 14" 14.5"  
15" 15.5" 16" 16.5"

Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than Front Seat Height 'C'.

**TWFTR3** **Angle Adjustable Footrest**

Standard

**TWFTR24** **Flip-Back Depth and Angle Adjustable Footrest** Not available when both 13"-15" seat width and abduction frame are ordered.

\$ 236.00

### High Mount Footrests

**Seat to Footrest 'E':** 5" 5.5" 6" 6.5" 7" 7.5" 8" 8.5" 9" 9.5" 10" 10.5"  
11" 11.5" 12" 12.5" 13"

Measure from front edge of seat sling to top rear of footrest.

**TWFTR4** **Angle Adjustable High Mount Footrest**

\$ 160.00

**TWFTR27** **Flip-Back Depth and Angle Adjustable High Mount Footrest** Not available when both 13"-15" seat width and abduction frame are ordered.

\$ 340.00

Seat Depth 'B'	3" Front Wheel		4" Front Wheel		5" Front Wheel	
	Standard	High-Mount	Standard	High-Mount	Standard	High-Mount
8"-11"	9"-12"	5"-8"	9"-12"	5"-8"		
9"-12"	10"-13"	6"-9"	10"-13"	6"-9"		
10"-13"	9"-13.5"	5"-9.5"	9"-13.5"	5"-9.5"	9"-14"	5"-10"
11"-14"	10"-14.5"	6"-10.5"	10"-14.5"	7"-11.5"	10"-15"	6"-11"
12"-15"	11"-15"	7"-11.5"	11"-15"	8"-12.5"	11"-15.5"	7"-12"
13"-16"	12"-16"	8"-12.5"	12"-16"	7.5"-12"	12"-16.5"	8"-13"
	Footrest Height Range					

**Footrest Width 'I'** | Measure from inside of front frame tube to inside of opposite front frame tube.

- TWFSW1 **No Taper** 2" narrower than seat width. Standard
- TWFSW4 **2" Taper** 4" narrower than seat width. Available with 12"-15" or 13"-16" seat depths only. Add 1.5" to minimum seat to footrest 'E' ranges. No charge

Seat Width 'A'	No Taper		2" Taper	
	8"	6"		
9"	7"			
10"	8"			
11"	9"			
12"	10"	8"		
13"	11"	9"		
14"	12"	10"		
15"	13"	11"		
	Footrest Width			

## Rear Wheel Positioning

**Rear Wheel Spacing 'J'** | Measure from outside of back post to inside of rear tire. 3/4" rear wheel spacing not available with 0° or 2° camber. Not all rear wheel spacing is available with all accessories.

- TWRWS1 **Fixed Rear Wheel Spacing** Standard  
0.75" 1" 1.25" 1.5" 1.75" 2"
- TWCBR12 **Adjustable Rear Wheel Spacing Camber Plugs** No charge
- TWCBR9 **Additional Spacers to Increase Rear Wheel Spacing** No charge

## Camber 'K'

- TWCBR15 **Aluminum Camber Tube** Standard
- TWCBR1: 0° 2° 4° 6° 8° 12°

## Backrest

### Back Type

- TWSB1 **Folding Aluminum Adjustable Height Back** With TiShaft aluminum back release bar. Standard
- TWSB11 **Backrest Locks When Folded** Standard
- TWSB22 **Height Adjustable Backrest Rigidizer Bar** \$ 170.00

### Push Handles

- TWSB6 **Integrated Push Handles** No charge
- TWSB7 **Integrated Fold-Down Push Handles** Requires tension adjustable by straps back upholstery. \$ 126.00
- TWSB8 **Bolt-On Push Handles** Requires 12"-14" or 13"-15" seat width range. \$ 219.00
- TWSB9 **Bolt-On Height-Adjustable Push Handles** Requires 12"-14" or 13"-15" seat width range. \$ 449.00
- TWSB20 **Bolt-On Center Mount Height Adjustable Push Handle** \$ 203.00  
Small (9.25") Standard (14.25") Tall (19.25")
- TWSB21 **Stroller Handle** \$ 279.00

**Back Height 'F'** | Set at 0.5" increments.

Push Handles	None / Bolt-On	Folding Back Height			Specify
		Short	Medium	Tall	
	Integrated	8.5"-11.5"	12"-16.5"	16.5"-20.5"	Set at:
	Fold-Down	10.5"-12.5"	13.5"-16.5"	17.5"-20.5"	Set at:
		13"-15.5"	14.5"-19"	18.5"-21"	Set at:
		Short Medium Tall			Specify
		Folding Back Height			

**Back Angle 'G'** | Measure angle from front of back post to floor. 90° will be perpendicular to the floor. To approximate the angle from the back post to the seat, subtract 2° from measurement G for every 0.5" of seat dump.

TWSBA1:                    80°      83°      86°      89°      92°      95°      98°      101°

## Back Upholstery

TWBUP2	<b>Velcro Adjustable Back Upholstery</b> <i>Not available with fold-down push handles.</i>	Standard
TWBUP3	<b>Tension Adjustable Back Upholstery By Straps, Padded Nylon</b> <i>Required for fold-down push handles.</i>	\$ 454.00
TWBUP5	<b>Tension Adjustable Back Upholstery By Straps, AIR Breathable</b> <i>Required for fold-down push handles.</i>	\$ 437.00
TWUPC3	<b>Ultrasuede Backrest Stripe, Black</b> <i>Not available with AIR breathable upholstery.</i>	\$ 116.00
TWBUP0	<b>Omit Back Upholstery</b>	No charge

## Front Wheels

See front seat height ranges chart for limitations.

TFWW24	<b>3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire</b> Silver Hub    Black Hub	\$ 197.00
TFWW14	<b>4"x0.75" Plastic Wheel w/Poly Tire</b> <i>Not available with Slipstream fork.</i>	No charge
TFWW3	<b>4"x0.75" Light-Up Micro Caster</b> <i>Not available with Slipstream fork.</i>	\$ 88.00
TFWW26	<b>4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire</b>	\$ 104.00
TFWW27	<b>4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire</b> Silver Hub    Black Hub	\$ 208.00
TFWW6	<b>5"x1" Plastic Wheel w/Poly Tire</b> <i>Not available with 90° front angle or Slipstream fork.</i>	No charge
TFWW28	<b>5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire</b> <i>Not available with 90° front angle.</i>	\$ 104.00
TFWW30	<b>5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire</b> <i>Not available with 90° front angle.</i> Silver Hub    Black Hub	\$ 274.00

## Precision Lock Front Forks

TWFK1	<b>TiLite Standard Fork</b>	No charge
TWFK2	<b>TiLite Slipstream Single-Sided Fork</b>	\$ 203.00
TWGL1	<b>Glide Suspension Fork</b> <i>Black anodize only. Select rider weight.</i> 125lbs    125-180lbs    180-265lbs	\$ 596.00
TWTS1	<b>Titanium Fork Stem</b>	\$ 66.00

## Rear Wheels

See rear seat height ranges chart for limitations.

TWRW14	<b>TiLite SHADOW Rear Wheel</b> <i>Select size.</i> 18" (400)    20" (451)    22" (501)    24" (540)	Standard
TWRW7	<b>Spinergy SPOX Rear Wheel</b> <i>Select size and spoke color. Orange, pink, purple and green spokes only available in 20" size.</i> 20" (451)    22" (501)    24" (540) Black    Red    Blue    White    Orange    Pink    Purple    Green	\$ 1,427
TWRW8	<b>Spinergy LX Rear Wheel</b> <i>Select size and spoke color.</i> 22" (501)    24" (540) Black    Red    Blue    White    Orange    Pink    Purple    Green	\$ 1,547
TWRW9	<b>Spinergy Design Lab LX Rear Wheel</b> <i>Black spokes only. Select size and anodizing color. If same color is selected, anodized hubs and rims match colour anodized package.</i> 22" (501)    24" (540) Gold    Green    Blue    Red    Purple    Orange	\$ 2,382
TWRW16	<b>Spinergy Carbon Blade Rear Wheel</b> <i>Select size. 2x2 twill weave carbon fiber.</i> 24" (540)	\$ 2,003
TWRW17	<b>X-Core Mag Rear Wheel</b> <i>Select size. Not available with short tabs, quad release axles, spoke guards or omit handrims options. Requires aluminum silver anodized handrims.</i> 20" (451)    22" (501)    24" (540)	\$ 269

## Axles

TWAXL1	<b>Stainless Quick Release Axle</b>	Standard
TWAXL2	<b>Titanium Quick Release Axle</b>	\$ 318.00
TWAXL3	<b>Quad Release Axle</b> <i>Not available with Mag rear wheels.</i>	\$ 187.00

## Pneumatic Rear Tires

TWRTR1	<b>Treaded Tire</b>	Standard
TWRTR2	<b>Airless Inserts</b> <i>Not available with 18" rear wheels.</i>	\$ 99.00
TWRTR4	<b>MBL SpeedLite Tire 20"(451), 22"(501), 24"(540)</b>	No charge
TWRTR30	<b>Primo Sentinel Tire 24"(540)</b>	\$ 175.00
TWRTR8	<b>Schwalbe Marathon Plus Evolution Tire 24"(540)</b>	\$ 202
TWRTR32	<b>Primo C-764 Knobby Tire 24"(540)</b>	\$ 106.00

## Solid Rear Tires

TWRTR22	<b>Primo Orion Tire 18" (400), 20"(451), 22"(501), 24"(540)</b> <i>Not available with Spinergy Carbon Blade rear wheels.</i>	\$ 121.00
TWRTR23	<b>Primo Active Tire 24"(540)</b> <i>Not available with Spinergy Carbon Blade rear wheels.</i>	\$ 148.00

## Handrims

TWRIM1	<b>Aluminium Handrim, Silver Anodized</b>	Standard
TWRIM21	<b>Aluminium Handrim, Black Anodized</b> <i>Not available with 18" rear wheels.</i>	No charge
TWRIM12	<b>Titanium Handrim, Satin Finish</b> <i>Not available with 18", 20" or 22" rear wheels, or 24" knobby tires.</i>	\$ 536.00
TWRIM20	<b>Q-Grip Handrim</b> <i>Exclusive high-friction, anti-slip coating. Not available with 18" or 20" rear wheels.</i>	\$ 410.00
TWRIM5	<b>Plastic Coated Handrim, Black</b> <i>Not available with 18" rear wheels.</i>	\$ 175.00
TWRIM10	<b>The Surge Handrim</b> <i>Large profile. All-In-One oval with Gription™ strip. Not available with 18" or 20" rear wheels.</i>	\$ 476.00
TWRIM11	<b>The Surge LT Handrim</b> <i>Small profile. All-In-One oval with Gription™ strip. Not available with 18" or 20" rear wheels.</i>	\$ 514.00
TWRIM14	<b>The Natural-Fit Handrim, Standard Grip</b> <i>Original, large profile. Not available with 18" or 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.</i>	\$ 454.00
TWRIM15	<b>The Natural-Fit Handrim, Super Grip</b> <i>Original, large profile. Not available with 18" or 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.</i>	\$ 514.00
TWRIM16	<b>The Natural-Fit Handrim, Without Thumb Grip</b> <i>Original, large profile.</i>	\$ 415.00
TWRIM17	<b>The Natural-Fit LT Handrim, Standard Grip</b> <i>Reduced weight, small profile. Not available with 18" or 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.</i>	\$ 546.00
TWRIM18	<b>The Natural-Fit LT Handrim, Super Grip</b> <i>Reduced weight, small profile. Not available with 18" or 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.</i>	\$ 585.00
TWRIM19	<b>The Natural-Fit LT Handrim, Without Thumb Grip</b> <i>Reduced weight, small profile.</i>	\$ 503.00
TWRIM22	<b>Omit Handrims</b>	No charge

## Handrim Tab Length

TWRIM7	<b>Long Tabs</b>	Standard
TWRIM8	<b>Short Tabs</b> <i>Not available with Mag rear wheels, The Natural-Fit handrims with thumb grip or Natural-Fit handrims on Shadow wheels.</i>	\$ 39.00

## Wheel Locks

TWWLK7	<b>Compact Composite Push to Lock</b> <i>Available with 0.75"-1.5" rear wheel spacing only.</i>	No charge
TWWLK8	<b>Compact Composite Pull to Lock</b> <i>Available with 0.75"-1.5" rear wheel spacing only.</i>	No charge
TWWLK10	<b>Extension Handles for Compact Composite Push/Pull to Lock</b>	\$ 94.00
TWWLK11	<b>MBL Composite Scissor Lock</b>	\$ 225.00
TWWLK12	<b>MBL Performance Aluminum Scissor Lock</b>	\$ 245
TWWLK3	<b>Uni-Lock</b>	\$ 181.00
	Push Pull	
TWWLK4	<b>Extension Handles for Uni-Lock</b>	\$ 94.00

## Seat Upholstery

TWSUP1	<b>Tension Adjustable Bolt-On Upholstery, Black</b>	Standard
TWSUP0	<b>Omit Seat Upholstery</b>	No charge

## Seat Cushions

TWCSH2	<b>Comfort Curve Cushion</b> <i>Stretch-Air cover.</i>	\$ 150.00
TWCSH23	<b>Comfort Embrace Zero Elevation Cushion</b> <i>Comfort-Tek cover.</i>	\$ 522.00
TWCSH11	<b>ROHO Low Profile Single Compartment Cushion</b>	\$ 711.00

## Side Guards

TWSDG13	<b>Side Guards Custom Cut to Profile of Rear Wheels</b> <i>Available on aluminum and carbon fiber non-fendered side guards.</i>	\$ 345.00
CMD007	<b>Custom Offset Profiled Side Guards (Removable or Fold-Down)</b> <i>Please specify desired offset to nearest 1/2" and side guards type in the Notes section of this order form. Side guard material must be selected. Some limitations may apply.</i>	CMD

### Aluminum

TWSDG2	<b>Aluminum Rigid Side Guards Removable</b>	\$ 318.00
	Extra Small    Small    Standard	
TWSDG3	<b>Aluminum Side Guards, Quad Removable w/Hole</b>	\$ 345.00
	Small    Standard	

### Carbon Fiber

TWSDG7	<b>Carbon Rigid Side Guards, Removable</b>	\$ 372.00
	Small    Standard	
TWSDG8	<b>Carbon Side Guards, Quad Removable w/Hole</b>	\$ 410.00
	Small    Standard	
TWSDG12	<b>Carbon Side Guards, Fendered Removable</b>	\$ 650.00

## Armrests

TWARM1	<b>Swing Away Tubular Armrests</b>	\$ 284.00
TWARM4	<b>Flip Back Tubular Armrests</b> <i>Not available with short seat back height or velcro adjustable back upholstery.</i>	\$ 345.00
TWARM2	<b>Removable Height Adjustable Desk Armrests with Rigid Side Guard</b> <i>10" armpad. Requires 10" minimum seat depth.</i>	\$ 394.00
TWARM5	<b>Armrests Covered in Black Ultrasuede</b> <i>Not available with desk arms.</i>	\$ 126.00

## Positioning

### Upper Body Positioning

TWUBP1	<b>Bodypoint PivotFit Dynamic Shoulder Harness</b> <i>Cam buckle mounts, rear pull.</i>	\$ 296.00
	Extra Small    Small	

### Pelvic Positioning

TWBLT1	<b>Velcro Adjustable Belt, Black Nylon Webbing</b>	\$ 61.00
TWBLT2	<b>Bodypoint 2" Belt, Push-Button Buckle (Auto Style)</b>	\$ 104.00
TWBLT4	<b>Bodypoint 1" Padded Belt, Push-Button Buckle (Auto Style)</b>	\$ 187.00
TWBLT5	<b>Bodypoint 1.5" Padded Belt, Push-Button Buckle (Auto Style)</b>	\$ 187.00
TWBLT8	<b>Q'Straint WC-19 Wheelchair Anchored Pelvic Belt</b> <i>Recommended for use with transit tie down.</i>	\$ 219.00
TWBLT10	<b>Bodypoint Optional Band Clamp Mounting Hardware</b>	\$ 22.00

### Lower Body Positioning

TWCLF15	<b>Bodypoint Ankle Huggers</b>	\$ 318.00
	Extra Small    Small	

## Calf Straps

TWCLF2	<b>Velcro Adjustable</b> Small (9"-12" footrest widths)    Medium (11"-15" footrest widths)	\$ 50.00
TWCLF3	<b>Bodypoint Padded Calf Strap, Velcro Adjustable</b> Medium (11"-15" footrest widths)	\$ 88.00
TWCLF1	<b>Buckle Adjustable Calf Strap, Black Nylon Webbing</b>	\$ 34.00

## Residual Limb Supports

TWRSL1	<b>Residual Limb Support, 8"W x 10"L, Single</b> Left    Right	\$ 563.00
TWRSL2	<b>Residual Limb Support, 8"W x 10"L, Pair</b>	\$ 1,126.00

## Anti-Tips

TWTIP1	<b>Rear Anti-Tips, Aluminum Black</b> <i>Seat widths 8"-10", 9"-11" and 10"-12" will ship with Center Mount Single Non-Swivel Caster Anti-Tip at no additional cost.</i>	\$ 214.00
TWTIP2	<b>Rear Anti-Tips, User Friendly Flip-Up</b> <i>Requires 12"-14" or 13"-15" seat widths.</i>	\$ 394.00
TWTIP5	<b>Center Mount Swivel Caster Anti-Tip</b>	\$ 279.00
TWTIP6	<b>Center Mount Non-Swivel Caster Anti-Tip</b>	\$ 279.00

## Accessories

TWGRD1	<b>Clear Lexan Spoke Guards</b> <i>Not available with Mag rear wheels.</i>	\$ 148.00
TWBCK3	<b>Permobil Medical Necessities Backpack</b>	\$ 104.00
TWPCH3	<b>Permobil Medical Necessities Seat Pouch</b>	\$ 66.00
TWUEST1	<b>Upper Extremity Support Tray</b> <i>Requires armrests.</i> 12"-14" seat widths    13"-15" seat widths	\$ 318.00
TWMPCT2	<b>Neoprene Impact Guards with Ultrasuede Stripe, Black</b>	\$ 116.00
TWCH1	<b>Crutch Holder</b>	\$ 160.00

## General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

## Notes

# TiFit Reference Diagrams

