

Order form

Custom ACTA-BACK® Deep

HCPCS code: E2617 / Custom back

Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process.

Customer information

Date: _____

Company name/ACCT #: _____

Client name: _____

Client weight: _____

Client height: _____

Client reference: _____

W/C manufacturer: _____


Model #: _____ Serial #: _____

Seat Frame Width: _____

Seat Frame Depth: _____

Seat Rail Tube Size: _____

Seating system measurements



NOTE: Seating system measurements do not affect any selections made throughout the order form.

A - Back support width	
B - Back support length	
C - Width between lateral trunk supports	
D - Width between lateral thigh supports	
E - Seat width	
F - Seat depth	
G - Effective seat depth	
H - Lateral trunk support height	

P.O. number: _____

Requested by: _____

Phone: _____ Fax: _____

Email: _____

Ship to: _____

Additional information:


1 Custom ACTA-BACK Deep size & cover material

Instructions: select the desired back support size by clicking on the corresponding box below or filling in the specific width and length measurements in the "custom back part number" box.

1.1 Back support size

HCPCS code: E2617

LATERAL WIDTH ADJUSTMENT (BASED ON 1 1/2" THICK FOAM SELECTION)



*Length (L) refers to the actual size dimension of the support from bottom to top edge.

Length*	22"	C-D-12W22H	C-D-14W22H	C-D-16W22H	C-D-18W22H	C-D-20W22H	C-D-22W22H	C-D-24W22H	C-D-26W22H	C-D-28W22H
	20"	C-D-12W20H	C-D-14W20H	C-D-16W20H	C-D-18W20H	C-D-20W20H	C-D-22W20H	C-D-24W20H	C-D-26W20H	C-D-28W20H
	18"	C-D-12W18H	C-D-14W18H	C-D-16W18H	C-D-18W18H	C-D-20W18H	C-D-22W18H	C-D-24W18H	C-D-26W18H	C-D-28W18H
	16"	C-D-12W16H	C-D-14W16H	C-D-16W16H	C-D-18W16H	C-D-20W16H	C-D-22W16H	C-D-24W16H	C-D-26W16H	C-D-28W16H
	14"	C-D-12W14H	C-D-14W14H	C-D-16W14H	C-D-18W14H	C-D-20W14H	C-D-22W14H	C-D-24W14H	C-D-26W14H	C-D-28W14H
	12"	C-D-12W12H	C-D-14W12H	C-D-16W12H	C-D-18W12H	C-D-20W12H	C-D-22W12H	C-D-24W12H	C-D-26W12H	C-D-28W12H
		12"	14"	16"	18"	20"	22"	24"	26"	28"
		\$1348 ⁰⁰					\$1481 ⁰⁰		\$2186 ⁰⁰	

Frame width

Custom back support part number:

C - D - _ W _ H

NOTE: Desired back support length is "H" dimension in part number above.


12-20" Width / 12-18" Length = \$1348.00 base price
22-24" Width / 12-18" Length = \$1481.00 base price
12-24" Width / 20-22" Length = \$1521.00 base price
26-28" Width / 12-22" Length = \$2186.00 base price
*If desired size is not listed above, a 60% up charge will be applied for materials.

If an odd width is needed, go 1" wider and choose the next even width.

1.2 Cover material options

REQUIRED

MSRP \$0.00




COVER FABRIC	PART NUMBER
STRETCH-AIR®	(CD-STRETCH-AIR)
COMFORT-TEK®	(CD-COMFORT-TEK)

Cover appearance may vary with hardware selection.

1.3 Military branch patch

OPTIONAL

MSRP \$0.00



2 Foam options

- Select only one type of foam per layer
- A minimum of one layer must be chosen if Foam in Place is not selected in Section 2.2.



Indicates standard ACTA-BACK foam configuration.
May use as a reference when choosing foam layers.




2.1 Custom foam selection

Rear layer
REQUIRED
H.R. MSRP \$0.00
Visco MSRP \$166.00

Middle layer
OPTIONAL
H.R. MSRP \$0.00
Visco MSRP \$166.00

Topper layer
OPTIONAL
H.R. MSRP \$0.00
Visco MSRP \$166.00

NOTE: 1/4" Crosslink foundation layer will come standard on all back supports. For other foam combinations, call Customer Support 800.736.0925.



MSRP \$0.00	Foam type	1/4"	1/2"	3/4"	1"	1 1/4"	1/4"	1/2"	3/4"	1"	1 1/4"	1/4"	1/2"	3/4"	1"	1 1/4"
	H.R. Medium	(CD-Rear-HRM-0.25)	(CD-Rear-HRM-0.50)	(CD-Rear-HRM-0.75)	(CD-Rear-HRM-1.00)	(CD-Rear-HRM-1.25)	(CD-Mid-HRM-0.25)	(CD-Mid-HRM-0.50)	(CD-Mid-HRM-0.75)	(CD-Mid-HRM-1.00)	(CD-Mid-HRM-1.25)	(CD-Top-HRM-0.25)	(CD-Top-HRM-0.50)	(CD-Top-HRM-0.75)	(CD-Top-HRM-1.00)	(CD-Top-HRM-1.25)
MSRP \$166.00	H.R. Soft	(CD-Rear-HRS-0.25)	(CD-Rear-HRS-0.50)	(CD-Rear-HRS-0.75)	(CD-Rear-HRS-1.00)	(CD-Rear-HRS-1.25)	(CD-Mid-HRS-0.25)	(CD-Mid-HRS-0.50)	(CD-Mid-HRS-0.75)	(CD-Mid-HRS-1.00)	(CD-Mid-HRS-1.25)	(CD-Top-HRS-0.25)	(CD-Top-HRS-0.50)	(CD-Top-HRS-0.75)	(CD-Top-HRS-1.00)	(CD-Top-HRS-1.25)
	Visco Medium		(CD-Rear-VIS-M-0.50)		(CD-Rear-VIS-M-1.00)			(CD-Mid-VIS-M-0.50)		(CD-Mid-VIS-M-1.00)			(CD-Top-VIS-M-0.50)		(CD-Top-VIS-M-1.00)	
	Visco Soft		(CD-Rear-VIS-S-0.50)		(CD-Rear-VIS-S-1.00)			(CD-Mid-VIS-S-0.50)		(CD-Mid-VIS-S-1.00)			(CD-Top-VIS-S-0.50)		(CD-Top-VIS-S-1.00)	
	Visco X-Soft		(CD-Rear-VIS-XS-0.50)		(CD-Rear-VIS-XS-1.00)			(CD-Mid-VIS-XS-0.50)		(CD-Mid-VIS-XS-1.00)			(CD-Top-VIS-XS-0.50)		(CD-Top-VIS-XS-1.00)	
	Visco XX-Soft		(CD-Rear-VIS-XXS-0.50)		(CD-Rear-VIS-XXS-1.00)			(CD-Mid-VIS-XXS-0.50)		(CD-Mid-VIS-XXS-1.00)			(CD-Top-VIS-XXS-0.50)		(CD-Top-VIS-XXS-1.00)	


3 Back support options

- A standard alignment has been optimized based on back support width if you choose not to select a specific width between the stays.
- Select the spacing between the moldable aluminum stays.
- Custom options are available by selecting your desired spacing below.

3.1 Moldable aluminum stays

NOTE:

- This back support comes with two patented aluminum stays that can be molded to support spinal curvature.
- Stays are offset from the spine and are held in place by pockets sewn into the rear of the cover.
- The stays are accessible while the user is seated against the back support through the zipper located at the top of the back support.



Space between stays

MSRP \$0.00


A = 1" (CD-STAY-1.0)
A = 1 1/2" (CD-STAY-1.5)
A = 2" (standard alignment) (CD-STAY-2.0)
A = 2 1/2" (CD-STAY-2.5)
A = 3" (CD-STAY-3.0)
A = 3 1/2" (CD-STAY-3.5)
A = 4" (CD-STAY-4.0)
A = 4 1/2" (CD-STAY-4.5)
Specify custom spacing _____ (CD-STAY-OTHER)
None (CD-STAY-NONE)

3.2 Hardware options

4 Point hardware is not available on 12"W or narrower back supports and may not fit all smaller back support lengths.

NOTE: Quick release and fixed hardware are adjustable for height, depth, width and angle of recline. Fixed requires tools for removal. 4-Point is designed to provide extra support for users with tone and for back supports 24" wide or greater. Compass 4 hardware may be installed on 3/4", 7/8", 1" and 1 1/8" tubing. Power mount hardware is height adjustable up to 8" in 1/2" increments but may be limited by the chair.


Compass® 4 Quick release (2 points of attachment) MSRP \$0.00



Quick release (2 point) (CD-Q-COMPASS)


Compass bundle
(CD-COMPASS-BUNDLE)

Compass 4 Quick release (4 points of attachment) MSRP \$0.00




Quick release (4 point) (CD-Q-COMPASS4)

Compass 4 Fixed (2 points of attachment) MSRP \$0.00



Fixed (2 point) (CD-NQ-COMPASS)

Compass 4 Fixed (4 points of attachment) MSRP \$0.00



Fixed (4 point) (CD-NQ-COMPASS4)

*Power mount hardware required for all power chairs listed below.

Compass Power Mount Hardware MSRP \$618.00

(Not available on back supports smaller than 12"L or 12" W, or sizes 12Wx16L & 12Wx18L)

*Using Compass Power Mount on Rovi requires mounting to the RECESSED PLANAR INTERFACE PLATE available when ordering the Motion Concepts Seating System. Drilling is required for mounting.

Quantum w/ Tru Balance 3
(CD-PMB-QUANTUM)

Permobil w/ Corpus 3G
(CD-PMB-PERMOBIL)

Quickie w/ 3.7 Power
(CD-PMB-QUICKIE)

Avid Rehab w/ Power T & R
(CD-PMB-AVID)

Rovi w/ Motion Concepts*
(CD-PMB-ROVI)

3.3 Color *OPTIONAL* MSRP \$0.00

NOTE: Accent colors can be applied via colored piping that is stitched to the flaps of the back support. If no colors are selected, our standard silver reflective piping will be used.



Color location

Placement of the piping may vary by back support size and mounting hardware.

Red (CD-RED)	Blue (CD-BLUE)
Purple (CD-PURPLE)	Green (CD-GREEN)
Pink (CD-PINK)	Remove (CD-REMOVE)
Orange (CD-ORANGE)	

3.4 Embroidery *OPTIONAL* MSRP \$0.00

Fill out sections A, B, and C.

A.

Thread color options

Red (CD-EMB-RED)	Blue (CD-EMB-BLUE)
Purple (CD-EMB-PURPLE)	Green (CD-EMB-GREEN)
Pink (CD-EMB-PINK)	Silver (CD-EMB-SILVER)

B.

Font/style options

Helvetica (CD-EMB-Helvetica)	Twilight (CD-EMB-Twilight)
Garamond (CD-EMB-Garamond)	Script (CD-EMB-Script)

C.

Word one

(Input is case sensitive. Please enter exactly as you want the words to look. Example: Marcus or marcus or MARCUS)


- AND/OR -

Word two

(Input is case sensitive. Please enter exactly as you want the words to look. Example: Marcus or marcus or MARCUS)

3.5 Privacy flap *OPTIONAL* MSRP \$44.00

NOTE: This flap connects with hook & loop attachment between the seat and back support for added privacy.



Privacy flap

Privacy flap (CD-PF-CV)

4 Thoracic lateral support

4.1 Single lateral depth MSRP \$0.00


Select depth of lateral desired from the options below.

*Length (L) refers to the actual size dimension of the support from bottom to top edge.

NOTE: Laterals are mounted at a fixed height position. Laterals are width and angle adjustable.

Lateral pad height changes with different depth selections.


If different wing height location is needed, please make comments in section 6.1.



L	Lateral depth	R
(CD-DEEP-LAT-L-OMIT)	Omit lateral	(CD-DEEP-LAT-R-OMIT)
(CD-DEEP-LAT-L-3)	3"D = (5"L)	(CD-DEEP-LAT-R-3)
(CD-DEEP-LAT-L-5)	5"D = (6"L)	(CD-DEEP-LAT-R-5)
(CD-DEEP-LAT-L-7)	7"D = (7"L)	(CD-DEEP-LAT-R-7)

Typical wing height mounting chart

If different wing height location is needed, please make comments in section 6.1.



5 BodiLink Head Support hardware & pads

Length (L) refers to the actual size dimension of the support from bottom to top edge.

5.1 Premium pad

HCPCS code: E0955

MSRP \$269.00

W		Pad Firmness	
L	6W x 3L	Soft	BL-HSP1-6W3L
		Extra soft	BL-HSP2-6W3L
L	8W x 4L	Soft	BL-HSP1-8W4L
		Extra soft	BL-HSP2-8W4L
L	10W x 5L with adjustable wings	Soft	BL-HSP1-10W5L
		Extra soft	BL-HSP2-10W5L
L	12W x 5L with adjustable wings	Soft	BL-HSP1-12W5L
		Extra soft	BL-HSP2-12W5L
L	14W x 5L with adjustable wings	Soft	BL-HSP1-14W5L
		Extra soft	BL-HSP2-14W5L
L	18W x 5L with adjustable wings	Soft	BL-HSP1-18W5L
		Extra soft	BL-HSP2-18W5L
W			



Premium head support pads angles can be formed by adjusting the plates in the pad.

5.2 Premium pad Cover material

GlideWear

For skin protection & shear reduction. Not fluid resistant.

COMFORT-TEK



For fluid protection & an easily cleaned surface.

STRETCH-AIR



For client comfort.

	GlideWear™ (GLIDEWEAR) (add \$30.00)
	COMFORT-TEK (COMFORT-TEK)
	STRETCH-AIR (STRETCH-AIR)

5.3 PT Hardware

Range of motion

PT MODEL	TOTAL # OF LINKS	TOTAL LINK LENGTH	MAX ANTERIOR REACH	MAX POSTERIOR REACH	MAX LATERAL REACH
MID 6	2	6"	5.5"	1"	N/A
LAT 6	3	6"	5.5"	1"	6"
MID 8	2	8"	7.5"	1"	N/A
LAT 8	3	8"	7.5"	1"	8"
MAX 8	4	8"	7.5"	3"	8"

A. PT Hardware, Removable


HCPCS E1028



Shape	MSRP	Part number
MID 6	\$235.00	(BL-HS-PT1RMMID6)
LAT 6	\$273.00	(BL-HS-PT1RMLAT6)
MID 8	\$235.00	(BL-HS-PT1RMMID8)
LAT 8	\$273.00	(BL-HS-PT1RMLAT8)
MAX 8*	\$344.00	(BL-HS-PT1RMMAX8)

*Provides the greatest amount of articulation and range of motion.

PT Hardware, Removable, Stem length option




Stem length	MSRP	Part number
Standard 7" stem	\$0.00	(HS-STEM-7)
**3" stem	\$0.00	(HS-STEM-3)

**3"stem is recommended for use with power mounting where there may be special constraints with power chair functions.

B. PT Hardware, Flip-down
HCPCS E1028

Flip-Down Hardware has 7" Stem



Shape	MSRP	Part number
MID 6	\$287.00	(BL-HS-PT1FDMID6)
LAT 6	\$326.00	(BL-HS-PT1FDLAT6)
MID 8	\$287.00	(BL-HS-PT1FDMID8)
LAT 8	\$326.00	(BL-HS-PT1FDLAT8)
MAX 8*	\$375.00	(BL-HS-PT1FDMAX8)

*Provides the greatest amount of articulation and range of motion.

5.4 Head support mounting option



Mounting block	MSRP
Posterior A mounting block	\$0.00

Part number
(HS-MB-POSTA)

5.5 Extra hardware

	Link length	MSRP	Extra part	
	•Extra 2" Link	\$65.00	(BL-HS-PT1LINK2)	




*Extra 3" Link	\$65.00	(BL-HS-PT1LINK3)	
----------------	---------	------------------	--

	Link length	MSRP	Extra part	*Extra 4" Link	\$65.00	(BL-HS-PT1LINK4)
Link length	MSRP	Extra part				
*Extra 4" Link	\$65.00	(BL-HS-PT1LINK4)				



Tool-less wing lock	MSRP	Extra part		Extra part		Extra part		Extra part	
1 Lock	\$24.00/each	(BL-HS-52WINGA1)		(BL-HS-52WINGA1)		(BL-HS-52WINGA1)		(BL-HS-52WINGA1)	



Additional	MSRP	Extra part	
1/4" Spacer plate	\$43.00/each	(BL-HS-SPCRA)	

*For PT hardware only

The following are trademarks and registered trademarks of Permobil: ACTA-BACK®, BodiLink®, COMFORT-TEK®, Compass®, and STRETCH-AIR®. GlideWear™ is a trademark of MIPS.

6 Additional information

6.1 Comments

INSTRUCTIONS: Please reference any section or subsection of the form when commenting.

Section #	Page #	Comments

INSTRUCTIONS: If preferred, please sketch ideas in the spaces provided below.

Submit by email

.....