

# FreeWheel

Mobility. Freedom.  
Independence.



permobil

## What is the FreeWheel?

The FreeWheel is a lightweight, durable third-wheel that instantly turns your wheelchair into an all-terrain chair.

It securely attaches to your solid footrest and safely lifts the front casters off the ground.

Whether it's participating in a race, at the park watching your child play sports, or camping and hitting the trails, the FreeWheel is the one piece of adaptive equipment you don't want to be without.



## Over curbs, on trails, through grass and sand.



You no longer need to spend time looking at the ground, worrying about every bump, stick or rock, or wasting energy wheeling over grass.

With the FreeWheel you can look at the world around you and experience life to the fullest.



## Folding Wheelchair? No problem.

The ADAPTOR for Folding Wheelchairs securely clamps on to the horizontal side tubing of your chair. Once set up for the individual chair, it comes on and off in seconds, making it easy to use.

## Stories from around the world.



Super easy to put on, the FreeWheel transforms my chair into an all terrain chair in seconds so I can go where I want to. The FreeWheel is awesome!

Mal, Australia

My goal is to use my FreeWheel to explore all the wonders of the world!

Catherine, Canada



Rough terrain, mounds and bumps, Freewheel is now the reason I live life on the edge.

Paul, United Kingdom



## Innovating for Individuals

Permobil Australia

8 Viewbrook Close  
Seven Hills NSW 2147  
Australia

Phone | 1300 845 483 or 02 8838 3100  
Email | [info.au@permobil.com](mailto:info.au@permobil.com)  
[www.permobil.com/en-au](http://www.permobil.com/en-au)

**permobil**