



PDG Fuze T50 - AADL W653



Dealer Name*: _____ Account #*: _____

City, Province, Post Code: _____

Phone: _____ Email: _____ ATP: _____

Rep Name: _____ P.O. Number: _____ Quote Number: _____

Client Reference*: _____ Client Weight: _____

Special considerations that need to be addressed here (e.g., diagnosis): _____ *Required field

October 2024, Revised January 2026. Prices are in Canadian Dollars (CAD)

Wheelchair & Frame

Fuze T50 50° Manual Tilt Wheelchair Also available: the narrow Fuze T50n with an overall chair width of only seat width + 4". See separate order form. W653 \$ 2,685.00

Frame Type

- Fixed Width Frame, 300lbs Weight Capacity** Select a fixed seat width below. No charge
- Fixed Width Frame, 400lbs Weight Capacity** Not compatible with transit tie down system. Requires flat free rear tires or heavy duty spoke rear wheels. Select a fixed seat width below. Prior Approval \$ 300.00
- Adjustable Width Frame, 250lbs Weight Capacity** Select an adjustable seat width below. \$ 200.00

Frame Options

Transit Tie Down System (TTDS) Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details. Tested to maximum weight capacity of 300lbs. No charge

Seat Measurements

Seat Width | Overall chair width is +7.5" with 12/16" rear wheels and +8.5" with 20/22/24" rear wheels.

| | | | |
|----------------------------|---------------------------------|------|-----------|
| Fixed: | 15" 16" 17" 18" 19" | | No charge |
| | 20" 21" 22" | W111 | \$ 250.00 |
| | 23" 24" | W111 | \$ 500.00 |
| Adjustable 15"-20": | 15" 16" 17" 18" 19" | | No charge |
| | 20" | W111 | \$ 250.00 |

Seat Depth

| | | | |
|----------------------|---------------------------------|------|-----------|
| Select Depth: | 15" 16" 17" 18" 19" | | No charge |
| | 20" 21" 22" | W111 | \$ 250.00 |

21"-22" seat depths are not adjustable and require straight stroller handle back. Not compatible with reclining or fold-down back.

Frame Finishes

Select Color: High Gloss Black Hyper Blue

Armrests

Armrest Type

No Armrests

No charge

Single Post, Height Adjustable Armrests

No charge

10-14" Range 12"-16" Range

Armrest Position: Rear Front

If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.

Armpads

Regular Armpads

No charge

Desk Length 2" x 10" Full Length 2" x 14"

Seat Height

Seat to floor height is adjustable within the range of the chosen axle plate (Low/Standard/High). Chairs with swing-away front rigging are set with 3° of pre-tilt to avoid interference. This increases the actual seat to floor height at the front of the seat by up to 1" compared to the selected value. Easy-Out and standard plug-in front rigging are not affected.

RW

| | 5" Casters | | | | | 6" Casters & 6" Forks | | | | | 6" Casters & 8" Forks | | | | | 8" Casters & 8" Forks | | | | |
|-----|----------------------|-----|-----|------|------|-----------------------|-----|------|------|------|-----------------------|------|------|------|------|-----------------------|------|------|------|------|
| 24" | | | Low | Low | Low | | Low | Low | Low | Std | Low | Low | Low | Std | Std | Low | Low | Std | Std | Std |
| 22" | | Low | Low | Low | Std | Low | Low | Low | Std | Std | Low | Low | Std | Std | Std | Low | Std | Std | Std | |
| 20" | Low | Low | Low | Std | Std | Low | Low | Std | Std | Std | Low | Std | Std | Std | | Std | Std | Std | | High |
| 16" | | Std | Std | Std | | Std | Std | Std | | High | Std | Std | | High | High | Std | | High | High | High |
| 12" | Std | Std | | High | High | Std | | High | High | High | | High | High | High | | High | High | High | High | High |
| | 13" | 14" | 15" | 16" | 17" | 14" | 15" | 16" | 17" | 18" | 15" | 16" | 17" | 18" | 19" | 16" | 17" | 18" | 19" | 20" |
| | Seat to Floor Height | | | | | | | | | | | | | | | | | | | |

Seat Pan

Seat Pan

No Seat Pan

No charge

Solid Aluminium Seat Pan

No charge

Front Casters

Urethane Casters

5" x 1.25" Caster

No charge

6" x 1.25" Caster

No charge

8" x 1.25" Caster

No charge

Rear Wheels

Composite Rear Wheel

No charge

12" 16" 20" 22" 24"

Heavy Duty Spoke Wheel

\$ 350.00

24"

Rear Wheel Position: Rear (More Stable) Mid (Default) Front (Less Stable)

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

Axles

Quick Release Rear Axles

No charge

Fixed Rear Axles

No charge

Rear Tires

Urethane Tires

No charge

Handrims | Handrims are only available with 20", 22" and 24" rear wheels.

| | |
|---|-----------|
| No Handrims | No charge |
| Standard Aluminium Handrims | No charge |
| Plastic Coated Handrims <i>Not compatible with drum brakes.</i> | \$ 175.00 |

Wheel Locks

| | |
|------------------------------|-----------|
| User Wheel Locks | No charge |
| Push to Lock Pull to Lock | |

Wheel Lock Options

| | |
|--|-----------|
| Wheel Lock Extensions <i>Not compatible with attendant wheel lock.</i> | No charge |
|--|-----------|

Anti-Tippers

| | |
|--|-----------|
| Rear Anti-Tippers <i>Required safety option.</i> | No charge |
|--|-----------|

Backrest

| | |
|---------------------|--|
| Back Post Style: | Stroller Handle Back |
| | <i>Cane handle back not compatible with power tilt with manual override.</i> |
| Back Post Type: | Straight 8° Bend |
| Back Post Height: | Short (17") Tall (23") |
| | <i>Measured from seat pan to handle.</i> |
| Folding Style: | Fold-Down |
| | <i>Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.</i> |
| Seat to Back Angle: | Set at 90° Set at 98° Set at 106° Set at 114° Set at 122° |

Backrest Options

| | |
|---|-----------|
| Angle/Height Adjustable Stroller Handle: <i>Not compatible with cane handle back.</i> | No charge |
| Button on Outside Button on Inside | |

Tilt Activation

Attendant Tilt Only

| | |
|--------------------------------------|-----------|
| Attendant Tilt Only, Right Back Post | No charge |
| Attendant Tilt Only, Left Back Post | No charge |

Tilt Limit | *Optional. If no tilt limit is selected the full tilt range will be available.*

| | |
|-------------|--|
| Tilt Limit: | Set at 20° Set at 30° Set at 40° |
|-------------|--|

Front Rigging

| | |
|------------------|-----------|
| No Front Rigging | No charge |
|------------------|-----------|

Easy-Out Front Rigging | *Effortless insertion and removal. No hanger brackets remain on the frame.*

| | |
|-----------|----------|
| Fixed 70° | \$ 95.00 |
|-----------|----------|

Swing Away Front Rigging | *Legacy design. Not compatible with Easy-Out.*

| | |
|---|----------------|
| Swing Away 70° | No charge |
| Swing Away 90° <i>Includes 11"-16" extension tubes. Not compatible with 8" casters.</i> | \$ 0.00 |
| Swing Away Elevating Legrests (ELR) <i>Includes calf pads.</i> | W106 \$ 375.00 |
| 1" Offset Hangers for Swing Away Front Rigging | \$ 95.00 |

Plug In Front Rigging | *Legacy design. Not compatible with Easy-Out.*

| | |
|-------------|-----------|
| Plug In 70° | No charge |
|-------------|-----------|

Footplates

- Composite Footplates \$ 0.00
- Aluminium Angle Adjustable Footplates W107 \$ 165.00

Extension Tubes | Consider the height of your seat cushion.

- Standard Extension Tubes, Short No charge
- Standard Extension Tubes, Long \$ 75.00

Footplate Height:

Specify:

| | | | | | |
|----------------|-------|-------------------|--------------------------|----------------|------------------|
| Ext | Short | 11" - 16" | 12" - 16" | 4" - 9" | 13" - 16" |
| | Long | 14" - 22" | 14" - 19" | 12" - 16" | 14" - 17" |
| | | Standard Easy-Out | 60°, 70°, ELR Swing-Away | 90° Swing-Away | Standard Plug-In |
| Legrest | | | | | |

Front Rigging Accessories

- Length Adjustable Heel Loops No charge
- Standard Calf Strap Basic Seating \$ 50.00

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes