

permobil

ROHO SEATING + POSITIONING




ROHO® AGILITY® BACKS THE NEXT GENERATION

AGILITY BACKS ARE EVEN BETTER NOW WITH:

- New, more comfortable, air and foam combinations
- Larger shell cutouts to create a whole new look
- Hardware that is smaller and easier to install
- An air, foam and cover ventilation system to improve breathability

NOW WITH A FIVE-YEAR WARRANTY



Wherever you're headed, the next generation of ROHO AGILITY is ready. Our full line of wheelchair backs will support your every move.

PERFORMANCE THAT LASTS

Selecting a backrest is an important decision. You need a back that performs optimally - time after time. We know you depend on us and we are here to deliver exciting new enhancements and a [market-changing, five-year warranty](#).*

Rest assured, the new ROHO AGILITY Back Systems will deliver. With the same benefits you have come to expect, AGILITY Backs still offer different DRY FLOATATION® and AIR FLOATATION™ inserts to provide superior comfort, protection and postural support.

AGILITY MEDICARE CODES

- E2611 Foam only
- E2613 Foam with 1" lumbar support pad
- E2615 Foam with optional air insert and 1" lumbar support pad
- E2620 Foam with optional air insert
- E2617 Custom AGILITY — custom fabricated wheelchair back cushion, any size, including any type mounting hardware

SPECIFICATIONS

AGILITY MAX CONTOUR

Depth	7"
Width	14", 15", 16", 17", 18", 19", 20"
Height	16", 18", 20"
Hardware	Fixed, Quick-Release, Direct Mount (Permobil only)

AGILITY MID CONTOUR

Depth	4"
Width	14", 15", 16", 17", 18", 19", 20"
Height	10", 13", 16", 18", 20"
Hardware	Fixed, Quick-Release, Direct Mount (Permobil only)

AGILITY MINIMUM CONTOUR

Depth	3" (2.5" with foam only configuration)
Width	14", 15", 16", 17", 18", 19", 20"
Height	10", 13", 16", 18"
Hardware	Fixed, Quick-Release, Direct Mount (Permobil only)

AGILITY CUSTOM BACKS

Completely customizable DRY FLOATATION insert.

Custom cell heights	1.25"
	2.25"
	3.25"
	4.25"

We also offer a full range of accessories, including covers, laterals and headrests, to complement your ROHO AGILITY Back System.

*The market-changing, five-year limited warranty applies to hardware and shell replacement parts only. Products may be covered by one or more U.S. and foreign patents and trademarks. The following are trademarks and registered trademarks of Permobil: ROHO®, DRY FLOATATION®, AGILITY®, and AIR FLOATATION™. ©2017 Permobil | Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. For coverage criteria in other countries, consult your local ROHO distributor.