

Titanium high performance wheelchairs

TiLite TR | TiLite ZR



permobil

TiLite TR

The best from every angle.

The sweeping curves of the titanium dual-tube frame combine striking design, enhanced performance, and optimal functionality. As our flagship rigid chair, the TR's efficiency and elegance move you to the pinnacle of performance.



TiLite ZR

Head turning style.

The ZR's titanium mono-tube frame means less truly is more. This minimalist aesthetic turns heads while the intentional frame design provides you the durability and vibration-dampening ride performance you expect from TiLite.




Tattoo options

- 1 Mosaic
- 2 Flower garden
- 3 Need for speed




Frame colours

High gloss paint




Metallic paint




Pearlescent paint




Matte paint




Candy paint



Matte textured paint




Natural Titanium



Component colours

Anodised



Specifications

TiLite TR

Frame features	Titanium Dual-Tube Rigid, TiFit, Limited adjustability
Total chair width*	22"
Seat width	12" – 20"
Front seat width (seat taper)	0", 1", 2"
Seat depth	12" – 20"
Front seat height	15" – 21"
Rear seat height	12.5" – 20"
Front frame angle	70°, 75°, 80°, 85°, or 90°
Camber	0°, 2°, 4°, 6° or 8°

Centre of gravity adjustment	5"
Folding seat back	Adj. height 9.5" – 22" Adj. angle 80° – 101°
Fixed seat back	Fixed height 5" – 20" Fixed angle 80° – 101°
Ergonomic seat	5" – 9"
Transit options	-
Transport weight*	4.2 kg
Weight limit	120 kg
Warranty	Five years on frame One year on TiLite components

TiLite ZR

Frame features	Titanium Mono-Tube Rigid, TiFit, Limited adjustability
Total chair width*	23"
Seat width	10" – 20"
Front seat width (seat taper)	0", 1.5"
Seat depth	10" – 20"
Front seat height	16" – 21"
Rear seat height	12.5" – 20"
Front frame angle	70°, 75°, 80°, 85°, or 90°
Camber	0°, 2°, 4°, 6° or 8°

Centre of gravity adjustment	2"
Folding seat back	Adj. height 9.5" – 22" Adj. angle 80° – 101°
Fixed seat back	Fixed height 8" – 20" Fixed angle 80° – 101°
Ergonomic seat	5" – 9"
Transit options	Transport bracket
Transport weight*	4.2 kg
Weight limit	120 kg Standard 136 kg Heavy duty
Warranty	Five years on frame One year on TiLite components

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.3 kg.