

Order form

Custom ACTA-Relief

HCPCS code: E2617 / Custom back

Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process.

Customer information

Date: _____

Company name/ACCT #: _____

Client name: _____

Client weight: _____

Client height: _____

Client reference: _____

W/C manufacturer: _____


Model #: _____ Serial #: _____

Seat Frame Width: _____

Seat Frame Depth: _____

Seat Rail Tube Size: _____

Seating system measurements



For best results, do not fill out in your browser. Interactive form should be completed using Adobe Reader after saving to your local drive. Then email or print and fax to Customer Support at orders.comfort@permobil.com or 877.522.8563.

DO NOT SEND PROTECTED HEALTH INFORMATION. IT IS NOT NEEDED TO MAKE THE PRODUCT YOU ARE REQUESTING.

P.O. number: _____

Requested by: _____

Phone: _____ Fax: _____

Email: _____

Ship to: _____

Additional information:

NOTE: Seating system measurements do not affect any selections made throughout the order form.

A - Back support width	
B - Back support length	
C - Width between lateral trunk supports	
D - Width between lateral thigh supports	
E - Seat width	
F - Seat depth	
G - Effective seat depth	
H - Lateral trunk support height	

ACTA-RELIEF-CD-ORDER REV111422

1 Back support size & cover material

INSTRUCTIONS: Select the desired back support size by clicking on the corresponding box below or filling in the specific width and length measurements in the "Custom back part number" box. ACTA-RELIEF IS LIMITED TO ONLY THE SIZES LISTED.

1.1 Back support size

HCPCS code: E2617



NOTE:

12"L Acta-Relief back supports are not available with 4 PT hardware when using lateral trunk supports.

10"L Acta-Relief back supports are not available with 4 PT or Power Mount hardware, lateral trunk supports or a head support.

Appearance and Boa® closure system may vary with size and hardware selection.

Sizes

*Length (L) refers to the actual size dimension of the support from bottom to top edge.

Length*	22"	C-R-14W22H	C-R-16W22H	C-R-18W22H	C-R-20W22H
	20"	C-R-14W20H	C-R-16W20H	C-R-18W20H	C-R-20W20H
	18"	C-R-14W18H	C-R-16W18H	C-R-18W18H	C-R-20W18H
	16"	C-R-14W16H	C-R-16W16H	C-R-18W16H	C-R-20W16H
	14"	C-R-14W14H	C-R-16W14H	C-R-18W14H	C-R-20W14H
	12"	C-R-14W12H	C-R-16W12H	C-R-18W12H	C-R-20W12H
	10"	C-R-14W10H	C-R-16W10H	C-R-18W10H	C-R-20W10H
		14"	16"	18"	20"

\$2,133.00

Frame width

Custom back support part number:

C - R - _ W _ H

NOTE: Desired back support length is "H" dimension in part number above
Acta-Relief is limited to only the sizes listed in chart above.

1.2 Cover material options

REQUIRED

MSRP \$0.00

STRETCH-AIR®

Stretch-Air has comfortable stretch, air transmission and heat dissipation. It is not fluid resistant.

- OR -

COMFORT-TEK®

Comfort-Tek is a smooth, easily cleaned surface. The fluid protection aids in infection control.


COVER FABRIC	PART NUMBER	
STRETCH-AIR®	(CD-STRETCH-AIR)	
COMFORT-TEK®	(CD-COMFORT-TEK)	

Cover appearance may vary with hardware selection.

1.3 Military branch patch

OPTIONAL

MSRP \$0.00



	Air Force (CD-AIR-FORCE-PATCH)
	Army (CD-ARMY-PATCH)
	Marines (CD-MARINES-PATCH)
	Navy (CD-NAVY-PATCH)

2 Foam options

- Select only one type of foam per layer
 - A minimum of one layer must be chosen if Foam in Place is not selected in Section 2.2.
- Indicates standard ACTA-BACK foam configuration.
May use as a reference when choosing foam layers.

2.1 Custom foam selection

NOTE: 1/4" Crosslink foundation layer will come standard on all back supports. For other foam combinations, call Customer Support 800.736.0925.

Rear layer
REQUIRED
H.R. MSRP \$0.00
Visco MSRP \$166.00

Middle layer
OPTIONAL
H.R. MSRP \$0.00
Visco MSRP \$166.00

Topper layer
OPTIONAL
H.R. MSRP \$0.00
Visco MSRP \$166.00

Foam type	1/4"	1/2"	3/4"	1"	1 1/4"
H.R. Medium	(CD-Rear-HRM-0.25)	(CD-Rear-HRM-0.50)	(CD-Rear-HRM-0.75)	(CD-Rear-HRM-1.00)	(CD-Rear-HRM-1.25)
H.R. Soft	(CD-Rear-HRS-0.25)	(CD-Rear-HRS-0.50)	(CD-Rear-HRS-0.75)	(CD-Rear-HRS-1.00)	(CD-Rear-HRS-1.25)
Visco Medium		(CD-Rear-VIS-M-0.50)		(CD-Rear-VIS-M-1.00)	
Visco Soft		(CD-Rear-VIS-S-0.50)		(CD-Rear-VIS-S-1.00)	
Visco X-Soft		(CD-Rear-VIS-XS-0.50)		(CD-Rear-VIS-XS-1.00)	
Visco XX-Soft		(CD-Rear-VIS-XXS-0.50)		(CD-Rear-VIS-XXS-1.00)	

Foam type	1/4"	1/2"	3/4"	1"	1 1/4"
(CD-Mid-HRM-0.25)	(CD-Mid-HRM-0.25)	(CD-Mid-HRM-0.50)	(CD-Mid-HRM-0.75)	(CD-Mid-HRM-1.00)	(CD-Mid-HRM-1.25)
(CD-Mid-HRS-0.25)	(CD-Mid-HRS-0.25)	(CD-Mid-HRS-0.50)	(CD-Mid-HRS-0.75)	(CD-Mid-HRS-1.00)	(CD-Mid-HRS-1.25)
		(CD-Mid-VIS-M-0.50)		(CD-Mid-VIS-M-1.00)	
		(CD-Mid-VIS-S-0.50)		(CD-Mid-VIS-S-1.00)	
		(CD-Mid-VIS-XS-0.50)		(CD-Mid-VIS-XS-1.00)	
		(CD-Mid-VIS-XXS-0.50)		(CD-Mid-VIS-XXS-1.00)	

Foam type	1/4"	1/2"	3/4"	1"	1 1/4"
(CD-Top-HRM-0.25)	(CD-Top-HRM-0.25)	(CD-Top-HRM-0.50)	(CD-Top-HRM-0.75)	(CD-Top-HRM-1.00)	(CD-Top-HRM-1.25)
(CD-Top-HRS-0.25)	(CD-Top-HRS-0.25)	(CD-Top-HRS-0.50)	(CD-Top-HRS-0.75)	(CD-Top-HRS-1.00)	(CD-Top-HRS-1.25)
		(CD-Top-VIS-M-0.50)		(CD-Top-VIS-M-1.00)	
		(CD-Top-VIS-S-0.50)		(CD-Top-VIS-S-1.00)	
		(CD-Top-VIS-XS-0.50)		(CD-Top-VIS-XS-1.00)	
		(CD-Top-VIS-XXS-0.50)		(CD-Top-VIS-XXS-1.00)	


- OR -

2.2 Foam in place kit

- Select extra soft, soft or medium density memory foam bag pour.
- Select the appropriate pour size based on your choice of a 3" or 4" deep back support border.
- Contact Customer Support with any questions 800.736.0925.

KIT INCLUDES: (Note: Price of Visco layers will be added to the order separately.)

- Back cushion cover with an optional 3" or 4" deep back, with a stretchy, pleated contact surface to allow for foam growth.
- Two aluminum stays sewn into the rear of the cover which can be molded to support spinal curvature prior to pouring FIP Chemicals.
- 1" of Visco Medium foam is cut to fit the shape of the back support and is placed loosely inside the cover. This will give the user some padding, prior to pouring the foam.
- 1/2" of Visco Soft foam is cut to fit the shape of the back support and is placed inside the front of the cover. This will give the final mold a plush feel for the user and will also smooth out any potential inconsistencies after pouring the foam.
- Two plastic bags are included to contain the foam chemicals as they expand.



NOTE: BAG POUR SIZE RECOMMENDATION CHART

Border thickness	Frame width	Recommended bag size
3"	10"W-16"W	1 bag
3"	17"W-22"W	1 1/2 bag
4"	10"W-14"W	1 bag
4"	15"W-20"W	1 1/2 bag

Extra soft memory foam

1 Bag FIP MSRP: \$371.00
1 1/2 Bag FIP MSRP: \$494.00

3" BACK BORDER	POUR SIZE	4" BACK BORDER
(CD-FIP3-1-XS)	1 BAG	(CD-FIP4-1-XS)
(CD-FIP3-1.5-XS)	1 1/2 BAG	(CD-FIP4-1.5-XS)

Soft memory foam

1 Bag FIP MSRP: \$371.00
1 1/2 Bag FIP MSRP: \$494.00

3" BACK BORDER	POUR SIZE	4" BACK BORDER
(CD-FIP3-1-S)	1 BAG	(CD-FIP4-1-S)
(CD-FIP3-1.5-S)	1 1/2 BAG	(CD-FIP4-1.5-S)

Medium memory foam

1 Bag FIP MSRP: \$371.00
1 1/2 Bag FIP MSRP: \$494.00

3" BACK BORDER	POUR SIZE	4" BACK BORDER
(CD-FIP3-1-M)	1 BAG	(CD-FIP4-1-M)
(CD-FIP3-1.5-M)	1 1/2 BAG	(CD-FIP4-1.5-M)

3 Back support options

NOTE: Quick release and Fixed hardware are adjustable for height, depth, width and angle of recline. Fixed requires tools for removal. 4-Point is designed to provide extra support for users with tone. Compass 4 hardware may be installed on 3/4", 7/8", 1" and 1 1/8" tubing. Power mount bracket is height adjustable up to 8" in 1/2" increments but may be limited by the chair.

3.1 Hardware options

4 point hardware may not fit all smaller back support lengths.

Compass® 4 Quick release (2 points of attachment) MSRP \$0.00	Compass 4 Quick release (4 points of attachment) MSRP \$0.00	Compass 4 Fixed (2 points of attachment) MSRP \$0.00	Compass 4 Fixed (4 points of attachment) MSRP \$0.00	Compass 4 Bundle (2 PT, 4 PT, Fixed & Quick release) MSRP \$577.00
Quick Release (2 point) (CD-Q-COMPASS)	Quick Release (4 point) (CD-Q-COMPASS4)	Fixed (2 point) (CD-NQ-COMPASS)	Fixed (4 point) (CD-NQ-COMPASS4)	Bundle (CD-COMPASS-BUNDLE)

Compass Power Mount Hardware


MSRP \$618.00

Not available on 10"L back supports.

Head support bracket not required for Acta-Relief when using Compass Power Mount hardware.

Quantum w/ Tru Balance 3 (CD-PMB-QUANTUM)
Permobil w/ Corpus 3G (CD-PMB-PERMOBIL)
Quickie w/ 3.7 Power (CD-PMB-QUICKIE)
Avid Rehab w/ Power T & R (CD-PMB-AVID)
Rovi w/ Motion Concepts* (CD-PMB-ROVI)


*Using Compass Power Mount on Rovi requires mounting to the RECESSED PLANAR INTERFACE PLATE available when ordering the Motion Concepts Seating System. Drilling is required for mounting.



*Power mount bracket required for all power chairs listed below.

3.2 Color *OPTIONAL* MSRP \$0.00

NOTE: Accent colors can be applied via colored piping that is stitched to the flaps of the back support. If no colors are selected, our standard silver reflective piping will be used.



Placement of the piping may vary by back support size and mounting hardware.

Red (CD-RED)	Blue (CD-BLUE)
Purple (CD-PURPLE)	Green (CD-GREEN)
Pink (CD-PINK)	Remove (CD-REMOVE)
Orange (CD-ORANGE)	

3.3 Embroidery *OPTIONAL* MSRP \$0.00

Fill out sections A, B, and C.

A. Thread color options

Red (CD-EMB-RED)	Blue (CD-EMB-BLUE)
Purple (CD-EMB-PURPLE)	Green (CD-EMB-GREEN)
Pink (CD-EMB-PINK)	Silver (CD-EMB-SILVER)

B. Font/style options

Helvetica (CD-EMB-Helvetica)	Twilight (CD-EMB-Twilight)
Garamond (CD-EMB-Garamond)	Script (CD-EMB-Script)

C. Word one

(Input is case sensitive. Please enter exactly as you want the words to look. Example: Marcus or marcus or MARCUS)


- AND/OR -

Word two

(Input is case sensitive. Please enter exactly as you want the words to look. Example: Marcus or marcus or MARCUS)

3.4 Privacy flap *OPTIONAL* MSRP \$44.00

NOTE: This flap connects with hook & loop attachment between the seat and back support for added privacy.



Privacy flap (CD-PF-CV)

4 BodiLink® Lateral Trunk Support pads & hardware

4.1 Premium pad size & shape

HCPCS code: E0956 MSRP \$101.00/ea.

• You may make one selection per column in section 4.1.

• Follow each column down to sections 4.2 and 4.3 and select a style and cover for each pad selected in section 4.1.


*Length (L) refers to the dimension of the support from top to bottom edge. Depth (D) refers to the dimension of the support from anterior to posterior edge.

	Shape	Size	Left	*Extra left	Right	*Extra right
Supportive		3"L x 4"D	BL-LTSP2S-3L4D-LH	BL-LTSP2S-3L4D-LH	BL-LTSP2S-3L4D-RH	BL-LTSP2S-3L4D-RH
		3"L x 5.5"D	BL-LTSP2S-3L5D-LH	BL-LTSP2S-3L5D-LH	BL-LTSP2S-3L5D-RH	BL-LTSP2S-3L5D-RH
		4"L x 5"D	BL-LTSP2S-4L5D-LH	BL-LTSP2S-4L5D-LH	BL-LTSP2S-4L5D-RH	BL-LTSP2S-4L5D-RH
		4"L x 6.5"D	BL-LTSP2S-4L6D-LH	BL-LTSP2S-4L6D-LH	BL-LTSP2S-4L6D-RH	BL-LTSP2S-4L6D-RH
		5"L x 6.5"D	BL-LTSP2S-5L6D-LH	BL-LTSP2S-5L6D-LH	BL-LTSP2S-5L6D-RH	BL-LTSP2S-5L6D-RH
		6"L x 6.5"D	BL-LTSP2S-6L6D-LH	BL-LTSP2S-6L6D-LH	BL-LTSP2S-6L6D-RH	BL-LTSP2S-6L6D-RH
Active		3"L x 4"D	BL-LTSP2A-3L4D-LH	BL-LTSP2A-3L4D-LH	BL-LTSP2A-3L4D-RH	BL-LTSP2A-3L4D-RH
		3"L x 5.5"D	BL-LTSP2A-3L5D-LH	BL-LTSP2A-3L5D-LH	BL-LTSP2A-3L5D-RH	BL-LTSP2A-3L5D-RH
		4"L x 5"D	BL-LTSP2A-4L5D-LH	BL-LTSP2A-4L5D-LH	BL-LTSP2A-4L5D-RH	BL-LTSP2A-4L5D-RH
		4"L x 6.5"D	BL-LTSP2A-4L6D-LH	BL-LTSP2A-4L6D-LH	BL-LTSP2A-4L6D-RH	BL-LTSP2A-4L6D-RH
		5"L x 6.5"D	BL-LTSP2A-5L6D-LH	BL-LTSP2A-5L6D-LH	BL-LTSP2A-5L6D-RH	BL-LTSP2A-5L6D-RH
		6"L x 6.5"D	BL-LTSP2A-6L6D-LH	BL-LTSP2A-6L6D-LH	BL-LTSP2A-6L6D-RH	BL-LTSP2A-6L6D-RH


**"EXTRA" parts may not fit depending on back support, pad, and hardware size and style combinations. Please contact customer support with questions.

4.2 Premium pad style

MSRP \$0.00/ea.




Or



Style options	Left	*Extra left	Right	*Extra right
Moldable	LTS-P2-MOLD	LTS-P2-MOLD	LTS-P2-MOLD	LTS-P2-MOLD
Rigid	LTS-P2-RIGID	LTS-P2-RIGID	LTS-P2-RIGID	LTS-P2-RIGID

4.3 Cover


Cover options	MSRP	Left	*Extra left	Right	*Extra right
No cover	\$0.00/ea.	LTS-P2-NO-CV	LTS-P2-NO-CV	LTS-P2-NO-CV	LTS-P2-NO-CV
COMFORT-TEK	\$0.00/ea.	LTS-P2-COMFORT	LTS-P2-COMFORT	LTS-P2-COMFORT	LTS-P2-COMFORT
STRETCH-AIR	\$0.00/ea.	LTS-P2-STRETCH	LTS-P2-STRETCH	LTS-P2-STRETCH	LTS-P2-STRETCH
GlideWear™	\$15.00/ea.	LTS-P2-GLIDE	LTS-P2-GLIDE	LTS-P2-GLIDE	LTS-P2-GLIDE
COMFORT-TEK with foam	\$15.00/ea.	LTS-P2-COMFORT-FM	LTS-P2-COMFORT-FM	LTS-P2-COMFORT-FM	LTS-P2-COMFORT-FM
STRETCH-AIR with foam	\$15.00/ea.	LTS-P2-STRETCH-FM	LTS-P2-STRETCH-FM	LTS-P2-STRETCH-FM	LTS-P2-STRETCH-FM
GlideWear with foam	\$31.00/ea.	LTS-P2-GLIDE-FM	LTS-P2-GLIDE-FM	LTS-P2-GLIDE-FM	LTS-P2-GLIDE-FM



4.4 GT Mounting hardware

- You may make one selection per column in section 4.4.
- Extra links and long mounting bracket (section 4.5) only for GT2 style hardware.

Both fixed and swing-away hardware styles are adjustable.



Slot mounting	MSRP	Left	*Extra left	Right	*Extra right
Swing-away HCPCS Code: E1028	\$339.00/ea.	BL-LTS-GT2SAK-LH	BL-LTS-GT2SAK-LH	BL-LTS-GT2SAK-RH	BL-LTS-GT2SAK-RH
Fixed	\$202.00/ea.	BL-LTS-GT2FXK-LH	BL-LTS-GT2FXK-LH	BL-LTS-GT2FXK-RH	BL-LTS-GT2FXK-RH

4.5 GT2 Style additional options



Extra hardware links	Left	Additional left	Right	Additional right
MSRP \$43.00/link	BL-LTS-GT2LINK	BL-LTS-GT2LINK	BL-LTS-GT2LINK	BL-LTS-GT2LINK
Long mounting bracket	Left	*Extra left	Right	*Extra right
For slot mount MSRP \$43.00/ea	BL-LTS-GT2LSLOT-LH	BL-LTS-GT2LSLOT-LH	BL-LTS-GT2LSLOT-RH	BL-LTS-GT2LSLOT-RH
Hardware cover	Left	*Extra left	Right	*Extra right
MSRP \$31.00/ea	BL-LTS-HCOVER	BL-LTS-HCOVER	BL-LTS-HCOVER	BL-LTS-HCOVER

5 BodiLink Head Support hardware & pads

Length (L) refers to the actual size dimension of the support from bottom to top edge.

5.1 Premium pad size

HCPCS code: E0955 MSRP \$269.00

W	Pad firmness		
6W x 3L	Soft	BL-HSP1-6W3L	
	Extra soft	BL-HSP2-6W3L	
8W x 4L	Soft	BL-HSP1-8W4L	
	Extra soft	BL-HSP2-8W4L	
10W x 5L with adjustable wings	Soft	BL-HSP1-10W5L	
	Extra soft	BL-HSP2-10W5L	
12W x 5L with adjustable wings	Soft	BL-HSP1-12W5L	
	Extra soft	BL-HSP2-12W5L	
14W x 5L with adjustable wings	Soft	BL-HSP1-14W5L	
	Extra soft	BL-HSP2-14W5L	
18W x 5L with adjustable wings	Soft	BL-HSP1-18W5L	
	Extra soft	BL-HSP2-18W5L	



Premium head support pads angles can be formed by adjusting the plates in the pad.

5.2 Premium pad cover material

GlideWear	For skin protection & shear reduction. Not fluid resistant.
COMFORT-TEK	For fluid protection & an easily cleaned surface.
STRETCH-AIR	For client comfort.

GlideWear (GLIDEWEAR) (add \$30.00)
COMFORT-TEK (COMFORT-TEK)
STRETCH-AIR (STRETCH-AIR)

5.3 PT Hardware

- You may make one selection per column.
- Choose PT hardware type in 5.3 A or B.
- If choosing A, select a stem size immediately after.

Range of motion

PT MODEL	TOTAL # OF LINKS	TOTAL LINK LENGTH	MAX ANTERIOR REACH	MAX POSTERIOR REACH	MAX LATERAL REACH
MID 6	2	6"	5.5"	1"	N/A
LAT 6	3	6"	5.5"	1"	6"
MID 8	2	8"	7.5"	1"	N/A
LAT 8	3	8"	7.5"	1"	8"
MAX 8	4	8"	7.5"	3"	8"

A. PT Hardware, Removable

HCPCS E1028



Shape	MSRP	Part number
MID 6	\$254.00	(BL-HS-PT1RMMID6)
LAT 6	\$295.00	(BL-HS-PT1RMLAT6)
MID 8	\$254.00	(BL-HS-PT1RMMID8)
LAT 8	\$295.00	(BL-HS-PT1RMLAT8)
MAX 8*	\$372.00	(BL-HS-PT1RMMAX8)

*Provides the greatest amount of articulation and range of motion.

PT Hardware, Removable, Stem length option




Stem length	MSRP	Part number
Standard 7" stem	\$0.00	(HS-STEM-7)
**3" stem	\$0.00	(HS-STEM-3)

**3" Stem is recommended for use with power mounting where there may be special constraints with power chair functions.

B. PT Hardware, Flip-down

HCPCS E1028


Flip-Down Hardware has 7" Stem



Shape	MSRP	Part number
MID 6	\$310.00	(BL-HS-PT1FDMID6)
LAT 6	\$353.00	(BL-HS-PT1FDLAT6)
MID 8	\$310.00	(BL-HS-PT1FDMID8)
LAT 8	\$353.00	(BL-HS-PT1FDLAT8)
MAX 8*	\$405.00	(BL-HS-PT1FDMAX8)

*Provides the greatest amount of articulation and range of motion.

5.4 Head support mounting option



BRACKET WITH MOUNTING BLOCK-
Acta-Relief
18W Shown Above
For Non-Power
Mount Only



MOUNTING BLOCK-POSTERIOR A
For Acta-Relief with
Power Mount Only

Mounting block or bracket with mounting block	MSRP
•Acta-Relief 14W Bracket with mounting block	\$84.00
•Acta-Relief 16W Bracket with mounting block	\$84.00
•Acta-Relief 18W Bracket with mounting block	\$84.00
•Acta-Relief 20W bracket with mounting block	\$84.00
Posterior A mounting block	\$0.00

Part number
(HS-MB-AR14W)
(HS-MB-AR16W)
(HS-MB-AR18W)
(HS-MB-AR20W)
(HS-MB-POSTA)

*For Acta-Relief, includes mounting bracket and POSTA block.
If you already have an Acta-Relief head support bracket, order POSTA Block.

5.5 Extra hardware

**For PT hardware only



Link length	MSRP
**Extra 2" link	\$65.00

Extra part
(BL-HS-PT1LINK2)



**Extra 3" link	\$65.00
-----------------	---------

(BL-HS-PT1LINK3)



**Extra 4" link	\$65.00
-----------------	---------

(BL-HS-PT1LINK4)




Tool-less wing lock	MSRP
**1 lock	\$25.00/each

Extra part
(BL-HS-52WINGA1)

Extra part
(BL-HS-52WINGA1)

Extra part
(BL-HS-52WINGA1)

Extra part
(BL-HS-52WINGA1)



BRACKET-
ACTA-RELIEF
18W Shown Above

Bracket size	MSRP
*Acta-Relief 14W bracket	\$173.00
*Acta-Relief 16W bracket	\$173.00
*Acta-Relief 18W bracket	\$173.00
*Acta-Relief 20W bracket	\$173.00

Part number
(BL-HS-AR14W)
(BL-HS-AR16W)
(BL-HS-AR18W)
(BL-HS-AR20W)



Additional	MSRP
1/4" Spacer plate	\$43.00/each

EXTRA PART
(BL-HS-SPCRA)

*Bracket only, does not include mounting block.

The following are trademarks and registered trademarks of Permobil: ACTA-BACK®, BodiLink®, COMFORT-TEK®, Compass®, and STRETCH-AIR®. GlideWear™ is a trademark of MIPS.

6 Additional information

6.1 Comments

INSTRUCTIONS: Please reference any section or subsection of the form when commenting.

Section # Page # Comments

Submit by email

INSTRUCTIONS: If preferred, please sketch ideas in the spaces provided below.

.....