

permobil

AERO T

Pure geometry.

Our most advanced aluminium chair, the Aero T's dual-tube frame minimises weight while maximising performance. Add to that full adjustability and TiFit, and you've got something special.



AERO T



TiFit. Personal Ergonomics.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do. TiFit...not standard.

Custom Configuration.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimise your riding experience. This includes choices of wheels, handrims, tyres, wheel locks, colour, tattoos, and more. TiLite offers more options for the individual.

Technical Specifications

Frame Features	Aluminium Dual-Tube Rigid TiFit Full Adjustability	Centre of Gravity Adjustment	5"
Total Chair Width*	23"	Front Seat Width (Seat Taper)	0", 1", 2", 3"
Seat Width	12" – 20"	Transit Options	Transport Bracket
Seat Depth	12" – 20"	Transport Weight*	5.4 kg.
Front Seat Height	15" – 21"	Weight Limit	120.2 kg Standard 158.7 kg Heavy Duty
Rear Seat Height	13" – 21"	Warranty	Five Years on Frame One Year on TiLite Components
Front Frame Angle	70°, 75°, 80°, 85°, or 90°	* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.27 kg.	
Seat Back Type	Folding Seat Back Adj. Height 8.5" – 21" Adj. Angle 80° – 101°		
Camber	0°, 2°, 4°, 6° or 8°		

Colour Options

High Gloss Paint



TiLite Red Burgundy Neon Pink Super Wet Black

Matte Paint



Deep Blue Ice Titanium Grey Black Lava

Candy Paint



Ruby Red Azure Blue Electric Plum

Pearlescent Paint



Midnight Blue Pearl White Pearl

Metallic Paint



Salsa Red Metallic Tangerine Metallic Acid Green Metallic Ocean Blue Metallic

Anodised Colours



Red Pepper Burnt Orange Sapphire Ultraviolet Fuchsia Black