

Power wheelchairs

Wheelchair matrix



permobil

Front wheel

Mid wheel



Rehab Series	F3 Corpus	F5 Corpus	F5 Corpus VS
Max user weight	300 lb	300 lb	300 lb
Max speed	6 mph 5 mph <i>without anti-tips</i>	7.5 mph 6 mph <i>without anti-tips</i>	7.5 mph
Range	16 mi (Grp 34s) 20 mi (Grp 24s)*	20 mi*	20 mi*
Base width	24"	25.5"	25.5"
Base length	40" <i>with anti-tips</i> 36" <i>without anti-tips</i>	42" <i>with anti-tips</i> 38" <i>without anti-tips</i>	43"
Base turning radius	26.5"	30"	30"
Ground clearance	3"	3"	3"
Weight w/ batteries	386 lb (Gr 34), 407 lb (Gr 24)	412 lb (Grp 24s)	432 lb (Grp 24s)
Seat-to-floor height	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"
Seat elevator	Optional 12"	Optional 14"	14"
Power tilt options	0° – 50° <i>Posterior</i> 5°, 10°, 20° & 30° <i>Anterior</i>	0° – 50° <i>Posterior</i> 5°, 10°, 20° & 45° <i>Anterior</i>	0° – 50° <i>Posterior</i> 45° <i>Anterior</i>
Recline options	85° – 120° (manual) 85° – 180° (power)	85° – 120° (manual) 85° – 180° (power)	85° – 180° (power)
Legrest options	85° – 180°	85° – 180°	85° – 180°
Seat depth	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)
Seat width	17", 19", 21", 23"	17", 19", 21", 23"	17", 19", 21", 23"
Backrest height	20", 23" – 28" (by 1" incr.)	20", 23" – 28" (by 1" incr.)	23" – 28" (by 1" incr.)



M1	M3 Corpus	M5 Corpus	M300 HD
300 lb	300 lb	300 lb	450 lb
5 mph	6 mph	7.5 mph	5 mph
16 mi*	16 mi (Grp 34s) 20 mi (Grp 24s)*	20 mi*	13 mi (Grp 34s) 20 mi (Grp 24s)*
24"	24"	25"	26"
35.5"	35.5"	37.5"	38"
20.5"	20.5"	21.5"	22"
3"	3"	3"	3"
319 lb (Gr 34)	385 lb (Grp 34s) 405 lb (Grp 24s)	423 lb (Grp 24s)	425 lb (Gr 34) 447 lb (Gr 24)
16", 17", 18"	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	18", 19", 20" (elevator) 18", 19", 20" (fixed)
—	Optional 12"	Optional 14"	Optional 8"
0° – 50° <i>Posterior</i>	0° – 50° <i>Posterior</i> 5°, 10° & 20° <i>Anterior</i>	0° – 50° <i>Posterior</i> 5°, 10° & 20° <i>Anterior</i>	0° – 45° <i>Posterior</i>
—	85° – 120° (manual) 85° – 180° (power)	85° – 120° (manual) 85° – 180° (power)	85° – 120° (manual) 85° – 150° (power)
85° – 140° 65° – 80° (adjustable swing-away)	85° – 180°	85° – 180°	85° – 175°
14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	19" – 25" (by 1" incr.)
14" – 22" (by 2" incr.)	17", 19", 21", 23"	17", 19", 21", 23"	19", 21", 23", 25"
18" – 29"	20", 23" – 28" (by 1" incr.)	20", 23" – 28" (by 1" incr.)	21", 24" – 29" (by 1" incr.)

*All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories. | *Actual driving range will vary based on driving and battery conditions.*



	Front wheel			Mid wheel	Rear wheel
					
Pediatric Series	Explorer Mini	Koala	K300 PS Jr.	M300 PS Jr.	K450 MX
Intended Age Range	12 months – 3 years	—	—	—	—
Max user weight & height	35 lb <i>Weight</i> 39" <i>Height</i>	80 lb <i>Weight</i>	165 lb <i>Weight</i>	165 lb <i>Weight</i>	125 lb <i>Weight</i>
Max speed	1.5 mph	4 mph	5 mph 6.5 mph	6 mph 7.5 mph	5 mph
Range	3.5 mi*	15 mi*	16 mi*	13 mi (Grp 34s) 20 mi (Grp 24s)*	17 mi*
Base width	19"	23"	24"	24"	25"
Base length	25"	32"	40"	36.5"	35"
Min turning radius	—	21"	25"	20"	28"
Ground clearance	1"	2.5"	3"	3"	2.5"
Weight w/ batteries	52 lb	226 lb (M 40)	272 lb (Gr 34)	313 lb (Gr 34) 335 lb (Gr 24)	337 lb (Gr 34)
Seat-to-floor height	12.1" – 21.7"	15" (elevator)	16.5" (elevator) 16.5", 17.25", 18.5" (fixed)	17" (elevator) 17", 17.75", 18.5" (fixed)	3" – 26" (elevator)
Seat elevator	—	8"	Optional 8"	Optional 8"	8"
Power tilt options	—	0° – 25°	0° – 45° (power or manual)	0° – 45° (power or manual)	0° – 45°
Recline options	—	85° – 120° (manual)	85° – 120° (manual)	85° – 120° (manual)	85° – 105° (manual)
Legrest options	—	90° – 150° (manual)	70° (swing away) 90 – 150° (center mount)	70° or 90° (swing away) 90° – 150° (center mount)	90° – 125° (manual)
Seat depth	—	10", 12", 14"	10" – 18"	10" – 18"	10" – 18" (by 1" incr.)
Seat width	—	10", 12", 14"	11" – 16" (by 1" incr.)	11" – 16" (by 1" incr.)	10", 12", 14", 16"
Backrest height	Total chair height 29" – 37"	14", 16", 20"	16.5" – 22.5" (by 1" incr.)	16.5" – 22.5" (by 1" incr.)	14", 16", 20" (MX) 16.5" – 22.5" (PS Jr)

*All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.
Actual driving range will vary based on driving and battery conditions.

Unmatched power positioning

Corpus is the most ergonomic and capable wheelchair seating system in the world. It is thoughtfully designed and configured to empower functional movements while supporting your optimal body position.



MyPermobil connectivity comes standard on all new M & F-Series Corpus & Corpus VS models. Download the app to get started today.

