

Power wheelchairs

# Koala



permobil



# Good things come in small packages

Designed to grow with a child for years to come, the Permobil Koala Mini-Flex is the ideal pediatric power wheelchair for older toddlers and young children.

Tilt feature offers both improved stability for going down ramps, and a chance for kids to rest after a long day of play

Parents can have extra peace of mind with optional features like the front bumper, lights, and a remote stop system for their little explorers



## Specifications

<b>Maximum user weight</b>	80 lb
<b>Maximum speed</b>	4 mph
<b>Range</b>	15 mi*
<b>Base width</b>	23"
<b>Base length</b>	32"
<b>Minimum turning radius</b>	21"
<b>Ground clearance</b>	2.5"
<b>Recommended battery type</b>	M40
<b>Weight including batteries</b>	226 lb
<b>Drive electronics</b>	R-net

<b>Seat-to-floor height</b>	15" (elevator)
<b>Optional seat elevator</b>	8"
<b>Power tilt</b>	0° – 25°
<b>Manual recline</b>	85° – 120°
<b>Manual legrest elevation</b>	90° – 150°
<b>Seat depth</b>	10", 12", 14"
<b>Seat width</b>	10", 12", 14"
<b>Backrest height</b>	14", 16", 20"

\*Actual driving range will vary based on driving and battery conditions.

Specifications may vary based upon configuration. Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions.

## Color options

