

Fitting / Order form

Prices valid from: 11/5/2020

M1

Midwheel drive chair with independent wheel suspension.
As standard equipped with airfilled tires, 60Ah maintenance free batteries and charger 24V/10A.
The seating system has an adjustable seat depth and comes with a Comfort Curve seat cushion. The armrests are independent height and angle adjustable.



Dealer information

Customer (name / code)	/
Contact	
Address	
Post code / City	/
Phone / Fax	/

	Article	Position	Price	Funding	Additional	☒ Required choice
	code	Left Right	excl VAT incl VAT	code	Info	☐ / ☐ Optional choice
						☐ / ☒ Default required option

Chassis

M1	108006-99-0	☒	4400.00			
----	-------------	-------------------------------------	---------	--	--	--

Drive wheels

WHEELS 3.00x8" AIR FILLED TIRES BLACK M1	107975-99-0	☒	0.00			
WHEELS 3.00x8" SOLID TIRES BLACK M1	107974-99-0	C	200.00			



Suspension setting

SUSPENSION ADJUSTMENT 0-110KG	107766-99-0	☒	0.00			
SUSPENSION ADJUSTMENT 111-150KG	107765-99-0	☐	0.00			
SUSPENSION ADJUSTMENT 126-150KG	107826-99-0	☐	0.00			

Seat height

SEAT HEIGHT STANDARD WITHOUT SEAT TILT	108026-99-0	☒	0.00			
SEAT HEIGHT (STF) +25MM WITHOUT SEAT TILT	108025-99-0	C	0.00			
SEAT HEIGHT (STF) +50MM WITHOUT SEAT TILT	108024-99-0	C	0.00			

Seat sizes

		Depth													
Size (mm)		360	285	410	435	460	485	510	535	560					
Width	360	☐	☐	☐	☐	☐	☐	☐	☐	☐					
	410	☐	☐	☐	☐	☐	☐	☐	☐	☐					
	460	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	510	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐			



Seat cushion

SEAT CUSHION CURVE	10xxxx-99-0	☒	0.00			
WITHOUT SEAT CUSHION	107250-99-0	C	0.00			




Mounting armrests

WIDTH BETWEEN ARMRESTS: NARROW	108283-99-0	☒	0.00			
WIDTH BETWEEN ARMRESTS: WIDE	108282-99-0	C	0.00			
ARMREST HIGH MOUNT 249-367MM	108034-99-0	☒	0.00			
ARMREST LOW MOUNT 178-308MM	108035-99-0	C	0.00			



Armrests

	Left				Right						
	Vinyl FOC	Textile FOC	Gel £122		Vinyl FOC	Textile FOC	Gel £122				
50mm x 200mm	☐	☐	☐		☐	☐	☐		0.00		
50mm x 350mm	☐	☐	☐		☐	☐	☐		0.00		
Corpus 100mm x 260mm	☐	☐	☐		☐	☐	☐		70.00		
Corpus 100mm x 355mm	☐	☐	☐		☐	☐	☐		70.00		
Corpus 100mm x 410mm	☐	☐	☐		☐	☐	☐		70.00		
Corpus 100mm x 460mm	☐	☐	☐		☐	☐	☐		70.00		



Legrest

LEGREST SWING AWAY SHORT	108318-99-0	☒	0.00			
LEGREST SWING AWAY	108158-99-0	C	0.00			



Notes

Accessories seat								
Headrest								
HEADREST ORIGINAL 340x130 LONG HOLDER, M1	108903-99-0	C	495.00					
HEADREST ORIGINAL 340x130 STD HOLDER, M1	108904-99-0	C	445.00					
HEADREST ORIGINAL 250x125 LONG HOLDER, M1	108902-99-0	C	445.00					
HEADREST ORIGINAL 250x125 STD HOLDER, M1	108901-99-0	C	395.00					
HEADREST ORIGINAL 180x110 LONG HOLDER, M1	108900-99-0	C	0.00					
HEADREST ORIGINAL 180x110 STD HOLDER, M1	108899-99-0	C	395.00					
Seatbelts								
SAFETY BELT,BUCKLE	100518-99-0	C	50.00					
SEATBELT, BODYPOINT 2- POINT, LARGE, PADDED	104516-99-0	C	85.00					
SEATBELT, BODYPOINT 2- POINT, MEDIUM, PADDED	104520-99-0	C	85.00					
BP HIP BELT, 2P C PULL LATCH 46CM PAD	108267-99-0	C	85.00					
BP HIP BELT, 2P C PULL LATCH 62CM PAD	108266-99-0	C	85.00					
BP HIP BELT, 4P R PULL PUSHBUTT. 46CM PAD	108265-99-0	C	115.00					
Joystick holder								
PANELHOLDER FIXED HIGH PS ARMREST LEFT	108210-99-0	O	0.00					
PANELHOLDER FIXED HIGH PS ARMREST RIGHT	108211-99-0	O	0.00					
PANELHOLDER SWING AWAY L=335MM	103473-99-0	O	0.00					
PANELHOLDER PARALLEL MOTION UNITRACK RIGHT	108037-99-0	O	150.00					
PANELHOLDER PARALLEL MOTION UNITRACK LEFT	108038-99-0	O	150.00					
Height adjustable holder								
HEIGHT ADJUSTABLE RIGHT JOYSTICK HOLDER	107960-99-0	O	110.00					
HEIGHT ADJUSTABLE LEFT JOYSTICK HOLDER	107963-99-0	O	110.00					
HEIGHT ADJUSTABLE JOYSTICK HOLDER	108005-99-0	O	110.00					
Alternative Joystick knobs								
JOYSTICK KNOB CHIN CUP	100960-99-0	C	29.00					
JOYSTICK KNOB T-GRIP	100961-99-0	C	29.00					
JOYSTICK KNOB BALL BIG	100962-99-0	C	29.00					
JOYSTICK KNOB SOFT BALL	100963-99-0	C	29.00					
JOYSTICK KNOB MUSHROOM	100964-99-0	C	29.00					
JOYSTICK KNOB STICK	100965-99-0	C	29.00					
JOYSTICK KNOB BALL SMALL	101165-99-0	C	29.00					
JOYSTICK KNOB STANDARD	101166-99-0	C	29.00					
JOYSTICK KNOB SOFT BALL BLACK (CONE)	105442-99-0	C	32.00					
JOYSTICK KNOB BODYPOINT 7.5CM U-SHAPED	105846-99-0	C	49.00					
JOYSTICK KNOB BODYPOINT 10CM U-SHAPED	105847-99-0	C	49.00					
JOYSTICK KNOB BODYPOINT RUBBER DOME	105848-99-0	C	45.00					
Notes				Total				