

# SmartDrive

Featuring the  
PushTracker E2

permobil



Rugged, durable  
and waterproof

Touch screen technology makes the  
SmartDrive as flexible, accessible  
and convenient as possible

## PushTracker E2

SmartDrive gives manual wheelchair users more control over their power-assist experience than ever before. The PushTracker E2 makes using your SmartDrive even easier, with gesture-activated controls and activity tracking technology.

8x more reliable Bluetooth® connection

## Specifications

<b>Dimensions</b>	1.85 x 2.06 x 0.5 in (46.9 x 52.2 x 12.9 mm)	<b>Contains Bluetooth® LE</b>	FCC ID 2AP42-WG12026 IC ID 24006-WG12026
<b>Colors</b>	Black	<b>Water resistance</b>	5 ATM (Swim-ready, up to 50 m)
<b>Operating system</b>	Wear OS by Google™	<b>NFC payment</b>	No
<b>Phone compatibility</b>	Android™, iPhone™	<b>Band width</b>	0.87 in (22 mm)
<b>Nominal operating voltage</b>	3.9 V	<b>Band material</b>	Silicone (interchangeable) - Latex-Free and completely Biocompatible
<b>Electric charge</b>	415 mA-h	<b>Case material</b>	Polycarbonate
<b>Wireless operating band</b>	Bluetooth / BLE 2.4 GHz (2.402 to 2.480 GHz) WiFi 802.11 b/g/n 2.4 GHz (2.412 to 2.472 GHz) GPS / GLONASS / Beidou 1575 MHz (1.559 to 1.610 GHz)	<b>Display</b>	AMOLED (400 x 400 pixel)
		<b>App version(s)</b>	Check Google Play™ Store

*\*All figures shown are for illustration purposes only and actual products may vary. Wear OS by Google, iPhone, Android, Bluetooth and Google Play are trademarks or registered trademarks of their respective owners.*