Power wheelchairs M300 HD



permobil



Rugged style, refined details

M300 HD takes all the things you love about Permobil's midwheel drive wheelchairs and amplifies them. Thanks to the stronger actuators, more powerful motors, and a reinforced seating system, it is built to withstand the rigors of your active heavy duty lifestyle. You'll benefit from the fully adjustable legrests, and can take your chair to the next level with the optional 8" seat elevator and up to four power seat functions.

Specifications

Maximum user weight	450 lb
Maximum speed	5 mph
Range	13 mi (Grp 34) 20 mi (Grp 24) [°]
Base width	26"
Base length	38"
Minimum turning radius	22"
Ground clearance	3"
Recommended battery type	Group 34 & Group 24
Weight including batteries	425 lb with Group 34 447 lb with Group 24
Drive electronics	R-Net 120A

Seat-to-floor height	18", 19", 20"
Optional seat elevator	8"
Power tilt	0° – 45°
Recline options	85° – 120° (manual) 85° – 150° (power)
Legrest elevation	85° - 175°
Backrest height	21", 24" – 29" (by 1" incr.)
Seat depth	19" – 25" (by 1" incr.)
Seat width	19", 21", 23", 25"

*Actual driving range will vary based on driving and battery conditions.

Specifications may vary based upon configuration. Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions.

Color options



permobil