



Table of contents

| Price List Item | Page |
|--------------------------------|---------|
| Cushion Sizing Chart | 2 |
| ROHO® Seat Cushions | 3 - 10 |
| Select Series | 3 - 4 |
| Single Compartment Series | 4 - 5 |
| Dual Compartment Series | 5 - 6 |
| Hybrid Series | 6 - 8 |
| Smart Series | 9 - 10 |
| Other Cushioning Products | 10 |
| Specialty Application Products | 11 |
| Seating Accessories and Covers | 12 |
| ROHO AGILITY® Back Supports | 13 - 17 |
| AGILITY Standard Series | 13 - 17 |
| AGILITY Active Series | 18 - 20 |
| Terms and Conditions | 21 |

ROHO® Cushion Sizing Chart

ROHO DRY FLOATATION® cushion sizing chart for:

- ROHO LOW PROFILE® Single or Dual Compartment Cushion
- ROHO MID PROFILE™ Single or Dual Compartment Cushion
- ROHO HIGH PROFILE® Single or Dual Compartment Cushion
- ROHO ENHANCER® Cushion
- ROHO MID PROFILE Single Compartment with Sensor Ready® Technology
- ROHO HIGH PROFILE Single Compartment with Sensor Ready Technology
- ROHO QUADTRO SELECT® LOW PROFILE Cushion
- ROHO QUADTRO SELECT MID PROFILE Cushion
- ROHO QUADTRO SELECT HIGH PROFILE Cushion
- ROHO CONTOUR SELECT® Cushion

Item numbers for cushions listed for use with the sizing chart include the number of cells used for width and depth. Refer to the model number sections in this price list to determine the correct item number to use. Example: ROHO QUADTRO SELECT HIGH PROFILE cushion - require a cushion to fit a chair size of 20" (51cm) x 18" (45.5cm), select QS1110C.

| Chair Size (Inches) | Chair Size (cm) | Number of Cells |
|------------------------|--------------------|--------------------|
| 10 | 25.5 | 5 |
| 11 - 12 | 28 - 30.5 | 6 |
| 13 | 33 | 7 |
| 14 - 15 | 35.5 - 38 | 8 |
| 16 - 17 | 40.5 - 43 | 9 |
| 18 - 19 | 45.5 - 48 | 10 |
| 20 - 21 | 51 - 53.5 | 11 |
| 22 | 56 | 12 |
| 23 - 24 | 58.5 - 61 | 13 |
| 25 | 63.5 | 14 |
| 26 - 27 | 66 - 68.5 | 15 |

Note: Availability of cushion sizes varies by model.

Select the correct size ROHO Cushion for your wheelchair

Step 1

Measure the wheelchair seat width and seat depth.

Note: Wheelchair manufacturers' published seat dimensions may vary.

Step 2

Use your wheelchair seat dimensions and the above table to help you determine the number of cells in cushion width and in depth.

ROHO seat cushions

- The following cushions are available in both custom modifications and custom sizes.
- Unless otherwise noted, all cushions include cover, pump, repair kit, and operating instructions.
- Class A1 & A2 products are eligible for contracted discounts as per provider purchasing contracts.
- Contact ROHO Customer Support for additional information.

Select Series Cushions

| Suggested Retail | Item | | | | | HCPCS* |
|----------------------|---|---|-----------|-----------|-----------|---|
| \$534.00 | QUADTRO SELECT HIGH PROFILE Cushion: <u>CLASS: A1</u> | | | | | E2624 |
| | QS66C | QS87C | QS910C | QS1011C | QS1111C |  |
| | QS67C | QS88C | QS911C | QS1012C | QS1112C | |
| | QS77C | QS89C | QS912C | QS1013C | QS1113C | |
| | QS78C | QS810C | QS108C | QS118C | QS1114C | |
| | QS79C | QS98C | QS109C | QS119C | | |
| | QS710C | QS99C | QS1010C | QS1110C | | |
| | | | | | | |
| \$594.00 | QUADTRO SELECT HIGH PROFILE Cushion: <u>CLASS: A1</u> | | | | | E2625 |
| | QS128C | QS1212C | QS139C | QS1313C | QS1410C | |
| | QS129C | QS1213C | QS1310C | QS1314C | QS1411C | |
| | QS1210C | QS1214C | QS1311C | QS148C | QS1412C | |
| | QS1211C | QS138C | QS1312C | QS149C | QS1413C | |
| \$534.00 \$624.00 | QUADTRO SELECT HIGH PROFILE Cushion: <u>CLASS: A2</u> | | | | | E2624 |
| | QSSMC* | Non Stock Small Sizes (11 cells wide or less) | | | | E2625 |
| | QSLGC* | Non Stock Large Sizes (12 cells wide or greater) | | | | |
| \$556.00 | QUADTRO SELECT HIGH PROFILE Cushion - Heavy Duty Cover: <u>CLASS: A1</u> | | | | | Non-Coded |
| | QS710HD* | QS910HD | QS109HD | QS1011HD | QS1111HD | |
| | QS99HD | QS911HD | QS1010HD | QS1110HD | | |
| \$616.00 | QUADTRO SELECT HIGH PROFILE Cushion - Heavy Duty Cover: <u>CLASS: A1</u> | | | | | Non-Coded |
| | QS139HD* | QS1312HD* | | | | |
| \$556.00 \$650.00 | QUADTRO SELECT HIGH PROFILE Cushion - Heavy Duty Cover: <u>CLASS: A2</u> | | | | | Non-Coded |
| | QSSMHD* | Non Stock Small Sizes (11 cells wide or less) | | | | |
| | QSLGHD* | Non Stock Large Sizes (12 cells wide or greater) | | | | |
| \$534.00 | QUADTRO SELECT MID PROFILE Cushion: <u>CLASS: A1</u> | | | | | E2624 |
| | QS77MPC | QS88MPC | QS99MPC | QS109MPC | | |
| | QS78MPC | QS89MPC | QS910MPC | QS1010MPC | | |
| \$534.00 \$624.00 | QUADTRO SELECT MID PROFILE Cushion: <u>CLASS: A2</u> | | | | | E2624 |
| | QSSMMP* | Non Stock Mid Profile Small Sizes (11 cells wide or less) | | | | E2625 |
| | QSLGMP* | Non Stock Mid Profile Large Sizes (12 cells wide or greater) | | | | |
| \$534.00 | QUADTRO SELECT LOW PROFILE Cushion: <u>CLASS: A1</u> | | | | | E2624 |
| | QS56LPC* | QS79LPC | QS98LPC | QS109LPC | QS1111LPC | |
| | QS66LPC | QS87LPC | QS99LPC | QS1010LPC | | |
| | QS67LPC | QS88LPC | QS910LPC | QS1011LPC | | |
| | QS77LPC | QS89LPC | QS911LPC | QS119LPC | | |
| | QS78LPC | QS810LPC | QS108LPC | QS1110LPC | | |
| | | | | | | |
| \$594.00 | QUADTRO SELECT LOW PROFILE Cushion: <u>CLASS: A1</u> | | | | | E2625 |
| | QS1210LPC | QS1211LPC | QS1212LPC | QS1310LPC | QS1311LPC | |
| \$534.00 \$624.00 | QUADTRO SELECT LOW PROFILE Cushion: <u>CLASS: A2</u> | | | | | E2624 |
| | QSSMLPC* | Non Stock Small Sizes (11 cells wide or less) | | | | E2625 |
| | QSLGLPC* | Non Stock Large Sizes (12 cells wide or greater) | | | | |
| \$738.00 | QUADTRO SELECT Custom Cushion: <u>CLASS: A2</u> | | | | | E2609 |
| | QSCUSTC* | Custom Cushion - with cover between 36 and 160 cells. Smaller than 36 cells or greater than 160 cells priced upon request. | | | | |


+ Decisions about how to code a product are the responsibility of the provider/supplier. Since Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. For coverage criteria in other countries, consult your local Permobil Seating & Positioning distributor.

*May take up to three weeks for shipment.

ROHO Seat Cushions

- The following cushions are available in both custom modifications and custom sizes.
- Unless otherwise noted, all cushions include cover, pump, repair kit, and operating instructions.
- Class A1 & A2 products are eligible for contracted discounts as per provider purchasing contracts.
- Contact ROHO Customer Support for additional information.

Select Series Cushions Continued

| Suggested Retail | Item | | | | | HCPCS+ |
|------------------|---|--|---------|---------|---------|---|
| \$534.00 | CONTOUR SELECT® Cushion: <u>CLASS: A1</u> | | | | | E2624 |
| | CS78C | CS810C | CS910C | CS1010C | CS1110C |  |
| | CS88C | CS98C | CS911C | CS1011C | CS1111C | |
| | CS89C | CS99C | CS109C | CS119C | | |
| \$594.00 | CONTOUR SELECT® Cushion: <u>CLASS: A1</u> | | | | | E2625 |
| | CS1210C | CS1211C | CS1212C | CS1310C | CS1311C | |
| \$534.00 | CONTOUR SELECT® Cushion: <u>CLASS: A2</u> | | | | | E2624 |
| | CSSMC* | Non Stock Small Sizes (11 cells wide or less) | | | | |
| \$624.00 | CSLGC* | Non Stock Large Sizes (12 cells wide or greater) | | | | E2625 |
| \$738.00 | CONTOUR SELECT® Cushion: <u>CLASS: A2</u> | | | | | E2609 |
| | CSCUSTC* Custom Cushion (multiple valves, deleted cells, cutouts, etc.) with cover between 36 and 160 cells. Smaller than 36 cells or greater than 160 cells priced upon request. | | | | | |

Single Compartment Series Cushions

| | | | | | | |
|----------|---|----------|----------|----------|---------|---|
| \$534.00 | HIGH PROFILE® Single Compartment Cushion: <u>CLASS: A1</u> | | | | | E2622 |
| | 1R66C | 1R88C | 1R910C | 1R1011C | 1R1111C |  |
| | 1R67C | 1R89C | 1R911C | 1R1012C | 1R1112C | |
| | 1R77C | 1R810C | 1R912C | 1R1013C | 1R1113C | |
| | 1R78C | 1R811C | 1R108C | 1R118C | 1R1114C | |
| | 1R79C | 1R98C | 1R109C | 1R119C | | |
| | 1R87C | 1R99C | 1R1010C | 1R1110C | | |
| \$594.00 | HIGH PROFILE® Single Compartment Cushion: <u>CLASS: A1</u> | | | | | E2623 |
| | 1R128C | 1R1213C | 1R1312C | 1R1410C | 1R1511C | |
| | 1R129C | 1R1214C | 1R1313C | 1R1411C | | |
| | 1R1210C | 1R138C | 1R1314C | 1R1412C | | |
| | 1R1211C | 1R1310C | 1R148C | 1R1413C | | |
| | 1R1212C | 1R1311C | 1R149C | 1R1414C | | |
| | | | | | | |
| \$534.00 | HIGH PROFILE® Single Compartment Cushion: <u>CLASS: A2</u> | | | | | E2622 |
| | 1RSMC* Non Stock Small Size (11 cells wide or less) | | | | | E2623 |
| \$624.00 | 1RLGC* Non Stock Large Sizes (12 cells wide or greater) | | | | | |
| \$556.00 | HIGH PROFILE Single Compartment Cushion - Heavy Duty: <u>CLASS: A1</u> | | | | | Non-Coded |
| | 1R99HD | 1R911HD | 1R1010HD | 1R1110HD | | |
| | 1R910HD | 1R109HD | 1R1011HD | 1R1111HD | | |
| \$616.00 | HIGH PROFILE Single Compartment Cushion - Heavy Duty: <u>CLASS: A1</u> | | | | | Non-Coded |
| | 1R1210HD | 1R1211HD | | | | |
| \$556.00 | HIGH PROFILE Single Compartment Cushion - Heavy Duty: <u>CLASS: A2</u> | | | | | Non-Coded |
| | 1RSMHD* Non Stock Small Sizes (11 cells wide or less) | | | | | |
| \$650.00 | 1RLGHD* Non Stock Large Sizes (12 cells wide or greater) | | | | | |



+ Decisions about how to code a product are the responsibility of the provider/supplier. Since Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. For coverage criteria in other countries, consult your local Permobil Seating & Positioning distributor.

*May take up to three weeks for shipment.



ROHO Seat Cushions

- The following cushions are available in both custom modifications and custom sizes.
- Unless otherwise noted, all cushions include cover, pump, repair kit, and operating instructions.
- Class A1 & A2 products are eligible for contracted discounts as per provider purchasing contracts.
- Contact ROHO Customer Support for additional information.

Single Compartment Series Cushions Continued

| Suggested Retail | Item | HCPCS ⁺ |
|----------------------|---|---|
| \$534.00 | MID PROFILE Single Compartment Cushion: <u>CLASS: A1</u> 1R77MPC 1R88MPC 1R99MPC 1R109MPC 1R1011MPC 1R1111MPC 1R78MPC 1R89MPC 1R910MPC 1R1010MPC 1R1110MPC | E2622  |
| \$534.00 \$624.00 | MID PROFILE Single Compartment Cushion: <u>CLASS: A2</u> 1RSMMPC* Non Stock Sizes (11 cells wide or less) 1RLGMPC* Non Stock Large Sizes (12 cells wide or greater) | Non-Coded |
| \$534.00 | LOW PROFILE Single Compartment Cushion: <u>CLASS: A1</u> 1R66LPC 1R79LPC 1R810LPC 1R911LPC 1R1011LPC 1R67LPC 1R87LPC 1R98LPC 1R108LPC 1R119LPC 1R77LPC 1R88LPC 1R99LPC 1R109LPC 1R1110LPC 1R78LPC 1R89LPC 1R910LPC 1R1010LPC 1R1111LPC | E2622  |
| \$594.00 | LOW PROFILE Single Compartment Cushion: <u>CLASS: A1</u> 1R1210LPC 1R1310LPC | E2623 |
| \$534.00 \$624.00 | LOW PROFILE Single Compartment Cushion: <u>CLASS: A2</u> 1RSMLPC* Non Stock Small Size (11 cells wide or less) 1RLGLPC* Non Stock Large Size (12 cells wide or greater) | E2622 E2623 |


Dual Compartment Series Cushions

| | | |
|--|---|--|
| \$559.00 | HIGH PROFILE Dual Compartment Cushion: <u>CLASS: A2</u> 2R88C 2R99C 2R109C 2R1012C 2R89C 2R910C 2R1010C 2R119C 2R98C 2R911C 2R1011C 2R1110C | E2624  |
| \$559.00 | LOW PROFILE Dual Compartment Cushion: <u>CLASS: A2</u> 2R88LPC 2R99LPC 2R109LPC 2R1011LPC 2R119LPC 2R89LPC 2R910LPC 2R1010LPC 2R1012LPC 2R1110LPC | E2624  |
| \$559.00 \$559.00 \$660.00 \$660.00 | Non-Stock Dual Compartment Cushion: <u>CLASS: A2</u> 2RSMC* HIGH PROFILE, Non Stock Small Sizes (11 cells wide or less) 2RSMLPC* LOW PROFILE, Non Stock Small Sizes (11 cells wide or less) 2RLGC* HIGH PROFILE, Non Stock Large Sizes (12 cells wide or greater) 2RLGLPC* LOW PROFILE, Non Stock Large Sizes (12 cells wide or greater) | E2624 E2624 E2625 E2625 |
| \$738.00 | Custom Dual Compartment Cushion: <u>CLASS: A2</u> CUSTC* Custom Cushion (multiple valves, deleted cells, cutouts, etc.) with cover between 36 and 160 cells. Smaller than 36 cells or greater than 160 cells priced upon request. | E2609 |

ROHO Seat Cushions


- The following cushions are available in both custom modifications and custom sizes.
- Unless otherwise noted, all cushions include cover, pump, repair kit, and operating instructions.
- Class A1 & A2 products are eligible for contracted discounts as per provider purchasing contracts.
- Contact ROHO Customer Support for additional information.

Dual Compartment Series Cushions Continued

| Suggested Retail | Item | | | | | HCPCS* |
|------------------|---|--|---------|----------|----------|---|
| \$534.00 | ENHANCER Cushion: <u>CLASS: A1</u> | | | | | E2624 |
| | ENH88C | ENH98C | ENH910C | ENH1010C | ENH1110C |  |
| | ENH89C | ENH99C | ENH109C | ENH1011C | | |
| \$534.00 | ENHANCER Cushion: <u>CLASS: A2</u> | | | | | E2624 |
| | ENHSMC* | Non Stock Small Sizes (11 cells wide or less) | | | | E2638 |
| \$624.00 | ENHLGC* | Non Stock Large Sizes (12 cells wide or greater) | | | | |
| \$738.00 | ENHANCER Cushion: <u>CLASS: A2</u> | | | | | E2609 |
| | ENHCUSTC* | Custom Cushion (multiple valves, deleted cells, cutouts, etc.) with cover between 36-160 cells. Smaller than 36 cells or greater than 160 cells priced upon request. | | | | |

Hybrid Series Cushions

| | | | | | | | | |
|----------|---|-------------|------------------------------------|-------------|-------------|-------------|------------------------------------|-------------|
| \$571.00 | Hybrid Select Cushion: <u>CLASS: A1</u> | | | | | | | |
| | Item Number | | Fits Chair Size (width x depth) | | Item Number | | Fits Chair Size (width x depth) | |
| | HS1414C | | 14" x 14" / 36 cm x 36 cm | | HS1717C | | 17" x 17" / 43 cm x 43 cm | |
| | HS1416C | | 14" x 16" / 36 cm x 41 cm | | HS1718C | | 17" x 18" / 43 cm x 46 cm | |
| | HS1515C | | 15" x 15" / 38 cm x 38 cm | | HS1720C | | 17" x 20" / 43 cm x 51 cm | |
| | HS1516C | | 15" x 16" / 38 cm x 41 cm | | HS1816C | | 18" x 16" / 46 cm x 41 cm | |
| | HS1517C | | 15" x 17" / 38 cm x 43 cm | | HS1817C | | 18" x 17" / 46 cm x 43 cm | |
| | HS1615C | | 16" x 15" / 41 cm x 38 cm | | HS1818C | | 18" x 18" / 46 cm x 46 cm | |
| | HS1616C | | 16" x 16" / 41 cm x 41 cm | | HS1820C | | 18" x 20" / 46 cm x 51 cm | |
| | HS1617C | | 16" x 17" / 41 cm x 43 cm | | HS2016C | | 20" x 16" / 51 cm x 41 cm | |
| | HS1618C | | 16" x 18" / 41 cm x 46 cm | | HS2017C | | 20" x 17" / 51 cm x 43 cm | |
| | HS1620C | | 16" x 20" / 41 cm x 51 cm | | HS2018C | | 20" x 18" / 51 cm x 46 cm | |
| | HS1716C | | 17" x 16" / 43 cm x 41 cm | | HS2020C | | 20" x 20" / 51 cm x 51 cm | |
| \$151.00 | Hybrid Select Replacement Foam Base: <u>CLASS: A1</u> | | | | | | | |
| | Item Number | Description | Item Number | Description | Item Number | Description | Item Number | Description |
| | 46400 | 14 x 14 | 46410 | 16 x 16 | 46417 | 17 x 27 | 46423 | 20 x 16 |
| | 46402 | 14 x 16 | 46411 | 16 x 17 | 46418 | 17 x 20 | 46424 | 20 x 17 |
| | 46405 | 15 x 15 | 46412 | 16 x 18 | 46419 | 18 x 16 | 46425 | 20 x 18 |
| | 46406 | 15 x 16 | 46413 | 16 x 20 | 46420 | 18 x 17 | 46426 | 20 x 20 |
| | 46407 | 15 x 17 | 46415 | 17 x 16 | 46421 | 18 x 18 | 46425 | 20 x 18 |
| | 46409 | 16 x 15 | 46416 | 17 x 17 | 46422 | 18 x 20 | | |
| | | | | | | | | |




E2624

+ Decisions about how to code a product are the responsibility of the provider/supplier. Since Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. For coverage criteria in other countries, consult your local Permobil Seating & Positioning distributor.
*May take up to three weeks for shipment.

ROHO Seat Cushions

- Unless otherwise noted, all cushions include cover, pump, repair kit, and operating instructions.
- Class A1 & A2 products are eligible for contracted discounts as per provider purchasing contracts.
- Contact ROHO Customer Support for additional information.

Hybrid Select Series Cushions Continued

| Suggested Retail | Item | | | | HCPCS* | |
|--|---|---------------------------------------|--------------|---------------------------------------|---|-------------|
| \$322.00 | Hybrid Select Replacement Air Pad: <u>CLASS: A1</u> | | | | | |
| | Item Number | Description | Item Number | Description | Item Number | Description |
| | HSAP14 | 14 Widths | HSAP16 | 20 Widths | HSAP18 | 18 Widths |
| | HSAP15 | 15 Widths | HSAP17 | 17 Widths | HSAP20 | 20 Widths |
| \$217.00 | Hybrid Select Replacement IT Air Pad: <u>CLASS: A1</u> | | | | | |
| | Item Number | Description | Item Number | Description | | |
| | HSIAP14-16 | IT Air Pad Hybrid Select 14-16 Widths | HSIAP17-20 | IT Air Pad Hybrid Select 17-20 Widths | | |
| \$552.00 | Hybrid Elite SR® Single Compartment Cushion: <u>CLASS: A1</u> | | | | | |
| | Item Number | Fits Chair Size (width x depth) | Item Number | Fits Chair Size (width x depth) |  E2622 | |
| | 1RHE1414C-SR | 14" x 14" / 36 cm x 36 cm | 1RHE1716C-SR | 17" x 16" / 43 cm x 41 cm | | |
| | 1RHE1415C-SR | 14" x 15" / 36 cm x 38 cm | 1RHE1717C-SR | 17" x 17" / 43 cm x 43 cm | | |
| | 1RHE1416C-SR | 14" x 16" / 36 cm x 41 cm | 1RHE1718C-SR | 17" x 18" / 43 cm x 46 cm | | |
| | 1RHE1515C-SR | 15" x 15" / 38 cm x 38 cm | 1RHE1720C-SR | 17" x 20" / 43 cm x 51 cm | | |
| | 1RHE1516C-SR | 15" x 16" / 38 cm x 41 cm | 1RHE1816C-SR | 18" x 16" / 46 cm x 41 cm | | |
| | 1RHE1517C-SR | 15" x 17" / 38 cm x 43 cm | 1RHE1817C-SR | 18" x 17" / 46 cm x 43 cm | | |
| | 1RHE1615C-SR | 16" x 15" / 41 cm x 38 cm | 1RHE1818C-SR | 18" x 18" / 46 cm x 46 cm | | |
| | 1RHE1616C-SR | 16" x 16" / 41 cm x 41 cm | 1RHE1820C-SR | 18" x 20" / 46 cm x 51 cm | | |
| | 1RHE1617C-SR | 16" x 17" / 41 cm x 43 cm | 1RHE2017C-SR | 20" x 17" / 51 cm x 43 cm | | |
| | 1RHE1618C-SR | 16" x 18" / 41 cm x 46 cm | 1RHE2018C-SR | 20" x 18" / 51 cm x 46 cm | | |
| | 1RHE1620C-SR | 16" x 20" / 41 cm x 51 cm | 1RHE2020C-SR | 20" x 20" / 51 cm x 51 cm | | |
| \$552.00 | Hybrid Elite SR® Single Compartment Cushion: <u>CLASS: A1</u> | | | | | |
| | Item Number | Fits Chair Size (width x depth) | | | | |
| | 1RHE2218C-SR | 22" x 18" / 56 cm x 46 cm | | | | |
| | 1RHE2220C-SR | 22" x 20" / 56 cm x 51 cm | | | | |
| | 1RHE2222C-SR | 22" x 22" / 56 cm x 58 cm | | | | |
| | 1RHE2418C-SR | 24" x 18" / 61 cm x 46 cm | | | | |
| | 1RHE2420C-SR | 24" x 20" / 61 cm x 51 cm | | | | |
| ^Hybrid Elite SR Single Compartment Cushions are also a part of the Smart Series Cushions. | | | | | | |




E2622

E2623

ROHO Seat Cushions

- Unless otherwise noted, all cushions include cover, pump, repair kit, and operating instructions.
- Class A1 & A2 products are eligible for contracted discounts as per provider purchasing contracts.
- Contact ROHO Customer Support for additional information.

Hybrid Series Cushions Continued



| Suggested Retail | Item | | | | HCPCS ⁺ |
|--|---|---|--------------------|--|--|
| \$552.00 | Hybrid Elite® Dual Compartment Cushion: <u>CLASS: A1</u> | | | |  E2622 |
| | Item Number | Fits Chair Size (width x depth) | Item Number | Fits Chair Size (width x depth) | |
| | 2RHE1414C | 14" x 14" / 36 cm x 36 cm | 2RHE1716C | 17" x 16" / 43 cm x 41 cm | |
| | 2RHE1415C | 14" x 15" / 36 cm x 38 cm | 2RHE1717C | 17" x 17" / 43 cm x 43 cm | |
| | 2RHE1416C | 14" x 16" / 36 cm x 41 cm | 2RHE1718C | 17" x 18" / 43 cm x 46 cm | |
| | 2RHE1515C | 15" x 15" / 38 cm x 38 cm | 2RHE1720C | 17" x 20" / 43 cm x 51 cm | |
| | 2RHE1516C | 15" x 16" / 38 cm x 41 cm | 2RHE1816C | 18" x 16" / 46 cm x 41 cm | |
| | 2RHE1517C | 15" x 17" / 38 cm x 43 cm | 2RHE1817C | 18" x 17" / 46 cm x 43 cm | |
| | 2RHE1615C | 16" x 15" / 41 cm x 38 cm | 2RHE1818C | 18" x 18" / 46 cm x 46 cm | |
| | 2RHE1616C | 16" x 16" / 41 cm x 41 cm | 2RHE1820C | 18" x 20" / 46 cm x 51 cm | |
| | 2RHE1617C | 16" x 17" / 41 cm x 43 cm | 2RHE2017C | 20" x 17" / 51 cm x 43 cm | |
| | 2RHE1618C | 16" x 18" / 41 cm x 46 cm | 2RHE2018C | 20" x 18" / 51 cm x 46 cm | |
| | 2RHE1620C | 16" x 20" / 41 cm x 51 cm | 2RHE2020C | 20" x 20" / 51 cm x 51 cm | |
| \$611.00 | Hybrid Elite® Dual Compartment Cushion: <u>CLASS: A1</u> | | | | E2623 |
| | Item Number | Fits Chair Size (width x depth) | Item Number | Fits Chair Size (width x depth) | |
| | 2RHE2218C | 22" x 18" / 56 cm x 46 cm | 2RHE2418C | 24" x 18" / 61 cm x 46 cm | |
| | 2RHE2220C | 22" x 20" / 56 cm x 51 cm | 2RHE2420C | 24" x 20" / 61 cm x 51 cm | |
| | 2RHE2222C | 22" x 22" / 56 cm x 58 cm | | | |
| ^Hybrid Elite Dual Compartment Cushions are also a part of the Dual Compartment Series Cushions | | | | | |
| ROHO Hybrid Elite Cushion Heavy Duty Cover Options: <u>CLASS: A1</u> | | | | | |
| \$647.00 | 1RHEHD-SR* | Hybrid Elite Single Valve Sensor Ready with Heavy Duty Cover, ≥ 20" | | | Non-Coded |
| \$647.00 | 2RHEHD* | Hybrid Elite Dual Valve with Heavy Duty Cover, 20" and under | | | |
| \$814.00 | 1RHEHDLG-SR* | Hybrid Elite Single Valve Sensor Ready with Heavy Duty Cover, ≤ 22" | | | |
| \$814.00 | 2RHEHDLG* | Hybrid Elite Dual Valve with Heavy Duty Cover, 22" and over | | | |

⁺ Decisions about how to code a product are the responsibility of the provider/supplier. Since Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. For coverage criteria in other countries, consult your local Permobil Seating & Positioning distributor.
 *May take up to three weeks for shipment.

ROHO Seat Cushions

- Unless otherwise noted, all cushions include cover, pump, repair kit, and operating instructions.
- Class A1 & A2 products are eligible for contracted discounts as per provider purchasing contracts.
- Contact ROHO Customer Support for additional information.

Smart Series Cushions

| Suggested Retail | Item | | | | HCPCS* |
|------------------|---|----------------|----------------|----------------|---|
| \$534.00 | HIGH PROFILE Single Compartment Cushion with Sensor Ready Technology: CLASS: A1 refer to page 2 for sizes (width & depth) | | | | E2622 |
| | 1R66H-CA-SR | 1R88H-CA-SR | 1R911H-CA-SR | 1R119H-CA-SR |  |
| | 1R67H-CA-SR | 1R89H-CA-SR | 1R108H-CA-SR | 1R1110H-CA-SR | |
| | 1R77H-CA-SR | 1R810H-CA-SR | 1R109H-CA-SR | 1R1111H-CA-SR | |
| | 1R78H-CA-SR | 1R98H-CA-SR | 1R1010H-CA-SR | 1R1112H-CA-SR | |
| | 1R79H-CA-SR | 1R99H-CA-SR | 1R1011H-CA-SR | | |
| | 1R87H-CA-SR | 1R910H-CA-SR | 1R1012H-CA-SR | | |
| \$594.00 | HIGH PROFILE Single Compartment Cushion with Sensor Ready Technology: CLASS: A1 | | | | E2622 |
| | 1R129H-CA-SR | 1R1210H-CA-SR | 1R1211H-CA-SR | 1R1212H-CA-SR | |
| \$594.00 | HIGH PROFILE Single Compartment Cushion with Sensor Ready Technology: CLASS: A1 | | | | E2623 |
| | 1R1310H-CA-SR | 1R1311H-CA-SR | 1R1411H-CA-SR | 1R1511H-CA-SR | |
| \$556.00 | HIGH PROFILE Single Compartment Cushion with Sensor Ready Technology with Heavy Duty Cover: CLASS: A1 refer to page 2 for sizes (width & depth) | | | | E2622 |
| | 1R66H-CH-SR* | 1R88H-CH-SR* | 1R108H-CH-SR* | 1R119H-CH-SR* | |
| | 1R67H-CH-SR* | 1R89H-CH-SR* | 1R109H-CH-SR* | 1R1110H-CH-SR* | |
| | 1R77H-CH-SR* | 1R810H-CH-SR* | 1R1010H-CH-SR* | 1R1111H-CH-SR* | |
| | 1R78H-CH-SR* | 1R98H-CH-SR* | 1R1011H-CH-SR* | 1R1112H-CH-SR* | |
| | 1R79H-CH-SR* | 1R99H-CH-SR* | 1R1012H-CH-SR* | | |
| | 1R87H-CH-SR* | 1R911H-CH-SR* | | | |
| | | | | | |
| \$616.00 | HIGH PROFILE Single Compartment Cushion with Sensor Ready Technology with Heavy Duty Cover: CLASS: A1 | | | | E2622 |
| | 1R129H-CH-SR* | 1R1210H-CH-SR* | 1R1211H-CH-SR* | 1R1212H-CH-SR* | |
| \$616.00 | HIGH PROFILE Single Compartment Cushion with Sensor Ready Technology with Heavy Duty Cover: CLASS: A1 | | | | E2623 |
| | 1R1310H-CH-SR* | 1R1311H-CH-SR* | 1R1411H-CH-SR | 1R1511H-CH-SR | |
| \$534.00 | MID PROFILE Single Compartment Cushion with Sensor Ready Technology: CLASS: A1 | | | | E2622 |
| | 1R77M-CA-SR | 1R810M-CA-SR | 1R108M-CA-SR | 1R119M-CA-SR |  |
| | 1R78M-CA-SR | 1R98M-CA-SR | 1R109M-CA-SR | 1R1110M-CA-SR | |
| | 1R87M-CA-SR | 1R99M-CA-SR | 1R1010M-CA-SR | 1R1111M-CA-SR | |
| | 1R88M-CA-SR | 1R910M-CA-SR | 1R1011M-CA-SR | 1R1112M-CA-SR | |
| | 1R89M-CA-SR | 1R911M-CA-SR | 1R1012M-CA-SR | | |
| | | | | | |
| | | | | | |
| | | | | | |
| \$594.00 | MID PROFILE Single Compartment Cushion with Sensor Ready Technology: CLASS: A1 | | | | E2622 |
| | 1R1210M-CA-SR | 1R1211M-CA-SR | 1R1212M-CA-SR | | |
| \$594.00 | MID PROFILE Single Compartment Cushion with Sensor Ready Technology: CLASS: A1 | | | | E2623 |
| | 1R1310M-CA-SR | 1R1311M-CA-SR | | | |


+ Decisions about how to code a product are the responsibility of the provider/supplier. Since Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. For coverage criteria in other countries, consult your local Permobil Seating & Positioning distributor.

*May take up to three weeks for shipment.



ROHO Seat Cushions

- Unless otherwise noted, all cushions include cover, pump, repair kit, and operating instructions.
- Class A1 & A2 products are eligible for contracted discounts as per provider purchasing contracts.
- Contact ROHO Customer Support for additional information.

Smart Series Cushions Continued

| Suggested Retail | Item | HCPCS ⁺ |
|------------------|---|--|
| \$556.00 | MID PROFILE Single Compartment Cushion with Sensor Ready Technology with Heavy Duty Cover: CLASS: A1 | |
| | 1R77M-CH-SR* 1R89M-CH-SR* 1R910M-CH-SR* 1R1010M-CH-SR* 1R1110M-CH-SR* | |
| | 1R78M-CH-SR* 1R810M-CH-SR* 1R911M-CH-SR* 1R1011M-CH-SR* 1R1111M-CH-SR* | E2622 |
| | 1R87M-CH-SR* 1R98M-CH-SR* 1R108M-CH-SR* 1R1012M-CH-SR* 1R1112M-CH-SR* | |
| | 1R88M-CH-SR* 1R99M-CH-SR* 1R109M-CH-SR* 1R119M-CH-SR* | |
| \$616.00 | MID PROFILE Single Compartment Cushion with Sensor Ready Technology with Heavy Duty Cover: CLASS: A1 | |
| | 1R1210M-CH-SR* 1R1211M-CH-SR* 1R1212M-CH-SR* | E2622 |
| | 1R1310M-CH-SR* 1R1311M-CH-SR* | E2623 |
| \$187.00 | Smart Check[®] by ROHO: CLASS: A1 | Non-Coded |
| | SMARTCHECK-01 Set up and push to check device for 1R HIGH PROFILE, MID PROFILE and HYBRID ELITE Sensor Ready cushions |  |

Other Cushioning Products

| | | | |
|----------|---|---------------------------------------|---|
| \$155.00 | MOSAIC[®] Cushion: CLASS: B | |  |
| | Item Number | Standard Sizes (width x depth) | |
| | MOS1616CA | 16" x 16" / 41 cm x 41 cm | |
| | MOS1618CA | 16" x 18" / 41 cm x 46 cm | |
| | MOS1717CA | 17" x 17" / 44 cm x 44 cm | |
| | The above listed MOSAIC Cushions include a standard cover. | | |
| | • Weight limit: 315 lbs. / 143 kg | | |
| \$184.00 | MOSAIC[®] Cushion with Heavy Duty Cover: CLASS: B | | |
| | Item Number | Standard Sizes (width x depth) | |
| | MOS1616HD | 16" x 16" / 41 cm x 41 cm | |
| | MOS1618HD | 16" x 18" / 41 cm x 46 cm | |
| | MOS1717HD | 17" x 17" / 44 cm x 44 cm | |
| | • Weight limit: 315 lbs. / 143 kg | | |
| \$244.00 | ROHO AirLITE[®] Cushion: CLASS: B | |  |
| | Item Number | Standard Sizes (width x depth) | |
| | AL1313HD | 13" x 13" / 33 cm x 33 cm | |
| | AL1315HD | 13" x 15" / 33 cm x 38 cm | |
| | AL1414HD | 14" x 14" / 36 cm x 36 cm | |
| | AL1416HD | 14" x 16" / 36 cm x 41 cm | |
| | AL1515HD | 15" x 15" / 38 cm x 38 cm | |
| | AL1517HD | 15" x 17" / 38 cm x 43 cm | |
| | AL1616HD | 16" x 16" / 41 cm x 41 cm | |
| | AL1618HD | 16" x 18" / 41 cm x 46 cm | |
| | AL1717HD | 17" x 17" / 43 cm x 43 cm | |
| | AL1718HD | 17" x 18" / 43 cm x 46 cm | |
| | AL1720HD | 17" x 20" / 43 cm x 51 cm | |
| | AL1816HD | 18" x 16" / 46 cm x 41 cm | |
| | AL1818HD | 18" x 18" / 46 cm x 46 cm | |
| | AL1820HD | 18" x 20" / 46 cm x 51 cm | |
| | AL1917HD | 19" x 17" / 48 cm x 43 cm | |
| | AL1919HD | 19" x 19" / 48 cm x 48 cm | |
| | AL2017HD | 20" x 17" / 51 cm x 43 cm | |
| | AL2018HD | 20" x 18" / 51 cm x 46 cm | |
| | • Weight limit: 300 lbs. / 136 kg | | |
| | • Does not include repair kit or pump. | | |

⁺ Decisions about how to code a product are the responsibility of the provider/supplier. Since Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. For coverage criteria in other countries, consult your local Permobil Seating & Positioning distributor.
^{*}May take up to three weeks for shipment.

Speciality Application Products

- Unless otherwise noted, all cushions include cover, pump, repair kit, and operating instructions.
- Contact ROHO Customer Support for additional information.

| Suggested Retail | Item | | |
|------------------|--|--|---|
| \$1,262.00 | Personal Recliner Cushioning System: <u>CLASS: A3</u> | |  |
| | Item Number | Product Description | |
| | RECLINERCUSTC* | Custom Recliner Wheelchair cushion, back and seat protection Made to fit per order • Priced upon request if exceeds 320 cells overall | |
| \$373.00 | MINI-MAX® Cushion: <u>CLASS: A3</u> | |  |
| | Item Number | Standard Sizes (width x depth) | |
| | MINI1515C | 15" x 15" / 38 cm x 38 cm | |
| | MINI1616C | 16" x 16" / 41 cm x 41 cm | |
| | MINI1816C | 18" x 16" / 46 cm x 41 cm | |
| | MINI1818C | 18" x 18" / 46 cm x 46 cm | |
| | MINI2018C | 18" x 20" / 46 cm x 51 cm | |
| | MINI2020C | 20" x 20" / 51 cm x 51 cm | |
| \$242.00 | PACK-IT® Cushion: <u>CLASS: A3</u> | |  |
| | Item Number | Standard Sizes (width x depth) | |
| | PACKIT | 16.25" x 9.5" / 41.5 cm x 24 cm | |
| \$135.00 | The ADAPTOR® Pad: <u>CLASS: A3</u> | |  |
| | Item Number | Standard Size (width x depth) | |
| | ADAPTOR | 9.75" X 13.5" / 25 cm x 34 cm • Self-contained air cells can be cut in sections to fit various surface shapes • Repair kit, pump, and cover not included | |
| \$308.00 | Toilet Seat Cushion: <u>CLASS: A3</u> | |  |
| | Item Number | Standard Size | |
| | SEATTOILET | • Fits all standard toilet seats • Does not include cover | |
| \$460.00 | Shower / Commode Cushion: <u>CLASS: A3</u> (LOW PROFILE® cells only) | |  |
| | Item Number | Standard Sizes (width x depth) | |
| | SEATCMD88LP | 14.75" x 14.75" / 37.5 cm x 37.5 cm | |
| | SEATCMD89LP | 14.75" x 16.5" / 37.5 cm x 42 cm | |
| | SEATCMD99LP | 16.5" x 16.5" / 42 cm x 42 cm | |
| | SEATCMD910LP | 16.5" x 18.25" / 42 cm x 46.5 cm | |
| | SEATCMD1010LP | 18.25" X 18.25" / 46.5 cm x 46.5 cm | |
| | SEATCMDCUST | Custom sizes available • Does not include cover | |
| \$116.00 | LTV® Seat Cushion: <u>CLASS: H</u> (non-medical cushion used for additional comfort) | |  |
| | Item Number | Standard Sizes (width x depth) | |
| | LTVCL | 18" x 16" / 45.5 cm x 40.5 cm - Black Quilted Fabric Cover | |
| \$175.00 | LTVBLK | 18" x 16" / 45.5 cm x 40.5 cm - Black Quilted ULTRALEATHER® Cover | |

* May take up to three weeks for shipment.

Seating Accessories and Covers

| Suggested Retail | Item | | | | | |
|--|--|---|--|-------------|-----------------------------------|--|
| \$67.00 | ROHO Contour Base (used with sling seat): CLASS: C | | | | | |
| | Item Number | Standard Sizes (width x depth) | | | | |
| | BASE1313 | 13.25" x 13" / 33.5 cm x 33 cm | | | | |
| | BASE1515 | 15.25" x 15" / 38.5 cm x 38.5 cm | | | | |
| | BASE1715 | 17" x 15.25" / 43 cm x 38.5 cm | | | | |
| | BASE1717 | 17" x 17.5" / 43 cm x 44.5 cm | | | | |
| \$56.00 (\$85.00: SIMTO) | ROHO Planar Solid Seat Insert (used with any seat): CLASS: C | | | | | |
| | Item Number | Standard Sizes (width x depth) | | Item Number | Standard Sizes (width x depth) | |
| | SI1212 | 11.5" x 11.5" / 29 cm x 29 cm | | SI1717 | 16.5" x 16.5" / 42 cm x 42 cm | |
| | SI1214 | 11.5" x 13.5" / 29 cm x 34.5 cm | | SI1718 | 16.5" x 17.5" / 42 cm x 44.5 cm | |
| | SI1414 | 13.5" x 13.5" / 34.5 cm x 34.5 cm | | SI1818 | 17.5" x 17.5" / 44.5 cm x 44.5 cm | |
| | SI1415 | 13.5" x 14.5" / 34.5 cm x 37 cm | | SI1820 | 17.5" x 19.5" / 44.5 cm x 49.5 cm | |
| | SI1416 | 13.5" x 15.5" / 34.5 cm x 39.5 cm | | SI1822 | 17.5" x 21.5" / 44.5 cm x 54.5 cm | |
| | SI1515 | 14.5" x 14.5" / 37 cm x 37 cm | | SI1824 | 17.5" x 23.5" / 44.5 cm x 59.5 cm | |
| | SI1516 | 14.5" x 15.5" / 37 cm x 39.5 cm | | SI2020 | 19.5" x 19.5" / 49.5 cm x 49.5 cm | |
| | SI1517 | 14.5" x 16.5" / 38 cm x 42 cm | | SI2022 | 19.5" x 21.5" / 49.5 cm x 54.5 cm | |
| | SI1616 | 15.5" x 15.5" / 39.5 cm x 39.5 cm | | SI2024 | 19.5" x 23.5" / 49.5 cm x 59.5 cm | |
| | SI1617 | 15.5" x 16.5" / 39.5 cm x 42 cm | | SI2424 | 24" x 24" / 61 cm x 61 cm | |
| | SI1618 | 15.5" x 17.5" / 39.5 cm x 44.5 cm | | SIMTO* | (Indicate measurements needed) | |
| | SI1620 | 15.5" x 19.5" / 39.5 cm x 49.5 cm | | | | |
| | • No weight limit; however the Solid Seat Insert must be properly sized to the cushion and the individual. | | | | | |
| \$61.00 \$61.00 \$83.00 \$83.00 \$69.00 \$83.00 \$69.00 \$72.00 | Covers: CLASS: C | | | | | |
| | Item Number | Product Description | | | | |
| | COV-SSA (size) | Select Series Cushion Cover - Please specify size. | | | | |
| | COV-A (model/size) | Standard Cushion Cover - Please specify cushion model and size. | | | | |
| | COV-ACUST (model/size) | Custom Cushion Cover - Please specify cushion model and size. | | | | |
| | COV-AMTO | Non-Standard Stock Cushion Cover - Please specify cushion model and size. | | | | |
| | COV-HE (size) | Hybrid Elite Cushion Cover - Please specify cushion model and size. | | | | |
| | COV-HD (model/size) | Heavy duty, fluid-resistant cover - Available to fit most standard-sized ROHO DRY FLOATATION, Hybrid Elite seat cushions and MOSAIC cushions. | | | | |
| | COV-HDAL (size) | AirLITE Heavy Duty Cushion Cover - Please specify size. | | | | |
| | COV-HS (size) | Hybrid Select Cover - Please specify size | | | | |
| \$10.00 \$14.00 \$64.00 \$64.00 | Miscellaneous Items: CLASS: C | | | | | |
| | Item Number | Product Description | | | | |
| | HP45 | ROHO Inflation Pump | | | | |
| | HP20 | ROHO Inflation Pump for support surface products | | | | |
| | Miscellaneous Modifications: CLASS: C | | | | | |
| | | Hook & Loop Attachment | | | | |
| | Drainage Holes | | | | | |
| \$754.00 | Custom cushions for back supports | | | | | |
| | Item Number | Product Description | | | | |
| | CUSTBKC | Custom neoprene back with cover | | | | |

* May take up to three weeks for shipment.

+ Repairs are only available for neoprene custom products and mattress overlays.

ROHO AGILITY Standard Series Back Support Systems

- All AGILITY Standard Series Back Supports include attachment hardware, cover, pump, and operating instructions.
- Contact Customer Support for additional information.
- AG6 = Max Contour
- -1 = Quick Release Hardware, -2 = Fixed Hardware, -P = Direct Mount Hardware

ROHO AGILITY Standard Series Sizing Matrix

A – Back Shell Length


A² – Back Shell and Foam Length

B – Back Cane Width Range
(there will be a 1.5cm loss of seat depth if an AGILITY Aluminum is fitted to the smallest wheelchair width within each size range)

C – Inner Width

D – Flat Center Section of Shell

E – Length of Wings



AGILITY Standard Series Max Contour Sizing Matrix

| | | | | | | | | | | | | |
|----------|------|-----|------|-----|-----------|-----------|----|-------|------|-------|------|----|
| AG6-1416 | 40.5 | 16" | 48.5 | 19" | 33-35.5 | 13" - 14" | 28 | 11" | 15 | 6" | 20.5 | 8" |
| AG6-1418 | 45.5 | 18" | 53.5 | 21" | 33-35.5 | 13" - 14" | 28 | 11" | 15 | 6" | 20.5 | 8" |
| AG6-1420 | 51 | 20" | 58.5 | 23" | 33-35.5 | 13" - 14" | 28 | 11" | 15 | 6" | 20.5 | 8" |
| AG6-1516 | 40.5 | 16" | 48.5 | 19" | 35.5-38 | 14" - 15" | 32 | 12.5" | 20.5 | 8" | 20.5 | 8" |
| AG6-1518 | 45.5 | 18" | 53.5 | 21" | 35.5-38 | 14" - 15" | 32 | 12.5" | 20.5 | 8" | 20.5 | 8" |
| AG6-1520 | 51 | 20" | 58.5 | 23" | 35.5-38 | 14" - 15" | 32 | 12.5" | 20.5 | 8" | 20.5 | 8" |
| AG6-1616 | 40.5 | 16" | 48.5 | 19" | 38-40.5 | 15" - 16" | 33 | 13" | 23 | 9" | 20.5 | 8" |
| AG6-1618 | 45.5 | 18" | 53.5 | 21" | 38-40.5 | 15" - 16" | 33 | 13" | 23 | 9" | 20.5 | 8" |
| AG6-1620 | 51 | 20" | 58.5 | 23" | 38-40.5 | 15" - 16" | 33 | 13" | 23 | 9" | 20.5 | 8" |
| AG6-1716 | 40.5 | 16" | 48.5 | 19" | 40.5-43 | 16" - 17" | 37 | 14.5" | 25.5 | 10" | 23 | 9" |
| AG6-1718 | 45.5 | 18" | 53.5 | 21" | 40.5-43 | 16" - 17" | 37 | 14.5" | 25.5 | 10" | 23 | 9" |
| AG6-1720 | 51 | 20" | 58.5 | 23" | 40.5-43 | 16" - 17" | 37 | 14.5" | 25.5 | 10" | 23 | 9" |
| AG6-1816 | 40.5 | 16" | 48.5 | 19" | 43-45.5 | 17" - 18" | 38 | 15" | 29 | 11.5" | 23 | 9" |
| AG6-1818 | 45.5 | 18" | 53.5 | 21" | 43-45.5 | 17" - 18" | 38 | 15" | 29 | 11.5" | 23 | 9" |
| AG6-1820 | 51 | 20" | 58.5 | 23" | 43-45.5 | 17" - 18" | 38 | 15" | 29 | 11.5" | 23 | 9" |
| AG6-1916 | 40.5 | 16" | 48.5 | 19" | 45.5-48.5 | 18" - 19" | 42 | 16.5" | 30.5 | 12" | 23 | 9" |
| AG6-1918 | 45.5 | 18" | 53.5 | 21" | 45.5-48.5 | 18" - 19" | 42 | 16.5" | 30.5 | 12" | 23 | 9" |
| AG6-1920 | 51 | 20" | 58.5 | 23" | 45.5-48.5 | 18" - 19" | 42 | 16.5" | 30.5 | 12" | 23 | 9" |
| AG6-2016 | 40.5 | 16" | 48.5 | 19" | 48.5-51 | 19" - 20" | 43 | 17" | 32 | 12.5" | 23 | 9" |
| AG6-2018 | 45.5 | 18" | 53.5 | 21" | 48.5-51 | 19" - 20" | 43 | 17" | 32 | 12.5" | 23 | 9" |
| AG6-2020 | 51 | 20" | 58.5 | 23" | 48.5-51 | 19" - 20" | 43 | 17" | 32 | 12.5" | 23 | 9" |

ROHO AGILITY Standard Series Back Support Systems

- All AGILITY Standard Series Back Supports include attachment hardware, cover, pump, and operating instructions.
- Refer to pages 12 & 13 for sizing and back dimension information.
- AG2 = Minimum Contour, AG3 = Mid Contour, AG6 = Max Contour
- -1 = Quick Release Hardware, -2 = Fixed Hardware, -P = Direct Mount Hardware
- Contact ROHO Customer Support for additional information.

AGILITY Standard Series - Max Contour



| | | | | | |
|----------|---|-------------|---------------------------|--------------|-----------------------|
| \$957.00 | AGILITY Max Contour Back System with Quick Release Hardware: | | | | Air Only E2620 |
| | <u>CLASS: M</u> | | | | |
| | Air Only | | Air and Lumbar Pad | | |
| | Item Number | | Item Number | | Air & Lumbar E2615 |
| | AG6-1416A-1 | AG6-1720A-1 | AG6-1416AL-1 | AG6-1720AL-1 | |
| | AG6-1418A-1 | AG6-1816A-1 | AG6-1418AL-1 | AG6-1816AL-1 | |
| | AG6-1420A-1 | AG6-1818A-1 | AG6-1420AL-1 | AG6-1818AL-1 | |
| | AG6-1516A-1 | AG6-1820A-1 | AG6-1516AL-1 | AG6-1820AL-1 | |
| | AG6-1518A-1 | AG6-1916A-1 | AG6-1518AL-1 | AG6-1916AL-1 | |
| | AG6-1520A-1 | AG6-1918A-1 | AG6-1520AL-1 | AG6-1918AL-1 | |
| | AG6-1616A-1 | AG6-1920A-1 | AG6-1616AL-1 | AG6-1920AL-1 | |
| | AG6-1618A-1 | AG6-2016A-1 | AG6-1618AL-1 | AG6-2016AL-1 | |
| | AG6-1620A-1 | AG6-2018A-1 | AG6-1620AL-1 | AG6-2018AL-1 | |
| | AG6-1716A-1 | AG6-2020A-1 | AG6-1716AL-1 | AG6-2020AL-1 | |
| | AG6-1718A-1 | | AG6-1718AL-1 | | |
| \$957.00 | AGILITY Max Contour Back System with Fixed Hardware: | | | | Air Only E2620 |
| | <u>CLASS: M</u> | | | | |
| | Air Only | | Air and Lumbar Pad | | |
| | Item Number | | Item Number | | Air & Lumbar E2615 |
| | AG6-1416A-2 | AG6-1816A-2 | AG6-1416AL-2 | AG6-1816AL-2 | |
| | AG6-1516A-2 | AG6-1916A-2 | AG6-1516AL-2 | AG6-1916AL-2 | |
| | AG6-1616A-2 | AG6-2016A-2 | AG6-1616AL-2 | AG6-2016AL-2 | |
| | AG6-1716A-2 | | AG6-1716AL-2 | | |




ROHO AGILITY Standard Series Back Support Systems

- All AGILITY Standard Series Back Supports include attachment hardware, cover, pump, and operating instructions.
- Refer to pages 12 & 13 for sizing and back dimension information.
- AG2 = Minimum Contour, AG3 = Mid Contour, AG6 = Max Contour
- -1 = Quick Release Hardware, -2 = Fixed Hardware, -P = Direct Mount Hardware
- Contact ROHO Customer Support for additional information.



AGILITY Standard Series - Max Contour Continued

| Suggested Retail | Item | | | | | HCPCS ⁺ |
|------------------|---|--------------|--------------|--------------|--------------|---|
| \$957.00 | AGILITY Max Contour Back System with Direct Mount Hardware with Air (Compatible with Permobil Power Wheelchairs with compatible Corpus Seating): CLASS: M | | | | | E2620 |
| | AG6-1416A-P | AG6-1520A-P | AG6-1716A-P | AG6-1818A-P | AG6-1920A-P |  |
| | AG6-1418A-P | AG6-1616A-P | AG6-1718A-P | AG6-1820A-P | AG6-2016A-P | |
| | AG6-1420A-P | AG6-1618A-P | AG6-1720A-P | AG6-1916A-P | AG6-2018A-P | |
| | AG6-1516A-P | AG6-1620A-P | AG6-1816A-P | AG6-1918A-P | AG6-2020A-P | |
| | AG6-1518A-P | | | | | |
| | | | | | | |
| \$957.00 | AGILITY Max Contour Back System with Direct Mount Hardware with Air and Lumbar (Compatible with Permobil Power Wheelchairs with compatible Corpus Seating): CLASS: M | | | | | E2615 |
| | AG6-1416AL-P | AG6-1520AL-P | AG6-1716AL-P | AG6-1818AL-P | AG6-1920AL-P |  |
| | AG6-1418AL-P | AG6-1616AL-P | AG6-1718AL-P | AG6-1820AL-P | AG6-2016AL-P | |
| | AG6-1420AL-P | AG6-1618AL-P | AG6-1720AL-P | AG6-1916AL-P | AG6-2018AL-P | |
| | AG6-1516AL-P | AG6-1620AL-P | AG6-1816AL-P | AG6-1918AL-P | AG6-2020AL-P | |
| | AG6-1518AL-P | | | | | |
| | | | | | | |

AGILITY Standard Series - Custom Back Support

| | | | | |
|--|----------------------------|-----------------|--------------------------------------|---|
| \$1,983.00 \$2,092.00 \$2,788.00 | AGILITY Custom Back | | | E2617 |
| | AG2CUSTC* | Minimum Contour | Contact Customer Support for details |  |
| | AG3CUSTC* | Mid Contour | Contact Customer Support for details | |
| | AG6CUSTC* | Max Contour | Contact Customer Support for details | |

AGILITY Standard Series - Laterals

| | | | | | |
|----------|---|-----------------|-----------------|--------------------|---|
| \$527.00 | AGILITY Back System Lateral Pads: CLASS: M | | | |  |
| | SWING-AWAY with Curved Pad | Pad Size | Pad Size | Offset Type | |
| | AG-LATSA-M-FL | Medium | 4 T x 4 L | Flush | |
| | AG-LATSA-M-OF | Medium | 4 T x 4 L | Offset | |
| | AG-LATSA-M-OF1 | Medium | 4 T x 4 L | 1" / 2.5 cm Offset | |
| | AG-LATSA-M-OF2 | Medium | 4 T x 4 L | 2" / 5 cm Offset | |
| | AG-LATSA-L-FL | Large | 4 T x 5 L | Flush | |
| | AG-LATSA-L-OF | Large | 4 T x 5 L | Offset | |
| | AG-LATSA-L-OF1 | Large | 4 T x 5 L | 1" / 2.5 cm Offset | |
| | AG-LATSA-L-OF2 | Large | 4 T x 5 L | 2" / 5 cm Offset | |
| | AG-LATSA-X-FL | X-Large | 5 T x 5 L | Flush | |
| | AG-LATSA-X-OF | X-Large | 5 T x 5 L | Offset | |
| | AG-LATSA-X-OF1 | X-Large | 5 T x 5 L | 1" / 2.5 cm Offset | |
| | AG-LATSA-X-OF2 | X-Large | 5 T x 5 L | 2" / 5 cm Offset | |
| \$236.00 | FIXED with Flat Pad | Pad Size | Pad Size | Offset Type | E0956 |
| | AG-LATFX-M-FL | Medium | 4 T x 4 L | Flush |  |
| | AG-LATFX-M-OF | Medium | 4 T x 4 L | Offset | |
| | AG-LATFX-L-FL | Large | 4 T x 5 L | Flush | |
| | AG-LATFX-L-OF | Large | 4 T x 5 L | Offset | |
| | AG-LATFX-X-FL | X-Large | 5 T x 5 L | Flush | |
| | AG-LATFX-X-OF | X-Large | 5 T x 5 L | Offset | |

+ Decisions about how to code a product are the responsibility of the provider/supplier. Since Medicare coding is subject to change, the provider should always confirm the HCPSC code and coverage criteria as part of the client assessment process. For coverage criteria in other countries, consult your local Permobil Seating & Positioning distributor.
*May take up to three weeks for shipment.

ROHO AGILITY Standard Series Back Support Systems

- All AGILITY Standard Series Back Supports include attachment hardware, cover, pump, and operating instructions.
- Refer to sizing chart on pages 12 & 13.
- AG2 = Minimum Contour, AG3 = Mid Contour, AG6 = Max Contour
- Contact ROHO Customer Support for additional information.

AGILITY Standard Series - Headrests

| Suggested Retail | Item |
|------------------|---|
| \$199.00 | Therafin Headrest for ROHO AGILITY Back Support Systems: <u>CLASS: M</u> |
| | AG-HRPAD-SM Small (5" Height x 8" Width) |
| | AG-HRPAD-MED Medium (5" Height x 10" Width) |
| \$244.00 | AG-HRPAD-LG Large (5" Height x 14" Width) |
| | Therafin Headrest for ROHO AGILITY Back Support Systems: <u>CLASS: M</u> |
| \$107.00 | AG-HRPAD-BRKT Multi-Axis Headrest Hardware |
| | AGILITY Headrest: <u>CLASS: C</u> |
| \$21.00 | AG-HR-BRKT Universal Headrest Adaptor Plate |
| \$21.00 | AG-HR-UT Headrest bracket and unitrack to fit 15"- 20" width |
| \$21.00 | AG-HR-UT14 Headrest bracket and unitrack to fit 14" width |




AGILITY Standard Series - Replacements

| | | | | | |
|---|--|--------------|--------------------|--|--------------------|
| Minimum and Mid Contours: \$75.00 Max Contour: \$87.00 | AGILITY Back Replacement Covers - Standard: <u>CLASS: C</u> | | | | |
| | Minimum Contour | | Mid Contour | | Max Contour |
| | COV-AG2-1410 | COV-AG2-1810 | COV-AG3-1410 | COV-AG3-1810 | COV-AG6-1416 |
| | COV-AG2-1413 | COV-AG2-1813 | COV-AG3-1413 | COV-AG3-1813 | COV-AG6-1418 |
| | COV-AG2-1416 | COV-AG2-1816 | COV-AG3-1416 | COV-AG3-1816 | COV-AG6-1420 |
| | COV-AG2-1418 | COV-AG2-1818 | COV-AG3-1418 | COV-AG3-1818 | COV-AG6-1616 |
| | COV-AG2-1513 | COV-AG2-1910 | COV-AG3-1420 | COV-AG3-1820 | COV-AG6-1618 |
| | COV-AG2-1610 | COV-AG2-1913 | COV-AG3-1510 | COV-AG3-1910 | COV-AG6-1620 |
| | COV-AG2-1613 | COV-AG2-1916 | COV-AG3-1513 | COV-AG3-1916 | COV-AG6-1816 |
| | COV-AG2-1616 | COV-AG2-1918 | COV-AG3-1518 | COV-AG3-1918 | COV-AG6-1818 |
| | COV-AG2-1618 | COV-AG2-2010 | COV-AG3-1610 | COV-AG3-1920 | COV-AG6-1820 |
| | COV-AG2-1710 | COV-AG2-2013 | COV-AG3-1613 | COV-AG3-2010 | COV-AG6-1920 |
| | COV-AG2-1713 | COV-AG2-2016 | COV-AG3-1616 | COV-AG3-2013 | COV-AG6-2016 |
| | COV-AG2-1716 | COV-AG2-2018 | COV-AG3-1618 | COV-AG3-2016 | COV-AG6-2018 |
| | COV-AG2-1718 | | COV-AG3-1620 | COV-AG3-2018 | COV-AG6-2020 |
| | | | COV-AG3-1710 | COV-AG3-2020 | |
| | | | COV-AG3-1713 | | |
| | AGILITY Minimum Contour Replacement Foam: <u>CLASS: C</u> | | | | |
| | Item Number For Foam With No Cut Out | | | Item Number For Foam With Cut Out | |
| | AG2FO-1410 | AG2FO-1710 | AG2FO-1910 | AG2FC-1410 | AG2FC-1710 |
| | AG2FO-1413 | AG2FO-1713 | AG2FO-1913 | AG2FC-1413 | AG2FC-1713 |
| | AG2FO-1416 | AG2FO-1716 | AG2FO-1916 | AG2FC-1416 | AG2FC-1716 |
| | AG2FO-1418 | AG2FO-1718 | AG2FO-1918 | AG2FC-1418 | AG2FC-1718 |
| | AG2FO-1513 | AG2FO-1810 | AG2FO-2010 | AG2FC-1513 | AG2FC-1810 |
| | AG2FO-1610 | AG2FO-1813 | AG2FO-2013 | AG2FC-1610 | AG2FC-1813 |
| | AG2FO-1613 | AG2FO-1816 | AG2FO-2016 | AG2FC-1613 | AG2FC-1816 |
| | AG2FO-1616 | AG2FO-1818 | AG2FO-2018 | AG2FC-1616 | AG2FC-1818 |
| | AG2FO-1618 | | | AG2FC-1618 | |

ROHO AGILITY Standard Series Back Support Systems

- All AGILITY Standard Series Back Supports include attachment hardware, cover, pump, and operating instructions.
- Refer to sizing chart on pages 12 & 13.
- AG2 = Minimum Contour, AG3 = Mid Contour, AG6 = Max Contour
- Contact ROHO Customer Support for additional information.

AGILITY Standard Series - Replacements Continued

| Suggested Retail | Item | | | | | | |
|------------------|---|---|------------|------------|------------|--|---|
| \$81.00 | AGILITY Mid Contour Replacement Foam with Cut Out: CLASS: C | | | | | | |
| | AG3FC-1410 | AG3FC-1510 | AG3FC-1616 | AG3FC-1810 | AG3FC-1910 | AG3FC-2010 | |
| | AG3FC-1413 | AG3FC-1513 | AG3FC-1618 | AG3FC-1813 | AG3FC-1916 | AG3FC-2013 | |
| | AG3FC-1416 | AG3FC-1518 | AG3FC-1620 | AG3FC-1816 | AG3FC-1916 | AG3FC-2016 | |
| | AG3FC-1418 | AG3FC-1610 | AG3FC-1710 | AG3FC-1818 | AG3FC-1918 | AG3FC-2018 | |
| | AG3FC-1420 | AG3FC-1613 | AG3FC-1713 | AG3FC-1820 | AG3FC-1920 | AG3FC-2020 | |
| \$76.00 | ROHO AGILITY Max Contour Replacement Foam with Cut Out: CLASS: C | | | | | | |
| | AG6FC-1416 | AG6FC-1616 | AG6FC-1620 | AG6FC-1818 | AG6FC-1920 | AG6FC-2018 | |
| | AG6FC-1418 | AG6FC-1618 | AG6FC-1816 | AG6FC-1820 | AG6FC-2016 | AG6FC-2020 | |
| | AG6FC-1420 | | | | | | |
| \$232.00 | AGILITY Back Support System Replacement Components: CLASS: C | | | | | | |
| | AG-HDWR-QR | QUICK RELEASE CANE CLAMP, AGILITY BACK PAIR | | | |  Quick Release | |
| | AG-HDWR-FX | FIXED CANE CLAMP, AGILITY BACK PAIR | | | | | |
| | AG-HDWR-DM14 | DIRECT MOUNT HDWR AGILITY 14W | | | | |  Fixed |
| | AG-HDWR-DM15 | DIRECT MOUNT HDWR AGILITY 15W | | | | | |
| | AG-HDWR-DM16 | DIRECT MOUNT HDWR AGILITY 16W | | | | | |
| | AG-HDWR-DM17 | DIRECT MOUNT HDWR AGILITY 17W | | | | | |
| | AG-HDWR-DM18 | DIRECT MOUNT HDWR AGILITY 18W | | | | |  Direct Mount |
| | AG-HDWR-DM19 | DIRECT MOUNT HDWR AGILITY 19W | | | | | |
| | AG-HDWR-DM20 | DIRECT MOUNT HDWR AGILITY 20W | | | | | |
| \$23.00 | AG-LUMBAR | LUMBAR PAD, AGILITY BACK | | | | | |
| \$12.00 | AG-WRENCH | ALLEN WRENCH, 4 MM AGILITY BACK | | | | | |

Back Support Accessories

ROHO Privacy Shield

| | Item Number | Sizes (width x height) | Fits Wheelchair Sizes (width x depth) |
|---------|-------------|---------------------------|--|
| \$65.00 | PS14* | 10" x 14" | 13" - 15" |
| | PS16* | 13" x 14" | 15.5" - 17" |
| | PS18* | 15" x 14" | 17.5" - 19" |
| | PS20* | 17" x 14" | 19.5" - 21" |

*May take up to three weeks for shipment.

ROHO AGILITY Active Series Back Support Systems

- All AGILITY Active Series Carbon Back supports include attachment hardware, a cover, pump, and operating instructions.
- Contact ROHO Customer Support for additional information.
- AG2 = Minimum Contour, AG3 = Mid Contour

AGILITY Active Series: Carbon Sizing Matrix

A¹ – Back Shell LengthE¹ – Length of Wings, ShellA² – Foam LengthE² – Length of Wings, Foam


B – Outside Cane to Outside Cane

F – Shell Depth

C – Inner Foam Width

G – Inner Foam Depth

D – Flat Center Section



Minimum Contour Sizing Matrix

| Number | A ¹ | A ² | B | C | D | E ¹ | E ² | F | G |
|---------------|----------------|----------------|---------------|---------------|---------------|----------------|----------------|-------------|-------------|
| AG2-1410A-2CB | 10" (25.5 cm) | 12" (30.5 cm) | 14" (35.5 cm) | 11" (28 cm) | 6" (15 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1413A-2CB | 13" (33 cm) | 15" (38 cm) | 14" (35.5 cm) | 11" (28 cm) | 6" (15 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1416A-2CB | 16" (40.5 cm) | 18" (45.5 cm) | 14" (35.5 cm) | 11" (28 cm) | 6" (15 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1510A-2CB | 10" (25.5 cm) | 12" (30.5 cm) | 15" (38 cm) | 12" (30.5 cm) | 7" (18 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1513A-2CB | 13" (33 cm) | 15" (38 cm) | 15" (38 cm) | 12" (30.5 cm) | 7" (18 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1516A-2CB | 16" (40.5 cm) | 18" (45.5 cm) | 15" (38 cm) | 12" (30.5 cm) | 7" (18 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1610A-2CB | 10" (25.5 cm) | 12" (30.5 cm) | 16" (40.5 cm) | 13" (33 cm) | 8" (20.5 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1613A-2CB | 13" (33 cm) | 15" (38 cm) | 16" (40.5 cm) | 13" (33 cm) | 8" (20.5 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1616A-2CB | 16" (40.5 cm) | 18" (45.5 cm) | 16" (40.5 cm) | 13" (33 cm) | 8" (20.5 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1710A-2CB | 10" (25.5 cm) | 12" (30.5 cm) | 17" (43 cm) | 14" (35.5 cm) | 9" (23 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1713A-2CB | 13" (33 cm) | 15" (38 cm) | 17" (43 cm) | 14" (35.5 cm) | 9" (23 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1716A-2CB | 16" (40.5 cm) | 18" (45.5 cm) | 17" (43 cm) | 14" (35.5 cm) | 9" (23 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1810A-2CB | 10" (25.5 cm) | 12" (30.5 cm) | 18" (45.5 cm) | 15" (38 cm) | 10" (25.5 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1813A-2CB | 13" (33 cm) | 15" (38 cm) | 18" (45.5 cm) | 15" (38 cm) | 10" (25.5 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1816A-2CB | 16" (40.5 cm) | 18" (45.5 cm) | 18" (45.5 cm) | 15" (38 cm) | 10" (25.5 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |

Mid Contour Sizing Matrix

| Number | A ¹ | A ² | B | C | D | E ¹ | E ² | F | G |
|---------------|----------------|----------------|---------------|---------------|---------------|----------------|----------------|------------|------------|
| AG3-1410A-2CB | 10" (25.5 cm) | 12" (30.5 cm) | 14" (35.5 cm) | 12" (30.5 cm) | 6" (15 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1413A-2CB | 13" (33 cm) | 15" (38 cm) | 14" (35.5 cm) | 12" (30.5 cm) | 6" (15 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1416A-2CB | 16" (40.5 cm) | 18" (45.5 cm) | 14" (35.5 cm) | 12" (30.5 cm) | 6" (15 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1510A-2CB | 10" (25.5 cm) | 12" (30.5 cm) | 15" (38 cm) | 13" (33 cm) | 7" (18 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1513A-2CB | 13" (33 cm) | 15" (38 cm) | 15" (38 cm) | 13" (33 cm) | 7" (18 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1516A-2CB | 16" (40.5 cm) | 18" (45.5 cm) | 15" (38 cm) | 13" (33 cm) | 7" (18 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1610A-2CB | 10" (25.5 cm) | 12" (30.5 cm) | 16" (40.5 cm) | 14" (35.5 cm) | 8" (20.5 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1613A-2CB | 13" (33 cm) | 15" (38 cm) | 16" (40.5 cm) | 14" (35.5 cm) | 8" (20.5 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1616A-2CB | 16" (40.5 cm) | 18" (45.5 cm) | 16" (40.5 cm) | 14" (35.5 cm) | 8" (20.5 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1710A-2CB | 10" (25.5 cm) | 12" (30.5 cm) | 17" (43 cm) | 15" (38 cm) | 9" (23 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1713A-2CB | 13" (33 cm) | 15" (38 cm) | 17" (43 cm) | 15" (38 cm) | 9" (23 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1716A-2CB | 16" (40.5 cm) | 18" (45.5 cm) | 17" (43 cm) | 15" (38 cm) | 9" (23 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1810A-2CB | 10" (25.5 cm) | 12" (30.5 cm) | 18" (45.5 cm) | 16" (40.5 cm) | 10" (25.5 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1813A-2CB | 13" (33 cm) | 15" (38 cm) | 18" (45.5 cm) | 16" (40.5 cm) | 10" (25.5 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1816A-2CB | 16" (40.5 cm) | 18" (45.5 cm) | 18" (45.5 cm) | 16" (40.5 cm) | 10" (25.5 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |

ROHO AGILITY Active Series Back Support Systems

- All AGILITY Active Series Aluminum Back supports include attachment hardware, a cover, and operating instructions.
- Contact ROHO Customer Support for additional information.
- AG2 = Minimum Contour, AG3 = Mid Contour

AGILITY Active Series: Aluminum Sizing Matrix



A¹ – Back Shell LengthE¹ – Length of Wings, ShellA² – Foam LengthE² – Length of Wings, Foam

B – Outside Cane to Outside Cane F – Shell Depth

C – Inner Foam Width

G – Inner Foam Depth

D – Flat Center Section



Minimum Contour Sizing Matrix

| Number | A ¹ | A ² | B | C | D | E ¹ | E ² | F | G |
|--------------|----------------|----------------|---------------|---------------|---------------|----------------|----------------|-------------|-------------|
| AG2-1408F-2A | 8" (20.5 cm) | 9" (23 cm) | 14" (35.5 cm) | 11" (28 cm) | 6" (15 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1410F-2A | 10" (25.5 cm) | 12" (30.5 cm) | 14" (35.5 cm) | 11" (28 cm) | 6" (15 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1508F-2A | 8" (20.5 cm) | 9" (23 cm) | 15" (38 cm) | 12" (30.5 cm) | 7" (18 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1510F-2A | 10" (25.5 cm) | 12" (30.5 cm) | 15" (38 cm) | 12" (30.5 cm) | 7" (18 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1608F-2A | 8" (20.5 cm) | 9" (23 cm) | 16" (40.5 cm) | 13" (33 cm) | 8" (20.5 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1610F-2A | 10" (25.5 cm) | 12" (30.5 cm) | 16" (40.5 cm) | 13" (33 cm) | 8" (20.5 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1708F-2A | 8" (20.5 cm) | 9" (23 cm) | 17" (43 cm) | 14" (35.5 cm) | 9" (23 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1710F-2A | 10" (25.5 cm) | 12" (30.5 cm) | 17" (43 cm) | 14" (35.5 cm) | 9" (23 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1808F-2A | 8" (20.5 cm) | 9" (23 cm) | 18" (45.5 cm) | 15" (38 cm) | 10" (25.5 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |
| AG2-1810F-2A | 10" (25.5 cm) | 12" (30.5 cm) | 18" (45.5 cm) | 15" (38 cm) | 10" (25.5 cm) | 4" (10 cm) | 5.5" (14 cm) | 3" (7.5 cm) | 3" (7.5 cm) |


Mid Contour Sizing Matrix

| Number | A ¹ | A ² | B | C | D | E ¹ | E ² | F | G |
|--------------|----------------|----------------|---------------|---------------|---------------|----------------|----------------|------------|------------|
| AG3-1408F-2A | 8" (20.5 cm) | 9" (23 cm) | 14" (35.5 cm) | 12" (30.5 cm) | 6" (15 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1410F-2A | 10" (25.5 cm) | 12" (30.5 cm) | 14" (35.5 cm) | 12" (30.5 cm) | 6" (15 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1508F-2A | 8" (20.5 cm) | 9" (23 cm) | 15" (38 cm) | 13" (33 cm) | 7" (18 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1510F-2A | 10" (25.5 cm) | 12" (30.5 cm) | 15" (38 cm) | 13" (33 cm) | 7" (18 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1608F-2A | 8" (20.5 cm) | 9" (23 cm) | 16" (40.5 cm) | 14" (35.5 cm) | 8" (20.5 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1610F-2A | 10" (25.5 cm) | 12" (30.5 cm) | 16" (40.5 cm) | 14" (35.5 cm) | 8" (20.5 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1708F-2A | 8" (20.5 cm) | 9" (23 cm) | 17" (43 cm) | 15" (38 cm) | 9" (23 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1710F-2A | 10" (25.5 cm) | 12" (30.5 cm) | 17" (43 cm) | 15" (38 cm) | 9" (23 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1808F-2A | 8" (20.5 cm) | 9" (23 cm) | 18" (45.5 cm) | 16" (40.5 cm) | 10" (25.5 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |
| AG3-1810F-2A | 10" (25.5 cm) | 12" (30.5 cm) | 18" (45.5 cm) | 16" (40.5 cm) | 10" (25.5 cm) | 5.25" (13 cm) | 6.75" (17 cm) | 4" (10 cm) | 4" (10 cm) |


ROHO AGILITY Active Series Back Support Systems

- All AGILITY Active Series Back supports include attachment hardware, a cover, and operating instructions. All AGILITY Active Series Carbon back supports also include a pump.
- Refer to page 20-21 for sizing and back dimension information.
- Contact ROHO Customer Support for additional information.
- AG2 = Minimum Contour, AG3 = Mid Contour

AGILITY Active Carbon Back Support

| Suggested Retail | Item | | | | | HCPCS ⁺ |
|------------------|--|---------------|---------------|---------------|---------------|---|
| \$1,202.00 | AGILITY Active Carbon Minimum Contour Back System with Air Only: CLASS: M | | | | | E2620 plus Carbon Upgrade |
| | AG2-1410A-2CB | AG2-1510A-2CB | AG2-1610A-2CB | AG2-1710A-2CB | AG2-1810A-2CB |  |
| | AG2-1413A-2CB | AG2-1513A-2CB | AG2-1613A-2CB | AG2-1713A-2CB | AG2-1813A-2CB | |
| | AG2-1416A-2CB | AG2-1516A-2CB | AG2-1616A-2CB | AG2-1716A-2CB | AG2-1816A-2CB | |
| \$1,202.00 | AGILITY Active Carbon Mid Contour Back Support System with Air Only: CLASS: M | | | | | E2620 plus Carbon Upgrade |
| | AG3-1410A-2CB | AG3-1510A-2CB | AG3-1610A-2CB | AG3-1710A-2CB | AG3-1810A-2CB | |
| | AG3-1413A-2CB | AG3-1513A-2CB | AG3-1613A-2CB | AG3-1713A-2CB | AG3-1813A-2CB | |
| | AG3-1416A-2CB | AG3-1516A-2CB | AG3-1616A-2CB | AG3-1716A-2CB | AG3-1816A-2CB | |

AGILITY Active Aluminum Back Support

| | | | | | | |
|----------|---|--------------|--------------|--------------|--------------|--|
| \$738.00 | AGILITY Active Aluminum Minimum Contour Back System with Foam Only: CLASS: M | | | | | E2620 |
| | AG2-1408F-2A | AG2-1508F-2A | AG2-1608F-2A | AG2-1708F-2A | AG2-1808F-2A |  |
| | AG2-1410F-2A | AG2-1510F-2A | AG2-1610F-2A | AG2-1710F-2A | AG2-1810F-2A | |
| \$738.00 | AGILITY Active Aluminum Mid Contour Back System with Foam Only: CLASS: M | | | | | E2620 |
| | AG3-1408F-2A | AG3-1508F-2A | AG3-1608F-2A | AG3-1708F-2A | AG3-1808F-2A | |
| | AG3-1410F-2A | AG3-1510F-2A | AG3-1610F-2A | AG3-1710F-2A | AG3-1810F-2A | |

AGILITY Active Series Replacements (Non-Coded)

| | | | | | | |
|----------|--|----------------|----------------|----------------|----------------|----------------|
| \$69.00 | AGILITY Active Carbon Replacement Foam: CLASS: C | | | | | |
| | AG2FC-1410CB | AG2FC-1610CB | AG2FC-1810CB | AG3FC-1410CB | AG3FC-1610CB | AG3FC-1810CB |
| | AG2FC-1413CB | AG2FC-1613CB | AG2FC-1813CB | AG3FC-1413CB | AG3FC-1613CB | AG3FC-1813CB |
| | AG2FC-1416CB | AG2FC-1616CB | AG2FC-1816CB | AG3FC-1416CB | AG3FC-1616CB | AG3FC-1816CB |
| | AG2FC-1510CB | AG2FC-1710CB | | AG3FC-1510CB | AG3FC-1710CB | |
| | AG2FC-1513CB | AG2FC-1713CB | | AG3FC-1513CB | AG3FC-1713CB | |
| | AG2FC-1516CB | AG2FC-1716CB | | AG3FC-1516CB | AG3FC-1716CB | |
| \$69.00 | AGILITY Active Aluminum Replacement Foam: CLASS: C | | | | | |
| | AG2FO-1408AC | AG2FO-1608AC | AG2FO-1808AC | AG3FO-1408AC | AG3FO-1608AC | AG3FO-1808AC |
| | AG2FO-1410AC | AG2FO-1610AC | AG2FO-1810AC | AG3FO-1410AC | AG3FO-1610AC | AG3FO-1810AC |
| | AG2FO-1508AC | AG2FO-1708AC | | AG3FO-1508AC | AG3FO-1708AC | |
| | AG2FO-1510AC | AG2FO-1710AC | | AG3FO-1510AC | AG3FO-1710AC | |
| \$103.00 | AGILITY Active Series Back Support Replacement Covers: CLASS: C | | | | | |
| | COV-AG2-1408CB | COV-AG2-1516CB | COV-AG2-1713CB | COV-AG3-1408CB | COV-AG3-1516CB | COV-AG3-1713CB |
| | COV-AG2-1410CB | COV-AG2-1608CB | COV-AG2-1716CB | COV-AG3-1410CB | COV-AG3-1608CB | COV-AG3-1716CB |
| | COV-AG2-1413CB | COV-AG2-1610CB | COV-AG2-1808CB | COV-AG3-1413CB | COV-AG3-1610CB | COV-AG3-1808CB |
| | COV-AG2-1416CB | COV-AG2-1613CB | COV-AG2-1810CB | COV-AG3-1416CB | COV-AG3-1613CB | COV-AG3-1810CB |
| | COV-AG2-1508CB | COV-AG2-1616CB | COV-AG2-1813CB | COV-AG3-1508CB | COV-AG3-1616CB | COV-AG3-1813CB |
| | COV-AG2-1510CB | COV-AG2-1708CB | COV-AG2-1816CB | COV-AG3-1510CB | COV-AG3-1708CB | COV-AG3-1816CB |
| | COV-AG2-1513CB | COV-AG2-1710CB | | COV-AG3-1513CB | COV-AG3-1710CB | |
| | | | | | | |

Terms and Conditions

1. Definitions:

- i. Buyer" shall mean the purchaser of any ROHO® Products.
- ii. Intellectual Property" means all copyrights, trademarks, trade names, patents, trade and corporate names, trade secrets, know how, formulas, technical information and other intellectual property rights concerning the Products.
- iii. "Person" shall mean an individual, association, joint venture, partnership, corporation, limited liability company, organization or other entity.
- iv. "Product" shall mean any goods sold by ROHO.
- v. "ROHO" shall mean ROHO, Inc., an Illinois corporation, all of its subsidiaries, its successors and assigns.

2. To the extent that the terms of any written agreement, signed by both ROHO and Buyer, are inconsistent with the provisions appearing hereon, the provisions of such executed agreement shall control. Except as provided herein, the provisions appearing hereon may not be altered or amended in any way except in writing signed by both ROHO and Buyer. ROHO's acceptance of any purchase order for Product submitted by Buyer to ROHO is expressly limited to Buyer's acceptance of all of the provisions hereof and any provisions contained in any document heretofore or hereafter submitted by Buyer to ROHO, which are in addition to or different from the provisions hereof, shall not become a part of the agreement between ROHO and Buyer unless said provisions are expressly agreed to in writing by ROHO; nor shall ROHO's filling of any purchase order submitted by Buyer for Products be construed as an acceptance of any such additional or different provisions.

3. Buyer's right to market, promote, sell or otherwise distribute ROHO Products is non-exclusive. Buyer agrees not to market, promote sell, or otherwise distribute ROHO Products outside of the United States, its territories and possessions.

4. Buyer agrees, at all times, to comply with the Product operation manual and other documents and guidelines produced and distributed by ROHO and shall not use, sell or otherwise distribute any Product for any purpose or use for which Product has not been designed or contrary to manufacturer and/or ROHO specifications. Buyer shall deliver all Products sold with its applicable operation manual, other documents and guidelines and all applicable accessories. Buyer agrees that all individuals dealing with the Products on its behalf will be trained in the selection of the appropriate Product and the individual setup and/or adjustment of the Products. Buyer shall insure that its customers and/or the customers' caregiver fully understand how to perform said setup and/or adjustments, all instructions and all warnings associated with the use of the Product and the fact that the use of any Product for any purpose or in any manner other than which it has been designed or any use inconsistent with the operation manual is strictly prohibited.

Buyer shall inform all of its customers about ROHO's standard Product warranties. Buyer shall not make any representations or statements inconsistent with the Product's operation manual, documents, guidelines and ROHO's standard Product warranties.

5. Buyer shall indemnify and hold ROHO, its officers, directors, employees, successors and assigns, harmless from any liability for any claim or action made upon or brought against ROHO, etc., by any Person, of whatever form or nature, including attorneys fees and other costs of legal defense, which may be sustained or incurred as a result of any acts or omissions of Buyer, or any of its directors, officers, employees or agents, including but not limited to:

- i. Reach of any provision of these terms and conditions;
- ii. Negligence or other tortuous actions;
- iii. Representations or statements for Product usage or alterations not specifically authorized in writing by ROHO; and
- iv. Violations of any laws, orders, regulations or ordinances.

Buyer shall, at its own expense, maintain comprehensive general liability insurance for bodily injury or death and property damage. This requirement shall in no way limit the liability of Buyer to ROHO. Buyer is responsible for any claims of damage or injury attributable to Buyer's actions.

6. Buyer recognizes that the Intellectual Property is an asset owned by ROHO and remains the property of ROHO. Buyer shall not:

- i. Use or permit the use of the Intellectual Property for any purpose other than as provided herein;
- ii. Permit any other Person to use the Intellectual Property;
- iii. Transfer or encumber the Intellectual Property;
- iv. Make or allow to be made any changes or alterations in the Intellectual Property, operations manual or other Product documents and guidelines;
- v. Remove, alter or deface any patent numbers, trade names, trademarks, notices, serial numbers, labels, tags, warning labels or other identifying marks, symbols or legends impressed on or affixed to any Product or container or package;
- vi. Disclose to others any information it receives from ROHO which is not generally known and available in the industry; or
- vii. Misrepresent ROHO or the Products.

7. The agreement between ROHO and Buyer is an Illinois contract and shall be governed by Illinois laws. ROHO and Buyer consent to the jurisdiction and venue of the Court sitting in St. Clair County, Illinois and waive the defense of inconvenient forum.

8. Pricing subject to change at ROHO's discretion with a 30-day notice.



Permobil Seating and Positioning
1501 South 74th Street • Belleville, IL 62223
Toll Free: 800-736-0925 • Local: 618-277-9173 • Fax: 618-277-6518
permobil.com • orders.roho@permobil.com

© 2019 Permobil. Products may be covered by one or more U.S. and foreign patents and trademarks, including:
ADAPTOR®, AGILITY®, AIRLITE®, CONTOUR SELECT®, DRY FLOATATION®, ENHANCER®, HIGH PROFILE®, Hybrid Elite®, Hybrid Elite SR®, LOW
PROFILE®, LTV®, MID PROFILE™, MINI-MAX®, MOSAIC®, PACK-IT®, QUADTRO SELECT®, ROHO®, Sensor Ready®, and Smart Check®.

ULTRALEATHER® is a registered trademark of Ultrafabrics, LLC.

Product information is changed as needed; current product information is available at permobil.com.

ROHO, Inc. is part of Permobil.