

Cushions

Comparison matrix



permobil



Comfort Curve®



Comfort Elements® with Gel



Comfort Support Pro 7 Series



ROHO® MOSAIC®



Comfort Ascent



ROHO® AirLITE



Comfort M2 Zero-Elevation M2 Anti-Thrust



Comfort Embrace™ Zero-Elevation Embrace™ Anti-Thrust

HPCS	General Use E2601/E2602	General Use E2601/E2602	General Use E2602	General Use E2601	Skin Protection E2603/E2604	Positioning E2605	Skin Protection and Positioning E2607/E2608	Skin Protection and Positioning E2607/E2608
Widths	8" - 24" (1" increments)	16"-24" (1" increments)	22"-30" (2" increments)	16" - 20" (1" increments)	16"- 24" (1" increments)	13" - 20" (1" increments)	10"-24" (1" increments)	10"-24" (1"increments)
Depths	8"-20" (1" increments)	16"-20" (1" increments)	18" and 20"	16" - 18" (1" increments)	16"-20" (1" increments)	13" - 20" (1" increments)	10"- 24" (1" increments)	10" - 24" (1" increments)
Cushion Thickness Dimensions	3" 2.25" inferior thigh	2" or 3"	3.5"	3"	4" 3" inferior thigh	2" - 3.75" contoured	Zero-elevation: ≤ 16" x 14", 3" ≥ 16" x 16", 4.25" Anti-thrust: ≤ 16" x 14", 2.75" ≥ 16" x 16", 3.75"	Zero-elevation: ≤ 16" x 14", 3" ≥ 16" x 16", 4" Anti-thrust: ≤ 16" x 14", 3" ≥ 16" x 16", 3"
Cushion Materials	Foam	Foam, QuadraGel® insert	Foam Optional QuadraGel® insert	Flowing Air	Foam	Foam, Air insert	Foam, QuadraGel® insert	Foam
Cover Options	COMFORT-TEK® STRETCH-AIR®	COMFORT-TEK® STRETCH-AIR®	Solace fabric with Kwik-Strap	Two-way stretch Heavy Duty cover	COMFORT-TEK® STRETCH-AIR®	Fluid-resistant, two-way stretch	COMFORT-TEK® STRETCH-AIR® COMFORT-TEK® with GlideWear® STRETCH-AIR® with GlideWear®	COMFORT-TEK® STRETCH-AIR® COMFORT-TEK® with GlideWear® STRETCH-AIR® with GlideWear®
Weight Capacity / Limit	12" wide and under - 100 lbs 14" wide - 150 lbs 16" wide - 250 lbs 18" wide - 300 lbs 20" wide - 350 lbs 22" wide and > - 400 lbs	16" wide - 250 lbs 18" wide - 300 lbs 20" wide - 350 lbs 22" wide and > - 400 lbs	750 lbs	315 lbs	16" wide - 250 lbs 18" wide - 300 lbs 20" wide - 350 lbs 22" wide and > - 400 lbs	300 lbs	14" wide and under- 200 lbs 16"wide - 250 lbs 18" wide - 300 lbs 20" wide - 350 lbs 22" wide and > - 400 lbs	14" wide and under- 200 lbs 16" wide - 250 lbs 18" wide - 300 lbs 20" wide - 350 lbs 22" wide and > - 400 lbs
Warranty	6 months cover 5 years cushion	6 months cover 5 years cushion	6 months cover 5 years cushion	6 months cover 12 months cushion	6 months cover 5 years cushion	6 months cover 24 months cushions	6 months cover 5 years cushion	6 months cover 5 years cushion
Add-ons	Incontinence liner	Incontinence liner	-	-	Incontinence liner Extra cover	-	GlideWear™ ATI Incontinence liner Extra cover Growth notches Rail cuts Colored piping	GlideWear™ ATI Incontinence liner Extra cover Growth notches Rail cuts Colored piping



Comfort
Saddle 7 Series



ROHO®
Single
Compartment



ROHO®
Hybrid Elite SR®
Single Compartment



ROHO®
HYBRID ELITE®
Dual Compartment



Vicair
Vicair® Adjuster O2



Comfort
Embrace with ATI



Comfort
M2 with ATI



Vicair
Vicair® Active O2

HPCS	Skin Protection and Positioning E2608	Adjustable Skin Protection E2622/E2623	Adjustable Skin Protection E2622/E2623	Adjustable Skin Protection E2622/E2623	Adjustable Skin Protection E2622/E2623	Adjustable Skin Protection and Positioning E2624/E2625	Adjustable Skin Protection and Positioning E2624/E2625	Adjustable Skin Protection and Positioning E2624/E2625
Widths	22" - 30" (2" increments)	10" - 27" (1" increments) 5 - 15 cells	14" - 18" (1" increments), 20", 22", and 24"	14" - 18" (1" increments), 20", 22", and 24"	14" - 20" (1" increments), 22" and 24"	10" - 24" (1" increments)	10" - 24" (1" increments)	14" - 24" (2" increments)
Depths	18" and 20"	10" - 25" (1" increments) 5 - 14 cells	14" - 18" (1" increments), 20" and 22"	14" - 18" (1" increments), 20" and 22"	14" - 22" (1" increments)	10" - 24" (1" increments)	10" - 24" (1" increments)	14" - 24" (2" increments)
Cushion Thickness Dimensions	4.5" 3.5" inferior thigh 5" abductor	HIGH PROFILE 4.25" cell height MID PROFILE 3.25" cell height LOW PROFILE 2.25" cell height	4.25"	4.25"	Standard: 3.93" / 10 cm Low: 2.36" / 6 cm	Zero-elevation: ≤ 16" x 14", 3" ≥ 16" x 16", 4" Anti-thrust: ≤ 16" x 14", 3" ≥ 16" x 16", 3"	Zero-elevation: ≤ 16" x 14", 3" ≥ 16" x 16", 4.25" Anti-thrust: ≤ 16" x 14", 2.75" ≥ 16" x 16", 3.75"	Standard: 3.93" / 10 cm Low: 2.36" / 6 cm
Cushion Materials	Foam Optional QuadraGel® insert	Flowing air	Foam Flowing air insert	Foam Flowing air insert	Non Flowing air	Foam	Foam, QuadraGel® insert	Non Flowing air
Cover Options	Solace fabric with Kwik-Strap	Two-way stretch Heavy duty cover (non coded)	Two-way stretch Heavy duty cover (non coded)	Two-way stretch Heavy duty cover (non coded)	Incotec Comfair Optional Kwik-Strap (limited sizes)	COMFORT-TEK® STRETCH-AIR® COMFORT-TEK® with GlideWear® STRETCH-AIR® with GlideWear®	COMFORT-TEK® STRETCH-AIR® COMFORT-TEK® with GlideWear® STRETCH-AIR® with GlideWear®	Active cover
Weight Capacity / Limit	750 lbs	No weight limit if appropriately sized to the individual	500 lbs	500 lbs	550 lbs	14" wide and under- 200 lbs 16" wide - 250 lbs 18" wide - 300 lbs 20" wide - 350 lbs 22" wide and > - 400 lbs	14" wide and under- 200 lbs 16" wide - 250 lbs 18" wide - 300 lbs 20" wide - 350 lbs 22" wide and > - 400 lbs	440 lbs
Warranty	6 months cover 5 years cushion	6 months cover 24 months cushion	6 months cover 24 months cushion	6 months cover 24 months cushion	24 months cover Lifetime on SmartCells™	6 months cover 5 years cushion	6 months cover 5 years cushion	24 months cover Lifetime on SmartCells™
Add-ons	-	Smart Check Privacy Shield Planar Solid Seat insert Contour base	Smart Check Privacy Shield	Privacy Shield	Extra Vicair SmartCells™	GlideWear™ ATI Incontinence liner Extra cover Growth notches Rail cuts Colored piping	GlideWear™ ATI Incontinence liner Extra cover Growth notches Rail cuts Colored piping	Extra Vicair SmartCells™



Vicair
Vicair® Vector O2



ROHO®
Dual Compartment



ROHO®
QUADTRO SELECT®



ROHO®
CONTOUR SELECT®



ROHO®
ENHANCER®



ROHO®
Hybrid Select



Comfort
Custom Inception

HCPCS	Adjustable Skin Protection and Positioning E2624/E2625	Adjustable Skin Protection and Positioning E2624/E2625	Adjustable Skin Protection and Positioning E2624/E2625	Adjustable Skin Protection and Positioning E2624/E2625	Adjustable Skin Protection and Positioning E2624/E2625	Adjustable Skin Protection and Positioning E2624	Custom Fabricated Seat Cushion E2609
Widths	14" - 20" (1" increments) 22" and 24"	10" - 27" (1" increments) 5 - 15 cells	10" - 27" (1" increments) 5 - 15 cells	10" - 27" (1" increments) 5 - 15 cells	10" - 27" (1" increments) 5 - 15 cells	14" - 18" (1" increments) and 20"	10" - 28" (1" increments)
Depths	14" - 22" (1" increments)	10" - 25" (1" increments) 5 - 14 cells	10" - 25" (1" increments) 5 - 14 cells	10" - 25" (1" increments) 5 - 14 cells	10" - 25" (1" increments) 5 - 14 cells	14" - 18" (1" increments) and 20"	10" - 28" (1" increments)
Cushion Thickness Dimensions	Standard: 3.93" / 10 cm Low: 2.36" / 6 cm	HIGH PROFILE 4.25" cell height LOW PROFILE 2.25" cell height	HIGH PROFILE 4.25" cell height MID PROFILE 3.25" cell height LOW PROFILE 2.25" cell height	2.25" - 4.25" due to contour	2.25" - 4.25" due to contour	5" 2.5" inferior thigh	Contact Customer Experience for additional information and sizes.
Cushion Materials	Non Flowing air	Flowing air	ISOFLO® Memory Control Flowing air	ISOFLO® Memory Control Flowing air	Flowing air	Foam ISOFLO® Memory Control Flowing air inserts	Foam ROHO® insert (Flowing air) QuadraGel® insert Ischial gel insert
Cover Options	Incotec Comfair Optional Kwik Strap (limited sizes)	Two-way stretch	Two-way stretch Heavy Duty	Two-way stretch	Two-way stretch	Two-way stretch	COMFORT-TEK® STRETCH-AIR® COMFORT-TEK® with GlideWear® STRETCH-AIR® with GlideWear®
Weight Capacity / Limit	550 lbs	No weight limit if appropriately sized to the individual	No weight limit if appropriately sized to the individual	No weight limit if appropriately sized to the individual	No weight limit if appropriately sized to the individual	500 lbs	No weight limit if appropriately sized to the individual
Warranty	24 months cover Lifetime on SmartCells™	6 months cover 24 months cushion	6 months cover 36 months cushion	6 months cover 36 months cushion	6 months cover 24 months cushion	6 months cover 36 months cushion	6 months cover 5 years cushion
Add-ons	Extra Vicair SmartCells™	Privacy Shield Planar Solid Seat insert Contour base	Privacy Shield Planar Solid Seat insert Contour base	Privacy Shield Planar Solid Seat insert Contour base	Privacy Shield Planar Solid Seat insert Contour base	Privacy Shield	GlideWear™ ATI Incontinence liner Extra cover Growth notches Rail cuts Coccyx cutout Leg length discrepancy Colored piping



ROHO®
Custom



Comfort
TruShape

HCPCS	Custom Fabricated Seat Cushion E2609	Custom Fabricated Seat Cushion E2609
Widths	10" - 27" (1" increments) 5 - 15 cells	Based on individual measurements
Depths	10" - 27" (1" increments) 5 - 15 cells	Based on individual measurements
Cushion Thickness Dimensions	2.25" - 4.25" cell height Contact Customer Experience for additional information.	Contact Customer Experience for additional information and sizes.
Cushion Materials	Flowing air ISOFOLO® Memory Control Single and Multiple Compartments	Foam ROHO® insert (Flowing air)
Cover Options	Two-way stretch Heavy duty cover	Neoprene Spacer fabric Breathe 4™ COMFORT-TEK® STRETCH-AIR®
Weight Capacity / Limit	No weight limit if appropriately sized to the individual	No weight limit if appropriately sized to the individual
Warranty	6 months cover 36 months cushion	6 months cover 5 years cushion
Add-ons	HD cover Planar Solid Seat insert Contour base	Double reinforcement Soft spots Soft spot overlay Ventilation holes Leg length discrepancy Undercuts Rail cuts Corner notches Extra cover



GlideWear™

GlideWear™ technology is designed to reduce the risk of pressure injury through the use of shear reduction forces allowing for skin protection. It offers 100% breathability for microclimate control and is latex free. Available on select Comfort cushions.

Permobil is the US distributor of Vicair® O2 cushions.
Vicair® is a registered trademark of Vicair BV.
GlideWear™ is a trademark of MIPS.

Adjustable Technology Insert

The Adjustable Technology Insert (ATI) is a moldable, high strength, light weight aluminum base that allows for individualized positioning and support.

Key benefits:

- Remoldable to accommodate the user's needs over time
- Designed to be placed beneath a compatible cushion, inside the cushion cover.
- Includes the metal insert, bending tool and operation manual



Smart Check™

The world's first personal electronic cushion monitoring system, Permobil's exclusive Smart Check device can help determine the proper inflation level for your ROHO® Sensor Ready cushion.

Key benefits:

- Stores your inflation range
- Press to check inflation level
- Real time feedback
- Ability to reset
- Detachable for storage



Privacy Shield

The ROHO® Privacy Shield is perfect for anyone who wants coverage between the bottom of the back support and the seat of their wheelchair. Designed to attach directly to the back canes of a wheelchair, it is easy to attach and detach for cleaning.