



# PDG Stellar GLT - AADL W321



Dealer Name\*: \_\_\_\_\_ Account #\*: \_\_\_\_\_

City, Province, Post Code: \_\_\_\_\_

Phone: \_\_\_\_\_ Email: \_\_\_\_\_ ATP: \_\_\_\_\_

Rep Name: \_\_\_\_\_ P.O. Number: \_\_\_\_\_ Quote Number: \_\_\_\_\_

Client Reference\*: \_\_\_\_\_ Client Weight: \_\_\_\_\_

Special considerations that need to be addressed here (e.g., diagnosis): \_\_\_\_\_ \*Required field

## October 2024, Revised January 2026. Prices are in Canadian Dollars (CAD)

### Wheelchair & Frame

Stellar GL 45° Manual Tilt Wheelchair Base Model	W321	\$ 2,685.00
Stellar GLT Lightweight Aluminum Frame Upgrade <i>Aluminum frame is not recommended for agitated users.</i>		\$ 300.00

#### Frame Type

Fixed Frame Width, 300lb Weight Capacity <i>Select a fixed seat width below.</i>	No charge
Adjustable Width Frame, 300lb Weight Capacity <i>Select an adjustable seat width below.</i>	\$ 200.00

#### User Weight | Required for setup of tilt uplift force.

<b>User Weight Range:</b>	Up to 75lbs	75lbs - 100lbs	100lbs - 115lbs	115lbs - 150lbs	150lbs - 200lbs	200lbs - 250lbs
	250lbs - 300lbs					

#### Frame Options

Transit Tie Down System (TTDS) <i>Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details. Not compatible with seat widths &gt; 22". Tested to maximum weight capacity of 300lbs.</i>	No charge
--	-----------

### Seat Measurements

#### Seat Width | Overall chair width is +9.5" with 12"/16" rear wheels, +10.5" with 20"/22"/24" rear wheels.

<b>Fixed:</b>	15"	16"	17"	18"	19"		No charge
	20"	21"	22"			W111	\$ 250.00
	23"	24"				W111	\$ 500.00
<b>Adjustable 15"-20":</b>	15"	16"	17"	18"	19"		No charge
	20"					W111	\$ 250.00

#### Seat Depth

<b>Select Depth:</b>	16"	17"	18"	19"		No charge
	20"	21"	22"		W111	\$ 250.00

*16"-20" seat depths are adjustable.*

### Frame Finishes

<b>Select Color:</b>	High Gloss Black	Hyper Blue
----------------------	------------------	------------

## Armrests

### Armrest Type

No Armrests No charge

Single Post, Height Adjustable Armrests No charge

8"-11" Range    10-14" Range    12"-16" Range

Armrest Position:    Rear    Front

*If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.*

### Armpads

Regular Armpads No charge

Desk Length 2" x 10"    Full Length 2" x 14"

### Clothing Guards

Clothing Guards, Single Post Armrests *Not compatible with 8"-11" height range.* \$ 125.00

### Armrest Extensions | *Not compatible with flip-up armrests or desk length armpads.*

Armrest Extensions - 3" Forward \$ 75.00

Armrest Extensions - 3" Rearward \$ 75.00

## Seat Height

### Seat Height - 6" Casters | *Seat height is adjustable within the range of the chosen lower frame type (Low/High)*

Low						High						
RW	24"											
	22"											
	20"											
	16"											
	12"											
	13"	14"	15"	16"	17"	18"	15"	16"	17"	18"	19"	20"
	Seat to Floor Height											

### Seat Height - 8" Casters | *Seat height is adjustable within the range of the chosen lower frame type (Low/High)*

Low				High					
RW	24"								
	16"								
	15"	16"	17"	18"	17"	18"	19"	20"	
	Seat to Floor Height								

## Seat Pan

No Seat Pan No charge

Solid Aluminium Seat Pan No charge

## Front Casters

### Urethane Casters

6" x 1.25" Caster No charge

8" x 1.25" Caster *Only compatible with 16" and 24" rear wheels.* No charge

## Rear Wheels

Composite Rear Wheel No charge

12"    16"    20"    22"    24"

Rear Wheel Position:    Rear    Front (Default)

*If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.*

### Axles

Quick Release Rear Axles No charge

Fixed Rear Axles No charge

## Rear Tires

Urethane Tires No charge

**Handrims** | *Handrims are only available with 20", 22" and 24" rear wheels.*

No Handrims No charge

Standard Aluminium Handrims No charge

Plastic Coated Handrims *Not compatible with drum brakes.* \$ 175.00

## Wheel Locks

User Wheel Locks No charge

Push to Lock Pull to Lock

## Wheel Lock Options

Wheel Lock Extensions *Not compatible with attendant wheel lock.* No charge

## Anti-Tippers

Rear Anti-Tippers *Required safety option.* No charge

## Backrest

Back Post Style: Stroller Handle Back

Back Post Height: Short (22") Tall (25")

*Measured from seat pan to handle.*

Folding Style: Fold-Down

*Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.*

Seat to Back Angle: Set at 90° Set at 98° Set at 106° Set at 114° Set at 122°

## Backrest Options

Angle/Height Adjustable Stroller Handle: No charge

Button on Outside Button on Inside

Spring Back Suspension *Not compatible with reclining back or dynamic rocker back.* \$ 105.00

## Tilt Activation

### Attendant Tilt Only

Attendant Tilt Only, Right Back Post No charge

Attendant Tilt Only, Left Back Post No charge

Two Trigger Setup *A trigger mounted on each side of the stroller handle. Reduces grip force required to activate tilt. Not compatible with reclining backrest or drum brakes.* No charge

**Tilt Limit** | *Optional. If no tilt limit is selected the full tilt range will be available.*

Tilt Limit: Set at 15° Set at 25° Set at 35°

**Anterior Pre-Tilt** | *Add up to 5° of Anterior Tilt to the tilt range.*

Anterior Pre-Tilt Set at 0° (Default) Set at 5°

## Front Rigging

No Front Rigging No charge

**Easy-Out Front Rigging** | *Effortless insertion and removal. No hanger brackets remain on the frame.*

Fixed 70° \$ 95.00

**Swing Away Front Rigging** | *Legacy design. Not compatible with Easy-Out.*

Swing Away 70° No charge

Swing Away Elevating Legrests (ELR) *Includes calf pads.* W106 \$ 375.00

1" Offset Hangers for Swing Away Front Rigging \$ 95.00

**Plug In Front Rigging** | *Legacy design. Not compatible with Easy-Out.*

Plug In 70° No charge

## Footplates

Composite Footplates		\$ 0.00
Aluminium Angle Adjustable Footplates	W107	\$ 165.00

### Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short	8"-12.5" for Easy-Out high mount; 14"-18" for Easy-Out; 13"-17" for swing away and plug in.	No charge
Standard Extension Tubes, Long	11"-16" for Easy-Out high mount; 18"-24" for Easy-Out; 15"-19" for swing away and plug in.	\$ 75.00

Footplate Height: Specify:

Ext	Short	11" - 16"	12" - 16"	13" - 16"
	Long	14" - 22"	14" - 19"	14" - 17"
		Standard Easy-Out	60°, 70°, ELR Swing-Away	Standard Plug-In
<b>Legrest</b>				

## Front Rigging Accessories

Length Adjustable Heel Loops		No charge
Standard Calf Strap	Basic Seating	\$ 50.00

## General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

## Notes