

TiLite Z Order Form

Rigid Frame

Effective June 2025 (New items in red)

Page 1 of 8

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:	User Weight:	Name:	Ship Via:
TiLite Sales Rep:		Address:	
TiLite Information:	Promo Code:	City:	State: Zip:
Ship Date:	Order #:	Phone:	Fax:

Please contact the Permobil Manual Customer Service Team at chairs.tilite@permobil.com - Phone: 800.736.0925 or Fax: 866.586.2413

FRAME STYLE

- ☐ Aluminum Rigid Frame (Full Adjustability).....\$2,750 K0005, E1235
- ☐ 250 lb. weight limit.....STD
- ☐ 275 lb. weight limit.....NCO
- ☐ Aluminum Rigid Frame with Transit Tie-Down.....\$2,995 K0005, E1235
- ☐ 250 lb. weight limit.....STD
- ☐ 275 lb. weight limit.....NCO

SMARTDRIVE POWER ASSIST SYSTEM



- ☐ Permobil Smart Drive MX2+* (N/A with rear seat heights less than 13.5").....\$7,529 E0986
- ☐ SwitchControl ButtonsNCO
- ☐ SwitchControl Button with Mono Jack Option.....NCO
- ☐ Buddy ButtonNCO
- ☐ Carrying Bag.....\$84
- ☐ Universal Padded Handle\$22
- *SmartDrive apps available on Apple App Store and Google Play Store for Apple Watch Series 5 or newer and Samsung Galaxy Watch 4 or newer.

TRANSIT OPTIONS

- ☐ Transit Tie-Down (Tested to WC-19. Consists of four Tie-Down Brackets and two Pelvic Belt Mounting Brackets). (Please refer to the user manual for important information regarding the Transit Tie-Down Options.).....\$363 E1023
- ☐ Transport Bracket Option (Consists of four Tie-Down Brackets) (Unoccupied use only).....\$379 E1022

FRAME FINISHES

- | | | |
|--|---|--|
| ☐ High Gloss Black | ☐ Red Pepper | ☐ Neon Pink |
| ☐ Ocean Blue Metallic | ☐ Tangerine Metallic | ☐ Acid Green Metallic |
| ☐ Midnight Blue Pearl | ☐ White Pearl | ☐ Black Lava |
| ☐ Titanium Grey | ☐ Pomegranate Red | ☐ Blue Mist |
| ☐ Indigo Purple | ☐ Hyper Blue | ☐ Flame Orange |
| ☐ OD Green | ☐ Silver Streak | ☐ Lavender Dusk |
| ☐ Black Galaxy | ☐ Seafoam Green | ☐ Daffodil Yellow.....STD |

COLOR ANODIZE PACKAGES

- ☐ Black.....STD
- ☐ Gold ☐ Green ☐ Purple ☐ Orange ☐ Blue ☐ Red.....\$195

TiFit: TAILORED INDIVIDUALIZED FRAME (Please see information on page 8 regarding frame measurements)

		Seat Width A									
Seat Depth B	12	13	14	15	16	17	18	19	20		
	12	13	14	15	16	17	18	19	20		
12	☐	☐									
13	☐	☐	☐	☐	☐						
14	☐	☐	☐	☐	☐	☐	☐	☐			
15	☐	☐	☐	☐	☐	☐	☐	☐	☐		
16	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	
17	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
18	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
19			☐	☐	☐	☐	☐	☐	☐	☐	☐
20			☐	☐	☐	☐	☐	☐	☐	☐	☐

Note: Configurations including 20" Seat Width\$559 E2201

Configurations including 20" Seat Depth\$675 E2203



Permobil.com
2701 West Court Street | Pasco, Washington 99301
509.352.2245 Tel | 866.586.2413 Toll-Free
chairs.tilite@permobil.com

TiLite® Z / June 2025

TiLite Z Order Form

Rigid Frame

Effective June 2025 (New items in red)

Page 2 of 8

CUSTOM FRAME DEPTH Measure from front of seat sling to beginning of frame bend. Frame bend will start 1" in front of seat sling unless Custom Frame Depth is specified, or the combination of specifications affect the rolling dynamics of the chair.

☐ 1" **STD**
☐ 0.5" ☐ 1.5" ☐ 2" **NCO**

FRONT SEAT HEIGHT c Measure from floor to top of seat tube at front edge of seat sling.

☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21" **STD**

REAR SEAT HEIGHT RANGES (Will be optimized to achieve desired rear seat height allowing +/- 1" adjustment)

Rear Wheel Size	Rear Seat Height d
20" (451)	13.0"-18.5"
22" (501)	14.0"-19.5"
24" (540)	15"-20.5"
25" (559)	15.5"-20.5"
26" (590)	16"-20.5"

REAR SEAT HEIGHT d Measure from floor to top of seat tube at back post.

☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5"
☐ 19" ☐ 19.5" ☐ 20" **STD**

FRONT ANGLE e Measure from frame front to floor.

☐ 75° ☐ 80° ☐ 85° (N/A with 6" Front Wheels.) ☐ 90° (N/A with 5" or 6" Front Wheels.) **STD**

CENTER OF GRAVITY f Measure from front of back post at the seat tube to the center of the rear axle.

☐ 1" ☐ 1.5" ☐ 2" ☐ 2.5" ☐ 3" ☐ 3.5" ☐ 4" ☐ 4.5" ☐ 5" **STD**
☐ Amputee Axle Plate (Center of Gravity positioned -1" behind back canes) **\$225**

E0959

FOOTREST OPTIONS

Footrest Width

SEAT WIDTH A (Outside Measurement)	No Taper	1.5" Taper (4" Narrower than Rear Seat Width)	3.5" Taper (6" Narrower than Rear Seat Width)
12"	☐ 9.5"	☐ 8"	
13"	☐ 10.5"	☐ 9"	
14"	☐ 11.5"	☐ 10"	
15"	☐ 12.5"	☐ 11"	☐ 9"
16"	☐ 13.5"	☐ 12"	☐ 10"
17"	☐ 14.5"	☐ 13"	☐ 11"
18"	☐ 15.5"	☐ 14"	☐ 12"
19"	☐ 16.5"	☐ 15"	☐ 13"
20"	☐ 17.5"	☐ 16"	☐ 14"

SEAT TO FOOTREST h Measure from front edge of seat sling to top rear of footrest

Standard Footrests Measurement must be at least 2" shorter than Front Seat Height

☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19"

☐ **Aluminum Open Loop Footrest** **STD**
☐ ABS Cover **\$30**
☐ Angle Adjustable Footrest **\$185**
☐ Flip-Back Angle Adjustable Footrest **\$244**

K0040

K0040

High Mount Footrest

☐ 7"-9" (with 4" risers) ☐ 9"-11" (with 2" risers) ☐ 11"-13" Set at _____ **\$155**

K0037, K0040

FOOTREST ACCESSORIES

☐ Safe Tread Non-Skid Adhesive Patch **\$20**

CALF STRAPS

Velcro Adjustable **\$44**

K0038

- ☐ Small (Recommended for 9" - 12" Footrest Widths.)
- ☐ Medium (Recommended for 11" - 15" Footrest Widths.)
- ☐ Large (Recommended for 14" - 17" Footrest Widths.)
- ☐ Extra Large (Recommended for 16" - 19" Footrest Widths.)

TiLite Z Order Form

Rigid Frame

Effective June 2025 (New items in red)

Page 3 of 8

Bodypoint® Padded Velcro Adjustable.....	\$109	K0038
☐ Small (Recommended for 8"- 11" Footrest Widths.)		
☐ Medium (Recommended for 11"- 15" Footrest Widths.)		
☐ Large (Recommended for 16"- 19" Footrest Widths.)		
Buckle Adjustable, Black Nylon Webbing	\$50	K0038

CLAMP ON RESIDUAL LIMB SUPPORT

☐ Left (8" width x 10" length).....	\$495	E1020
☐ Right (8" width x 10" length).....	\$495	E1020
☐ Pair (8" width x 10" length).....	\$895	E1020

SEAT BACK CONFIGURATIONS

SEAT BACK TYPE

☐ Folding Aluminum Adjustable Height Back.....	STD
---	-----

PUSH HANDLES

☐ Integrated Push Handles.....	STD
☐ Integrated Fold-Down Push Handles.....	\$125
☐ Bolt on Height Adjustable Push Handles (Requires 4" deep rigidizer bar.).....	\$363
☐ Omit Push Handles.....	NCO

SEAT BACK HEIGHT i Measure from top of seat tube at rear of frame to top of backrest upholstery.

Push Handles	Short	Medium	Tall	Set at: (Use .5" increments)
None / Bolt-On	9"-11.5"	12"-15.5"	16"-20.5"	Set at: _____
Integrated	10.0"-12.5"	13.5"-16.5"	17.5"-20.5"	Set at: _____
Integrated Fold Down	12" – 13.5"	14" – 17.5"	18" – 20.5"	Set at: _____

BACKREST OPTIONS

RIGIDIZER BAR

☐ 3" deep OR ☐ 4" deep (Standard is 3" deep rigidizer bar).....	NCO
☐ Locking Backrest (locks in up and down position).....	STD
☐ Non-Locking Backrest (locks in up position only).....	NCO
☐ Depth Adjustable Back/Frame (N/A with 19" or 20" Seat Depths. Allows 2" of rearward adjustment in 1" increments. Bolt On Height Adjustable Push Handles have a limited height range with Depth Adjustable Back).....	\$185

SEAT BACK ANGLE j Measure angle from front of back post to floor.

☐ 80° ☐ 83° ☐ 86° ☐ 89° ☐ 92° ☐ 95° ☐ 98° ☐ 101°	STD
--	-----

BACK UPHOLSTERY

☐ Standard Back Upholstery (Tension Adjustable. Includes Privacy Flap).....	STD	
☐ AIR Breathable Back Upholstery (Tension Adjustable. Includes Privacy Flap).....	\$369	E2611
☐ Removable Lumbar Contour Pack ☐ Removable Lateral Contour Pack		
☐ Omit Back Upholstery.....	NCO	
4" Military Service Patch.....	\$115	
☐ Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard ☐ Purple Heart		

PERMOBIL BACK SUPPORTS (Please see information on page 8 regarding product shipment)

ROHO® AGILITY® ACTIVE ALUMINUM BACK SUPPORT

Shell Height ☐ 8" ☐ 10"

Contour and Insert

☐ Minimum Contour AGILITY® Active Aluminum with Foam	\$738	E2620
☐ Mid Contour AGILITY® Active Aluminum with Foam	\$738	E2620
☐ ROHO Privacy Shield	\$65	

ROHO® AGILITY® ACTIVE CARBON

Shell Height ☐ 10" ☐ 13" ☐ 16"

Contour and Insert

☐ Minimum Contour, Foam.....	\$1,202	E2620
☐ Mid Contour, Foam.....	\$1,202	E2620
☐ ROHO Privacy Shield.....	\$65	

TiLite Z Order Form

Rigid Frame

Effective June 2025 (New items in red)

Page 4 of 8

COMFORT ACTA-BACK® (Backs over 20" wide or tall subject to increase pricing)

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"\$623 E2616/E2614
☐ Comfort-Tek Cover ☐ Stretch-Air Cover.....NCO
☐ Comfort Privacy Flap.....\$44

COMFORT ACTA-EMBRACE® (Backs over 20" wide or tall subject to increase pricing)

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"\$649 E2613/E2614
☐ Comfort-Tek Cover ☐ Stretch-Air Cover.....NCO
☐ Comfort Privacy Flap.....\$44

COMFORT ACTA-RELIEF (Backs over 20" wide or tall subject to increase pricing)

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"\$886 E2615
☐ Comfort-Tek Cover ☐ Stretch-Air Cover.....NCO
☐ Comfort Privacy Shield.....\$44

SEAT UPHOLSTERY AND SEATING

☐ Tension Adjustable Seat Upholstery.....STD
☐ Omit Seat Upholstery.....NCO

PERMOBIL SEAT CUSHIONS (Please see information on page 8 regarding product shipment)

Comfort Company Cushions (Cushions over 20" wide subject to increase pricing)

☐ Comfort Curve Cushion Comfort-Tek cover.....\$122 E2601, E2602
☐ Comfort Acta-Embrace™ Zero Elevation Cushion with Comfort-Tek cover.....\$428 E2607, E2608
☐ Comfort M2 Zero Elevation Cushion with Comfort-Tek cover\$474 E2607, E2608

ROHO® Cushions (Cushions over 20" wide subject to increase pricing)

☐ ROHO Single Compartment ☐ Low Profile ☐ Mid Profile ☐ High Profile\$534 E2622, E2623
☐ ROHO Quadro Select ☐ Low Profile ☐ Mid Profile ☐ High Profile\$534 E2624, E2625
☐ ROHO Hybrid Elite ☐ Single Compartment ☐ Dual Compartment.....\$552 E2622, E2623
☐ ROHO Hybrid Select.....\$571 E2624

CUSTOM CONFIGURATION OPTIONS

FRONT WHEELS

3" Casters

☐ 3"x. 1.4" LiteSpeed-7 Billet Aluminum Wheel w/ Soft Roll Tire ☐ Silver Hub **OR** ☐ Black Hub.....\$205

4" Casters

☐ 4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire.....STD
☐ 4"x.75" Plastic Wheel w/Poly Tire (N/A with Slipstream fork).....\$65
☐ 4"x.75" Light-Up Micro-Caster (N/A with Slipstream fork).....\$125
☐ 4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire
☐ Silver Hub **OR** ☐ Black Hub.....\$217

5" Casters

☐ 5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire (N/A with 90° Front Angle **E**).....STD
☐ 5"x1" Plastic Wheel w/Poly Tire (N/A with 90° Front Angle **L** or Slipstream fork).....\$65
☐ 5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire
☐ Silver Hub **OR** ☐ Black Hub.....\$233

6" Casters (N/A with 85° or 90° Front Angle **E** or Slipstream fork)

☐ 6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire (N/A with 85° or 90° Front Angle **E** or Slipstream fork).....STD
☐ 6"x1" Plastic Wheel w/Poly Tire.....\$65
☐ 6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire (N/A with 85° or 90° Front Angle **E** or Slipstream fork)
☐ Silver Hub **OR** ☐ Black Hub.....\$244

FRONT FORKS

☐ TiLite Standard Fork.....STD
☐ TiLite Slipstream Single-Sided Fork.....\$185
☐ Glide Suspension Fork (Max user weight 265 lbs, N/A with ano).....\$526 E1015
Rider's Weight: ☐ 125 lbs. ☐ 125-180 lbs. ☐ 180-265 lbs.

TiLite Z Order Form

Rigid Frame

Effective June 2025 (New items in red)

Page 5 of 8

REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height D limitations)

STANDARD REAR WHEELS

	20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
TiLite Wheels	☐ STD	☐ STD	☐ STD	☐ STD	☐ STD

When ordering 24" wheel with treaded pneumatic or treaded tire with airless inserts and aluminum black or silver anodized handrims and long tabs, the chair will ship with TiLite standard rear wheel. Otherwise, the chair will ship TiLite Shadow rear wheel.

ALTERNATIVE REAR WHEELS

	20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
X-Core Mag®	☐ \$225	☐ \$225	☐ \$225		
Spinergy Spox ☐ Black ☐ Red ☐ Blue ☐ White	☐ \$1,195	☐ \$1,195	☐ \$1,195	☐ \$1,195	☐ \$1,195
Spinergy LX ☐ Black ☐ Red ☐ Blue ☐ White ☐ Bright Orange ☐ Flamingo Pink ☐ Purple		☐ \$1,295	☐ \$1,295	☐ \$1,295	☐ \$1,295
Spinergy Design Lab LX (Available with black spokes only) ☐ Gold Anodize ☐ Green Anodize ☐ Blue Anodize ☐ Red Anodize ☐ Purple Anodize ☐ Orange Anodize		☐ \$1,995	☐ \$1,995	☐ \$1,995	☐ \$1,995

X-Core Mags requires Aluminum Silver Anodized Handrims with long tabs. N/A with Heavy Duty Frame, Quad Release Axles, Spoke Guards

PNEUMATIC REAR TIRES

☐ Treaded (N/A with 25"(559) rear wheels).....	STD	
☐ Treaded with Airless Inserts (N/A with 25"(559) rear wheels).....	\$145	E2213
☐ MBL SpeedLite.....	NCO	
☐ Primo Sentinel.....	\$159	
☐ Schwalbe® Marathon Plus Evolution (N/A with Mag or 22" (501) wheels).....	\$169	
☐ Primo C-763 Knobby (Available on 24"(540) rear wheels.).....	\$129	

SOLID REAR TIRES

☐ Primo Active (Available on 24"(540) and 25"(559) rear wheels.).....	\$160
☐ Primo Orion (Replaces Primo Express) (N/A with 26"(590) Rear Wheels).....	\$160

CAMBER κ (Rear wheel spacing will be set to minimum value allowed by your configuration)

☐ 0° ☐ 2° ☐ 4° ☐ 6°	STD
---	-----

AXLES

☐ Stainless Quick Release.....	STD
☐ Quad Release (N/A with Mag or 24"(540)/ 25" (559) Spinergy Spox or Spinergy LX RearWheels).....	\$169

HANDRIMS

☐ Aluminum – ☐ Silver Ano OR ☐ Black Ano	STD
☐ Plastic Coated – Black.....	\$174
☐ Omit Handrims.....	NCO

The Natural-Fit® – Original, Larger Profile (N/A with 20"(451) Rear Wheels.)

☐ Standard Grip (Requires Long Tabs).....	\$428
☐ Super Grip (Requires Long Tabs).....	\$466
☐ Without Thumb Grip.....	\$390

The Natural-Fit® LT – Reduced Weight, Smaller Profile (N/A with 20"(451) Rear Wheels.)

☐ Standard Grip (Requires Long Tabs).....	\$450
☐ Super Grip (Requires Long Tabs).....	\$514
☐ Without Thumb Grip.....	\$423

The Surge (All-In-One Oval with Gription™ Strip)

☐ The Surge (Larger Profile).....	\$495
☐ The Surge LT (Smaller Profile)	\$466

Q Grip (High-Friction, Anti-Slip Coating)

☐ Q-Grip (N/A with 20" (451) Rear Wheels)	\$495
--	-------

TAB LENGTH

☐ Long Tabs.....	STD
☐ Short Tabs (N/A with Mag Rear Wheels.).....	\$39

SPOKE GUARDS

☐ Clear Lexan (N/A with Mag Rear Wheels.).....	\$179	K0065
---	-------	-------

TiLite Z Order Form
Rigid Frame

Effective June 2025 (New items in red)

Page 6 of 8

WHEEL LOCKS

Aluminum Wheel Locks

- ☐ Push to LockSTD
- ☐ Pull to LockSTD
- ☐ Extension Handles for Aluminum Push/Pull to Lock\$75 E0961

Composite Wheel Locks

- ☐ Compact Composite Push to Lock\$129
- ☐ Compact Composite Pull to Lock\$129
- ☐ Extension Handles for Compact Composite Push/Pull to Lock\$93 E0961
- ☐ Compact Composite Scissor Lock\$158

Additional Wheel Locks

- ☐ Grade Aids, Push to Lock (Requires pneumatic tires)\$395 E0974

ALUMINUM SIDE GUARDS

- Rigid Removable ☐ Small ☐ Standard\$217
- Quad Removable w/ Hole ☐ Small ☐ Standard\$233
- Fold Down Side Guards ☐ Small ☐ Standard\$363

ARMRESTS

- ☐ OmitNCO

Standard Armrest

- ☐ Fixed Height Non-Removable ArmrestSTD
- ☐ Short OR ☐ Tall

Swing Away Tubular

- ☐ Height Adjustable\$244 E0973
- Range: ☐ 6.5" to 9.5" OR ☐ 9.5" to 12.5"
- ☐ Height Adjustable Flip Back (Height Adjustable from 8.5" to 11")\$342 E0973
- ☐ Height Adjustable with Sideguards (Side guard will be NCO)\$425 E0973
- Range: ☐ 6.5" to 9.5" OR ☐ 9.5" to 12.5"

T-Style Armrests (Requires minimum seat depth of 15")

- ☐ Removeable Height Adjustable T Arm\$375 E0973
- ☐ 8" to 12" OR ☐ 10" to 14"
- ☐ Desk Length Pad (10") OR ☐ Full Length Pad (14")

POSITIONING BELTS

- ☐ Velcro Adjustable, Black Nylon Webbing\$61 E0978
- ☐ Bodypoint® 2" Push-Button Buckle (Auto-Style)\$103 E0978
- ☐ Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style, N/A with 19" -20" Rear Widths)\$185 E0978
- ☐ Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style, N/A with 19" -20" Rear Widths)\$196 E0978
- ☐ Bodypoint® EvoFlex® Pelvic Stabilizer Push Button\$201 E0978
- ☐ Q'Strain® WC-19 Wheelchair-Anchored Pelvic Belt (Recommended for use with Transit Tie Down.)\$217 E0978
- ☐ Bodypoint® Optional Band Clamp Mounting Hardware\$70

ACCESSORIES

ANTI-TIPS

- ☐ Rear Aluminum - Black\$174 E0971
- ☐ User Friendly Flip-up\$363 E0971

FUNCTIONAL ACCESSORIES

- ☐ Upper Extremity Support Tray
- ☐ 12" to 14" OR ☐ 13" to 15" Seat Widths A (Requires T Style Armrests.)\$315 E0950
- ☐ Luggage Carriers\$158
- ☐ Crutch Holder - ☐ Left OR ☐ Right\$244
- E2207

SOFT GOOD ACCESSORIES

- ☐ Permobil Medical Necessities Bag\$103
- ☐ Permobil Medical Necessities Seat Pouch\$66
- ☐ Bodypoint® Mobility Bag\$109
- ☐ Neoprene Impact Guards, Black\$195

TiLite Z Order Form

Rigid Frame

Effective June 2025 (New items in red)

Page 7 of 8

GENERAL

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at permobilus.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Straint® is a registered trademark of Q'Straint Inc.

NOTES




QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR

Measurement Guide

- A** Seat Width - Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.
- B** Seat Depth - Measure from front of back post to the front edge of seat sling. Frame bend will start .5" in front of seat sling unless Custom Frame Depth is specified. Custom Frame Depth is measured from front of seat sling to beginning of frame bend.
- C** Front Seat Height - Measure from floor to top of seat tube at front edge of seat sling.
- D** Rear Seat Height - Measure from floor to top of seat tube at back post.
- E** Front Angle - Measure from frame front to floor.
- F** Center of Gravity - Measure from front of back post at the seat tube to the center of the rear axle.
- G** Footrest Width - Measure from inside of front frame tube to inside of opposite front frame tube.
- H** Seat to Footrest - Measure from front edge of seat sling to top rear of footrest.
- I** Seat Back Height - Measure from top of back post to top of seat tube at rear of frame.
- J** Seat Back Angle - Measure angle from front of back post to floor.
- K** Camber – this will affect overall width of chair.



Seating and Positioning Shipping Information

The following Seating and Positioning products will ship with the chair. Cushions are stocked with Comfort Tek Cover and backs are stocked with Stretch Air Cover. All other Seating and Positioning products will be drop-shipped from the manufacturer. (Width x Depth)

Curve Cushion	Elements Cushion	Acta Embrace Cushion	M2 Cushion Zero Elevation	ROHO Single High Profile Cushion	ROHO Quadtro Cushion
18" x 18"	18" x 18"	18" x 18"	18" x 18"	18" x 18"	18" x 18"
16" x 16"	16" x 18"	16" x 18"	18" x 16"	16" x 18"	16" x 16"
16" x 18"	18" x 16"	16" x 16"	20" x 16"	16" x 16"	20" x 16"
			20" x 20"		20" x 18"
			16" x 16"		

Acta Back	Acta Embrace Back	Acta Relief Back	ROHO Agility Mid Active	ROHO Agility Min Active
18" x 12"	18" x 16"	18" x 16"	18" x 10"	18" x 10"
18" x 16"	16" x 16"	16" x 16"	16" x 10"	16" x 10"
16" x 16"		18" x 12"		
16" x 12"		16" x 12"		