## Vicair® Active O2 Low

HCPCS CODING: E2624/E2625 ADJUSTABLE SKIN PROTECTION & POSITIONING



For best results, do not fill out in your browser. Interactive form should be completed using Adobe Reader <u>after</u> saving to your local drive. Then email or print and fax to Customer Support at orders.comfort@permobil.com or 1.406.522.8563.

\*DO NOT SEND PROTECTED HEALTH INFORMATION. IT IS NOT NEEDED TO MAKE THE PRODUCT YOU ARE REQUESTING.\*



## SIZE & STYLE

Choose the size needed for the Active O2 Low.

#### **ACTIVE 02 LOW**

= Nonstock items available at extended lead times and a \$40.00 upcharge

		<b>~</b>	<b></b> \$470	6.00 —	<b></b>	<b>←</b> \$52	5.00 ->
	22"			OCTL-1822	OCTL-2022	OCTL-2222	OCTL-2422
I	20"		OCTL-1620	OCTL-1820	OCTL-2020	OCTL-2220	OCTL-2420
DEPTH	18"	OCTL-1418	OCTL-1618	OCTL-1818	OCTL-2018	OCTL-2218	
Δ	16"	OCTL-1416	OCTL-1616	OCTL-1816	OCTL-2016		
	14"	OCTL-1414					
		14"	16"	18"	20"	22"	24"
				WIE	OTH		

## 2 ACTIVE O2 ACCESSORIES



# REMOVABLE SOLID SEAT PAN & HARDWARE KIT

HCPCS Code: E2231 MSRP \$387.00

The solid seat pan kit includes a slotted aluminum pan and attaching hardware to accommodate 7/8" or 1" tubing. The kit also comes with two different cross bars to accommodate various wheelchair frame types. This will fit both folding and non-folding seat rails. Lateral thigh support hardware can be easily attached.





EXTRA VICAIR® AIR CELLS MSRP \$43.00 / 50 air cell pack



QUESTIONS? PLEASE CONTACT CUSTOMER SUPPORT FOR ASSISTANCE. | 800.736.0925

# BODILINK® LATERAL PELVIC/THIGH SUPPORT HARDWARE & PADS

Please select sizes in the appropriate fabric & style charts below. Contact Customer Support for custom sizes that are not listed below 800.736.0925.



#### PREMIUM PAD SIZE & SHAPE

HCPCS Code: E0953

MSRP \$62.00

- You may make one selection per column (Left, Extra Left, Right, and Extra Right) in section 3.1.
- · Make a fabric selection in 3.2.

\*Length (L) refers to the actual size dimension of the support from top to bottom edge. Depth (D) refers to the actual size dimension from anterior to posterior edge.



LEFT		
BL-LPTSP2Z-3L4D-LH		
BL-LPTSP2Z-3L8D-LH		
BL-LPTSP2Z-3L12D-LH		
BL-LPTSP2Z-5L6D-LH		

EXTRA LEFT	
BL-LPTSP2Z-3L4D-LH	
BL-LPTSP2Z-3L8D-LH	
BL-LPTSP2Z-3L12D-LH	
BL-LPTSP2Z-5L6D-LH	

RIGHT
BL-LPTSP2Z-3L4D-RH
BL-LPTSP2Z-3L8D-RH
BL-LPTSP2Z-3L12D-RH
BL-LPTSP2Z-5L6D-RH

ı	EXTRA RIGHT
	BL-LPTSP2Z-3L4D-RH
	BL-LPTSP2Z-3L8D-RH
	BL-LPTSP2Z-3L12D-RH
	BL-LPTSP2Z-5L6D-RH



### PREMIUM PAD COVER

WITH FOAM INSERT

COVER OPTIONS	MSRP
COMFORT-TEK®	\$0.00/ea.
STRETCH-AIR®	\$0.00/ea.
GLIDEWEAR®	\$13.00/ea.

LEFT		
LPTS-P2-COMFORT-FM		
LPTS-P2-STRETCH-FM		
LPTS-P2-GLIDE-FM		

EXTRA LEFT
LPTS-P2-COMFORT-FM
LPTS-P2-STRETCH-FM
LPTS-P2-GLIDE-FM

RIGHT
LPTS-P2-COMFORT-FM
LPTS-P2-STRETCH-FM
LPTS-P2-GLIDE-FM

ı	EXTRA RIGHT		
	LPTS-P2-COMFORT-FM		
	LPTS-P2-STRETCH-FM		
	LPTS-P2-GLIDE-FM		





#### **STRETCH-AIR®** FOR USER COMFORT & HEAT DISSIPATION



FOR SKIN PROTECTION & SHEAR REDUCTION. NOT FLUID RESISTANT.



#### **BODILINK® LATERAL PELVIC/** THIGH SUPPORT HARDWARE

- You may make one selection per column (Left, Extra Left, Right, and Extra Right) in section 3.3 from group A, B, or C.
- TT long extention arms (section 3.4) will only work with TT style hardware.




### TT HARDWARE, SLOT MOUNT

TT HARDWARE SLOT MOUNT- FIXED









HARDWARE SIZE	MAX CUSHION THICKNESS CLEARANCE	RECOMMENDED CUSHION THICKNESS
SMALL	3.0"	2.0"- 3.0"
MEDIUM	4.5"	3.0"- 4.5"
LARGE	6.5"	4.0"- 6.5"

	SIZE	MSRP		
	SMALL	\$141.00		
FIXED	MEDIUM	\$141.00		
_	LARGE	\$141.00		
OR				
BLE (	SMALL	\$248.00		
REMOVABLE E1028	MEDIUM	\$248.00		
REM	LARGE	\$248.00		

LEFT	
	BL-LPTS-TT1FXSL1-LH
	BL-LPTS-TT1FXSL2-LH
	BL-LPTS-TT1FXSL3-LH
	BL-LPTS-TT1RMSL1-LH
	BL-LPTS-TT1RMSL2-LH
	BL-LPTS-TT1RMSL3-LH

BL-LPTS-TT1FXSL1-LH
BL-LPTS-TT1FXSL2-LH
BL-LPTS-TT1FXSL3-LH
BL-LPTS-TT1RMSL1-LH
BL-LPTS-TT1RMSL2-LH
BL-LPTS-TT1RMSL3-LH

**EXTRA LEFT** 

RIGHT	
BL-LPTS-TT1FXSL1-RH	
BL-LPTS-TT1FXSL2-RH	
BL-LPTS-TT1FXSL3-RH	
BL-LPTS-TT1RMSL1-RH	
BL-LPTS-TT1RMSL2-RH	
BL-LPTS-TT1RMSL3-RH	

	BL-LPTS-TT1FXSL1-RH
	BL-LPTS-TT1FXSL2-RH
	BL-LPTS-TT1FXSL3-RH
ľ	
	BL-LPTS-TT1RMSL1-RH
	BL-LPTS-TT1RMSL2-RH
	BL-LPTS-TT1RMSL3-RH

**EXTRA RIGHT** 

VICAIR-ACTIVE O2 LOW-ORDER REV111920


#### B. TT HARDWARE, POWER MOUNT











HARDWARE SIZE	MAX CUSHION THICKNESS CLEARANCE	RECOMMENDED CUSHION THICKNESS
SMALL	3.5"	2.0"- 3.5"
MEDIUM	5.5"	3.0"- 5.5"
LARGE	7.5"	4.0"- 7.5"

	SIZE	MSRP	
FIXED	SMALL	\$185.00	
	MEDIUM	\$185.00	
	LARGE	\$185.00	
OR			
REMOVABLE E1028	SMALL	\$285.00	
	MEDIUM	\$285.00	
	LARGE	\$285.00	

	LEFT
	BL-LPTS-TT1FXPWL1-LH
	BL-LPTS-TT1FXPWL2-LH
	BL-LPTS-TT1FXPWL3-LH
•	
	BL-LPTS-TT1RMPWL1-LH
	BL-LPTS-TT1RMPWL2-LH
	BL-LPTS-TT1RMPWL3-LH

EXTRA LEFT		
	BL-LPTS-TT1FXPWL1-LH	
	BL-LPTS-TT1FXPWL2-LH	
	BL-LPTS-TT1FXPWL3-LH	
•		
	BL-LPTS-TT1RMPWL1-LH	
	BL-LPTS-TT1RMPWL2-LH	
	BL-LPTS-TT1RMPWL3-LH	

LEFT	RIGHT	
T1FXPWL1-LH	BL-LPTS-TT1FXPWL1-RH	
T1FXPWL2-LH	BL-LPTS-TT1FXPWL2-RH	
T1FXPWL3-LH	BL-LPTS-TT1FXPWL3-RH	
T1RMPWL1-LH	BL-LPTS-TT1RMPWL1-RH	
T1RMPWL2-LH	BL-LPTS-TT1RMPWL2-RH	
T1RMPWL3-LH	BL-LPTS-TT1RMPWL3-RH	

EXTRA RIGHT		
	BL-LPTS-TT1FXPWL1-RH	
	BL-LPTS-TT1FXPWL2-RH	
	BL-LPTS-TT1FXPWL3-RH	
•		
	BL-LPTS-TT1RMPWL1-RH	
	BL-LPTS-TT1RMPWL2-RH	
	BL-LPTS-TT1RMPWL3-RH	

#### TT HARDWARE, POWER MOUNT OPTION

**NOTE:** If choosing the Power Mount Hardware, choose the appropriate power mount type below.

POWER MOUNT	MSRP
QUANTUM	add \$0.00
PERMOBIL	add \$0.00
QUICKIE	add \$0.00
ROVI	add \$0.00
AVID REHAB	add \$0.00

PAI	RT NUMBER
	LPTS-TT1-PW1
	LPTS-TT1-PW2
	LPTS-TT1-PW3
	LPTS-TT1-PW4
	LPTS-TT1-PW5



### TT LONG EXTENSION ARM 4"

MSRP \$37.00/arm



LEFT	
	BL-LPTS-TT1LARM4

EXTRA LEFT	
	BL-LPTS-TT1LARM4

RIGHT	
	BL-LPTS-TT1LARM4

EXTRA RIGHT		
	BL-LPTS-TT1LARM4	

# BODILINK® MEDIAL KNEE/THIGH SUPPORT HARDWARE & PADS



#### **PAD SIZE & SHAPE**

HCPCS CODING: E0957 MSRP \$130.00

	STYLE	SIZE	PART NUMBER	
<b>^</b>		2.5"W x 3.5"D	BL-MKTSP1-3W4D	
w	WEDGE	3.5"W x 5"D	BL-MKTSP1-4W5D	
<b>V</b>		4"W x 6"D	BL-MKTSP1-4W6D	
$\uparrow$	w °°	2.5"W x 3.5"D	BL-MKTSP2-3W4D	
w		3.5"W x 5"D	BL-MKTSP2-4W5D	
<b>V</b>		5"W x 6"D	BL-MKTSP2-5W6D	
,	OVAL D		·	

- You may make one selection per column.
- Choose between wedge or oval in 4.1.
- Choose a cover material in 4.2.
- Choose hardware size in 4.3.



### **PAD COVER MATERIAL**

GlideWear

FOR SKIN PROTECTION & SHEAR REDUCTION. NOT FLUID RESISTANT.

COMFORT-TEK®



FOR FLUID PROTECTION & AN EASILY CLEANED SURFACE.

STRETCH-AIR®




FOR USER COMFORT & HEAT DISSIPATION.

GlideWear® (MKTSP-GLIDE) (add \$15.00)
Comfort-Tek® (MKTSP-COMFORT)
Stretch-Air® (MKTSP-STRETCH)


VICAIR-ACTIVE O2 LOW-ORDER REV111920



**HCPCS CODING:** E1028 MSRP \$180.00




SIZE	PART NUMBER	
SMALL	BL-MKTS-ST1SASL1	
MEDIUM	BL-MKTS-ST1SASL2	
LARGE	BL-MKTS-ST1SASL3	





Inferior Thigh Thickness



SMALL	MEDIUM	LARGE

SIZE	RECOMMENDED INFERIOR THIGH THICKNESS
Small	1.00"- 3.00"
Medium	2.00"- 5.00"
Large	3.00"- 7.00"

VICAIR-ACTIVE O2 LOW-ORDER REV111920

SUBMIT BY EMAIL