



TiLite 2GX - Order Form



Dealer Name*:		Account #*:
City, State, Zip:		
Phone:	Email:	ATP:
Rep Name:	P.O. Number:	Quote Number:
Client Reference*:	Client Weight:	
Special considerations that need to be addressed here (e.g., diagnosis):		

*Required field

Frame Style

- ☐ G2FS1 **2GX Titanium Folding Manual Wheelchair** Full adjustability, 265lb weight limit. Light weight and high strength. Not available with frame transit tie-down option. K0005 \$ 4,692.00

Frame Options

- ☐ XFS1 **Swing Away Front Frame** Standard
- ☐ XFS2 **Fixed Front Frame** No charge
- ☐ XFS3 **Hemi Height Frame** Available with a front seat height of 13.5"-20". Not with depth adjustable back/frame, 12"-13" seat depth or heavy duty frame. No charge
- ☐ XHDR1 **Heavy Duty Frame** 350lb weight limit. Includes tension adjustable by straps back upholstery. Available with 18"-22" seat widths and 16"-20" seat depths. Requires LiteSpeed billet aluminum front wheels and 0° camber. Not available with 90° front end or single sided forks. No charge
- ☐ XPWR1 **Power Adaptable Frame** Requires LiteSpeed billet aluminum front wheels and 0° camber. Not available with single-sided forks. No charge
- ☐ XRFR1 **Reinforced Frame** No charge

Transit Options

- ☐ XWTB **Transport Bracket Option** Consists of four tie-down brackets. Unoccupied use only. Not available when swing away frame and heavy duty frame are both selected. Not with swing away tubular armrests or swing away flip back tubular armrests. Requires a minimum seat depth of 14", or 15" with hemi height frame. E1022 \$ 379.00

SmartDrive Power Assist System

- ☐ XMAXSD1 **Permobil SmartDrive MX2+** Not available with rear seat heights less than 13.5". Maximum user weight capacity of 331 lbs. E0986 \$ 7,529.00
- ☐ XMAXSD18 **Switch Control Buttons** No charge
- ☐ XMAXSD8 **SwitchControl Button with Mono Jack Option** No charge
- ☐ XMAXSD19 **Buddy Button** For use with SwitchControl button with mono jack option. No charge
- ☐ XMAXSD4 **Carrying Bag** \$ 84.00
- ☐ XMAXSD10 **Universal Padded Handle** \$ 22.00

Frame Finishes

☐ XFF1 **Satin Finish** No charge

☐ XSCBT1 Buffing Pads for Satin Frame \$ 22.00

☐ XFF **Powder Coated Frame** *Select below.* No charge

Select color:

☐ XFF8 High Gloss Black	☐ XFF10 Red Pepper	☐ XFF29 Neon Pink	☐ XFF14 Ocean Blue Metallic
☐ XFF15 Tangerine Metallic	☐ XFF30 Acid Green Metallic	☐ XFF17 Midnight Blue Pearl	☐ XFF18 White Pearl
☐ XFF19 Black Lava	☐ XFF27 Titanium Grey	☐ XFF59 Pomegranate Red	☐ XFF62 Blue Mist
☐ XFF57 Indigo Purple	☐ XFF40 Hyper Blue	☐ XFF41 Flame Orange	☐ XFF42 OD Green
☐ XFF60 Silver Streak	☐ XFF63 Lavendar Dusk	☐ XFF58 Black Galaxy	☐ XFF61 Seafoam Green
☐ XFF64 Daffodil Yellow			

☐ XFF3 **Tattoos** *Tattoos on front frame tubes. Finished in satin. Tattoos will be trimmed and adjusted to fit your frame.* \$ 363.00

☐ XFF65 Liquid ☐ XFF66 Bloom ☐ XFF67 Rebel ☐ XFF68 Vector



Liquid



Bloom



Rebel



Vector

Frame Measurements

Measurements shown below are our most popular selections. Not seeing what you want? Contact customer service or view our Creative Mobility Division (CMD) catalog [here](#) to view additional options. Please refer to Reference Diagrams at end of form for dimensional information.

Seat Width 'A' | *Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.*

XSW1: ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" Standard

XSW2: ☐ 20" ☐ 21" ☐ 22" E2201 \$ 536.00

21"-22" seat widths only available with heavy duty frame.

Seat Depth 'B' | *Measure from 0.5" in front of back posts to the front edge of seat sling. 12"-14" seat depths not with 25"-26" rear wheels.*

XSD1: ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" Standard

XSD2: ☐ 20" E2203 \$ 568.00

XSD3: ☐ Custom Frame Depth. Specify: No charge

Frame depth longer than seat depth. Use 1" increments with a maximum of 2". A custom frame depth with a maximum of 20" is \$546.

Front Seat Height 'C' | *Measure from floor to top of seat tube at front edge of seat sling. 13.5"-14" front seat height requires hemi height frame. 20.5"-21.5" front seat height not with hemi height frame.*

XFSH1: ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21"

Front Wheel	8" Front Wheel	19.5"-21"	18.5"-19.5"
	6" Front Wheel	17.5"-20"	16.5"-19"
	6" Front Wheel w/Glide	19"	17.5"-18"
	5" Front Wheel	16.5"-19.5"	16"-18.5"
	5" Front Wheel w/Glide	17.5"-18.5"	16.5"-18"
	4" Front Wheel	15.5"-19"	14.5"-18"
	4" Front Wheel w/Glide	17"-18"	15.5"-17"
	3" Front Wheel	14.5"-18.5"	13.5"-17.5"
	3" Front Wheel w/Glide	16.5"-17.5"	15"-17"

Std Front Seat Height

Hemi Front Seat Height

Front Seat Height Range

Rear Seat Height 'D' | *Measure from floor to top of seat tube at back post.*

XRSH1: ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20"

Rear Wheel	26"	15"-20"
	25"	14"-19"
	24"	14"-19"
	22"	13"-18"

Rear Seat Height Range

Back Supports

Comfort Back Supports | Backs over 20" wide or tall subject to increased pricing.

☐ SBRS19	Comfort Unite 1 Back Replaces Comfort Acta-Back. Inc 2-point quick release Compass hardware. Select shell height and cover.	E2613	\$ 712
	☐ 12" ☐ 14" ☐ 16" ☐ 18"		No charge
	☐ 20"		\$ 100
	☐ Comfort-Tek Cover ☐ Diamond Spacer Cover		
	☐ Privacy Shield		\$ 52
☐ SBRS21	Comfort Unite 1 5" Wing Back Replaces Comfort Acta-Back Deep w/5" Wing. Inc 2-point quick release Compass hardware. Select shell height and cover.	E2620	\$ 845
	☐ 12" ☐ 14" ☐ 16" ☐ 18"		No charge
	☐ 20"		\$ 100
	☐ Comfort-Tek Cover ☐ Diamond Spacer Cover		
	☐ Privacy Shield		\$ 52
☐ SBRS20	Comfort Unite 2 Back Replaces Comfort Acta-Back Deep w/3" Wing. Inc 2-point quick release Compass hardware. Select shell height and cover.	E2615	\$ 820
	☐ 12" ☐ 14" ☐ 16" ☐ 18"		No charge
	☐ 20"		\$ 100
	☐ Comfort-Tek Cover ☐ Diamond Spacer Cover		
	☐ Privacy Shield		\$ 52
☐ SBRS22	Comfort Embrace Back Replaces Comfort Acta-Embrace Back. Inc 2-point quick release Compass hardware. Select shell height.	E2613	\$ 700
	☐ 10" ☐ 12" ☐ 14" ☐ 16"		
	☐ Comfort-Tek Cover ☐ Diamond Spacer Cover		
	☐ Privacy Shield		\$ 52
☐ SBRS23	Comfort Relief Back Replaces Comfort Acta-Relief Back. Inc 2-point quick release Compass hardware. Select shell height.	E2615	\$ 886
	☐ 10" ☐ 12" ☐ 14" ☐ 16"		
	☐ Comfort-Tek Cover ☐ Diamond Spacer Cover		
	☐ Privacy Shield		\$ 52

Backrest

Back Type

☐ XSB13	Fixed Angle, Height Adjustable Back	Standard
☐ XLB1	Folding Stabilizer Bar Required with transit tie-down option unless depth adjustable backrest is selected. Includes tension adjustable by strap back upholstery.	\$ 423.00
☐ XSB14	Fixed Angle, Height Adjustable Back w/8° Bend Not available with heavy duty frame.	No charge
☐ XSB18	Folding, Angle and Height Adjustable Back Not available with heavy duty frame.	\$ 250.00
☐ XSB16	Depth Adjustable Back Frame Not available with hemi height frame, heavy duty frame, 19"-20" seat depths, swing away tubular or swing away flip back armrests. Seat depths 12"-17" allow 3" of rearward adjustment. Seat depth of 18" allows 2" of rearward adjustment. Requires tension adjustable seat upholstery.	\$ 244.00

Push Handles

☐ XSB6	Integrated Push Handles	No charge
☐ XSB7	Integrated Fold-Down Push Handles Requires tension adjustable by straps back upholstery. Not available with heavy duty frame.	\$ 125.00
☐ XSB9	Bolt-On Height Adjustable Push Handles Requires tension adjustable by straps back upholstery. Not available with heavy duty frame, folding angle and height adjustable back, swing away tubular or swing away flip back tubular armrests.	\$ 363.00

Fixed Angle Adjustable Back Height 'E' | 11"-15" back height not with heavy duty frame, swing away tubular armrests 12"-15", bolt-on height adjustable push handles or folding stabilizer bar. Set at 1" increments.

Push Handles	None	☐ 11"-15"	☐ 16"-20"	Set at:
	Integrated	☐ 11"-15"	☐ 16"-20"	Set at:
	Fold-Down	☐ 11"-15"	☐ 16"-20"	Set at:
	Bolt-On Height Adj	N/A	☐ 16"-20"	Set at:
		Short	Tall	Specify
Fixed Angle Adj Back Height				

Folding Angle Adjustable Back Height 'E' | Not with folding stabilizer bar or heavy duty frame. Set at 1" increments.

Push Handles	None	☐ 11"-13"	☐ 13.5"-16.5"	☐ 16.5"-19.5"	Set at:
	Integrated	☐ 12.5"-14.5"	☐ 13.5"-16.5"	☐ 17"-19.5"	Set at:
	Fold-Down	☐ 15"-17"	☐ 16"-19"	☐ 18.5"-20.5"	Set at:
		Short	Medium	Tall	Specify
Folding Angle Adj Back Height					

Front Rigging

Front Angle 'H' | Front angle is measured from inside of front frame tube to top of bottom frame tube.

- ☐ XFA3 **Swing Away Front End** 85° not with 6" front wheels. 90° requires aluminum footplate and not with heavy duty frame or 5"-6" front wheels. No charge
- ☐ 70° ☐ 80° ☐ 85° ☐ 90° ☐ XFA0 Omit Hangers
- ☐ XFA6 **Fixed Front End** 85° not with 6" front wheels. Standard
- ☐ 70° ☐ 80° ☐ 85°

Select one seat to footrest length and one footrest type.

Standard Footrests | For front angles 70°-85°. Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than front seat height. 70° fixed front has a minimum seat to footrest of 14.5".

Seat to Footrest 'E': ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20"

- ☐ XFTR5 **Flip-Up Footrest** No charge
- ☐ XFTR9 **Depth Adjustable Flip-Up Footrest** K0045 \$ 120.00
- ☐ XFTR10 **Depth and Angle Adjustable Flip-Up Footrest** Available for 15"-20" seat to footrest only. K0040, K0045 \$ 293.00
- ☐ XFTR11 **Angle Adjustable Flip-Up Footrest** Available for 15"-20" seat to footrest only. K0040, K0045 \$ 185.00
- ☐ XHLP1 **Heel Loops** Not with seat width ≤ 13" or one piece footrests. E0951 \$ 55.00

High Mount Footrests | Not with luggage carriers or impact guards.

Seat to Footrest 'E': ☐ 7" ☐ 7.5" ☐ 8" ☐ 8.5" ☐ 9" ☐ 9.5" ☐ 10" ☐ 10.5" ☐ 11" ☐ 11.5" ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5"

- ☐ XFTR20 **High Mount Flip-Up** K0037 \$ 202.00
- ☐ XFTR15 **Angle Adjustable High Mount Flip-Up Footrest** K0040, K0037 \$ 352.00
- ☐ XFTR16 **Depth and Angle Adjustable High Mount Flip-Up Footrest** K0040, K0037 \$ 406.00
- ☐ XHLP1 **Heel Loops** Not with seat width ≤ 13" or one piece footrests. E0951 \$ 55.00

90° Footrests

Seat to Footrest 'E': ☐ 7" ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"

- ☐ XFTR5 **Flip-Up Footrest** No charge
- ☐ XFTR9 **Depth Adjustable Flip-Up Footrest** K0045 \$ 120.00
- ☐ XFTR10 **Depth and Angle Adjustable Flip-Up Footrest** Available for 15"-20" seat to footrest only. K0040, K0045 \$ 293.00
- ☐ XHLP1 **Heel Loops** Not with seat width ≤ 13" or one piece footrests. E0951 \$ 55.00

Elevating Legrests

Seat to Footrest 'E': ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐ 20" ☐ 21"

- ☐ XFA5 **Depth Adjustable Flip-Up Footrest, Single** E0990 \$ 185.00
- ☐ Left ☐ Right
- ☐ XFTR17 **Depth Adjustable Flip-Up Footrest, Pair** \$ 363.00
- ☐ XHLP1 **Heel Loops** Not with seat width ≤ 13" or one piece footrests. E0951 \$ 55.00

Camber Options

Camber 'K'

XCBR8: ☐ 0° ☐ 2° ☐ 4° ☐ 6°

E2601 Standard

Front Wheels

See front seat height ranges chart for front seat height 'C' limitations.

☐	XFW24	3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$ 206.00
		☐ Silver Hub ☐ Black Hub	
☐	XFW14	4"x0.75" Plastic Wheel w/Poly Tire <i>Not available with power adaptable frame or heavy duty frame.</i>	No charge
☐	XFW3	4"x0.75" Light-Up Micro Caster <i>Not available with power adaptable frame or heavy duty frame.</i>	\$ 88.00
☐	XFW26	4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with power assist, power add-on system or heavy duty frame.</i>	No charge
☐	XFW27	4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$ 217.00
		☐ Silver Hub ☐ Black Hub	
☐	XFW6	5"x1" Plastic Wheel w/Poly Tire <i>Not available with power adaptable frame, heavy duty frame or with 90° front angle.</i>	No charge
☐	XFW28	5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with power adaptable frame, heavy duty frame or with 90° front angle.</i>	No charge
☐	XFW30	5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not available with 90° front angle.</i>	\$ 233.00
		☐ Silver Hub ☐ Black Hub	
☐	XFW9	6"x1" Plastic Wheel w/Poly Tire <i>Not available with power adaptable frame, heavy duty frame or 85°-90° front angle.</i>	No charge
☐	XFW32	6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with power adaptable frame, heavy duty frame or 85°-90° front angle.</i>	No charge
☐	XFW33	6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not available with 85°-90° front angle.</i>	\$ 244.00
		☐ Silver Hub ☐ Black Hub	
☐	XFW16	8"x1" Plastic Wheel w/Poly Tire <i>Not available with power assist frame, power add-on frame, heavy duty frame or 12"-14" seat depths. Only with 70° front angle.</i>	No charge

Front Forks

☐	XFK1	TiLite Standard Fork	No charge
☐	XGL1	Glide Suspension Fork <i>Black anodize only. Not available with 8" front wheels or heavy duty frame. Select rider weight.</i>	E1015 \$ 526.00
		☐ 125lbs ☐ 125-180lbs ☐ 180-265lbs	

Rear Wheels

See rear seat height ranges chart for rear seat height 'D' limitations.

☐	XOAD1	One Arm Drive System <i>Available with 24" mag wheels and standard aluminum handrim only. Rear seat height, camber or axle restrictions may apply.</i>	E0958 \$ 1,298.00
		☐ Install on Left ☐ Install on Right	
☐	XRW14	TiLite SHADOW Rear Wheel <i>Select size. 25" and 26" wheels require 15" minimum seat depth.</i>	Standard
		☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)	
☐	XRW7	Spinergy SPOX Rear Wheel <i>Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth.</i>	\$ 1,195
		☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)	
		☐ Black ☐ Red ☐ Blue ☐ White ☐ Green	
☐	XRW8	Spinergy LX Rear Wheel <i>Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth. Only available with heavy duty frame when 0° or 2° camber is selected.</i>	\$ 1,295
		☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)	
		☐ Black ☐ Red ☐ Blue ☐ White ☐ Orange ☐ Pink ☐ Purple ☐ Green	
☐	ZRW9	Spinergy Design Lab LX Rear Wheel <i>Black spokes only. Select size and anodizing color. 25" and 26" wheels require 15" minimum seat depth.</i>	\$ 1,995
		☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)	
		☐ Gold ☐ Green ☐ Blue ☐ Red ☐ Purple ☐ Orange	
☐	XRW16	Spinergy Carbon Blade Rear Wheel <i>Select size. Not available with heavy duty frame. 25" wheels require 15" minimum seat depth.</i>	\$ 1,795
		☐ 24" (540) ☐ 25" (559)	
☐	XRW17	X-Core Mag Rear Wheel <i>Select size. Not available with heavy duty frame, short tabs, quad release axles, spoke guards or omit handrims options. Requires aluminum silver anodized handrims.</i>	\$ 225
		☐ 22" (501) ☐ 24" (540)	

Axles

☐	XAXL1	Stainless Quick Release Axle		Standard
☐	XAXL2	Titanium Quick Release Axle		\$ 277.00
☐	XAXL3	Quad Release Axle <i>Not available with Mag or 26" (590) Spinergy Spox or LX rear wheels.</i>		\$ 169.00

Pneumatic Rear Tires

☐	XRTR1	Treaded Tire 22"(501), 24"(540), 26"(590) <i>Not available with 25" rear wheels.</i>		Standard
☐	XRTR2	Airless Inserts <i>Not available with Spinergy Carbon Blade rear wheels.</i>	E2213	\$ 98.00
☐	XRTR31	MBL SpeedLite Tire 22"(501), 24"(540), 25"(559), 26"(590)		No charge
☐	XRTR30	Primo Sentinel Tire 24"(540), 25"(559), 26"(590)		\$ 136.00
☐	XRTR8	Schwalbe Marathon Plus Evolution Tire 24"(540), 25"(559), 26"(590)		\$ 169
☐	XRTR32	Primo C-764 Knobby Tire 24"(540)		\$ 104.00

Solid Rear Tires

☐	XRTR22	Primo Orion Tire 22"(501), 24"(540), 25"(559) <i>Not available with Spinergy Carbon Blade rear wheels.</i>		\$ 120.00
☐	XRTR23	Primo Active Tire 24"(540), 25"(559) <i>Not available with Spinergy Carbon Blade rear wheels.</i>	E2213	\$ 147.00

Handrims

☐	XRIM1	Aluminum Handrim, Silver Anodized		Standard
☐	XRIM21	Aluminum Handrim, Black Anodized		No charge
☐	XRIM12	Titanium Handrim, Satin Finish <i>Not with 22"(501) rear wheels or 24"(540) knobby tires.</i>		\$ 509.00
☐	XRIM20	Q-Grip Handrim <i>Exclusive high-friction, anti-slip coating.</i>		\$ 358.00
☐	XRIM5	Plastic Coated Handrim, Black		\$ 174.00
☐	XRIM10	The Surge Handrim <i>Large profile. All-In-One oval with Gription strip.</i>		\$ 433.00
☐	XRIM11	The Surge LT Handrim <i>Small profile. All-In-One oval with Gription strip.</i>		\$ 466.00
☐	XRIM14	The Natural-Fit Handrim, Standard Grip <i>Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 428.00
☐	XRIM15	The Natural-Fit Handrim, Super Grip <i>Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 466.00
☐	XRIM16	The Natural-Fit Handrim, Without Thumb Grip <i>Original, large profile.</i>		\$ 390.00
☐	XRIM17	The Natural-Fit LT Handrim, Standard Grip <i>Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 440.00
☐	XRIM18	The Natural-Fit LT Handrim, Super Grip <i>Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 514.00
☐	XRIM19	The Natural-Fit LT Handrim, Without Thumb Grip <i>Reduced weight, small profile.</i>		\$ 423.00
☐	XRIM22	Omit Handrims		No charge

Handrim Tab Length

☐	XRIM7	Long Tabs		Standard
☐	XRIM8	Short Tabs <i>Not available with Mag rear wheels or The Natural-Fit handrims with thumb grip.</i>		\$ 28.00

Wheel Locks

☐	XWLK7	Compact Composite Wheel Lock <i>Depending on configuration there may be compatibility issues with 12"-13" seat depths. Uni-Lock will be substituted at no charge.</i>		No charge
☐	XWLK8	Compact Composite Pull to Lock <i>Depending on configuration there may be compatibility issues with 12"-13" seat depths. Uni-Lock will be substituted at no charge.</i>		No charge
☐	XWLK10	Extension Handles for Compact Composite Push/Pull to Lock	E0961	\$ 93.00
☐	XWLK11	MBL Composite Scissor Lock <i>Not available with swing away frame.</i>		\$ 158.00
☐	XWLK12	MBL Performance Aluminum Scissor Lock <i>Not available with swing away frame.</i>		\$ 175
☐	XWLK3	Uni-Lock		\$ 179.00
		☐ Push to Lock ☐ Pull to Lock		
☐	XWLK4	Extension Handles for Uni-Lock	E0961	\$ 93.00
☐	XWLK15	Grade Aids, Push to Lock <i>Not available with solid rear tires.</i>	E0974	\$ 190.00

Back Upholstery

☐ XBUP1	Padded Nylon Back Upholstery <i>Black only. Not available with heavy duty frame, folding back, folding stabilizer bar, depth adjustable frame, fold-down push handles or ultrasuede backrest stripe.</i>		Standard
☐ XBUP3	Tension Adjustable Back Upholstery By Straps, Padded Nylon <i>Included with heavy duty frame. Required for aluminum fendered side guards and fold-down push handles.</i>	E2611	\$ 369.00
	☐ XUPC3 Ultrasuede Backrest Stripe, Black		\$ 115.00
☐ XBUP5	Tension Adjustable Back Upholstery By Straps, AIR Breathable <i>Included with heavy duty frame. Required for aluminum fendered side guards and fold-down push handles. Not with ultrasuede backrest stripe.</i>	E2611	\$ 369.00
☐ XBUP0	Omit Back Upholstery		No charge

Seat Upholstery

☐ XSUP4	Tension Adjustable Seat Upholstery, Black		Standard
☐ XSUP0	Omit Seat Upholstery		No charge

Seat Cushions

☐ XCSH2	Comfort Curve Cushion <i>Stretch-Air cover.</i>	E2601	\$ 122.00
☐ XCSH23	Comfort Acta-Embrace Comfort-Tek Zero Elevation Cushion	E2607	\$ 428.00
☐ XCSH15	ROHO Hybrid Elite SR Single Compartment Cushion <i>Fits seat widths 14"-19"</i>	E2622	\$ 552.00
☐ XCSH14	ROHO Quadtro Select Cushion <i>Fits seat widths 14"-19".</i>	E2624	\$ 534.00
	☐ Low Profile ☐ Mid Profile ☐ High Profile		
☐ XCSH16	ROHO Hybrid Select Cushion <i>Fits seat widths 14"-18", and 20"</i>	E2624	\$ 571.00
☐ XCSH17	Vicar Active O2 Low <i>Fits even seat widths 14"-20". Non-stock upcharges may apply.</i>	E2624	\$ 542.00

Side Guards

☐ XSDG2	Aluminum Side Guards, Rigid Removable	E0956, E1034	\$ 217.00
	☐ Small ☐ Standard		
☐ XSDG3	Aluminum Side Guards, Quad Removable w/Hole	E0956, E1034	\$ 233.00
	☐ Small ☐ Standard		
☐ XSDG7	Carbon Side Guards, Rigid Removable <i>2x2 Twill Weave</i>	E0956, E1034	\$ 336.00
	☐ Small ☐ Standard		
☐ XSDG8	Carbon Side Guards, Quad Removable w/Hole <i>2x2 Twill Weave</i>	E0956, E1034	\$ 347.00
	☐ Small ☐ Standard		
☐ XSDG6	Aluminum Side Guards, Fendered Adjustable <i>Requires tension adjustable by straps back upholstery. Not available with angle adjustable fold-down back.</i>	E0956, E1034	\$ 363.00
☐ CMD007	Custom Offset Profiled Side Guards (Removable or Fold-Down) <i>Please specify desired offset to nearest 0.5" and side guards type in the Notes section of this order form. Side guard material must be selected. Some limitations may apply.</i>		CMD

Armrests

☐ XARM1	Swing Away Tubular Armrests <i>Select range. Not available with heavy duty frame, folding stabilizer bar with short or medium back height ranges, or depth adjustable back.</i>	E0973	\$ 244.00
	☐ 9"-12" ☐ 12"-15"		
☐ XARM4	Swing Away Flip-Back Tubular Armrests <i>Not available with heavy duty frame, folding stabilizer bar or depth adjustable back.</i>	E0973	\$ 342.00
☐ XARM2	Removable Height Adjustable Desk Arm with Rigid Side Guard <i>Select range and pad size. Requires 14" minimum seat depth.</i>	E0973	\$ 363.00
	☐ 9"-12" ☐ 12"-15"		
	☐ Desk Length (10") ☐ Full Length (14")		
☐ XARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop <i>Select range and pad size. Requires 16" minimum seat depth.</i>	E0973	\$ 406.00
	☐ 9"-12" ☐ 12"-15"		
	☐ Desk Length (10") ☐ Full Length (14")		
☐ XARM5	Armrests Covered in Black Ultrasuede <i>Not available with desk arms.</i>		\$ 120.00

Positioning Belts

☐	XBLT1	Velcro Adjustable Belt, Black Nylon Webbing	E0978	\$ 61.00
☐	XBLT2	Bodypoint 2" Belt, Push-Button Buckle (Auto Style)	E0978	\$ 103.00
☐	XBLT4	Bodypoint 1" Padded Belt, Push-Button Buckle (Auto Style) <i>Not available with 19"-22" rear seat widths.</i>	E0978	\$ 196.00
☐	XBLT5	Bodypoint 1.5" Padded Belt, Push-Button Buckle (Auto Style) <i>Not available with 19"-22" rear seat widths.</i>	E0978	\$ 185.00
☐	XBLT10	Bodypoint Optional Band Clamp Mounting Hardware	K0108	\$ 22.00
☐	XBLT16	Bodypoint EvoFlex Pelvic Stabilizer, Push Button	E0978	\$ 201.00

Calf Straps

☐	XCLF2	Velcro Adjustable ☐ Small (9"-12" Footrest Widths) ☐ Medium (11"-15" Footrest Widths) ☐ Large (14"-17" Footrest Widths) ☐ Xtra Large (16"-19" Footrest Widths)	K0038	\$ 44.00
☐	XCLF3	Bodypoint® Padded Velcro Adjustable ☐ Medium (12"-15" Footrest Widths) ☐ Large (16"-20" Footrest Widths)	K0038	\$ 109.00
☐	XCLF1	Buckle Adjustable, Black Nylon Webbing	K0038	\$ 34.00

Residual Limb Supports

☐	XRSL1	Residual Limb Support, 8"W x 10"L, Single ☐ Left ☐ Right	E1020	\$ 428.00
☐	XRSL2	Residual Limb Support, 8"W x 10"L, Pair	E1020	\$ 855.00

Anti-Tips

☐	XTIP1	Rear Anti-Tips, Aluminum Black	E0971	\$ 174.00
☐	XTIP2	Rear Anti-Tips, User Friendly Flip-Up	E0971	\$ 363.00
☐	XTIP9	Tip Assist	E0971	\$ 115.00

Accessories

☐	XGRD1	Clear Spoke Guards <i>Not available with mag rear wheels.</i>	K0065	\$ 147.00
☐	XBCK3	Permobil Medical Necessities Backpack		\$ 103.00
☐	XPCH3	Permobil Medical Necessities Seat Pouch		\$ 66.00
☐	XPCH4	Bodypoint Mobility Bag		\$ 109.00
☐	XUEST1	Upper Extremity Support Tray <i>Requires armrests.</i> ☐ 12"-14" Seat Widths ☐ 13"-15" Seat Widths	E0950	\$ 315.00
☐	XLUG1	Luggage Carriers <i>Not available with high mount swing away footrests or 1-piece flip-up footrest.</i>		\$ 158.00
☐	XMPCT2	Black Neoprene Impact Guards with Ultrasuede Stripe <i>Not available with high mount footrests, 90° hangers or elevating legrests.</i>		\$ 115.00
☐	XCH1	Crutch Holder	E2207	\$ 244.00
☐	XAFP1	4" Military Service Patch ☐ Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard		\$ 115.00

General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

TiFit Reference Diagrams

