


		Permobil M3 Tilt - Corpus® 6 km/h								
	Permobil Account #:			Permobil M3- (706826)						
	Permobil PO #			Standard Options						
	Permobil Delivery Address	1	☐	● R-Net EL LED	● 65 Ah AGM Gel Batteries					
		2	☐	● MyPermobil Fleet Management	● Corpus 3G® Contour Back					
		3	☐	● Solid Drive Wheels & Castors	● Manual Adjust Backrest 135°					
4		☐	● Fixed Seat Height	● Manual Legrest Mid Mount 150°						
RRP	£8665		● 50° Posterior Tilt	● Footplate Narrow 16cm						
Motability Weekly Rental	£51.96		● 420 mm - 570 mm Seat Width	● Bodilink® Headrest Medium						
Max User Weight	150kg	Client Weight			kg					
			● 370 mm - 570 mm Seath Depth	● Reflector Kit						
			● Stretch Air Ergo Seat Cushion	● Charger Voltpro AGM/Gel IP44						
Configurable Options										
Accent Colour										
	Aqua	☐		Sky	☐					
	Violet	☐		Flamingo	☐					
	Canary	☐		Spring	☐					
	Arctic	☐		Sunset	☐					
	Midnight	☐								
Seat Size										
Size (mm)		Depth								
		370	395	420	445	470	495	520	545	570
Width	420	☐	☐	☐	☐	☐	☐	☐	☐	☐
	470	☐	☐	☐	☐	☐	☐	☐	☐	☐
	520	☐	☐	☐	☐	☐	☐	☐	☐	☐
	570	☐	☐	☐	☐	☐	☐	☐	☐	☐
Corpus 3G Backrest										
Size (mm)		Height								
		470	520	545	570	595	620	645	670	
Width	360	☐	☐	☐	☐	☐	☐	☐	☐	
	410	☐	☐	☐	☐	☐	☐	☐	☐	
	460	☐	☐	☐	☐	☐	☐	☐	☐	
	510	☐	☐	☐	☐	☐	☐	☐	☐	
Bodilink® Thigh Guides										
		Left	Right							
Small Thigh Support Stretch Air 90x110 (B)		☐	☐							
Medium Thigh Support Stretch Air 200x90 (C)		☐	☐							
X-Large Thigh Support Stretch Air 300x90 (G)		☐	☐							
Without Thigh Guides		☐	☐							
Joystick Holder										
		Left	Right							
435mm Swing Away Panelholder		☐	☐							
Parallel Motion Unitrack Panelholder		☐	☐							
				e-mail				Customer Support_UK@permobil.com		
				Motability CRN						
				Seat Assembly						
				Without Anterior/Negative tilt				☐		
				5° Active Seating				☐		
				Ergo Supports Back						
				Large - (7cm deep)				☐		
				Small - (5cm deep)				☐		
				Armrest Height						
				230mm Armrest Height (Pos 1)				☐		
				255mm Armrest Height (Pos 2)				☐		
				280mm Armrest Height (Pos 3)				☐		
				305mm Armrest Height (Pos 4)				☐		
				330mm Armrest Height (Pos 5)				☐		
				Armrest Style						
				Left				Right		
				Vinyl	Textile			Vinyl	Textile	
				260	☐	☐			☐	☐
				335	☐	☐			☐	☐
				410	☐	☐			☐	☐
				460	☐	☐			☐	☐

permobil		Permobil M3 Tilt - ROHO® 6 km/h																																													
		Permobil Account #:				Permobil M3- (706825)																																									
		Permobil PO #				Standard Options																																									
		Permobil Delivery Address		1 ☐		● R-Net EL LED		● 65 Ah AGM Gel Batteries																																							
				2 ☐		● MyPermobil Fleet Management		● ROHO® Max Contour Agility Back																																							
		3 ☐		● Solid Drive Wheels & Castors		● Manual Adjust Backrest 135°																																									
		4 ☐		● Fixed Seat Height		● Manual Legrest Mid Mount 150°																																									
RRP		£8920		● 50° Posterior Tilt		● Footplate Narrow 16cm																																									
Motability Weekly Rental		£53.10		● 420 mm - 570 mm Seat Width		● Bodilink® Headrest Medium																																									
Max User Weight		150kg		Client Weight		_____ kg		● 370 mm - 570 mm Seath Depth																																							
						● Stretch Air Ergo Seat Cushion		● Charger Voltpro AGM/Gel IP44																																							
Configurable Options																																															
Accent Colour																																															
	Aqua	☐		Sky	☐		Violet	☐		Flamingo	☐		Chilli	☐																																	
	Canary	☐		Spring	☐		Arctic	☐		Sunset	☐		Midnight	☐																																	
Seat Size										Seat Assembly																																					
		Depth										Without Anterior/Negative tilt ☐																																			
Size (mm)		370	395	420	445	470	495	520	545	570	5° Active Seating ☐																																				
Width	420	☐	☐	☐	☐	☐	☐	☐	☐	☐																																					
	470	☐	☐	☐	☐	☐	☐	☐	☐	☐																																					
	520	☐	☐	☐	☐	☐	☐	☐	☐	☐																																					
	570	☐	☐	☐	☐	☐	☐	☐	☐	☐																																					
ROHO® Agility Max Contour										Height Agility Lateral supports																																					
		Height			Mounting Height					Low Agility Lateral (Standard) ☐																																					
Size (mm)		405	455	510	(mm)	405	455	510	High Agility Lateral ☐																																						
Width	355	☐	☐	☐	470	☐																																									
	380	☐	☐	☐	495	☐																																									
	405	☐	☐	☐	520	☐	☐																																								
	430	☐	☐	☐	545	☐	☐	☐																																							
	455	☐	☐	☐	570	☐	☐	☐																																							
	485	☐	☐	☐	595	☐	☐	☐																																							
	510	☐	☐	☐	620		☐	☐																																							
										645			☐																																		
Bodilink® Thigh Guides										Armrest Height																																					
					Left	Right	230mm Armrest Height (Pos 1) ☐																																								
Small Thigh Support Stretch Air 90x110 (B)					☐	☐	255mm Armrest Height (Pos 2) ☐																																								
Medium Thigh Support Stretch Air 200x90 (C)					☐	☐	280mm Armrest Height (Pos 3) ☐																																								
X-Large Thigh Support Stretch Air 300x90 (G)					☐	☐	305mm Armrest Height (Pos 4) ☐																																								
Without Thigh Guides					☐	☐	330mm Armrest Height (Pos 5) ☐																																								
Joystick Holder										Armrest Style																																					
					Left	Right	<table border="1"> <tr> <td></td> <td colspan="2">Left</td> <td></td> <td colspan="2">Right</td> </tr> <tr> <td></td> <td>Vinyl</td> <td>Textile</td> <td></td> <td>Vinyl</td> <td>Textile</td> </tr> <tr> <td>260</td> <td>☐</td> <td>☐</td> <td></td> <td>☐</td> <td>☐</td> </tr> <tr> <td>335</td> <td>☐</td> <td>☐</td> <td></td> <td>☐</td> <td>☐</td> </tr> <tr> <td>410</td> <td>☐</td> <td>☐</td> <td></td> <td>☐</td> <td>☐</td> </tr> <tr> <td>460</td> <td>☐</td> <td>☐</td> <td></td> <td>☐</td> <td>☐</td> </tr> </table>						Left			Right			Vinyl	Textile		Vinyl	Textile	260	☐	☐		☐	☐	335	☐	☐		☐	☐	410	☐	☐		☐	☐	460	☐	☐		☐	☐
	Left			Right																																											
	Vinyl	Textile		Vinyl	Textile																																										
260	☐	☐		☐	☐																																										
335	☐	☐		☐	☐																																										
410	☐	☐		☐	☐																																										
460	☐	☐		☐	☐																																										
435mm Swing Away Panelholder					☐	☐																																									
Parallel Motion Unitrack Panelholder					☐	☐																																									
					<table border="1"> <tr> <td>e-mail</td> <td colspan="5">Customer_Support_UK@permobil.com</td> </tr> <tr> <td colspan="2">Motability CRN</td> <td colspan="4"></td> </tr> </table>					e-mail	Customer_Support_UK@permobil.com					Motability CRN																															
e-mail	Customer_Support_UK@permobil.com																																														
Motability CRN																																															

		Permobil M3 Tilt - Corpus® 10 km/h														
		Permobil Account #:				Permobil M3- (706832)										
		Permobil PO #				Standard Options										
		Permobil Delivery Address		1 ☐		● R-Net EL LED		● 65 Ah AGM Gel Batteries								
				2 ☐		● MyPermobil Fleet Management		● Corpus 3G® Contour Back								
		3 ☐		● Solid Drive Wheels & Castors		● Manual Adjust Backrest 135°										
		4 ☐		● Fixed Seat Height		● Manual Legrest Mid Mount 150°										
RRP		£9132		● 50° Posterior Tilt		● Footplate Narrow 16cm										
Motability Weekly Rental		£54.04		● 420 mm - 570 mm Seat Width		● Bodilink® Headrest Medium										
Max User Weight		150kg		Client Weight		_____ kg		● 370 mm - 570 mm Seath Depth								
						● Stretch Air Ergo Seat Cushion		● Charger Voltpro AGM/Gel IP44								
Configurable Options																
Accent Colour																
	Aqua	☐		Sky	☐		Violet	☐		Flamingo	☐		Chilli	☐		
	Canary	☐		Spring	☐		Arctic	☐		Sunset	☐		Midnight	☐		
Seat Size										Seat Assembly						
		Depth										Without Anterior/Negative tilt ☐				
Size (mm)		370	395	420	445	470	495	520	545	570	5° Active Seating ☐					
Width	420	☐	☐	☐	☐	☐	☐	☐	☐	☐						
	470	☐	☐	☐	☐	☐	☐	☐	☐	☐						
	520	☐	☐	☐	☐	☐	☐	☐	☐	☐						
	570	☐	☐	☐	☐	☐	☐	☐	☐	☐						
Corpus 3G Backrest										Postural Support Belt						
		Height								2 Point Hip Belt M L XL						
Size (mm)		470	520	545	570	595	620	645	670	Bodypoint® padded ☐ ☐ ☐						
Width	360	☐	☐	☐	☐	☐	☐	☐	☐							
	410	☐	☐	☐	☐	☐	☐	☐	☐							
	460	☐	☐	☐	☐	☐	☐	☐	☐							
	510	☐	☐	☐	☐	☐	☐	☐	☐							
Bodilink® Thigh Guides										Ergo Supports Back						
				Left	Right	Large - (7cm deep) ☐										
Small Thigh Support Stretch Air 90x110 (B)				☐	☐	Small - (5cm deep) ☐										
Medium Thigh Support Stretch Air 200x90 (C)				☐	☐											
X-Large Thigh Support Stretch Air 300x90 (G)				☐	☐											
Without Thigh Guides				☐	☐											
Joystick Holder										Armrest Height						
				Left	Right	230mm Armrest Height (Pos 1) ☐										
435mm Swing Away Panelholder				☐	☐	255mm Armrest Height (Pos 2) ☐										
Parallel Motion Unitrack Panelholder				☐	☐	280mm Armrest Height (Pos 3) ☐										
						305mm Armrest Height (Pos 4) ☐										
						330mm Armrest Height (Pos 5) ☐										
Armrest Style										Customer Support						
				Left	Right	e-mail Customer_Support_UK@permobil.com										
				Vinyl	Textile	Motability CRN										

		Permobil M3 Tilt - ROHO® 10 km/h												
		Permobil Account #:				Permobil M3- (706831)								
		Permobil PO #				Standard Options								
		Permobil Delivery Address		1 ☐		● R-Net EL LED		● 65 Ah AGM Gel Batteries						
				2 ☐		● MyPermobil Fleet Management		● ROHO® Max Contour Agility Back						
		3 ☐		● Solid Drive Wheels & Castors		● Manual Adjust Backrest 135°								
		4 ☐		● Fixed Seat Height		● Manual Legrest Mid Mount 150°								
RRP		£9287		● 50° Posterior Tilt		● Footplate Narrow 16cm								
Motability Weekly Rental		£54.74		● 420 mm - 570 mm Seat Width		● Bodilink® Headrest Medium								
Max User Weight	150kg	Client Weight			kg	● 370 mm - 570 mm Seath Depth		● LED Lights						
						● Stretch Air Ergo Seat Cushion		● Charger Voltpro AGM/Gel IP44						
Configurable Options														
Accent Colour														
	Aqua	☐		Sky	☐		Violet	☐		Flamingo	☐		Chilli	☐
	Canary	☐		Spring	☐		Arctic	☐		Sunset	☐		Midnight	☐
Seat Size														
		Depth												
Size (mm)		370	395	420	445	470	495	520	545	570				
Width	420	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	470	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	520	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	570	☐	☐	☐	☐	☐	☐	☐	☐	☐				
ROHO® Agility Max Contour														
		Height			Mounting Height									
Size (mm)		405	455	510	(mm)	405	455	510						
Width	355	☐	☐	☐	470	☐								
	380	☐	☐	☐	495	☐								
	405	☐	☐	☐	520	☐	☐							
	430	☐	☐	☐	545	☐	☐	☐						
	455	☐	☐	☐	570	☐	☐	☐						
	485	☐	☐	☐	595	☐	☐	☐						
	510	☐	☐	☐	620		☐	☐						
				645			☐							
Bodilink® Thigh Guides														
				Left	Right									
Small Thigh Support Stretch Air 90x110 (B)				☐	☐									
Medium Thigh Support Stretch Air 200x90 (C)				☐	☐									
X-Large Thigh Support Stretch Air 300x90 (G)				☐	☐									
Without Thigh Guides				☐	☐									
Joystick Holder														
				Left	Right									
435mm Swing Away Panelholder				☐	☐									
Parallel Motion Unitrack Panelholder				☐	☐									
Seat Assembly														
Without Anterior/Negative tilt				☐										
5° Active Seating				☐										
Postural Support Belt														
2 Point Hip Belt		M	L	XL										
Bodypoint® padded		☐	☐	☐										
Height Agility Lateral supports														
Low Agility Lateral (Standard)				☐										
High Agility Lateral				☐										
Armrest Height														
230mm Armrest Height (Pos 1)				☐										
255mm Armrest Height (Pos 2)				☐										
280mm Armrest Height (Pos 3)				☐										
305mm Armrest Height (Pos 4)				☐										
330mm Armrest Height (Pos 5)				☐										
Armrest Style														
	Left			Right										
	Vinyl	Textile		Vinyl	Textile									
260	☐	☐		☐	☐									
335	☐	☐		☐	☐									
410	☐	☐		☐	☐									
460	☐	☐		☐	☐									
e-mail Customer_Support_UK@permobil.com														
Motability CRN														

permobil		Permobil M3 Lift & Tilt - Corpus® 6 km/h																																												
		Permobil Account #:				Permobil M3- (706824)																																								
		Permobil PO #				Standard Options <table border="1"> <tr> <td>● R-Net EL LED</td> <td>●65 Ah AGM Gel Batteries</td> </tr> <tr> <td>● MyPermobil Fleet Management</td> <td>●Corpus 3G® Contour Back</td> </tr> <tr> <td>●Solid Drive Wheels & Castors</td> <td>●Manual Adjust Backrest 135°</td> </tr> <tr> <td>●300mm Seat Elevation</td> <td>●Manual Legrest Mid Mount 150°</td> </tr> <tr> <td>●50° Posterior Tilt</td> <td>●Footplate Narrow 16cm</td> </tr> <tr> <td>●420 mm - 570 mm Seat Width</td> <td>●Bodilink® Headrest Medium</td> </tr> <tr> <td>●370 mm - 570 mm Seath Depth</td> <td>●Reflector Kit</td> </tr> <tr> <td>●Stretch Air Ergo Seat Cushion</td> <td>●Charger Voltpro AGM/Gel IP44</td> </tr> </table>				● R-Net EL LED	●65 Ah AGM Gel Batteries	● MyPermobil Fleet Management	●Corpus 3G® Contour Back	●Solid Drive Wheels & Castors	●Manual Adjust Backrest 135°	●300mm Seat Elevation	●Manual Legrest Mid Mount 150°	●50° Posterior Tilt	●Footplate Narrow 16cm	●420 mm - 570 mm Seat Width	●Bodilink® Headrest Medium	●370 mm - 570 mm Seath Depth	●Reflector Kit	●Stretch Air Ergo Seat Cushion	●Charger Voltpro AGM/Gel IP44																					
		● R-Net EL LED	●65 Ah AGM Gel Batteries																																											
		● MyPermobil Fleet Management	●Corpus 3G® Contour Back																																											
●Solid Drive Wheels & Castors	●Manual Adjust Backrest 135°																																													
●300mm Seat Elevation	●Manual Legrest Mid Mount 150°																																													
●50° Posterior Tilt	●Footplate Narrow 16cm																																													
●420 mm - 570 mm Seat Width	●Bodilink® Headrest Medium																																													
●370 mm - 570 mm Seath Depth	●Reflector Kit																																													
●Stretch Air Ergo Seat Cushion	●Charger Voltpro AGM/Gel IP44																																													
Permobil Delivery Address		1 ☐ 2 ☐ 3 ☐ 4 ☐																																												
RRP		£9333																																												
Motability Weekly Rental		£54.86																																												
Max User Weight	150kg	Client Weight			kg																																									
Configurable Options																																														
Accent Colour																																														
	Aqua	☐		Sky	☐		Violet	☐		Flamingo	☐		Chilli	☐																																
	Canary	☐		Spring	☐		Arctic	☐		Sunset	☐		Midnight	☐																																
Seat Size										Seat Assembly																																				
		Depth									Without Anterior/Negative tilt ☐																																			
Size (mm)		370	395	420	445	470	495	520	545	570	10° Active Seating ☐																																			
Width	420	☐	☐	☐	☐	☐	☐	☐	☐	☐																																				
	470	☐	☐	☐	☐	☐	☐	☐	☐	☐																																				
	520	☐	☐	☐	☐	☐	☐	☐	☐	☐																																				
	570	☐	☐	☐	☐	☐	☐	☐	☐	☐																																				
Corpus 3G Backrest										Ergo Supports Back																																				
		Height								Large - (7cm deep) ☐																																				
Size (mm)		470	520	545	570	595	620	645	670	Small - (5cm deep) ☐																																				
Width	360	☐	☐	☐	☐	☐	☐	☐	☐																																					
	410	☐	☐	☐	☐	☐	☐	☐	☐																																					
	460	☐	☐	☐	☐	☐	☐	☐	☐																																					
	510	☐	☐	☐	☐	☐	☐	☐	☐																																					
Bodilink® Thigh Guides										Armrest Height																																				
				Left	Right	230mm Armrest Height (Pos 1) ☐																																								
Small Thigh Support Stretch Air 90x110 (B)				☐	☐	255mm Armrest Height (Pos 2) ☐																																								
Medium Thigh Support Stretch Air 200x90 (C)				☐	☐	280mm Armrest Height (Pos 3) ☐																																								
X-Large Thigh Support Stretch Air 300x90 (G)				☐	☐	305mm Armrest Height (Pos 4) ☐																																								
Without Thigh Guides				☐	☐	330mm Armrest Height (Pos 5) ☐																																								
Joystick Holder										Armrest Style																																				
				Left	Right	<table border="1"> <tr> <td></td> <td colspan="2">Left</td> <td></td> <td colspan="2">Right</td> </tr> <tr> <td></td> <td>Vinyl</td> <td>Textile</td> <td></td> <td>Vinyl</td> <td>Textile</td> </tr> <tr> <td>260</td> <td>☐</td> <td>☐</td> <td></td> <td>☐</td> <td>☐</td> </tr> <tr> <td>335</td> <td>☐</td> <td>☐</td> <td></td> <td>☐</td> <td>☐</td> </tr> <tr> <td>410</td> <td>☐</td> <td>☐</td> <td></td> <td>☐</td> <td>☐</td> </tr> <tr> <td>460</td> <td>☐</td> <td>☐</td> <td></td> <td>☐</td> <td>☐</td> </tr> </table>						Left			Right			Vinyl	Textile		Vinyl	Textile	260	☐	☐		☐	☐	335	☐	☐		☐	☐	410	☐	☐		☐	☐	460	☐	☐		☐	☐
	Left			Right																																										
	Vinyl	Textile		Vinyl	Textile																																									
260	☐	☐		☐	☐																																									
335	☐	☐		☐	☐																																									
410	☐	☐		☐	☐																																									
460	☐	☐		☐	☐																																									
435mm Swing Away Panelholder				☐	☐																																									
Parallel Motion Unitrack Panelholder				☐	☐																																									
						e-mail Customer_Support_UK@permobil.com Motability CRN																																								

		Permobil M3 Lift & Tilt - ROHO® 6 km/h														
		Permobil Account #:				Permobil M3- (706823)										
		Permobil PO #				Standard Options										
		Permobil Delivery Address		1 ☐		● R-Net EL LED		● 65 Ah AGM Gel Batteries								
				2 ☐		● MyPermobil Fleet Management		● ROHO® Max Contour Agility Back								
		3 ☐		● Solid Drive Wheels & Castors		● Manual Adjust Backrest 135°										
		4 ☐		● 300mm Seat Elevation		● Manual Legrest Mid Mount 150°										
RRP		£9320		● 50° Posterior Tilt		● Footplate Narrow 16cm										
Motability Weekly Rental		£54.89		● 420 mm - 570 mm Seat Width		● Bodilink® Headrest Medium										
Max User Weight	150kg	Client Weight			kg	● 370 mm - 570 mm Seath Depth		● Reflector Kit								
						● Stretch Air Ergo Seat Cushion		● Charger Voltpro AGM/Gel IP44								
Configurable Options																
Accent Colour																
	Aqua	☐		Sky	☐		Violet	☐		Flamingo	☐		Chilli	☐		
	Canary	☐		Spring	☐		Arctic	☐		Sunset	☐		Midnight	☐		
Seat Size										Seat Assembly						
		Depth										Without Anterior/Negative tilt ☐				
Size (mm)		370	395	420	445	470	495	520	545	570	10° Active Seating ☐					
Width	420	☐	☐	☐	☐	☐	☐	☐	☐	☐						
	470	☐	☐	☐	☐	☐	☐	☐	☐	☐						
	520	☐	☐	☐	☐	☐	☐	☐	☐	☐						
	570	☐	☐	☐	☐	☐	☐	☐	☐	☐						
ROHO® Agility Max Contour										Height Agility Lateral supports						
		Height			Mounting Height					Low Agility Lateral (Standard) ☐						
Size (mm)		405	455	510	(mm)	405	455	510	High Agility Lateral ☐							
Width	355	☐	☐	☐	470	☐										
	380	☐	☐	☐	495	☐										
	405	☐	☐	☐	520	☐	☐									
	430	☐	☐	☐	545	☐	☐	☐								
	455	☐	☐	☐	570	☐	☐	☐								
	485	☐	☐	☐	595	☐	☐	☐								
	510	☐	☐	☐	620		☐	☐								
					645			☐								
Bodilink® Thigh Guides										Armrest Height						
				Left	Right	230mm Armrest Height (Pos 1) ☐										
Small Thigh Support Stretch Air 90x110 (B)				☐	☐	255mm Armrest Height (Pos 2) ☐										
Medium Thigh Support Stretch Air 200x90 (C)				☐	☐	280mm Armrest Height (Pos 3) ☐										
X-Large Thigh Support Stretch Air 300x90 (G)				☐	☐	305mm Armrest Height (Pos 4) ☐										
Without Thigh Guides				☐	☐	330mm Armrest Height (Pos 5) ☐										
Joystick Holder										Armrest Style						
				Left	Right	Left					Right					
						Vinyl					Textile					
435mm Swing Away Panelholder				☐	☐	260					☐					
Parallel Motion Unitrack Panelholder				☐	☐	335					☐					
						410					☐					
						460					☐					
						e-mail					Customer_Support_UK@permobil.com					
						Motability CRN										

		Permobil M3 Tilt & Recline - Corpus® 10 km/h																					
		Permobil Account #:				Permobil M3- (706820) Standard Options <table border="1"> <tr> <td>• R-Net EL LED</td> <td>• 65 Ah AGM Gel Batteries</td> </tr> <tr> <td>• MyPermobil Fleet Management</td> <td>• Corpus 3G® Contour Back</td> </tr> <tr> <td>• Solid Drive Wheels & Castors</td> <td>• Powered Adjust Backrest 1805°</td> </tr> <tr> <td>• Fixed Seat Height</td> <td>• Manual Legrest Mid Mount 150°</td> </tr> <tr> <td>• 50° Posterior Tilt</td> <td>• Footplate Narrow 16cm</td> </tr> <tr> <td>• 420 mm - 570 mm Seat Width</td> <td>• Bodilink® Headrest Medium</td> </tr> <tr> <td>• 370 mm - 570 mm Seath Depth</td> <td>• LED Lights</td> </tr> <tr> <td>• Stretch Air Ergo Seat Cushion</td> <td>• Charger Voltpro AGM/Gel IP44</td> </tr> </table>		• R-Net EL LED	• 65 Ah AGM Gel Batteries	• MyPermobil Fleet Management	• Corpus 3G® Contour Back	• Solid Drive Wheels & Castors	• Powered Adjust Backrest 1805°	• Fixed Seat Height	• Manual Legrest Mid Mount 150°	• 50° Posterior Tilt	• Footplate Narrow 16cm	• 420 mm - 570 mm Seat Width	• Bodilink® Headrest Medium	• 370 mm - 570 mm Seath Depth	• LED Lights	• Stretch Air Ergo Seat Cushion	• Charger Voltpro AGM/Gel IP44
		• R-Net EL LED	• 65 Ah AGM Gel Batteries																				
		• MyPermobil Fleet Management	• Corpus 3G® Contour Back																				
		• Solid Drive Wheels & Castors	• Powered Adjust Backrest 1805°																				
• Fixed Seat Height	• Manual Legrest Mid Mount 150°																						
• 50° Posterior Tilt	• Footplate Narrow 16cm																						
• 420 mm - 570 mm Seat Width	• Bodilink® Headrest Medium																						
• 370 mm - 570 mm Seath Depth	• LED Lights																						
• Stretch Air Ergo Seat Cushion	• Charger Voltpro AGM/Gel IP44																						
Permobil PO #																							
Permobil Delivery Address		1	☐																				
		2	☐																				
		3	☐																				
		4	☐																				
RRP		£9333																					
Motability Weekly Rental		£54.94																					
Max User Weight	150kg	Client Weight			kg																		

Configurable Options											
Accent Colour											
	Aqua	☐		Sky	☐		Violet	☐		Flamingo	☐
	Canary	☐		Spring	☐		Arctic	☐		Sunset	☐
										Chilli	☐
										Midnight	☐

Seat Size										
		Depth								
Size (mm)		370	395	420	445	470	495	520	545	570
Width	420	☐	☐	☐	☐	☐	☐	☐	☐	☐
	470	☐	☐	☐	☐	☐	☐	☐	☐	☐
	520	☐	☐	☐	☐	☐	☐	☐	☐	☐
	570	☐	☐	☐	☐	☐	☐	☐	☐	☐

Corpus 3G Backrest									
		Height							
Size (mm)		470	520	545	570	595	620	645	670
Width	360	☐	☐	☐	☐	☐	☐	☐	☐
	410	☐	☐	☐	☐	☐	☐	☐	☐
	460	☐	☐	☐	☐	☐	☐	☐	☐
	510	☐	☐	☐	☐	☐	☐	☐	☐

Seat Assembly			
Without Anterior/Negative tilt			☐
5° Active Seating			☐

Postural Support Belt			
2 Point Hip Belt	M	L	XL
Bodypoint® padded	☐	☐	☐

Ergo Supports Back	
Large - (7cm deep)	☐
Small - (5cm deep)	☐

Armrest Height	
230mm Armrest Height (Pos 1)	☐
255mm Armrest Height (Pos 2)	☐
280mm Armrest Height (Pos 3)	☐
305mm Armrest Height (Pos 4)	☐
330mm Armrest Height (Pos 5)	☐

Bodilink® Thigh Guides		
	Left	Right
Small Thigh Support Stretch Air 90x110 (B)	☐	☐
Medium Thigh Support Stretch Air 200x90 (C)	☐	☐
X-Large Thigh Support Stretch Air 300x90 (G)	☐	☐
Without Thigh Guides	☐	☐

Joystick Holder		
	Left	Right
435mm Swing Away Panelholder	☐	☐
Parallel Motion Unitrack Panelholder	☐	☐

Armrest Style						
	Left			Right		
	Vinyl	Textile		Vinyl	Textile	
	260	☐		☐	☐	☐
	335	☐		☐	☐	☐
	410	☐		☐	☐	☐
460	☐	☐	☐	☐	☐	

e-mail	Customer_Support_UK@permobil.com
Motability CRN	

		Permobil M3 Tilt, Recline & ELR - Corpus® 6 km/h																			
		Permobil Account #:		Permobil M3- (706821) Standard Options <table border="1"> <tr> <td>• R-Net EL LED</td> <td>•65 Ah AGM Gel Batteries</td> </tr> <tr> <td>• MyPermobil Fleet Management</td> <td>•Corpus 3G® Contour Back</td> </tr> <tr> <td>•Solid Drive Wheels & Castors</td> <td>•Powered Backrest 180°</td> </tr> <tr> <td>•Fixed Seat Height</td> <td>•Powered Legrest Mid Mount 180 °</td> </tr> <tr> <td>•50° Posterior Tilt</td> <td>•Footplate Narrow 16cm</td> </tr> <tr> <td>•420 mm - 570 mm Seat Width</td> <td>•Bodilink® Headrest Medium</td> </tr> <tr> <td>•370 mm - 570 mm Seath Depth</td> <td>•Reflector Kit</td> </tr> <tr> <td>•Stretch Air Ergo Seat Cushion</td> <td>•Charger Voltpro AGM/Gel IP44</td> </tr> </table>		• R-Net EL LED	•65 Ah AGM Gel Batteries	• MyPermobil Fleet Management	•Corpus 3G® Contour Back	•Solid Drive Wheels & Castors	•Powered Backrest 180°	•Fixed Seat Height	•Powered Legrest Mid Mount 180 °	•50° Posterior Tilt	•Footplate Narrow 16cm	•420 mm - 570 mm Seat Width	•Bodilink® Headrest Medium	•370 mm - 570 mm Seath Depth	•Reflector Kit	•Stretch Air Ergo Seat Cushion	•Charger Voltpro AGM/Gel IP44
		• R-Net EL LED	•65 Ah AGM Gel Batteries																		
		• MyPermobil Fleet Management	•Corpus 3G® Contour Back																		
		•Solid Drive Wheels & Castors	•Powered Backrest 180°																		
•Fixed Seat Height	•Powered Legrest Mid Mount 180 °																				
•50° Posterior Tilt	•Footplate Narrow 16cm																				
•420 mm - 570 mm Seat Width	•Bodilink® Headrest Medium																				
•370 mm - 570 mm Seath Depth	•Reflector Kit																				
•Stretch Air Ergo Seat Cushion	•Charger Voltpro AGM/Gel IP44																				
Permobil PO #																					
Permobil Delivery Address																					
RRP																					
Motability Weekly Rental		Client Weight		kg																	
Max User Weight		150kg		£9333																	
£54.94		1 ☐		2 ☐																	
3 ☐		4 ☐																			
Configurable Options																					
Accent Colour																					
Aqua ☐		Sky ☐		Violet ☐																	
Canary ☐		Spring ☐		Arctic ☐																	
Flamingo ☐		Chilli ☐		Sunset ☐																	
Midnight ☐																					
Seat Size																					
Size (mm)		Depth																			
370 395 420 445 470 495 520 545 570																					
Width		420 470 520 570																			
Corpus 3G Backrest																					
Size (mm)		Height																			
470 520 545 570 595 620 645 670																					
Width		360 410 460 510																			
Bodilink® Thigh Guides																					
Left		Right																			
Small Thigh Support Stretch Air 90x110 (B)																					
Medium Thigh Support Stretch Air 200x90 (C)																					
X-Large Thigh Support Stretch Air 300x90 (G)																					
Without Thigh Guides																					
Joystick Holder																					
Left		Right																			
435mm Swing Away Panelholder																					
Parallel Motion Unitrack Panelholder																					
Seat Assembly																					
Without Anterior/Negative tilt																					
5° Active Seating																					
Ergo Supports Back																					
Large - (7cm deep)																					
Small - (5cm deep)																					
Armrest Height																					
230mm Armrest Height (Pos 1)																					
255mm Armrest Height (Pos 2)																					
280mm Armrest Height (Pos 3)																					
305mm Armrest Height (Pos 4)																					
330mm Armrest Height (Pos 5)																					
Armrest Style																					
Left		Right																			
Vinyl		Textile																			
260																					
335																					
410																					
460																					
Customer Support																					
e-mail		Customer Support		UK@permobil.com																	
Motability CRN																					

		Permobil M3 Lift & Tilt - ROHO® 10 km/h												
		Permobil Account #:				Permobil M3- (706827)								
		Permobil PO #				Standard Options								
		Permobil Delivery Address		1 ☐		● R-Net EL LED		● 65 Ah AGM Gel Batteries						
				2 ☐		● MyPermobil Fleet Management		● ROHO® Max Contour Agility Back						
		3 ☐		● Solid Drive Wheels & Castors		● Manual Adjust Backrest 135°								
		4 ☐		● Fixed Seat Height		● Manual Legrest Mid Mount 150°								
RRP		£9333		● 50° Posterior Tilt		● Footplate Narrow 16cm								
Motability Weekly Rental		£54.94		● 420 mm - 570 mm Seat Width		● Bodilink® Headrest Medium								
Max User Weight	150kg	Client Weight			kg	● 370 mm - 570 mm Seath Depth		● LED Lights						
						● Stretch Air Ergo Seat Cushion		● Charger Voltpro AGM/Gel IP44						
Configurable Options														
Accent Colour														
	Aqua	☐		Sky	☐		Violet	☐		Flamingo	☐		Chilli	☐
	Canary	☐		Spring	☐		Arctic	☐		Sunset	☐		Midnight	☐
Seat Size														
		Depth												
Size (mm)		370	395	420	445	470	495	520	545	570				
Width	420	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	470	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	520	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	570	☐	☐	☐	☐	☐	☐	☐	☐	☐				
ROHO® Agility Max Contour														
		Height			Mounting Height									
Size (mm)		405	455	510	(mm)	405	455	510						
Width	355	☐	☐	☐	470	☐								
	380	☐	☐	☐	495	☐								
	405	☐	☐	☐	520	☐	☐							
	430	☐	☐	☐	545	☐	☐	☐						
	455	☐	☐	☐	570	☐	☐	☐						
	485	☐	☐	☐	595	☐	☐	☐						
	510	☐	☐	☐	620		☐	☐						
					645			☐						
Bodilink® Thigh Guides														
				Left	Right									
Small Thigh Support Stretch Air 90x110 (B)				☐	☐									
Medium Thigh Support Stretch Air 200x90 (C)				☐	☐									
X-Large Thigh Support Stretch Air 300x90 (G)				☐	☐									
Without Thigh Guides				☐	☐									
Joystick Holder														
				Left	Right									
435mm Swing Away Panelholder				☐	☐									
Parallel Motion Unitrack Panelholder				☐	☐									
Seat Assembly														
Without Anterior/Negative tilt				☐										
10° Active Seating				☐										
Postural Support Belt														
2 Point Hip Belt		M	L	XL										
Bodypoint® padded		☐	☐	☐										
Height Agility Lateral supports														
Low Agility Lateral (Standard)				☐										
High Agility Lateral				☐										
Armrest Height														
230mm Armrest Height (Pos 1)				☐										
255mm Armrest Height (Pos 2)				☐										
280mm Armrest Height (Pos 3)				☐										
305mm Armrest Height (Pos 4)				☐										
330mm Armrest Height (Pos 5)				☐										
Armrest Style														
		Left				Right								
		Vinyl	Textile			Vinyl		Textile						
260		☐	☐			☐		☐						
335		☐	☐			☐		☐						
410		☐	☐			☐		☐						
460		☐	☐			☐		☐						
e-mail		Customer_Support_UK@permobil.com												
Motability CRN														

		Permobil M3 Tilt & ELR - Corpus® 10 km/h												
		Permobil Account #:				Permobil M3- (706828)								
		Permobil PO #				Standard Options								
		Permobil Delivery Address		1 ☐		● R-Net EL LED		● 65 Ah AGM Gel Batteries						
				2 ☐		● MyPermobil Fleet Management		● Corpus 3G® Contour Back						
		3 ☐		● Solid Drive Wheels & Castors		● Manual Adjust Backrest 135°								
		4 ☐		● Fixed Seat Height		● Powered Legrest Mid Mount 180 °								
RRP		£9333		● 50° Posterior Tilt		● Footplate Narrow 16cm								
Motability Weekly Rental		£54.94		● 420 mm - 570 mm Seat Width		● Bodilink® Headrest Medium								
Max User Weight	150kg	Client Weight			kg	● 370 mm - 570 mm Seath Depth		● LED Lights						
						● Stretch Air Ergo Seat Cushion		● Charger Voltpro AGM/Gel IP44						
Configurable Options														
Accent Colour														
	Aqua	☐		Sky	☐		Violet	☐		Flamingo	☐		Chilli	☐
	Canary	☐		Spring	☐		Arctic	☐		Sunset	☐		Midnight	☐
Seat Size														
		Depth												
Size (mm)		370	395	420	445	470	495	520	545	570				
Width	420	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	470	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	520	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	570	☐	☐	☐	☐	☐	☐	☐	☐	☐				
Corpus 3G Backrest														
		Height												
Size (mm)		470	520	545	570	595	620	645	670					
Width	360	☐	☐	☐	☐	☐	☐	☐	☐					
	410	☐	☐	☐	☐	☐	☐	☐	☐					
	460	☐	☐	☐	☐	☐	☐	☐	☐					
	510	☐	☐	☐	☐	☐	☐	☐	☐					
Seat Assembly														
		Without Anterior/Negative tilt ☐												
		5° Active Seating ☐												
Postural Support Belt														
2 Point Hip Belt		M	L	XL										
Bodypoint® padded		☐	☐	☐										
Ergo Supports Back														
		Large - (7cm deep) ☐												
		Small - (5cm deep) ☐												
Armrest Height														
		230mm Armrest Height (Pos 1) ☐												
		255mm Armrest Height (Pos 2) ☐												
		280mm Armrest Height (Pos 3) ☐												
		305mm Armrest Height (Pos 4) ☐												
		330mm Armrest Height (Pos 5) ☐												
Bodilink® Thigh Guides														
		Left	Right											
Small Thigh Support Stretch Air 90x110 (B)		☐	☐											
Medium Thigh Support Stretch Air 200x90 (C)		☐	☐											
X-Large Thigh Support Stretch Air 300x90 (G)		☐	☐											
Without Thigh Guides		☐	☐											
Joystick Holder														
		Left	Right											
435mm Swing Away Panelholder		☐	☐											
Parallel Motion Unitrack Panelholder		☐	☐											
		Armrest Style												
		Left				Right								
		Vinyl	Textile			Vinyl	Textile							
260		☐	☐			☐	☐							
335		☐	☐			☐	☐							
410		☐	☐			☐	☐							
460		☐	☐			☐	☐							
e-mail		Customer_Support_UK@permobil.com												
Motability CRN														

		Permobil M3 Lift & Tilt - ROHO® 10 km/h												
		Permobil Account #:				Permobil M3- (706829)								
		Permobil PO #				Standard Options								
		Permobil Delivery Address		1 ☐		● R-Net EL LED		● 65 Ah AGM Gel Batteries						
				2 ☐		● MyPermobil Fleet Management		● ROHO® Max Contour Agility Back						
		3 ☐		● Solid Drive Wheels & Castors		● Manual Adjust Backrest 135°								
		4 ☐		● 300mm Seat Elevation		● Manual Legrest Mid Mount 150°								
RRP		£9333		● 50° Posterior Tilt		● Footplate Narrow 16cm								
Motability Weekly Rental		£54.94		● 420 mm - 570 mm Seat Width		● Bodilink® Headrest Medium								
Max User Weight	150kg	Client Weight			kg	● 370 mm - 570 mm Seath Depth		● LED Lights						
						● Stretch Air Ergo Seat Cushion		● Charger Voltpro AGM/Gel IP44						
Configurable Options														
Accent Colour														
	Aqua	☐		Sky	☐		Violet	☐		Flamingo	☐		Chilli	☐
	Canary	☐		Spring	☐		Arctic	☐		Sunset	☐		Midnight	☐
Seat Size														
		Depth												
Size (mm)		370	395	420	445	470	495	520	545	570				
Width	420	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	470	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	520	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	570	☐	☐	☐	☐	☐	☐	☐	☐	☐				
ROHO® Agility Max Contour														
		Height			Mounting Height									
Size (mm)		405	455	510	(mm)	405	455	510						
Width	355	☐	☐	☐	470	☐								
	380	☐	☐	☐	495	☐								
	405	☐	☐	☐	520	☐	☐							
	430	☐	☐	☐	545	☐	☐	☐						
	455	☐	☐	☐	570	☐	☐	☐						
	485	☐	☐	☐	595	☐	☐	☐						
510	☐	☐	☐	620		☐	☐							
					645			☐						
Bodilink® Thigh Guides														
				Left	Right									
Small Thigh Support Stretch Air 90x110 (B)				☐	☐									
Medium Thigh Support Stretch Air 200x90 (C)				☐	☐									
X-Large Thigh Support Stretch Air 300x90 (G)				☐	☐									
Without Thigh Guides				☐	☐									
Joystick Holder														
				Left	Right									
435mm Swing Away Panelholder				☐	☐									
Parallel Motion Unitrack Panelholder				☐	☐									
Seat Assembly														
Without Anterior/Negative tilt						☐								
10° Active Seating						☐								
Postural Support Belt														
2 Point Hip Belt		M	L	XL										
Bodypoint® padded		☐	☐	☐										
Height Agility Lateral supports														
Low Agility Lateral (Standard)						☐								
High Agility Lateral						☐								
Armrest Height														
230mm Armrest Height (Pos 1)						☐								
255mm Armrest Height (Pos 2)						☐								
280mm Armrest Height (Pos 3)						☐								
305mm Armrest Height (Pos 4)						☐								
330mm Armrest Height (Pos 5)						☐								
Armrest Style														
		Left				Right								
		Vinyl	Textile			Vinyl		Textile						
260	☐	☐			☐		☐							
335	☐	☐			☐		☐							
410	☐	☐			☐		☐							
460	☐	☐			☐		☐							
e-mail Customer_Support_UK@permobil.com														
Motability CRN														

		Permobil M3 Lift & Tilt - Corpus® 10 km/h												
		Permobil Account #:				Permobil M3- (706830)								
		Permobil PO #				Standard Options								
		Permobil Delivery Address		1 ☐		● R-Net EL LED		● 65 Ah AGM Gel Batteries						
				2 ☐		● MyPermobil Fleet Management		● Corpus 3G® Contour Back						
		3 ☐		● Solid Drive Wheels & Castors		● Manual Adjust Backrest 135°								
		4 ☐		● 300mm Seat Elevation		● Manual Legrest Mid Mount 150°								
RRP		£9333		● 50° Posterior Tilt		● Footplate Narrow 16cm								
Motability Weekly Rental		£54.94		● 420 mm - 570 mm Seat Width		● Bodilink® Headrest Medium								
Max User Weight	150kg	Client Weight			kg	● 370 mm - 570 mm Seath Depth		● LED Lights						
						● Stretch Air Ergo Seat Cushion		● Charger Voltpro AGM/Gel IP44						
Configurable Options														
Accent Colour														
	Aqua	☐		Sky	☐		Violet	☐		Flamingo	☐		Chilli	☐
	Canary	☐		Spring	☐		Arctic	☐		Sunset	☐		Midnight	☐
Seat Size														
		Depth												
Size (mm)		370	395	420	445	470	495	520	545	570				
Width	420	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	470	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	520	☐	☐	☐	☐	☐	☐	☐	☐	☐				
	570	☐	☐	☐	☐	☐	☐	☐	☐	☐				
Corpus 3G Backrest														
		Height												
Size (mm)		470	520	545	570	595	620	645	670					
Width	360	☐	☐	☐	☐	☐	☐	☐	☐					
	410	☐	☐	☐	☐	☐	☐	☐	☐					
	460	☐	☐	☐	☐	☐	☐	☐	☐					
	510	☐	☐	☐	☐	☐	☐	☐	☐					
Bodilink® Thigh Guides														
		Left	Right											
Small Thigh Support Stretch Air 90x110 (B)		☐	☐											
Medium Thigh Support Stretch Air 200x90 (C)		☐	☐											
X-Large Thigh Support Stretch Air 300x90 (G)		☐	☐											
Without Thigh Guides		☐	☐											
Joystick Holder														
		Left	Right											
435mm Swing Away Panelholder		☐	☐											
Parallel Motion Unitrack Panelholder		☐	☐											
Seat Assembly														
Without Anterior/Negative tilt		☐												
10° Active Seating		☐												
Postural Support Belt														
2 Point Hip Belt		M	L	XL										
Bodypoint® padded		☐	☐	☐										
Ergo Supports Back														
Large - (7cm deep)		☐												
Small - (5cm deep)		☐												
Armrest Height														
230mm Armrest Height (Pos 1)		☐												
255mm Armrest Height (Pos 2)		☐												
280mm Armrest Height (Pos 3)		☐												
305mm Armrest Height (Pos 4)		☐												
330mm Armrest Height (Pos 5)		☐												
Armrest Style														
	Left			Right										
	Vinyl	Textile		Vinyl	Textile									
260	☐	☐		☐	☐									
335	☐	☐		☐	☐									
410	☐	☐		☐	☐									
460	☐	☐		☐	☐									
e-mail		Customer_Support_UK@permobil.com												
Motability CRN														