

Customer and Order Information			
Purchase Order *	Tag		v.0101
Dealer/Provider Name *			
Purchaser Name and Contact Information *			
Shipping Address *			
User Information (if available)	Weight _____ lbs	Height _____ ft/in	*required

Order Notes/Instructions

Base Chair			
Bentley LT-R 20° Manual Tilt Wheelchair Reverse Wheel Layout enhances user access for self-propulsion. pdgmobility.com/blt-r			
<input checked="" type="checkbox"/> 2701000	Bentley 20° Manual Tilt Wheelchair	E1161	\$ 3,150
<input checked="" type="checkbox"/> 2701009	Lightweight Aluminum Frame and Reverse Wheel Layout Upgrade	K0108	\$ 300
Choose the standard Bentley with a steel frame for agitated or heavier users. See separate order form.			

Frame Type & Weight Capacity			
Frame Type All frames offer 300 lbs Weight Capacity.			
<input type="checkbox"/> 2702010	Fixed Width Frame	Select a Fixed Seat Width below.	
<input type="checkbox"/> 2702025	Adjustable Width Frame	Select an Adjustable Seat Width below.	\$ 100
User Weight Range Required for setup of tilt uplift force.			
<input type="checkbox"/> 2702430	Up to 75 lbs	<input type="checkbox"/> 2702431	75 lbs to 100 lbs
<input type="checkbox"/> 2702433	115 lbs to 150 lbs	<input type="checkbox"/> 2702434	150 lbs to 200 lbs
<input type="checkbox"/> 2702436	250 lbs to 300 lbs	<input type="checkbox"/> 2702432	100 lbs to 115 lbs
		<input type="checkbox"/> 2702435	200 lbs to 250 lbs

Frame Options	
<input type="checkbox"/> 2702300	Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.

Seat Width			Overall chair width is +10.5".
Fixed Seat Widths Order with Fixed Frame.			
<input type="checkbox"/> 2703015	15" Seat Width	<input type="checkbox"/> 2703016	16" Seat Width
<input type="checkbox"/> 2703018	18" Seat Width	<input type="checkbox"/> 2703019	19" Seat Width
<input type="checkbox"/> 2703020	20" Seat Width	<input type="checkbox"/> 2703021	21" Seat Width
<input type="checkbox"/> 2703023	23" Seat Width	<input type="checkbox"/> 2703022	22" Seat Width
<input type="checkbox"/> 2703024	24" Seat Width		E2201 \$ 375
			E2201 \$ 375
			E2202 \$ 575
Adjustable Seat Widths Order with Adjustable Frame.			
<input type="checkbox"/> 2703112	12"-16" Range - Set at 12"	<input type="checkbox"/> 2703113	12"-16" Range - Set at 13"
<input type="checkbox"/> 2703115	12"-16" Range - Set at 15"	<input type="checkbox"/> 2703116	12"-16" Range - Set at 16"
<input type="checkbox"/> 2703215	15"-20" Range - Set at 15"	<input type="checkbox"/> 2703216	15"-20" Range - Set at 16"
<input type="checkbox"/> 2703218	15"-20" Range - Set at 18"	<input type="checkbox"/> 2703217	15"-20" Range - Set at 17"
<input type="checkbox"/> 2703220	15"-20" Range - Set at 20"		E2201 \$ 375

Seat Depth

Standard Seat Depths (Adjustable)

<input type="checkbox"/> 2704016	16" Seat Depth	<input type="checkbox"/> 2704017	17" Seat Depth		
<input type="checkbox"/> 2704018	18" Seat Depth	<input type="checkbox"/> 2704019	19" Seat Depth		
<input type="checkbox"/> 2704020	20" Seat Depth			E2203	\$ 525

Custom Seat Depths

<input type="checkbox"/> 2704021	21" Seat Depth			E2203	\$ 525
<input type="checkbox"/> 2704022	22" Seat Depth			E2204	\$ 800

Frame Colour

PDG Premium Colours

View colour swatches on PDG's website.

pdgmobility.com/colours

<input type="checkbox"/> 2705001	Black Shimmer	<input type="checkbox"/> 2705015	Sky Blue	<input type="checkbox"/> 2705016	Poppy Red
<input type="checkbox"/> 2705003	Fine Texture Black	<input type="checkbox"/> 2705010	Lollypop Blue	<input type="checkbox"/> 2705011	Lollypop Red
<input type="checkbox"/> 2705017	Arctic Grey	<input type="checkbox"/> 2705018	Tropical Teal	<input type="checkbox"/> 2705019	Tangerine Orange
<input type="checkbox"/> 2705020	Olive Green	<input type="checkbox"/> 2705021	Lush Lilac	<input type="checkbox"/> 2705012	Bubble Gum Pink
<input type="checkbox"/> 2705022	Leaf Green	<input type="checkbox"/> 2705023	Boysenberry Purple	<input type="checkbox"/> 2705014	Daffodil Yellow
<input type="checkbox"/> 2705200	Custom Colour	Contact PDG for specifications.			\$ 300

Armrests

<input type="checkbox"/> 2706000	No Armrests				
	Height Adjustable Armrests (Single Post)			E0973	\$ 345
<input type="checkbox"/> 2706111	8"-11" Range	<input type="checkbox"/> 2706114	10"-14" Range	<input type="checkbox"/> 2706116	12"-16" Range
<input type="checkbox"/> 2706100	Fixed Height Armrests (12")	Full Length Armpad only.			
	Armrest Height	If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.			
	Seat Pan to top of Armrest:	_____	Does not include height of Armpad		
	Armrest Position	Indicate preference for Armrest position. If left blank, the default setting will be used.			
	In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.				
<input type="checkbox"/> 2706401	Rear	Default Position	<input type="checkbox"/> 2706403	Front	

Armpads

Regular Armpads

<input type="checkbox"/> 2707100	Desk Length - 2" x 10"	<input type="checkbox"/> 2707101	Full Length - 2" x 14"		
	Waterfall Armpads				\$ 115
<input type="checkbox"/> 2707105	Desk Length - 2" x 10"	<input type="checkbox"/> 2707106	Full Length - 2" x 14"		
	Wide Armpads	May interfere with wheels in some tilt configurations.			\$ 200
<input type="checkbox"/> 2707110	Desk Length - 3" x 10"	<input type="checkbox"/> 2707111	Full Length - 3" x 12"		
	Gel Armpads	May interfere with wheels in some tilt configurations.			\$ 135
<input type="checkbox"/> 2707120	Full Length - 2" x 12"	<input type="checkbox"/> 2707121	Full Length - 3.5" x 12"		
	Contour Arm Trough with Hand Pad	May interfere with wheels in some tilt configurations.			E2209 \$ 125
<input type="checkbox"/> 2707130	Right (4.5" x 18")	<input type="checkbox"/> 2707131	Left (4.5" x 18")		

Armrest Options

Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change.

<input type="checkbox"/> 2707210	1" Total Additional Width				\$ 95
<input type="checkbox"/> 2707211	2" Total Additional Width				\$ 195
<input type="checkbox"/> 2707310	Clothing Guards	Not compatible with 8"-11" Height range.			\$ 110
	Armrest Extensions	Not compatible with Desk Length Armpads.			\$ 75
<input type="checkbox"/> 2707400	3" Extension - Forward	<input type="checkbox"/> 2707401	3" Extension - Rearward		

Seat Pan

<input type="checkbox"/> 2708000	No Seat Pan	
<input type="checkbox"/> 2708001	Solid Aluminum Seat Pan	E2231 \$ 165

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low							High						
	RW	12"	13"	14"	15"	16"	17"	RW	14"	15"	16"	17"	18"	19"
20"	<input type="checkbox"/>	20"	<input type="checkbox"/>											
22"		<input type="checkbox"/>	22"		<input type="checkbox"/>									
24"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters	Low						High						
	RW	12"	13"	14"	15"	17"	RW	14"	15"	16"	17"	18"	19"
24"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

Urethane Casters 8" Casters are compatible only with 24" Large Wheels.		
<input type="checkbox"/> 2711061	6" x 1.25" Caster size	<input type="checkbox"/> 2711081 8" x 1.25" Caster size
<input type="checkbox"/> 2711063	6" x 1.50" Caster size	\$ 85
<input type="checkbox"/> 2711084	8" x 2.00" Caster size	\$ 135
Pneumatic Casters 8" Casters are compatible only with 24" Large Wheels.		
<input type="checkbox"/> 2711161	6" x 1.25" Caster size	E2214 & E2215 \$ 110
<input type="checkbox"/> 2711181	8" x 1.25" Caster size	E2214 & E2215 \$ 120
<input type="checkbox"/> 2711184	8" x 2.00" Caster size	E2214 & E2215 \$ 150
<input type="checkbox"/> 2711260	Soft-Roll Casters 6" x 1.50" Caster size. Urethane tire & Silver Anodized Aluminum hub.	\$ 275
<input type="checkbox"/> 2711400	Frog Legs Suspension Forks Only compatible with 24" Large Wheels. Requires 6" Soft-Roll Caster.	E1015 \$ 595

Large Wheels

Composite Wheels		
<input type="checkbox"/> 2712020	20" Wheel size	<input type="checkbox"/> 2712022 22" Wheel size
<input type="checkbox"/> 2712024	24" Wheel size	\$ 75
Spoke Wheels		
<input type="checkbox"/> 2712120	20" Wheel size	<input type="checkbox"/> 2712122 22" Wheel size
<input type="checkbox"/> 2712124	24" Wheel size	
 Large Wheel Position Indicate Large Wheel positioning preference. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.		
<input type="checkbox"/> 2712401	Rear Less stable	<input type="checkbox"/> 2712402 Mid Default Position
<input type="checkbox"/> 2712403	Front More stable. Not compatible with 20"-22" Seat Depth.	

Axle Type

<input type="checkbox"/> 2713000	Quick Release Axles
<input type="checkbox"/> 2713010	Fixed Axles

Large Tires

<input type="checkbox"/> 2714000	Urethane Tires	
<input type="checkbox"/> 2714010	Pneumatic Tires May result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212 \$ 60
<input type="checkbox"/> 2714020	Pneumatic Tires with Flat-Free Inserts May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	E2211 & E2213 \$ 95

Handrims

<input type="checkbox"/> 2715000	No Handrims	
<input type="checkbox"/> 2715010	Standard Aluminum Handrims	
<input type="checkbox"/> 2715020	Plastic Coated Handrims	\$ 150
<input type="checkbox"/> 2715030	Projection Handrims Composite Wheels only. 10 Projections, 30°, 1-1.8" Length.	\$ 495

Wheel Locks & Brakes

User Wheel Locks

- 2716010 Push to Lock 2716020 Pull to Lock

User-activated Unilateral Wheel Lock

 Pull to Lock.

\$ 795

- 1016310 Left Activation 1016311 Right Activation

Wheel Lock Options

- 2716100 **Wheel Lock Extensions** *E0961* \$ 80

- 2716200 **Hill Holders** Requires Pneumatic Large Tires with Flat-Free Inserts. *E0974* \$ 450

Anti-Tippers

Front Anti-Tippers

 Required Safety Option. Must choose one type.

- 2717010 Standard *E0971* \$ 125

- 2717015 Spring-loaded Dampens the forward-tip movement of the chair during transfers and allows for improved ground clearance. Not removable. *E0971* \$ 200

Back Posts

Back Post Style

- 2718010 Stroller Handle Back

Back Post Height

 Measured from the Seat Pan to the Handle.

- 2718425 Short (22") 2718431 Tall (25")

Back Post Folding

 Fold-down Back is not compatible with Reclining Backrest, Dynamic Rocker Back and will be substituted with a Non-folding Back.

- 2719220 Fold-down

Seat to Back Angle

 Seat Angles above 98° may cause interference at Seat Depths above 19".

- 2719090 Set at 90° 2719098 Set at 98° 2719106 Set at 106°

- 2719114 Set at 114° 2719122 Set at 122°

Back Options

Angle/Height Adjustable Stroller Handle

 Not compatible with Headrest Mounting Plate.

\$ 150

- 2719010 Button on Outside 2719011 Button on Inside

Endomorph Back

 Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".

\$ 395

Not compatible with Dynamic Rocker Back, Reclining Backrest. Requires Straight Back Post.

- 2718102 1" Offset per side 2718103 1.5" Offset per side 2718104 2" Offset per side

- 2719210 **10" Backrest Extension** Requires Angle/Height Adjustable Stroller Handle. *E0966* \$ 150

- 2719212 **Headrest Mounting Plate** \$ 95

Not compatible with Backrest Extension, Angle/Height Adjustable Stroller Handle. Not required if mounting to a solid Back.

- 2719230 **Spring Back Suspension** Not compatible with Reclining Back. \$ 105

- 2719240 **Seating Dynamics Dynamic Rocker Back Interface** *E2398* \$ 2,995

Not compatible with Fold-down Back, Endomorph Back, 30° Recline, Standard Upholstery. Finished Back Post Height is 4" taller than selected Back Post

 Choose the standard **Bentley** with a steel frame for highly agitated users. See separate order form.

Back Upholstery

- 2720000 **No Back Upholstery**

Standard Upholstery

E2611/E2612 \$ 195

- 2720418 Short (18" Height) 2720421 Regular (21" Height) 2720424 Tall (24" Height)

Tension Adjustable Upholstery

 Height of Upholstery is 16". Not available with Cane Handle Back.

- 2720450 For 14" to 22" Seat Width *E2611/E2612* \$ 295

- 2720451 For 23" to 24" Seat Width *E2611/E2612* \$ 375

See next page for more options >

Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.

22" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>								
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

25" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>								
Regular				<input type="checkbox"/>					
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>						

Tilt Activation

Attendant Tilt only Select trigger mounting location below.

- 2723011 Right Back Post 2723012 Left Back Post
- 2723013 Two Trigger Setup A Trigger mounted on each side of the Stroller Handle. Reduces grip force required to activate tilt. Not compatible with Reclining Backrest, Drum Brakes.

User Tilt only Full Length Arm Pads only. Select trigger mounting location below.

\$ 100

- 2723021 Right Armrest 2723022 Left Armrest

Attendant & User Tilt Not compatible with Tall (25") Back Post, Reclining Backrest. Full Length Arm Pads only. Select trigger mounting locations below.

\$ 195

- 2723031 Right Back Post 2723032 Left Back Post
- 2723041 Right Armrest 2723042 Left Armrest

Power Tilt 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt

Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both. Not available for 12" STF Height.

- 2723400 Standard Power Tilt E1002 \$ 2,295
- 2723410 Power Tilt with Manual Tilt Override Allows for instant manual release of tilt in case of power tilt failure. Minimum Seat Width is 16". E1002 \$ 2,795
- 2723421 Toggle Switch Installed on Right side 2723422 Installed on Left side
 - 2723431 Button Switch Installed on Right side 2723432 Installed on Left side
 - 2723440 Stealth Egg Switch Not installed
 - Optional for Egg Switch: 2723445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. E1028 \$ 625

Specify Headrest to be used: _____

Recline

- 2723310 **Reclining Backrest (Manual)** E1225 \$ 895
30° maximum recline. Not compatible with Backrest Extension, Fold-down Back, Tension Adjustable Upholstery, Dynamic Rocker Back, Endomorph Back. Requires longer Frame. Maximum Seat Depth is 19". Finished Back Post Height is 3" taller than selected Back Post Height.

Front Rigging

- 2725000 **No Front Rigging**
- Easy-Out Front Rigging** Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out

 - 2725410 Fixed 70° Not available for 12" STF Height. \$ 95
 - 2725420 Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°) 12" STF Height is only compatible with 0°-60°. K0108 \$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment.
Optional for Universal AA: 2725421 AA High-mount Option For shorter Lower Leg Lengths of 8"-12.5".
- 2725110 **Standard Plug In 70° Front Rigging** Legacy Design. Not compatible with Easy-Out.

Footplates

- 2726010 **Composite Footplates**
- 2726020 **Aluminum Angle Adjustable Footplates** K0040 \$ 200
- 2726030 **One Piece Fixed Footplate** Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. \$ 350
- 2726040 **One Piece Flip-up Footplate** Requires selection of Composite or Aluminum Angle Adjustable. \$ 395
- 2726050 **Interlocking Angle Adjustable Footplates** Available for Seat Widths up to 21" only. K0040 \$ 200

Footrest/Legrest Accessories

<input type="checkbox"/>	2727012	Calf Pads	Includes reinforced premium Calf Pad mount.		\$ 200
<input type="checkbox"/>	2727032	Gel Knee Pads	Compatible with Universal Angle Adjustable Front Rigging only.		\$ 200
<input type="checkbox"/>	2727040	Length Adjustable Heel Loops		E0951	\$ 50
<input type="checkbox"/>	2727050	Bodypoint Padded Calf Strap		K0038	\$ 195

Front Rigging Extension Tubes

Regular Extension Tubes No selection required for Easy-Out AA High-mount Option.

<input type="checkbox"/>	2728010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Standard Plug In Front Rigging	
<input type="checkbox"/>	2728020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Standard Plug In Front Rigging	\$ 75



Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

		Spoke Guards	Made from durable and lightweight aluminum.	K0065	\$ 195
<input type="checkbox"/>	2730020	Black		<input type="checkbox"/>	2730021 Frame Color
<input type="checkbox"/>	2730050	Protective Lower Padding			\$ 795
		Stealth Flip Down Abductor	Includes wedge-style pad. Requires Solid Seat Pan.	E0957 & E1028	\$ 650
<input type="checkbox"/>	2730060	Small Pad	3"x4"x3"x2" size	<input type="checkbox"/>	2730062 Large Pad 5"x5"x4"x3" size
<input type="checkbox"/>	2730080	Gel Hip Guides	Reduces available Seat Width by 2".	E0956 & E1028	\$ 495
<input type="checkbox"/>	2730110	IV Pole		E0776	\$ 295
<input type="checkbox"/>	2730120	O2 Holder		E2208	\$ 395
<input type="checkbox"/>	2730200	Q'Straint Lap Belt	Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.		\$ 245
<input type="checkbox"/>	2730222	Medical Necessities Backpack	Black with Permobil Logo.		\$ 95
<input type="checkbox"/>	2730169	Rubber Shroud for Caster Journals			\$ 50
<input type="checkbox"/>	2730170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.		\$ 10

On-Chair Seating & Positioning

Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not guarantee coverage or payment.

Please send your completed order form by email to sales.pdg@permobil.com.