

# BODILINK® HEAD SUPPORT ORDER FORM



PRICING IN CANADIAN DOLLARS

Company Name/ACCT #: \_\_\_\_\_

P.O. Number: \_\_\_\_\_

Requested By: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

Ship To: \_\_\_\_\_

Client Reference: \_\_\_\_\_

For best results, do not fill out in your browser. Interactive form should be completed using Adobe Reader after saving to your local drive. Then email or print and fax to Customer Support at sales.ca@permobil.com or 1.877.636.8944.



PT Max 8 Hardware  
Shown on ACTA-BACK®



GT Lateral Hardware  
Shown on ACTA-BACK®

## 1 HEAD SUPPORT PADS

- You may make one selection per column.
- Choose between premium (1.1) -or- basic (1.3).
- Choose a cover material in 1.2 or 1.4.

### 1.1 PREMIUM PAD SIZE MSRP \$247.00

Length (L) refers to the actual size dimension of the support from bottom to top edge.

W		Pad Firmness	
L	6W x 3L	Soft	BL-HSP1-6W3L
		Extra Soft	BL-HSP2-6W3L
L	8W x 4L	Soft	BL-HSP1-8W4L
		Extra Soft	BL-HSP2-8W4L
L	10W x 5L <i>with adjustable wings</i>	Soft	BL-HSP1-10W5L
		Extra Soft	BL-HSP2-10W5L
L	12W x 5L <i>with adjustable wings</i>	Soft	BL-HSP1-12W5L
		Extra Soft	BL-HSP2-12W5L
L	14W x 5L <i>with adjustable wings</i>	Soft	BL-HSP1-14W5L
		Extra Soft	BL-HSP2-14W5L
L	18W x 5L <i>with adjustable wings</i>	Soft	BL-HSP1-18W5L
		Extra Soft	BL-HSP2-18W5L



Premium head support pad angles can be formed by adjusting the plates in the pad.

### 1.2 PREMIUM PAD COVER MATERIAL

	FOR SKIN PROTECTION & SHEAR REDUCTION. NOT FLUID RESISTANT.
	FOR FLUID PROTECTION & AN EASILY CLEANED SURFACE.
	FOR PATIENT COMFORT.

	GlideWear™ (GLIDEWEAR) (add \$82.00)
	COMFORT-TEK (COMFORT-TEK)
	STRETCH-AIR (STRETCH-AIR)

### 1.3 BASIC PAD SIZE MSRP \$247.00

Length (L) refers to the actual size dimension of the support from bottom to top edge.

W		
L	6W x 5L	BL-HSP3-6W5L
L	8W x 5L	BL-HSP3-8W5L
L	10W x 6L	BL-HSP3-10W6L
L	12W x 6L	BL-HSP3-12W6L
L	14W x 6L	BL-HSP3-14W6L
L	18W x 6L	BL-HSP3-18W6L



Basic head support pads can be formed to any shape using the moldable steel insert inside the pad.

### 1.4 BASIC PAD COVER MATERIAL

	FOR FLUID PROTECTION & AN EASILY CLEANED SURFACE.
	FOR PATIENT COMFORT & HEAT DISSIPATION.

	COMFORT-TEK (COMFORT-TEK)
	STRETCH-AIR (STRETCH-AIR)

# 2 HEAD SUPPORT HARDWARE


- Make one selection per column in subsection A, B, or C.
- If A is chosen, choose a stem length option.
- Choose which mounting block option you require in 2.2.

## 2.1 HARDWARE

### PT HARDWARE RANGE OF MOTION


PT MODEL	TOTAL # OF LINKS	TOTAL LINK LENGTH	MAX ANTERIOR REACH	MAX POSTERIOR REACH	MAX LATERAL REACH
MID 6	2	6"	5.5"	1"	N/A
LAT 6	3	6"	5.5"	1"	6"
MID 8	2	8"	7.5"	1"	N/A
LAT 8	3	8"	7.5"	1"	8"
MAX 8	4	8"	7.5"	3"	8"

### A. PT HARDWARE, REMOVABLE



SHAPE	MSRP	PART NUMBER
MID 6	\$355.00	(BL-HS-PT1RMMID6)
LAT 6	\$465.00	(BL-HS-PT1RMLAT6)
MID 8*	\$355.00	(BL-HS-PT1RMMID8)
LAT 8*	\$465.00	(BL-HS-PT1RMLAT8)
MAX 8**	\$519.00	(BL-HS-PT1RMMAX8)

- \*Recommended for ACTA-Relief back supports. See options for mounting hardware selections.
- \*\*Provides the greatest amount of articulation and range of motion.




### PT HARDWARE, REMOVABLE, STEM LENGTH OPTION

STEM LENGTH	MSRP	PART NUMBER
STANDARD 7" STEM	\$0.00	(HS-STEM-7)
***3" STEM	\$0.00	(HS-STEM-3)

- \*\*\*3" Stem is recommended for use with power mounting where there may be special constraints with power chair functions.

### B. PT HARDWARE, FLIP-DOWN


Flip-Down Hardware has 7" Stem



SHAPE	MSRP	PART NUMBER
MID 6	\$399.00	(BL-HS-PT1FDMID6)
LAT 6	\$504.00	(BL-HS-PT1FDLAT6)
MID 8*	\$399.00	(BL-HS-PT1FDMID8)
LAT 8*	\$504.00	(BL-HS-PT1FDLAT8)
MAX 8**	\$562.00	(BL-HS-PT1FDMAX8)

- \*Recommended for ACTA-Relief back supports. See options for mounting hardware selections.
- \*\*Provides the greatest amount of articulation and range of motion.

### C. GT HARDWARE, REMOVABLE



### GT HARDWARE RANGE OF MOTION

GT MODEL	TOTAL # OF LINKS	TOTAL LINK LENGTH	MAX. ANTERIOR REACH	MAX. POSTERIOR REACH	MAX. LATERAL REACH
LAT	5	10.5"	10"	2"	7.5"
MID	3	8.5"	8"	2"	N/A

STYLE	MSRP	PART NUMBER
LATERAL	\$333.00	(BL-HS-GT3RMLAT)
MID-LINE	\$296.00	(BL-HS-GT3RMMID)

## 2.2 HEAD SUPPORT MOUNTING OPTION



**MOUNTING BLOCK-POSTERIOR A or B**


**A:** For ACTA-BACK®, ACTA-EMBRACE®, Visco-Back, and ACTA-Relief with Power Mount only

**B:** For Radius Back® or OEM wood/t-nut back supports




**MOUNTING BLOCK-ANTERIOR A**

For ROHO® AGILITY and other OEM back supports




**MOUNTING BLOCK-QUANTUM A**

For Quantum Tru-Comfort back supports




**MOUNTING BLOCK-PERMOBIL A**

For Permobil Corpus 3G back supports (23"-28"H)



**MOUNTING BLOCK-PERMOBIL B**

For Permobil Corpus 3G back supports (20"H)



**BRACKET WITH MOUNTING BLOCK-ACTA-RELIEF**


18W shown above for Non-Power Mount only

MOUNTING BLOCK OR BRACKET WITH MOUNTING BLOCK	MSRP	PART NUMBER
POSTERIOR A MOUNTING BLOCK	\$0.00	(HS-MB-POSTA)
POSTERIOR B MOUNTING BLOCK	\$0.00	(HS-MB-POSTB)
ANTERIOR A MOUNTING BLOCK	\$0.00	(HS-MB-ANTA)
QUANTUM A MOUNTING BLOCK	\$147.00	(HS-MB-PW1A)
PERMOBIL A MOUNTING BLOCK	\$147.00	(HS-MB-PW2A)
PERMOBIL B MOUNTING BLOCK	\$147.00	(HS-MB-PW2B)
*ACTA-RELIEF 14W BRACKET WITH MOUNTING BLOCK	\$153.00	(HS-MB-AR14W)
*ACTA-RELIEF 16W BRACKET WITH MOUNTING BLOCK	\$153.00	(HS-MB-AR16W)
*ACTA-RELIEF 18W BRACKET WITH MOUNTING BLOCK	\$153.00	(HS-MB-AR18W)
*ACTA-RELIEF 20W BRACKET WITH MOUNTING BLOCK	\$153.00	(HS-MB-AR20W)

\*For ACTA-Relief, includes mounting bracket and Posterior A block. If you already have an ACTA-Relief head support bracket, order Posterior A block.

## 3 HEAD SUPPORT ADDITIONAL MOUNTING HARDWARE

### 3.1 EXTRA LINKS & WING LOCKS FOR PT HARDWARE ONLY



LINK LENGTH	MSRP	EXTRA PART
EXTRA 2" LINK	\$108.00	(BL-HS-PT1LINK2)



EXTRA 3" LINK	\$108.00	(BL-HS-PT1LINK3)
---------------	----------	------------------



EXTRA 4" LINK	\$108.00	(BL-HS-PT1LINK4)
---------------	----------	------------------



TOOL-LESS WING LOCK	MSRP	EXTRA PART
1 LOCK	\$44.00/each	(BL-HS-52WINGA1)

EXTRA PART
(BL-HS-52WINGA1)

EXTRA PART
(BL-HS-52WINGA1)

EXTRA PART
(BL-HS-52WINGA1)

### 3.2 EXTRA MOUNTING HARDWARE



**MOUNTING BLOCK-POSTERIOR A or B**


**A:** For ACTA-BACK®, ACTA-EMBRACE®, Visco-Back, and ACTA-Relief with Power Mount only

**B:** For Radius Back® or OEM wood/t-nut back supports




**MOUNTING BLOCK-ANTERIOR A**

For ROHO® AGILITY and other OEM back supports




**MOUNTING BLOCK-QUANTUM A**

For Quantum Tru-Comfort back supports




**MOUNTING BLOCK-PERMOBIL A**

For Permobil Corpus 3G back supports (23"-28"H)



**MOUNTING BLOCK-PERMOBIL B**

For Permobil Corpus 3G back supports (20"H)



**BRACKET-ACTA-RELIEF**

18W shown above



MOUNTING BLOCK OR BRACKET	MSRP	PART NUMBER
POSTERIOR A MOUNTING BLOCK	\$113.00	(BL-HS-MBPOSTA)
POSTERIOR B MOUNTING BLOCK	\$113.00	(BL-HS-MBPOSTB)
ANTERIOR A MOUNTING BLOCK	\$113.00	(BL-HS-MBANTA)
QUANTUM A MOUNTING BLOCK	\$113.00	(BL-HS-MBPW1A)
PERMOBIL A MOUNTING BLOCK	\$113.00	(BL-HS-MBPW2A)
PERMOBIL B MOUNTING BLOCK	\$113.00	(BL-HS-MBPW2B)
*ACTA-RELIEF 14W BRACKET	\$306.00	(BL-HS-AR14W)
*ACTA-RELIEF 16W BRACKET	\$306.00	(BL-HS-AR16W)
*ACTA-RELIEF 18W BRACKET	\$306.00	(BL-HS-AR18W)
*ACTA-RELIEF 20W BRACKET	\$306.00	(BL-HS-AR20W)

\*Bracket only, does not include mounting block.

ADDITIONAL	MSRP	EXTRA PART
1/4" SPACER PLATE	\$82.00/each	(BL-HS-SPCRA)

The following are registered trademarks of Permobil Seating and Positioning: ACTA-BACK®, ACTA-EMBRACE®, AGILITY®, BodiLink®, COMFORT-TEK®, RADIUS BACK®, and STRETCH-AIR®. GlideWear™ is a trademark of Tamarack Habilitation Technologies, Inc.