

permobil
POWER WHEELCHAIRS

TOUGH ENOUGH



M300
CORPUS **HD**




RUGGED STYLE, REFINED DETAILS

M300 HD takes all the things you love about the Permobil's Mid Wheel Drive and amplifies them. Outfitted with stronger actuators, more powerful motors, and a reinforced seating system; then completed with elegant features like fully adjustable leg rests and up to four power seat functions, the M300 HD is tough enough.

COLOUR OPTIONS



FABRIC OPTIONS



SPECIFICATIONS

Maximum User Weight	204 kg	Seat-to-Floor Height	18", 19", 20" (45.7cm, 48.2cm, 50.8cm)
Maximum Speed	8 km/h	Optional Seat Elevator	8" (20.3cm)
Range	28 km*	Power Tilt	0° – 45° (power)
Base Width	26" (66cm)	Recline Options	85° – 120° (manual) 85° – 150° (power)
Base Length	38" (96.5cm)	Legrest Elevation	85° - 175°
Minimum Turning Radius	22" (56cm)	Backrest Height	21", 24" – 29" (by 1" incr.)
Ground Clearance	3" (7.6cm)	Seat Depth	19" – 25" (by 1" incr.)
Recommended Battery Type	Group 24 & Group 34	Seat Width	19", 21", 23", 25" (48.2cm, 53.3cm, 58.4cm, 63.5cm)
Weight Including Batteries	202 kg with Group 24 193 kg with Group 34	Armrest Pad Lengths	13", 16", 18" (33cm, 40.6cm, 45.7cm)
Drive Electronics	R-Net 120A	*Actual driving range will vary based on driving and battery conditions. †Refer to Owner Manual for safe driving practices.	
Maximum Slope	7.5° Up, 7.5° Down†		