

# M2 ATI ORDER FORM

permobil

PRICING IN CANADIAN DOLLARS

Company Name/ACCT #: \_\_\_\_\_

P.O. Number: \_\_\_\_\_

Requested By: \_\_\_\_\_


Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

Ship To: \_\_\_\_\_

Client Reference: \_\_\_\_\_

For best results, do not fill out in your browser. Interactive form should be completed using Adobe Reader [after](#) saving to your local drive. Then email or print and fax to Customer Support at sales.ca@permobil.com or 1.877.636.8944.



## 1 SIZE, STYLE & COVER

### STYLES

ZERO ELEVATION

ANTI-THRUST

DEPTH	24"								18 x 24	19 x 24	20 x 24	21 x 24	22 x 24	23 x 24	24 x 24	
	23"								18 x 23	19 x 23	20 x 23	21 x 23	22 x 23	23 x 23	24 x 23	
	22"						16 x 22	17 x 22	18 x 22	19 x 22	20 x 22	21 x 22	22 x 22	23 x 22	24 x 22	
	21"					15 x 21	16 x 21	17 x 21	18 x 21	19 x 21	20 x 21	21 x 21	22 x 21	23 x 21	24 x 21	
	20"				14 x 20	15 x 20	16 x 20	17 x 20	18 x 20	19 x 20	20 x 20	21 x 20	22 x 20	23 x 20	24 x 20	
	19"			14 x 19	15 x 19	16 x 19	17 x 19	18 x 19	19 x 19	20 x 19	21 x 19	22 x 19	23 x 19	24 x 19		
	18"			12 x 18	13 x 18	14 x 18	15 x 18	16 x 18	17 x 18	18 x 18	19 x 18	20 x 18	21 x 18	22 x 18	23 x 18	24 x 18
	17"	11 x 17	12 x 17	13 x 17	14 x 17	15 x 17	16 x 17	17 x 17	18 x 17	19 x 17	20 x 17	21 x 17	22 x 17	23 x 17	24 x 17	
	16"	10 x 16	11 x 16	12 x 16	13 x 16	14 x 16	15 x 16	16 x 16	17 x 16	18 x 16	19 x 16	20 x 16	21 x 16	22 x 16	23 x 16	24 x 16
	15"	10 x 15	11 x 15	12 x 15	13 x 15	14 x 15	15 x 15	16 x 15	17 x 15	18 x 15	19 x 15	20 x 15				
14"	10 x 14	11 x 14	12 x 14	13 x 14	14 x 14	15 x 14	16 x 14	17 x 14	18 x 14	19 x 14						
13"	10 x 13	11 x 13	12 x 13	13 x 13	14 x 13	15 x 13	16 x 13	17 x 13								
12"	10 x 12	11 x 12	12 x 12	13 x 12	14 x 12	15 x 12	16 x 12	17 x 12								
11"	10 x 11	11 x 11	12 x 11	13 x 11	14 x 11	15 x 11										
10"	10 x 10	11 x 10	12 x 10	13 x 10	14 x 10	15 x 10										
	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"	
	WIDTH															

### PRICING

CUSHION SIZE	COMFORT-TEK OR STRETCH-AIR	WITH GLIDEWEAR
Standard	\$735.00	\$762.00
Bariatric >20"	\$919.00	\$945.00

### M2 ATI PART NUMBERS

CUSHION STYLE	STRETCH-AIR	COMFORT-TEK
Zero Elevation	M2-SM-WWDD	M2-FM-WWDD
Anti-Thrust	M2-SAM-WWDD	M2-FAM-WWDD
	STRETCH-AIR + GLIDEWEAR	COMFORT-TEK + GLIDEWEAR
Zero Elevation	M2-SQM-WWDD	M2-FQM-WWDD
Anti-Thrust	M2-SQAM-WWDD	M2-FQAM-WWDD

WWDD to be replaced with sizes from table; WW - width in inches & DD - depth in inches

### YOUR ORDER PART #

Please summarize your cushion style, size, and cover options here.

M2-

Example: M2-FQM-1816 for a M2 ATI Zero Elevation cushion with a Comfort-Tek cover & added GlideWear.

### COVER & LINER OPTIONS



#### COMFORT-TEK®

For fluid protection & an easily cleaned surface.



#### STRETCH-AIR®

For patient comfort & heat dissipation.



#### GLIDEWEAR™

Promotes healing & reduces risk of skin breakdown\*

\*A GlideWear® Shear Reduction Panel is incorporated into the cover under the most at-risk areas to allow for micro-movements without the deformation of tissue.  
U.S. GLIDEWEAR® Patent No. 8,646,459



#### INCONTINENCE LINER

Liner provides extra incontinence protection for the cushion. Add at time of cushion purchase for only \$53.00.

(INCON-LINER)

### ACCESSORIES

#### EXTRA COMFORT-TEK OUTER COVER

##### Regular:

(M2--WWDD-CV)

One = \$125.00

Two = \$250.00

##### with GlideWear:

(+CC-GLIDEWEAR)

One = \$150.00

Two = \$300.00

#### EXTRA STRETCH-AIR OUTER COVER

##### Regular:

(M2--WWDD-CV)

One = \$125.00

Two = \$250.00

##### with GlideWear:

(+CC-GLIDEWEAR)

One = \$150.00

Two = \$300.00

#### EXTRA INCONTINENCE LINER

##### Off cushion:

(BINC-M2--WWDD-CV)

One = \$120.00

Two = \$240.00


## 2 M2 ACCESSORIES

\*All accessories are cosmetic changes to the cushion except for options in gray, which are add-on items.

### 2.1 GROWTH NOTCHES MSRP \$0.00

If growth notches are selected, a 1¼" wide cut will be taken from each side to allow the cushion to fit between the chair canes.

#### DEPTH OF NOTCHES




C = 1" D (GROWTH-NOTCH-1)

C = 2" D (GROWTH-NOTCH-2)

### 2.2 RAIL CUTS MSRP \$55.00

A standard 1 1/2" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.




B = 1" Rail Cut (RAIL-CUT-1.0)

B = 1½" Rail Cut (RAIL-CUT-1.5)

B = 2" Rail Cut (RAIL-CUT-2.0)

### 2.3 COLOR - COMES STANDARD SILVER REFLECTIVE PIPING MSRP \$55.00

Accent colors can be applied to the cushion via colored piping. If no colors are selected, our standard silver reflective piping is used.



Red (RED)

Blue (BLUE)

Purple (PURPLE)

Green (GREEN)


Pink (PINK)

Remove (REMOVE)

Orange (ORANGE)

### 2.4 REMOVABLE SOLID SEAT PAN & HARDWARE KIT MSRP \$657.00

The solid seat pan kit includes a slotted aluminum pan and attaching hardware to accommodate 7/8" or 1" tubing. The kit also comes with two different cross bars to accommodate various wheelchair frame types. This will fit both folding and non-folding seat rails. Lateral thigh support hardware can be easily attached.



Solid Seat Pan & Hardware Kit  
(RSSP-N-HARD)


# 3 BODILINK® LATERAL PELVIC/THIGH SUPPORT HARDWARE & PADS

Please select sizes in the appropriate fabric & style charts below. Contact Customer Support for custom sizes that are not listed below at 800.265.9830.

## 3.1 PREMIUM PAD MSRP \$142.00/ea.

- You may make one selection per column (Left, Extra Left, Right, and Extra Right) in section 3.1 for Premium Pads, or 3.3 for Basic Pads.
- Make a pad cover fabric selection in 3.2 or 3.4.

\*Length (L) refers to the actual size dimension of the support from top to bottom edge. Depth (D) refers to the actual size dimension from anterior to posterior edge.

STYLE	SIZE	LEFT	EXTRA LEFT	RIGHT	EXTRA RIGHT
	3.5"L x 4"D	BL-LPTSP2Z-3L4D-LH	BL-LPTSP2Z-3L4D-LH	BL-LPTSP2Z-3L4D-RH	BL-LPTSP2Z-3L4D-RH
	3.5"L x 8"D	BL-LPTSP2Z-3L8D-LH	BL-LPTSP2Z-3L8D-LH	BL-LPTSP2Z-3L8D-RH	BL-LPTSP2Z-3L8D-RH
	3.5"L x 12"D	BL-LPTSP2Z-3L12D-LH	BL-LPTSP2Z-3L12D-LH	BL-LPTSP2Z-3L12D-RH	BL-LPTSP2Z-3L12D-RH
	5.5"L x 6"D	BL-LPTSP2Z-5L6D-LH	BL-LPTSP2Z-5L6D-LH	BL-LPTSP2Z-5L6D-RH	BL-LPTSP2Z-5L6D-RH

## 3.2 PREMIUM PAD COVER WITH FOAM INSERT



**COMFORT-TEK**  
FOR FLUID PROTECTION & AN EASILY CLEANED SURFACE.



**STRETCH-AIR**  
FOR PATIENT COMFORT & HEAT DISSIPATION.




FOR SKIN PROTECTION & SHEAR REDUCTION.  
**NOT FLUID RESISTANT.**

COVER OPTIONS	MSRP	LEFT	EXTRA LEFT	RIGHT	EXTRA RIGHT
COMFORT-TEK	\$0.00/ea.	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM
STRETCH-AIR	\$0.00/ea.	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM
GLIDEWEAR®	\$34.00/ea.	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM

## 3.3 BASIC PAD MSRP \$132.00/ea.

\*Length (L) refers to the actual size dimension of the support from top to bottom edge. Depth (D) refers to the actual size dimension from anterior to posterior edge.

STYLE	SIZE	LEFT	EXTRA LEFT	RIGHT	EXTRA RIGHT
	3"L x 4"D	BL-LPTSP1Z-3L4D-LH	BL-LPTSP1Z-3L4D-LH	BL-LPTSP1Z-3L4D-RH	BL-LPTSP1Z-3L4D-RH
	3"L x 5"D	BL-LPTSP1Z-3L5D-LH	BL-LPTSP1Z-3L5D-LH	BL-LPTSP1Z-3L5D-RH	BL-LPTSP1Z-3L5D-RH
	4"L x 4"D	BL-LPTSP1Z-4L4D-LH	BL-LPTSP1Z-4L4D-LH	BL-LPTSP1Z-4L4D-RH	BL-LPTSP1Z-4L4D-RH
	4"L x 6"D	BL-LPTSP1Z-4L6D-LH	BL-LPTSP1Z-4L6D-LH	BL-LPTSP1Z-4L6D-RH	BL-LPTSP1Z-4L6D-RH
	4"L x 8"D	BL-LPTSP1Z-4L8D-LH	BL-LPTSP1Z-4L8D-LH	BL-LPTSP1Z-4L8D-RH	BL-LPTSP1Z-4L8D-RH
	4"L x 10"D	BL-LPTSP1Z-4L10D-LH	BL-LPTSP1Z-4L10D-LH	BL-LPTSP1Z-4L10D-RH	BL-LPTSP1Z-4L10D-RH
	4"L x 12" D	BL-LPTSP1Z-4L12D-LH	BL-LPTSP1Z-4L12D-LH	BL-LPTSP1Z-4L12D-RH	BL-LPTSP1Z-4L12D-RH
	4"L x 14" D	BL-LPTSP1Z-4L14D-LH	BL-LPTSP1Z-4L14D-LH	BL-LPTSP1Z-4L14D-RH	BL-LPTSP1Z-4L14D-RH
	4"L x 16" D	BL-LPTSP1Z-4L16D-LH	BL-LPTSP1Z-4L16D-LH	BL-LPTSP1Z-4L16D-RH	BL-LPTSP1Z-4L16D-RH
	5"L x 7" D	BL-LPTSP1Z-5L7D-LH	BL-LPTSP1Z-5L7D-LH	BL-LPTSP1Z-5L7D-RH	BL-LPTSP1Z-5L7D-RH

## 3.4 BASIC PAD COVER



**COMFORT-TEK**  
FOR FLUID PROTECTION & AN EASILY CLEANED SURFACE.



**STRETCH-AIR**  
FOR PATIENT COMFORT & HEAT DISSIPATION.




FOR SKIN PROTECTION & SHEAR REDUCTION.  
**NOT FLUID RESISTANT.**

COVER OPTIONS	MSRP	LEFT	EXTRA LEFT	RIGHT	EXTRA RIGHT
COMFORT-TEK	\$0.00/ea.	LPTS-P1-COMFORT	LPTS-P1-COMFORT	LPTS-P1-COMFORT	LPTS-P1-COMFORT
STRETCH-AIR	\$0.00/ea.	LPTS-P1-STRETCH	LPTS-P1-STRETCH	LPTS-P1-STRETCH	LPTS-P1-STRETCH
GLIDEWEAR®	\$34.00/ea.	LPTS-P1-GLIDE	LPTS-P1-GLIDE	LPTS-P1-GLIDE	LPTS-P1-GLIDE


# BODILINK® LATERAL PELVIC/THIGH SUPPORT HARDWARE

- You may make one selection per column (Left, Extra Left, Right, and Extra Right) in section 3.5 from group A, B, -OR- C.
- TT long extension arms will only work with TT style hardware.
- GT hardware extra links will only work with GT style hardware.

## A. TT HARDWARE, SLOT MOUNT




HARDWARE SIZE  
VS. CUSHION THICKNESS




HARDWARE SIZE	MAX CUSHION THICKNESS CLEARANCE	RECOMMENDED CUSHION THICKNESS
SMALL	3.0"	2.0"- 3.0"
MEDIUM	4.5"	3.0"- 4.5"
LARGE	6.5"	4.0"- 6.5"

	SIZE	MSRP	LEFT		EXTRA LEFT		RIGHT		EXTRA RIGHT	
FIXED	SMALL	\$207.00/ea.		BL-LPTS-TT1FXSL1-LH		BL-LPTS-TT1FXSL1-LH		BL-LPTS-TT1FXSL1-RH		BL-LPTS-TT1FXSL1-RH
	MEDIUM	\$207.00/ea.		BL-LPTS-TT1FXSL2-LH		BL-LPTS-TT1FXSL2-LH		BL-LPTS-TT1FXSL2-RH		BL-LPTS-TT1FXSL2-RH
	LARGE	\$207.00/ea.		BL-LPTS-TT1FXSL3-LH		BL-LPTS-TT1FXSL3-LH		BL-LPTS-TT1FXSL3-RH		BL-LPTS-TT1FXSL3-RH
OR REMOVABLE	SMALL	\$296.00/ea.		BL-LPTS-TT1RMSL1-LH		BL-LPTS-TT1RMSL1-LH		BL-LPTS-TT1RMSL1-RH		BL-LPTS-TT1RMSL1-RH
	MEDIUM	\$296.00/ea.		BL-LPTS-TT1RMSL2-LH		BL-LPTS-TT1RMSL2-LH		BL-LPTS-TT1RMSL2-RH		BL-LPTS-TT1RMSL2-RH
	LARGE	\$296.00/ea.		BL-LPTS-TT1RMSL3-LH		BL-LPTS-TT1RMSL3-LH		BL-LPTS-TT1RMSL3-RH		BL-LPTS-TT1RMSL3-RH

## B. TT HARDWARE, POWER MOUNT



HARDWARE SIZE  
VS. CUSHION THICKNESS



HARDWARE SIZE	MAX CUSHION THICKNESS CLEARANCE	RECOMMENDED CUSHION THICKNESS
SMALL	3.5"	2.0"- 3.5"
MEDIUM	5.5"	3.0"- 5.5"
LARGE	7.5"	4.0"- 7.5"

	SIZE	MSRP	LEFT		EXTRA LEFT		RIGHT		EXTRA RIGHT	
FIXED	SMALL	\$268.00/ea.		BL-LPTS-TT1FXPWL1-LH		BL-LPTS-TT1FXPWL1-LH		BL-LPTS-TT1FXPWL1-RH		BL-LPTS-TT1FXPWL1-RH
	MEDIUM	\$268.00/ea.		BL-LPTS-TT1FXPWL2-LH		BL-LPTS-TT1FXPWL2-LH		BL-LPTS-TT1FXPWL2-RH		BL-LPTS-TT1FXPWL2-RH
	LARGE	\$268.00/ea.		BL-LPTS-TT1FXPWL3-LH		BL-LPTS-TT1FXPWL3-LH		BL-LPTS-TT1FXPWL3-RH		BL-LPTS-TT1FXPWL3-RH
OR REMOVABLE	SMALL	\$350.00/ea.		BL-LPTS-TT1RMPWL1-LH		BL-LPTS-TT1RMPWL1-LH		BL-LPTS-TT1RMPWL1-RH		BL-LPTS-TT1RMPWL1-RH
	MEDIUM	\$350.00/ea.		BL-LPTS-TT1RMPWL2-LH		BL-LPTS-TT1RMPWL2-LH		BL-LPTS-TT1RMPWL2-RH		BL-LPTS-TT1RMPWL2-RH
	LARGE	\$350.00/ea.		BL-LPTS-TT1RMPWL3-LH		BL-LPTS-TT1RMPWL3-LH		BL-LPTS-TT1RMPWL3-RH		BL-LPTS-TT1RMPWL3-RH

## TT HARDWARE, POWER MOUNT OPTION

**NOTE:** If choosing the Power Mount Hardware, choose the appropriate power mount type below.

POWER MOUNT	MSRP	PART NUMBER
QUANTUM	add \$0.00	LPTS-TT1-PW1
PERMOBIL	add \$0.00	LPTS-TT1-PW2
QUICKIE	add \$0.00	LPTS-TT1-PW3
ROVI	add \$0.00	LPTS-TT1-PW4
AVID REHAB	add \$0.00	LPTS-TT1-PW5

## TT LONG EXTENSION ARM 4"


MSRP \$55.00/arm




LEFT		EXTRA LEFT		RIGHT		EXTRA RIGHT	
	BL-LPTS-TT1LARM4		BL-LPTS-TT1LARM4		BL-LPTS-TT1LARM4		BL-LPTS-TT1LARM4

## C. GT HARDWARE, ONLY SLOT MOUNT


GT HARDWARE  
SLOT MOUNT- FIXED



GT HARDWARE  
SLOT MOUNT-  
REMOVABLE



NO MEDIAL



1" MEDIAL




HARDWARE	MAX CUSHION THICKNESS CLEARANCE	RECOMMENDED CUSHION THICKNESS
NO MEDIAL	*5.0"	2.5"- 5.0"
1" MEDIAL	*4.0"	2.25"- 4.0"

\*If you need to increase the height or medial reach, add additional links below.

SIZE	MSRP	LEFT	EXTRA LEFT	RIGHT	EXTRA RIGHT
FIXED	\$207.00/ea.	BL-LPTS-GT2FXS-LH	BL-LPTS-GT2FXS-LH	BL-LPTS-GT2FXS-RH	BL-LPTS-GT2FXS-RH
REMOVABLE	\$247.00/ea.	BL-LPTS-GT2RMS-LH	BL-LPTS-GT2RMS-LH	BL-LPTS-GT2RMS-RH	BL-LPTS-GT2RMS-RH

## ADDITIONAL GT OPTIONS



### GT ADDITIONAL HARDWARE LINKS



MSRP \$55.00/link  
(2 Additional 1.25" Links Allowed  
Per Hardware Ordered)

LEFT	EXTRA LEFT	RIGHT	EXTRA RIGHT
BL-LPTS-GT2LINK	BL-LPTS-GT2LINK	BL-LPTS-GT2LINK	BL-LPTS-GT2LINK
BL-LPTS-GT2LINK	BL-LPTS-GT2LINK	BL-LPTS-GT2LINK	BL-LPTS-GT2LINK

## 4 BODILINK® MEDIAL KNEE/THIGH SUPPORT HARDWARE & PADS

### 4.1 PAD SIZE & SHAPE

MSRP \$193.00/ea.

STYLE	SIZE	PART NUMBER
 <b>WEDGE</b>	2.5"W x 3.5"D	BL-MKTSP1-3W4D
	3.5"W x 5"D	BL-MKTSP1-4W5D
	4"W x 6"D	BL-MKTSP1-4W6D
 <b>OVAL</b>	2.5"W x 3.5"D	BL-MKTSP2-3W4D
	3.5"W x 5"D	BL-MKTSP2-4W5D
	5"W x 6"D	BL-MKTSP2-5W6D

### 4.2 PAD COVER MATERIAL

- You may make one selection per column.
- Choose between wedge or oval in 4.1.
- Choose a cover material in 4.2.
- Choose hardware size in 4.3.

**GlideWear**  
INNOVATIVE SKIN PROTECTION

FOR SKIN PROTECTION  
& SHEAR REDUCTION.  
NOT FLUID RESISTANT.

**COMFORT-TEK**



FOR FLUID PROTECTION  
& AN EASILY CLEANED  
SURFACE.

**STRETCH-AIR**




FOR PATIENT COMFORT  
& HEAT DISSIPATION.

GlideWear™ (MKTSP-GLIDE) (add \$34.00)
COMFORT-TEK (MKTSP-COMFORT)
STRETCH-AIR (MKTSP-STRETCH)



### 4.3 MEDIAL KNEE/THIGH SUPPORT SWING-AWAY HARDWARE

MSRP \$266.00/ea.



SIZE	PART NUMBER
SMALL	BL-MKTS-ST1SASL1
MEDIUM	BL-MKTS-ST1SASL2
LARGE	BL-MKTS-ST1SASL3

\*Swing-Away Hardware only, not available in Fixed Hardware  
\*\*MKTS Hardware must be mounted to Removable Solid Seat Pan in section 2.4  
\*\*\*If Leg Length Discrepancy is needed on Removable Solid Seat Pan, Offset Orientation must be selected



SMALL

MEDIUM

LARGE

SIZE	RECOMMENDED INFERIOR THIGH THICKNESS
Small	1.00"- 3.00"
Medium	2.00"- 5.00"
Large	3.00"- 7.00"