

Aluminium adjustable wheelchairs

TiLite Aero T | TiLite Aero Z



permobil

TiLite Aero T

Our most advanced aluminium chair, the Aero T's dual-tube frame minimises weight while maximising performance.

Custom designed and built to your specifications, with available Ergo seating

Find your perfect ride balance with our adjustable centre of gravity

Adjustable seating angle to maximise comfort

Dual-tube geometry provides better load distribution and less frame flex

Adjustable castor fork angle to ensure a precise ride



TiLite Aero Z

Distinguished by a clean, minimal look, the mono-tube frame's 1.25" tubing gives the Aero Z strength and durability.

Custom designed and built to your specifications, ergo seating available

Adjustable centre of gravity to find that perfect ride balance

Adjustable seating angle to maximise comfort

1.25" diameter tubing to maximise performance and minimise visual impact

Adjustable castor fork angle to ensure a precision ride




Tattoo options

- 1 Mosaic
- 2 Flower Garden
- 3 Need for Speed




Frame colours


High gloss paint




Matte textured paint




Candy paint



Metallic paint




Matte paint



Pearlescent paint



Natural titanium



Component colours

Anodised



Specifications

TiLite Aero T

Frame features	Aluminium Dual-Tube Rigid, TiFit, full adjustability
Total chair width*	23"
Seat width	12" – 20"
Seat depth	12" – 20"
Front seat height	15" – 21"
Rear seat height	13" – 20"
Front frame angle	70°, 75°, 80°, 85°, or 90°
Camber	0°, 2°, 4°, or 6°
Centre of gravity adjustment	5"

Folding seat back	Adj. height 9.5" – 22" Adj. angle 80° – 101°
Front seat width (seat taper)	0", 1", 2"
Ergonomic seat	6", 7", 8"
Transit options	Q-Straint lock down approved
Transport weight*	5.4 kg
Weight limit	120 kg Standard 159 kg Heavy duty
Warranty	Five years on frame One year on TiLite components

TiLite Aero Z

Frame features	Aluminium Mono-Tube Rigid, TiFit, full adjustability
Total chair width*	23"
Seat width	10" – 20"
Seat depth	10" – 20"
Front seat height	16" – 21"
Rear seat height	12.5" – 20"
Front frame angle	70°, 75°, 80°, 85°, or 90°
Camber	0°, 2°, 4°, or 6°
Centre of gravity adjustment	5"

Folding seat back	Adj. height 9.5" – 22" Adj. angle 80° – 101°
Front seat width (seat taper)	0", 1.5"
Ergonomic seat	6", 7", 8"
Transit options	Transit Tie-Down Transport Bracket
Transport weight*	5.6 kg
Weight limit	120 kg
Warranty	Five years on frame One year on TiLite components

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.3kg.