

Combination Agreement

The below mentioned products are in principle appropriate to be combined with the below mentioned **max mobility** products provided that the combination is made expertly and according to the restrictions of both manufacturers. The restrictions can be asked at the manufacturers.

SUNRISE MEDICAL product	Max mobility product
Quickie Helium	SmartDrive MX2
Quickie Easy Life	
Quickie Argon 2	
Quickie Neon 2	
Quickie Xenon 2	
Quickie Youngster 3 (24 " wheels)	
Quickie Simba (24" wheels)	
Quickie Impera	
Quickie Emineo	

Please note that mounting and limitations has to be effected according to the Sunrise Medical and max mobility mounting instructions by authorized and skilled staff only.

Sunrise Medical exclude and does not allow mounting the Smart Drive MX2 to a wheelchair equipped with a carbon tube axle.

The maximum values and restrictions (e.g. user weight) mentioned in the user manual of both manufacturers must be observed.

Machining, bending, welding or bracing on any safety relevant components is not allowed.

The combination agreement has been made in two copies, one for each party.

Authorized representatives of the parties have executed this agreement.

SUNRISE MEDICAL Europe
Kahlbachring 2-4
69254 Malsch
Germany

max mobility, LLC
5425 Crossings Blvd.
Antioch, TN37013
USA

Malsch, 9th December 2016



Gerd Weis
(Senior test engineer)

Antioch, 12th December 2016



Ben Hemkens
(Engineering Manager)