



PDG Bentley LT-R - AADL W340



Dealer Name*: _____ Account #*: _____

City, Province, Post Code: _____

Phone: _____ Email: _____ ATP: _____

Rep Name: _____ P.O. Number: _____ Quote Number: _____

Client Reference*: _____ Client Weight: _____

Special considerations that need to be addressed here (e.g., diagnosis): _____ *Required field

October 2024, Revised January 2026. Prices are in Canadian Dollars (CAD)

Wheelchair & Frame

Bentley 20° Manual Tilt Wheelchair Base Model	W340	\$ 2,685.00
Bentley LT-R Lightweight Aluminum Frame and Reverse Wheel Layout Upgrade <i>Reverse wheel layout enhances user access for self-propulsion. Choose standard Bentley with steel frame for agitated or heavier users.</i>		\$ 300.00

Frame Type | All frames offer 300lbs weight capacity.

Fixed Width Frame <i>Select a fixed seat width below.</i>	No charge
Adjustable Width Frame <i>Select an adjustable seat width below.</i>	\$ 200.00

User Weight | Required for setup of tilt uplift force.

User Weight Range:	Up to 75lbs	75lbs - 100lbs	100lbs - 115lbs	115lbs - 150lbs	150lbs - 200lbs	200lbs - 250lbs
	250lbs - 300lbs					

Frame Options

Transit Tie Down System (TTDS) <i>Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details.</i>	No charge
---	-----------

Seat Measurements

Seat Width | Overall chair width is +10.5".

Fixed:	15"	16"	17"	18"	19"		No charge
	20"	21"	22"			W111	\$ 250.00
	23"	24"				W111	\$ 500.00
Adjustable 15"-20":	15"	16"	17"	18"	19"		No charge
	20"					W111	\$ 250.00

Seat Depth

Select Depth:	16"	17"	18"	19"		No charge
	20"	21"	22"		W111	\$ 250.00

16" - 20" seat depths are adjustable.

Frame Finishes

Select Color:	High Gloss Black	Hyper Blue
----------------------	------------------	------------

Armrests

Armrest Type

No Armrests No charge

Single Post, Height Adjustable Armrests No charge

8"-11" Range 10-14" Range 12"-16" Range

Armrest Position: Rear (Default) Front

If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.

Armpads

Regular Armpads No charge

Desk Length 2" x 10" Full Length 2" x 14"

Clothing Guards

Clothing Guards, Single Post Armrests *Not compatible with 8"-11" height range.* \$ 125.00

Armrest Extensions | *Not compatible with flip-up armrests or desk length armpads.*

Armrest Extensions - 3" Forward \$ 75.00

Armrest Extensions - 3" Rearward \$ 75.00

Seat Height

Seat Height - 6" Casters

RW

	Low						High					
20"												
22"												
24"												
	12"	13"	14"	15"	16"	17"	14"	15"	16"	17"	18"	19"
	Seat to Floor Height											

Seat Height - 8" Casters

	Low				High			
24"								
	14"	15"	16"	17"	16"	17"	18"	19"
	Seat to Floor Height							

Seat Pan

No Seat Pan No charge

Solid Aluminium Seat Pan No charge

Rear Casters

Urethane Casters

6" x 1.25" Caster No charge

8" x 1.25" Caster *Only compatible with 24" rear wheels.* No charge

Drive Wheels

Composite Drive Wheel No charge

20" 22" 24"

Drive Wheel Position: Rear (Less Stable) Mid (Default) Front (More Stable)

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference. Front position not compatible with 20"-22" seat depths.

Axles

Quick Release Axles No charge

Fixed Axles No charge

Tires

Urethane Tires No charge

Handrims

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	\$ 175.00

Wheel Locks

User Wheel Locks	No charge
Push to Lock Pull to Lock	

Wheel Lock Options

Wheel Lock Extensions	No charge
-----------------------	-----------

Anti-Tippers

Must select one option.

Front Anti-Tippers, Standard <i>Required safety option.</i>	No charge
---	-----------

Front Anti-Tippers, Spring Loaded <i>Required safety option. Dampens the forward-tip movement of the chair during transfers and allows for improved ground clearance. Not removable.</i>	\$ 200.00
--	-----------

Backrest

Stroller Handle Back

Back Post Height:	Short (22")	Tall (25")			
	<i>Measured from seat pan to handle.</i>				
Folding Style:	Fold-Down				
	<i>Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.</i>				
Seat to Back Angle:	Set at 90°	Set at 98°	Set at 106°	Set at 114°	Set at 122°
	<i>Back angles >98° in combination with seat depths >19" may limit range of tilt.</i>				

Backrest Options

Angle/Height Adjustable Stroller Handle: <i>Not compatible with headrest mounting plate.</i>	No charge
--	-----------

Button on Outside Button on Inside

Spring Back Suspension <i>Not compatible with reclining backrest or dynamic rocker back.</i>	\$ 105.00
--	-----------

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only, Right Back Post	No charge
--------------------------------------	-----------

Attendant Tilt Only, Left Back Post	No charge
-------------------------------------	-----------

Two Trigger Setup <i>A trigger mounted on each side of the stroller handle. Reduces grip force required to activate tilt. Not compatible with reclining backrest or drum brakes.</i>	No charge
--	-----------

Front Rigging

No Front Rigging	No charge
------------------	-----------

Easy-Out Front Rigging | *Effortless insertion and removal. No hanger brackets remain on the frame.*

Fixed 70°	\$ 95.00
-----------	----------

Plug In Front Rigging | *Legacy design. Not compatible with Easy-Out.*

Plug In 70°	No charge
-------------	-----------

Footplates

Composite Footplates		\$ 0.00
Aluminium Angle Adjustable Footplates	W107	\$ 165.00

Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short		No charge
Standard Extension Tubes, Long		\$ 75.00

Footplate Height: Specify:

Ext	Short	11" - 16"	12" - 16"	13" - 16"
	Long	14" - 22"	14" - 19"	14" - 17"
		Standard Easy-Out	60°, 70°, ELR Swing-Away	Standard Plug-In
Legrest				

Front Rigging Accessories

Length Adjustable Heel Loops		No charge
Standard Calf Strap	Basic Seating	\$ 50.00

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes