

BODILINK® LATERAL PELVIC/THIGH SUPPORT ORDER FORM

PRICING IN CANADIAN DOLLARS

permobil



Lateral Pelvic Thigh Support with TT Power Mount Removable Hardware



Removable Solid Seat Pan

Company Name/ACCT #: _____

P.O. Number: _____

Requested By: _____

Phone: _____ Fax: _____

Email: _____

Ship To: _____

Client Reference: _____

For best results, do not fill out in your browser. Interactive form should be completed using Adobe Reader after saving to your local drive. Then email or print and fax to Customer Support at sales.ca@permobil.com or 1.877.636.8944.

1 SUPPORT PADS

Please select sizes in the appropriate fabric & style charts below. Contact Customer Support for custom sizes that are not listed below at 800.265.9830.

1.1 PREMIUM PAD MSRP \$142.00/ea.

- You may make one selection per column (Left, Extra Left, Right, and Extra Right) in section 1.1 for Premium Pads, or 1.3 for Basic Pads.
- Make a pad cover fabric selection in 1.2 or 1.4.

*Length (L) refers to the actual size dimension of the support from top to bottom edge. Depth (D) refers to the actual size dimension from anterior to posterior edge.

STYLE	SIZE	LEFT	EXTRA LEFT	RIGHT	EXTRA RIGHT
	3.5"L x 4"D	BL-LPTSP2Z-3L4D-LH	BL-LPTSP2Z-3L4D-LH	BL-LPTSP2Z-3L4D-RH	BL-LPTSP2Z-3L4D-RH
	3.5"L x 8"D	BL-LPTSP2Z-3L8D-LH	BL-LPTSP2Z-3L8D-LH	BL-LPTSP2Z-3L8D-RH	BL-LPTSP2Z-3L8D-RH
	3.5"L x 12"D	BL-LPTSP2Z-3L12D-LH	BL-LPTSP2Z-3L12D-LH	BL-LPTSP2Z-3L12D-RH	BL-LPTSP2Z-3L12D-RH
	5.5"L x 6"D	BL-LPTSP2Z-5L6D-LH	BL-LPTSP2Z-5L6D-LH	BL-LPTSP2Z-5L6D-RH	BL-LPTSP2Z-5L6D-RH

1.2 PREMIUM PAD COVER WITH FOAM INSERT

COMFORT-TEK®
FOR FLUID PROTECTION & AN
EASILY CLEANED SURFACE.

STRETCH-AIR®
FOR PATIENT COMFORT &
HEAT DISSIPATION.

GlideWear
FOR SKIN PROTECTION
& SHEAR REDUCTION.
NOT FLUID RESISTANT.

COVER OPTIONS	MSRP	LEFT	EXTRA LEFT	RIGHT	EXTRA RIGHT
COMFORT-TEK	\$0.00/ea.	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM
STRETCH-AIR	\$0.00/ea.	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM
GLIDEWEAR™	\$34.00/ea.	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM

1.3 BASIC PAD

MSRP \$132.00/ea.

*Length (L) refers to the actual size dimension of the support from top to bottom edge. Depth (D) refers to the actual size dimension from anterior to posterior edge.

STYLE	SIZE	LEFT	EXTRA LEFT	RIGHT	EXTRA RIGHT
	3"L x 4"D	BL-LPTSP1Z-3L4D-LH	BL-LPTSP1Z-3L4D-LH	BL-LPTSP1Z-3L4D-RH	BL-LPTSP1Z-3L4D-RH
	3"L x 5"D	BL-LPTSP1Z-3L5D-LH	BL-LPTSP1Z-3L5D-LH	BL-LPTSP1Z-3L5D-RH	BL-LPTSP1Z-3L5D-RH
	4"L x 4"D	BL-LPTSP1Z-4L4D-LH	BL-LPTSP1Z-4L4D-LH	BL-LPTSP1Z-4L4D-RH	BL-LPTSP1Z-4L4D-RH
	4"L x 6"D	BL-LPTSP1Z-4L6D-LH	BL-LPTSP1Z-4L6D-LH	BL-LPTSP1Z-4L6D-RH	BL-LPTSP1Z-4L6D-RH
	4"L x 8"D	BL-LPTSP1Z-4L8D-LH	BL-LPTSP1Z-4L8D-LH	BL-LPTSP1Z-4L8D-RH	BL-LPTSP1Z-4L8D-RH
	4"L x 10"D	BL-LPTSP1Z-4L10D-LH	BL-LPTSP1Z-4L10D-LH	BL-LPTSP1Z-4L10D-RH	BL-LPTSP1Z-4L10D-RH
	4"L x 12"D	BL-LPTSP1Z-4L12D-LH	BL-LPTSP1Z-4L12D-LH	BL-LPTSP1Z-4L12D-RH	BL-LPTSP1Z-4L12D-RH
	4"L x 14"D	BL-LPTSP1Z-4L14D-LH	BL-LPTSP1Z-4L14D-LH	BL-LPTSP1Z-4L14D-RH	BL-LPTSP1Z-4L14D-RH
	4"L x 16"D	BL-LPTSP1Z-4L16D-LH	BL-LPTSP1Z-4L16D-LH	BL-LPTSP1Z-4L16D-RH	BL-LPTSP1Z-4L16D-RH
	5"L x 7"D	BL-LPTSP1Z-5L7D-LH	BL-LPTSP1Z-5L7D-LH	BL-LPTSP1Z-5L7D-RH	BL-LPTSP1Z-5L7D-RH

1.4 BASIC PAD COVER

COMFORT-TEK
FOR FLUID PROTECTION & AN EASILY CLEANED SURFACE.

STRETCH-AIR
FOR PATIENT COMFORT & HEAT DISSIPATION.

GlideWear
FOR SKIN PROTECTION & SHEAR REDUCTION. NOT FLUID RESISTANT.


COVER OPTIONS	MSRP	LEFT	EXTRA LEFT	RIGHT	EXTRA RIGHT
COMFORT-TEK	\$0.00/ea.	LPTS-P1-COMFORT	LPTS-P1-COMFORT	LPTS-P1-COMFORT	LPTS-P1-COMFORT
STRETCH-AIR	\$0.00/ea.	LPTS-P1-STRETCH	LPTS-P1-STRETCH	LPTS-P1-STRETCH	LPTS-P1-STRETCH
GLIDEWEAR	\$34.00/ea.	LPTS-P1-GLIDE	LPTS-P1-GLIDE	LPTS-P1-GLIDE	LPTS-P1-GLIDE

2 MOUNTING HARDWARE


2.1 BODILINK LATERAL PELVIC/THIGH SUPPORT HARDWARE

- You may make one selection per column (Left, Extra Left, Right, and Extra Right) in section 2.1 from group A, B, -OR- C.
- TT long extension arms will only work with TT style hardware.
- GT hardware extra links will only work with GT style hardware.

A. TT HARDWARE, SLOT MOUNT



HARDWARE SIZE VS. CUSHION THICKNESS




HARDWARE SIZE	MAX CUSHION THICKNESS CLEARANCE	RECOMMENDED CUSHION THICKNESS
SMALL	3.0"	2.0" - 3.0"
MEDIUM	4.5"	3.0" - 4.5"
LARGE	6.5"	4.0" - 6.5"


	SIZE	MSRP	LEFT	EXTRA LEFT	RIGHT	EXTRA RIGHT
FIXED	SMALL	\$207.00/ea.	BL-LPTS-TT1FXSL1-LH	BL-LPTS-TT1FXSL1-LH	BL-LPTS-TT1FXSL1-RH	BL-LPTS-TT1FXSL1-RH
	MEDIUM	\$207.00/ea.	BL-LPTS-TT1FXSL2-LH	BL-LPTS-TT1FXSL2-LH	BL-LPTS-TT1FXSL2-RH	BL-LPTS-TT1FXSL2-RH
	LARGE	\$207.00/ea.	BL-LPTS-TT1FXSL3-LH	BL-LPTS-TT1FXSL3-LH	BL-LPTS-TT1FXSL3-RH	BL-LPTS-TT1FXSL3-RH
OR REMOVABLE	SMALL	\$296.00/ea.	BL-LPTS-TT1RMSL1-LH	BL-LPTS-TT1RMSL1-LH	BL-LPTS-TT1RMSL1-RH	BL-LPTS-TT1RMSL1-RH
	MEDIUM	\$296.00/ea.	BL-LPTS-TT1RMSL2-LH	BL-LPTS-TT1RMSL2-LH	BL-LPTS-TT1RMSL2-RH	BL-LPTS-TT1RMSL2-RH
	LARGE	\$296.00/ea.	BL-LPTS-TT1RMSL3-LH	BL-LPTS-TT1RMSL3-LH	BL-LPTS-TT1RMSL3-RH	BL-LPTS-TT1RMSL3-RH

B. TT HARDWARE, POWER MOUNT


TT HARDWARE
POWER MOUNT- FIXED



TT HARDWARE
POWER MOUNT- REMOVABLE



HARDWARE SIZE
VS. CUSHION THICKNESS



HARDWARE SIZE	MAX CUSHION THICKNESS CLEARANCE	RECOMMENDED CUSHION THICKNESS
SMALL	3.5"	2.0"- 3.5"
MEDIUM	5.5"	3.0"- 5.5"
LARGE	7.5"	4.0"- 7.5"

	SIZE	MSRP	LEFT		EXTRA LEFT		RIGHT		EXTRA RIGHT	
FIXED	SMALL	\$268.00/ea.		BL-LPTS-TT1FXPWL1-LH		BL-LPTS-TT1FXPWL1-LH		BL-LPTS-TT1FXPWL1-RH		BL-LPTS-TT1FXPWL1-RH
	MEDIUM	\$268.00/ea.		BL-LPTS-TT1FXPWL2-LH		BL-LPTS-TT1FXPWL2-LH		BL-LPTS-TT1FXPWL2-RH		BL-LPTS-TT1FXPWL2-RH
	LARGE	\$268.00/ea.		BL-LPTS-TT1FXPWL3-LH		BL-LPTS-TT1FXPWL3-LH		BL-LPTS-TT1FXPWL3-RH		BL-LPTS-TT1FXPWL3-RH
OR REMOVABLE	SMALL	\$350.00/ea.		BL-LPTS-TT1RMPWL1-LH		BL-LPTS-TT1RMPWL1-LH		BL-LPTS-TT1RMPWL1-RH		BL-LPTS-TT1RMPWL1-RH
	MEDIUM	\$350.00/ea.		BL-LPTS-TT1RMPWL2-LH		BL-LPTS-TT1RMPWL2-LH		BL-LPTS-TT1RMPWL2-RH		BL-LPTS-TT1RMPWL2-RH
	LARGE	\$350.00/ea.		BL-LPTS-TT1RMPWL3-LH		BL-LPTS-TT1RMPWL3-LH		BL-LPTS-TT1RMPWL3-RH		BL-LPTS-TT1RMPWL3-RH


TT HARDWARE, POWER MOUNT OPTION

NOTE: If choosing the Power Mount Hardware, choose the appropriate power mount type below.

POWER MOUNT	MSRP	PART NUMBER
QUANTUM	add \$0.00	LPTS-TT1-PW1
PERMOBIL	add \$0.00	LPTS-TT1-PW2
QUICKIE	add \$0.00	LPTS-TT1-PW3
ROVI	add \$0.00	LPTS-TT1-PW4
AVID REHAB	add \$0.00	LPTS-TT1-PW5

TT LONG EXTENSION ARM 4"

MSRP \$55.00/arm




LEFT		EXTRA LEFT		RIGHT		EXTRA RIGHT	
	BL-LPTS-TT1LARM4		BL-LPTS-TT1LARM4		BL-LPTS-TT1LARM4		BL-LPTS-TT1LARM4


QUESTIONS? PLEASE CONTACT CUSTOMER SUPPORT FOR ASSISTANCE. | 800.265.9830

C. GT HARDWARE, ONLY SLOT MOUNT


GT HARDWARE
SLOT MOUNT- FIXED



GT HARDWARE
SLOT MOUNT- REMOVABLE



NO MEDIAL



1" MEDIAL




HARDWARE	MAX CUSHION THICKNESS CLEARANCE	RECOMMENDED CUSHION THICKNESS
NO MEDIAL	*5.0"	2.5"- 5.0"
1" MEDIAL	*4.0"	2.25"- 4.0"

*If you need to increase the height or medial reach, add additional links below.

SIZE	MSRP	LEFT		EXTRA LEFT		RIGHT		EXTRA RIGHT	
FIXED	\$207.00/ea.		BL-LPTS-GT2FXS-LH		BL-LPTS-GT2FXS-LH		BL-LPTS-GT2FXS-RH		BL-LPTS-GT2FXS-RH
REMOVABLE	\$247.00/ea.		BL-LPTS-GT2RMS-LH		BL-LPTS-GT2RMS-LH		BL-LPTS-GT2RMS-RH		BL-LPTS-GT2RMS-RH

ADDITIONAL GT OPTIONS



GT ADDITIONAL HARDWARE LINKS
MSRP \$55.00/link (2 additional 1.25" links allowed per hardware ordered)

LEFT		EXTRA LEFT		RIGHT		EXTRA RIGHT	
	BL-LPTS-GT2LINK		BL-LPTS-GT2LINK		BL-LPTS-GT2LINK		BL-LPTS-GT2LINK
	BL-LPTS-GT2LINK		BL-LPTS-GT2LINK		BL-LPTS-GT2LINK		BL-LPTS-GT2LINK

3 SOLID SEAT PANS

3.1


REMOVABLE SOLID SEAT PAN & HARDWARE KIT

MSRP \$657.00

The solid seat pan kit includes a slotted aluminum pan and attaching hardware to accommodate 7/8" or 1" tubing. The kit also comes with two different cross bars to accommodate various wheelchair frame types. This will fit both folding and non-folding seat rails. Lateral thigh support hardware can be easily attached.

Please fill in the below with dimensions. Example: RRSP-1-1416 for a 14"W x 16"D solid seat pan.

RSSP-1- _____



SIZE RANGE:
Width: 10" - 28"
Depth: 10" - 22"