

Order form

Custom ACTA-BACK®

HCPCS code: E2617 / Custom back

"Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process."

Customer information

Date: _____

Company name/ACCT #: _____

Client name: _____

Client weight: _____

Client height: _____

Client reference: _____

W/C manufacturer: _____


Model #: _____ Serial #: _____

Seat Frame Width: _____

Seat Frame Depth: _____

Seat Rail Tube Size: _____

Seating system measurements



NOTE: Seating system measurements do not affect any selections made throughout the order form.

A - Back support width	
B - Back support length	
C - Width between lateral trunk supports	
D - Width between lateral thigh supports	
E - Seat width	
F - Seat depth	
G - Effective seat depth	
H - Lateral trunk support height	

P.O. number: _____

Requested by: _____

Phone: _____ Fax: _____

Email: _____

Ship to: _____

Additional information:


1 Custom ACTA-BACK size & cover material

Instructions: select the desired back support size by clicking on the corresponding box below or filling in the specific width and length measurements in the "custom back part number" box.

1.1 Back support size

HCPCS code: E2617

Example of ACTA-BACK geometry



*Length (L) refers to the actual size dimension of the support from bottom to top edge.

Length*	\$1316 ⁰⁰	22"	C-A-10W22H	C-A-12W22H	C-A-14W22H	C-A-16W22H	C-A-18W22H	C-A-20W22H	C-A-22W22H	C-A-24W22H	C-A-26W22H	C-A-28W22H
		20"	C-A-10W20H	C-A-12W20H	C-A-14W20H	C-A-16W20H	C-A-18W20H	C-A-20W20H	C-A-22W20H	C-A-24W20H	C-A-26W20H	C-A-28W20H
		18"	C-A-10W18H	C-A-12W18H	C-A-14W18H	C-A-16W18H	C-A-18W18H	C-A-20W18H	C-A-22W18H	C-A-24W18H	C-A-26W18H	C-A-28W18H
		16"	C-A-10W16H	C-A-12W16H	C-A-14W16H	C-A-16W16H	C-A-18W16H	C-A-20W16H	C-A-22W16H	C-A-24W16H	C-A-26W16H	C-A-28W16H
		14"	C-A-10W14H	C-A-12W14H	C-A-14W14H	C-A-16W14H	C-A-18W14H	C-A-20W14H	C-A-22W14H	C-A-24W14H	C-A-26W14H	C-A-28W14H
		12"	C-A-10W12H	C-A-12W12H	C-A-14W12H	C-A-16W12H	C-A-18W12H	C-A-20W12H	C-A-22W12H	C-A-24W12H	C-A-26W12H	C-A-28W12H
		10"	C-A-10W10H	C-A-12W10H	C-A-14W10H	C-A-16W10H	C-A-18W10H	C-A-20W10H	C-A-22W10H	C-A-24W10H	C-A-26W10H	C-A-28W10H
			10"	12"	14"	16"	18"	20"	22"	24"	26"	28"
			\$1164 ⁰⁰						\$1278 ⁰⁰		\$1887 ⁰⁰	
			Frame width									

Custom back support part number:

C - A - _ W _ H

NOTE: Desired back support length is "H" dimension in part number above.

10-20" Width / 10-18" Length = \$1164 base price
22-24" Width / 10-18" Length = \$1278 base price
10-24" Width / 19-22" Length = \$1316 base price
26-28" Width / 10-22" Length = \$1887 base price
*If desired size is not listed above, a 60% up-charge will be applied for materials.

If an odd width is needed, go 1" wider and choose the next even width.

Please note that 10L sizes are not available with lateral trunk supports or a head support.


1.2 Cover material options

REQUIRED

MSRP \$0.00



- OR -




COVER FABRIC	PART NUMBER
STRETCH-AIR®	(CD-STRETCH-AIR)
COMFORT-TEK®	(CD-COMFORT-TEK)

Cover appearance may vary with hardware selection.

1.3 Military branch patch

OPTIONAL

MSRP \$0.00



Air Force	(CD-AIR-FORCE-PATCH)
Army	(CD-ARMY-PATCH)
Marines	(CD-MARINES-PATCH)
Navy	(CD-NAVY-PATCH)

2 Foam options

- Select only one type of foam per layer
- A minimum of one layer must be chosen if Foam in Place is not selected in Section 2.2.



Indicates standard ACTA-BACK foam configuration.
May use as a reference when choosing foam layers.




2.1 Custom foam selection

Rear layer
REQUIRED
H.R. MSRP \$0.00
Visco MSRP \$166.00

Middle layer
OPTIONAL
H.R. MSRP \$0.00
Visco MSRP \$166.00

Topper layer
OPTIONAL
H.R. MSRP \$0.00
Visco MSRP \$166.00

NOTE: 1/4" Crosslink foundation layer will come standard on all back supports. For other foam combinations, call Customer Support 800.736.0925.



Foam type	1/4"	1/2"	3/4"	1"	1 1/4"	1/4"	1/2"	3/4"	1"	1 1/4"	1/4"	1/2"	3/4"	1"	1 1/4"
	(CD-Rear-HRM-0.25)	(CD-Rear-HRM-0.50)	(CD-Rear-HRM-0.75)	(CD-Rear-HRM-1.00)	(CD-Rear-HRM-1.25)	(CD-Mid-HRM-0.25)	(CD-Mid-HRM-0.50)	(CD-Mid-HRM-0.75)	(CD-Mid-HRM-1.00)	(CD-Mid-HRM-1.25)	(CD-Top-HRM-0.25)	(CD-Top-HRM-0.50)	(CD-Top-HRM-0.75)	(CD-Top-HRM-1.00)	(CD-Top-HRM-1.25)
H.R. Medium	(CD-Rear-HRS-0.25)	(CD-Rear-HRS-0.50)	(CD-Rear-HRS-0.75)	(CD-Rear-HRS-1.00)	(CD-Rear-HRS-1.25)	(CD-Mid-HRS-0.25)	(CD-Mid-HRS-0.50)	(CD-Mid-HRS-0.75)	(CD-Mid-HRS-1.00)	(CD-Mid-HRS-1.25)	(CD-Top-HRS-0.25)	(CD-Top-HRS-0.50)	(CD-Top-HRS-0.75)	(CD-Top-HRS-1.00)	(CD-Top-HRS-1.25)
H.R. Soft															
Visco Medium		(CD-Rear-VIS-M-0.50)		(CD-Rear-VIS-M-1.00)			(CD-Mid-VIS-M-0.50)		(CD-Mid-VIS-M-1.00)			(CD-Top-VIS-M-0.50)		(CD-Top-VIS-M-1.00)	
Visco Soft		(CD-Rear-VIS-S-0.50)		(CD-Rear-VIS-S-1.00)			(CD-Mid-VIS-S-0.50)		(CD-Mid-VIS-S-1.00)			(CD-Top-VIS-S-0.50)		(CD-Top-VIS-S-1.00)	
Visco X-Soft		(CD-Rear-VIS-XS-0.50)		(CD-Rear-VIS-XS-1.00)			(CD-Mid-VIS-XS-0.50)		(CD-Mid-VIS-XS-1.00)			(CD-Top-VIS-XS-0.50)		(CD-Top-VIS-XS-1.00)	
Visco XX-Soft		(CD-Rear-VIS-XXS-0.50)		(CD-Rear-VIS-XXS-1.00)			(CD-Mid-VIS-XXS-0.50)		(CD-Mid-VIS-XXS-1.00)			(CD-Top-VIS-XXS-0.50)		(CD-Top-VIS-XXS-1.00)	


- OR -

2.2 Foam in place kit

- Select extra soft, soft or medium density memory foam bag pour.
- Select the appropriate pour size based on your choice of a 3" or 4" deep back support.
- Contact Customer Support with any questions 800.736.0925

KIT INCLUDES: (Note: Price of Visco layers will be added to the order separately.)

- Back cushion cover with an optional 3" or 4" deep back, with a stretchy, pleated contact surface to allow for foam growth.
- Two aluminum stays sewn into the rear of the cover which can be molded to support spinal curvature prior to pouring FIP Chemicals.
- 1" of Visco Medium foam is cut to fit the shape of the back support and is placed loosely inside the cover. This will give the user some padding, prior to pouring the foam.
- 1/2" of Visco Soft foam is cut to fit the shape of the back support and is placed inside the front of the cover. This will give the final mold a plush feel for the user and will also smooth out any potential inconsistencies after pouring the foam.
- Two plastic bags are included to contain the foam chemicals as they expand.



NOTE: BAG POUR SIZE RECOMMENDATION CHART

Border thickness	Frame width	Recommended bag size
3"	10"W-16"W	1 BAG
3"	17"W-22"W	1 1/2 BAG
4"	10"W-14"W	1 BAG
4"	15"W-20"W	1 1/2 BAG

Extra soft memory foam

1 Bag FIP MSRP: \$371.00
1 1/2 Bag FIP MSRP: \$494.00

3" Back border	Pour size	4" Back border
(CD-FIP3-1-XS)	1 BAG	(CD-FIP4-1-XS)
(CD-FIP3-1.5-XS)	1 1/2 BAG	(CD-FIP4-1.5-XS)

Soft memory foam

1 Bag FIP MSRP: \$371.00
1 1/2 Bag FIP MSRP: \$494.00

3" Back border	Pour size	4" Back border
(CD-FIP3-1-S)	1 BAG	(CD-FIP4-1-S)
(CD-FIP3-1.5-S)	1 1/2 BAG	(CD-FIP4-1.5-S)

Medium memory foam

1 Bag FIP MSRP: \$371.00
1 1/2 Bag FIP MSRP: \$494.00

3" Back border	Pour size	4" Back border
(CD-FIP3-1-M)	1 BAG	(CD-FIP4-1-M)
(CD-FIP3-1.5-M)	1 1/2 BAG	(CD-FIP4-1.5-M)


3 Back support options

- A standard alignment has been optimized based on back support width if you choose not to select a specific width between the stays.
- Select the spacing between the moldable aluminum stays.
- Custom options are available by selecting your desired spacing below.

3.1 Moldable aluminum stays

NOTE:

- This back support comes with two patented aluminum stays that can be molded to support spinal curvature.
- Stays are offset from the spine and are held in place by pockets sewn into the rear of the cover.
- The stays are accessible while the user is seated against the back support through the zipper located at the top of the back support.



Space between stays

MSRP \$0.00


A = 1" (CD-STAY-1.0)
A = 1 1/2" (CD-STAY-1.5)
A = 2" (standard alignment) (CD-STAY-2.0)
A = 2 1/2" (CD-STAY-2.5)
A = 3" (CD-STAY-3.0)
A = 3 1/2" (CD-STAY-3.5)
A = 4" (CD-STAY-4.0)
A = 4 1/2" (CD-STAY-4.5)
Specify custom spacing _____ (CD-STAY-OTHER)
None (CD-STAY-NONE)

3.2 Hardware options

4 Point hardware may not fit all smaller back support lengths.


NOTE: Quick release and fixed hardware are adjustable for height, depth, width and angle of recline. Fixed requires tools for removal. 4-Point is designed to provide extra support for users with tone and for back supports 24" wide or greater. Compass 4 hardware may be installed on 3/4", 7/8", 1" and 1 1/8" tubing. Power mount bracket is height adjustable up to 8" in 1/2" increments but may be limited by the chair. 4 point hardware is required on backs 24-28W.

Compass® 4 Quick release (2 points of attachment) MSRP \$0.00




Quick release (2 point) (CD-Q-COMPASS)

Compass 4 Quick release (4 points of attachment) MSRP \$0.00




Quick release (4 point) (CD-Q-COMPASS4)

Compass 4 Fixed (2 points of attachment) MSRP \$0.00



Fixed (2 point) (CD-NQ-COMPASS)

Compass 4 Fixed (4 points of attachment) MSRP \$0.00



Fixed (4 point) (CD-NQ-COMPASS4)

Compass Power Mount Hardware MSRP \$618.00

(Not available on back supports smaller than 12"L or 12" W, or sizes 12Wx16L & 12Wx18L)

Power mount bracket required for all power chairs listed below.

Quantum w/ Tru Balance 3 (CD-PMB-QUANTUM)
Permobil w/ Corpus 3G (CD-PMB-PERMOBIL)
Quickie w/ 3.7 Power (CD-PMB-QUICKIE)
Avid Rehab w/ Power T & R (CD-PMB-AVID)
Rovi w/ Motion Concepts* (CD-PMB-ROVI)

*Using Compass Power Mount on Rovi requires mounting to the RECESSED PLANAR INTERFACE PLATE available when ordering the Motion Concepts Seating System. Drilling is required for mounting.

3.3 Color *OPTIONAL* MSRP \$0.00

NOTE: Accent colors can be applied via colored piping that is stitched to the flaps of the back support. If no colors are selected, our standard silver reflective piping will be used.



Placement of the piping may vary by back support size and mounting hardware.

Red (CD-RED)	Blue (CD-BLUE)
Purple (CD-PURPLE)	Green (CD-GREEN)
Pink (CD-PINK)	Remove (CD-REMOVE)
Orange (CD-ORANGE)	

3.4 Embroidery *OPTIONAL* MSRP \$0.00

Fill out sections A, B, and C.

A.

Thread color options

Red (CD-EMB-RED)	Blue (CD-EMB-BLUE)
Purple (CD-EMB-PURPLE)	Green (CD-EMB-GREEN)
Pink (CD-EMB-PINK)	Silver (CD-EMB-SILVER)

B.

Font/style options

Helvetica (CD-EMB-Helvetica)	Twilight (CD-EMB-Twilight)
Garamond (CD-EMB-Garamond)	Script (CD-EMB-Script)

C.

Word one

(Input is case sensitive. Please enter exactly as you want the words to look. Example: Marcus or marcus or MARCUS)

- AND/OR -


Word two

(Input is case sensitive. Please enter exactly as you want the words to look. Example: Marcus or marcus or MARCUS)

Order Number: _____

3.5 Privacy flap *OPTIONAL* MSRP \$44.00

NOTE: This flap connects with hook & loop attachment between the seat and back support for added privacy.



Privacy flap (CD-PF-CV)

4 BodiLink® Lateral Trunk Support pads & hardware

4.1 Premium pad size & shape

HPCPS code: E0956 MSRP \$101.00/ea.

- You may make one selection per column in section 4.1.
- Follow each column down to sections 4.2 and 4.3 and select a style and cover for each pad selected in section 4.1.

*Length (L) refers to the dimension of the support from top to bottom edge. Depth (D) refers to the dimension of the support from anterior to posterior edge.

	Shape	Size	Left	*Extra left	Right	*Extra right
Supportive		3"L x 4"D	BL-LTSP2S-3L4D-LH	BL-LTSP2S-3L4D-LH	BL-LTSP2S-3L4D-RH	BL-LTSP2S-3L4D-RH
		3"L x 5"D	BL-LTSP2S-3L5D-LH	BL-LTSP2S-3L5D-LH	BL-LTSP2S-3L5D-RH	BL-LTSP2S-3L5D-RH
		4"L x 5"D	BL-LTSP2S-4L5D-LH	BL-LTSP2S-4L5D-LH	BL-LTSP2S-4L5D-RH	BL-LTSP2S-4L5D-RH
		4"L x 6.5"D	BL-LTSP2S-4L6D-LH	BL-LTSP2S-4L6D-LH	BL-LTSP2S-4L6D-RH	BL-LTSP2S-4L6D-RH
		5"L x 6.5"D	BL-LTSP2S-5L6D-LH	BL-LTSP2S-5L6D-LH	BL-LTSP2S-5L6D-RH	BL-LTSP2S-5L6D-RH
		6"L x 6.5"D	BL-LTSP2S-6L6D-LH	BL-LTSP2S-6L6D-LH	BL-LTSP2S-6L6D-RH	BL-LTSP2S-6L6D-RH
Active		3"L x 4"D	BL-LTSP2A-3L4D-LH	BL-LTSP2A-3L4D-LH	BL-LTSP2A-3L4D-RH	BL-LTSP2A-3L4D-RH
		3"L x 5"D	BL-LTSP2A-3L5D-LH	BL-LTSP2A-3L5D-LH	BL-LTSP2A-3L5D-RH	BL-LTSP2A-3L5D-RH
		4"L x 5"D	BL-LTSP2A-4L5D-LH	BL-LTSP2A-4L5D-LH	BL-LTSP2A-4L5D-RH	BL-LTSP2A-4L5D-RH
		4"L x 6.5"D	BL-LTSP2A-4L6D-LH	BL-LTSP2A-4L6D-LH	BL-LTSP2A-4L6D-RH	BL-LTSP2A-4L6D-RH
		5"L x 6.5"D	BL-LTSP2A-5L6D-LH	BL-LTSP2A-5L6D-LH	BL-LTSP2A-5L6D-RH	BL-LTSP2A-5L6D-RH
		6"L x 6.5"D	BL-LTSP2A-6L6D-LH	BL-LTSP2A-6L6D-LH	BL-LTSP2A-6L6D-RH	BL-LTSP2A-6L6D-RH


**EXTRA" parts may not fit depending on back support, pad, and hardware size and style combinations. Please contact customer support with questions.

4.2 Premium pad style

MSRP \$0.00/ea.

	Style options	Left	*Extra left	Right	*Extra right
Or	MOLDABLE	LTS-P2-MOLD	LTS-P2-MOLD	LTS-P2-MOLD	LTS-P2-MOLD
	RIGID	LTS-P2-RIGID	LTS-P2-RIGID	LTS-P2-RIGID	LTS-P2-RIGID


4.3 Premium pad cover




Cover options	MSRP	Left	*Extra left	Right	*Extra right
No cover	\$0.00/ea.	LTS-P2-NO-CV	LTS-P2-NO-CV	LTS-P2-NO-CV	LTS-P2-NO-CV
COMFORT-TEK	\$0.00/ea.	LTS-P2-COMFORT	LTS-P2-COMFORT	LTS-P2-COMFORT	LTS-P2-COMFORT
STRETCH-AIR	\$0.00/ea.	LTS-P2-STRETCH	LTS-P2-STRETCH	LTS-P2-STRETCH	LTS-P2-STRETCH
GlideWear	\$15.00/ea.	LTS-P2-GLIDE	LTS-P2-GLIDE	LTS-P2-GLIDE	LTS-P2-GLIDE
COMFORT-TEK with Foam	\$15.00/ea.	LTS-P2-COMFORT-FM	LTS-P2-COMFORT-FM	LTS-P2-COMFORT-FM	LTS-P2-COMFORT-FM
STRETCH-AIR with Foam	\$15.00/ea.	LTS-P2-STRETCH-FM	LTS-P2-STRETCH-FM	LTS-P2-STRETCH-FM	LTS-P2-STRETCH-FM
GlideWear with Foam	\$31.00/ea.	LTS-P2-GLIDE-FM	LTS-P2-GLIDE-FM	LTS-P2-GLIDE-FM	LTS-P2-GLIDE-FM

4.4 GT2 Mounting hardware

Both fixed and swing-away hardware styles are adjustable.




Swing-away shown



Slot mount



Fixed shown



****Track mount**
Swing-away shown

- You may make one selection per column in section 4.4.
- **TRACK MOUNT** is only available when using Compass 4 hardware. If using Compass Power Mount Hardware, select slot mounting lateral hardware.
- "Extra" track mount hardware does not come with an extra track
- Extra links and long mounting bracket (section 4.5) only for GT2 style hardware.

Slot mounting	MSRP	Left	*Extra left	Right	*Extra right
Swing-away HPCPS Code: E1028	\$339.00/ea.	BL-LTS-GT2SAS-LH	BL-LTS-GT2SAS-LH	BL-LTS-GT2SAS-RH	BL-LTS-GT2SAS-RH
Fixed	\$202.00/ea.	BL-LTS-GT2FXS-LH	BL-LTS-GT2FXS-LH	BL-LTS-GT2FXS-RH	BL-LTS-GT2FXS-RH

**Track mounting	MSRP	Left	*Extra left	Right	*Extra right
Swing-away HPCPS Code: E1028	\$339.00/ea.	BL-LTS-GT2SAT-LH	BL-LTS-GT2SATX-LH	BL-LTS-GT2SAT-RH	BL-LTS-GT2SATX-RH
Fixed	\$202.00/ea.	BL-LTS-GT2FXT-LH	BL-LTS-GT2FXTX-LH	BL-LTS-GT2FXT-RH	BL-LTS-GT2FXTX-RH

4.5 GT2 Additional options



1 1/4"

Extra hardware links	Left	Additional left	Right	Additional right
MSRP \$43.00/link	BL-LTS-GT2LINK	BL-LTS-GT2LINK	BL-LTS-GT2LINK	BL-LTS-GT2LINK



adds 1" of adjustment per bracket

Long mounting bracket	Left	*Extra left	Right	*Extra right
For Slot Mount MSRP \$43.00/ea	BL-LTS-GT2LSLOT-LH	BL-LTS-GT2LSLOT-LH	BL-LTS-GT2LSLOT-RH	BL-LTS-GT2LSLOT-RH



Hardware cover	Left	*Extra left	Right	*Extra right
MSRP \$31.00/ea	BL-LTS-HCOVER	BL-LTS-HCOVER	BL-LTS-HCOVER	BL-LTS-HCOVER

5 BodiLink Head Support hardware & pads

Length (L) refers to the actual size dimension of the support from bottom to top edge.

- You may make one selection per column.
- Choose pad firmness and size in 5.1.
- Choose a cover material in 5.2.

5.1 Premium pad size

HPCPS CODING: E0955 MSRP \$269.00

Pad firmness			
W	L	Soft	BL-HSP1-6W3L
6W x 3L		Extra soft	BL-HSP2-6W3L
W	L	Soft	BL-HSP1-8W4L
8W x 4L		Extra soft	BL-HSP2-8W4L
W	L	Soft	BL-HSP1-10W5L
10W x 5L with adjustable wings		Extra soft	BL-HSP2-10W5L
W	L	Soft	BL-HSP1-12W5L
12W x 5L with adjustable wings		Extra soft	BL-HSP2-12W5L
W	L	Soft	BL-HSP1-14W5L
14W x 5L with adjustable wings		Extra soft	BL-HSP2-14W5L
W	L	Soft	BL-HSP1-18W5L
18W x 5L with adjustable wings		Extra soft	BL-HSP2-18W5L



Premium head support pads angles can be formed by adjusting the plates in the pad.

5.2 Premium pad cover material



For skin protection & shear reduction. Not fluid resistant.



For fluid protection & an easily cleaned surface.



For client comfort.

GlideWear™ (GLIDEWEAR) (add \$30.00)
COMFORT-TEK (COMFORT-TEK)
STRETCH-AIR (STRETCH-AIR)

5.3 PT Hardware

- You may make one selection per column.
- Choose PT hardware type in 5.3 A or B.
- If choosing A, select a stem size immediately after.

Range of motion

PT MODEL	TOTAL # OF LINKS	TOTAL LINK LENGTH	MAX ANTERIOR REACH	MAX POSTERIOR REACH	MAX LATERAL REACH
MID 6	2	6"	5.5"	1"	N/A
LAT 6	3	6"	5.5"	1"	6"
MID 8	2	8"	7.5"	1"	N/A
LAT 8	3	8"	7.5"	1"	8"
MAX 8	4	8"	7.5"	3"	8"

A. PT Hardware, Removable

HCPCS E1028



Shape	MSRP	Part number
MID 6	\$254.00	(BL-HS-PT1RMMID6)
LAT 6	\$295.00	(BL-HS-PT1RMLAT6)
MID 8	\$254.00	(BL-HS-PT1RMMID8)
LAT 8	\$295.00	(BL-HS-PT1RMLAT8)
MAX 8*	\$372.00	(BL-HS-PT1RMMAX8)

*Provides the greatest amount of articulation and range of motion.

PT Hardware, Removable, Stem length option




Stem length	MSRP	Part number
Standard 7" stem	\$0.00	(HS-STEM-7)
**3" stem	\$0.00	(HS-STEM-3)

**3" stem is recommended for use with power mounting where there may be special constraints with power chair functions.

B. PT Hardware, Flip-down

HCPCS E1028

Flip-down hardware has 7" stem



Shape	MSRP	Part number
MID 6	\$310.00	(BL-HS-PT1FDMID6)
LAT 6	\$353.00	(BL-HS-PT1FDLAT6)
MID 8	\$310.00	(BL-HS-PT1FDMID8)
LAT 8	\$353.00	(BL-HS-PT1FDLAT8)
MAX 8*	\$405.00	(BL-HS-PT1FDMAX8)

*Provides the greatest amount of articulation and range of motion.

5.4 Head support mounting option



Mounting block or bracket	MSRP
Posterior A mounting block	\$0.00

Part number
(HS-MB-POSTA)

5.5 Extra hardware



Link length	MSRP	Extra part	
*Extra 2" Link	\$65.00	(BL-HS-PT1LINK2)	



*Extra 3" Link	\$65.00	(BL-HS-PT1LINK3)	
----------------	---------	------------------	--



*Extra 4" Link	\$65.00	(BL-HS-PT1LINK4)	
----------------	---------	------------------	--



Tool-less wing lock	MSRP	Extra part	Extra part	Extra part	Extra part
*1 Lock	\$25.00/each	(BL-HS-52WINGA1)	(BL-HS-52WINGA1)	(BL-HS-52WINGA1)	(BL-HS-52WINGA1)



Additional	MSRP	Extra part	
1/4" Spacer plate	\$43.00/each	(BL-HS-SPCRA)	

*For PT hardware only

The following are trademarks and registered trademarks of Permobil: ACTA-BACK®, BodiLink®, COMFORT-TEK®, Compass®, and STRETCH-AIR®. GlideWear™ is a trademark of MIPS.

6 Additional information

6.1 Comments

INSTRUCTIONS: Please reference any section or subsection of the form when commenting.

Section #	Page #	Comments

INSTRUCTIONS: If preferred, please sketch ideas in the spaces provided below.

Submit by email