



ROHO® ENHANCER® Cushion

The multiple cell heights of the ENHANCER Cushion, combined with the two independent compartment (one consisting of the ischial well and leg channeling cells, the other consisting of the perimeter and pommel cells), provide positioning channels for your legs, giving you side-to-side stability and optimum skin/soft tissue protection.

The ENHANCER Cushion is the only ROHO cushion that allows you to isolate the ischial well and increase the lateral support in the greater trochanter area.

- ✓ Two independent inflation compartments
- ✓ Positioning and stability
- ✓ Skin/soft tissue protection

Dimensions

Width: 14.75 in. to 20 in. /
37.5 cm to 50.5 cm

Depth: 14.75 in. to 18.25 in. /
37.5 cm to 46.5 cm

Height: 2.25 in. to 4.25 in. /
5.5 cm to 10.5 cm (Range
due to varying cell heights)

Weight: Approx. 3.0 lbs. /
1.4 kg

Weight limit

No weight limit.
However, the cushion
must be properly sized to
the individual.

Limited warranty

Cushion: 24 months
Cover: 6 months



Our evidence shows the ENHANCER is ideal for people who:

- Require increased sitting stability
- Require control or positioning of the pelvis
- Have issues with leg or thigh control, either abduction or adduction
- Need to offload the ITs
- Are at risk for skin/soft tissue breakdown or deep tissue injury