



M- & F-Series Corpus FAQs

150kg max user weight

What models are available with a maximum user weight of 150kg?

M3 Corpus, M5 Corpus, F3 Corpus and F5 Corpus

Why are we making this change to the maximum user weight now?

The M- & F-Series Corpus models have been available in the European market with a maximum user weight of 150kg, and have passed all applicable ISO testing at the 150kg weight capacity. By increasing the user weight in the Americas and APAC regions, we are aligning the product offering globally and in turn, are able to offer our industry-leading products to more consumers.

Are there any differences in the product offering for individuals over 150kg?

There are no differences in the F-Series Corpus models for individuals over 136kg. However, there are additional restrictions for the M-Series Corpus models as follows:

- ActiveReach will not be available on M3 or M5 for individuals greater than 136kg.
- There are different speed restriction conditions for the M3 and M5 for individuals greater than 136kg. Please speak with your Permobil sales manager for details.

Does this increased user weight apply to existing users?

No, as the wheelchair was sold prior to the official change, the product does not meet the necessary regulatory and/or technical requirements to promote the higher weight capacity. An upgrade package will not be available.