

permobil

# TiLite TWIST

Play. Giggle. Grow.

Designed from the ground up for kids, the TWIST encourages exploration, independence, and, yes, playfulness. There are no parts or growth kits to buy. Instead, the TWIST grows with them.



Other chairs say they can grow too, but parts are required. With the TWIST, growth is built in. Sized from 8" x 8" up to 15" x 16", the TWIST offers 2" of growth in seat width and 3" of growth in seat depth.



Custom configuration allows you to add options and components to your child's chair to meet their needs and to optimise their riding experience. This includes choices of wheels, handrims, tyres, wheel locks, colour, and more.

## Specifications

Frame features	Aluminium mono-tube rigid Growable Full adjustability	Camber	0°, 2°, 4°, 6°, 8°, or 12°
Total chair width*	19"	Centre of gravity adjustment	6"
Seat width	8" – 15" (2" of built-in Seat Width Adj.)	Transport options	Transit tie-down
Seat depth	8" – 16" (3" of built-in Seat Depth Adj.)	Transport weight*	5.4 kg
Front seat height	13" – 20" (up to 2.25" of built-in front seat height adj.)	Weight limit	74.8 kg
Rear seat height	12" – 19.5" (up to 3.5" of built-in rear seat height adj.)	Warranty	Five years on frame One year on TiLite components
Front frame angle	70°, 80°, or 90°	* Weights and widths based on 12" x 12" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weights. For reference, a pair of standard 22" rear wheels weighs 3.17 kg.	
Seat back type	Folding seat back Adj. height 8.5" – 21" Adj. angle 80° – 101°		

## Frame Colours

### High Gloss Paint




TiLite Red Neon Pink Super Wet Black

### Metallic Paint



Tangerine Metallic Acid Green Metallic Ocean Blue Metallic

### Pearlescent Paint



Midnight Blue Pearl White Pearl

### Anodised



Black Blue Red Gold

### Candy Paint




Ruby Red Azure Blue Electric Plum

### Matte Textured Paint



Flame Orange OD Green Hyper Blue

### Matte Paint



Titanium Grey Black Lava

### Component Colours



Purple Orange Green