

permobil | ROHO SEATING + POSITIONING



ROHO® AGILITY® Sizing Matrix

Table of Contents

AGILITY Aluminum Sizing Matrix


- 4 Minimum Contour
- 6 Mid Contour
- 8 Max Contour

AGILITY CARBON™ Sizing Matrix

- 10 Minimum Contour
- 10 Mid Contour



AGILITY Aluminum Sizing Matrix



- A** – Back Shell Height
- A²** – Back Shell and Foam Height
- B** – Back Cane Width Range*
- C** – Inner Width
- D** – Flat Center Section of Shell
- E** – Length of Wings


- There will be a 1.5cm loss of seat depth if an AGILITY Aluminum is fitted to the smallest wheelchair width within each size range.
- The measurements are for reference only. Permobil strongly recommends clinicians complete a full client assessment and equipment trial to ensure optimal fit and use of an AGILITY back.
- ©2017 ROHO, Inc. Please visit Permobilus.com to receive additional information regarding product description and intended use.
- The following are trademarks and registered trademarks of Permobil | ROHO Seating and Positioning.: AGILITY®, AGILITY CARBON™, and ROHO®.

Minimum Contour

HCPCS CODES: E2611, E2613, E2615, E2620

NUMBER	A		A ²		B		C		D		E	
AG2-1410	10"	25.5 cm	13"	33 cm	13" – 14"	33 – 35.5 cm	10.5"	26.5 cm	6"	15 cm	4"	10 cm
AG2-1413	13"	33 cm	16"	40.5 cm	13" – 14"	33 – 35.5 cm	10.5"	26.5 cm	6"	15 cm	4"	10 cm
AG2-1416	16"	40.5 cm	19"	48.5 cm	13" – 14"	33 – 35.5 cm	10.5"	26.5 cm	6"	15 cm	4"	10 cm
AG2-1418	18"	45.5 cm	21"	53.5 cm	13" – 14"	33 – 35.5 cm	10.5"	26.5 cm	6"	15 cm	4"	10 cm
AG2-1510	10"	25.5 cm	13"	33 cm	14" – 15"	35.5 – 38 cm	11.5"	29 cm	6.5"	16.5 cm	4"	10 cm
AG2-1513	13"	33 cm	16"	40.5 cm	14" – 15"	35.5 – 38 cm	11.5"	29 cm	6.5"	16.5 cm	4"	10 cm
AG2-1516	16"	40.5 cm	19"	48.5 cm	14" – 15"	35.5 – 38 cm	11.5"	29 cm	6.5"	16.5 cm	4"	10 cm
AG2-1518	18"	45.5 cm	21"	53.5 cm	14" – 15"	35.5 – 38 cm	11.5"	29 cm	6.5"	16.5 cm	4"	10 cm
AG2-1610	10"	25.5 cm	13"	33 cm	15" – 16"	38 – 40.5 cm	12"	30.5 cm	7"	18 cm	4"	10 cm
AG2-1613	13"	33 cm	16"	40.5 cm	15" – 16"	38 – 40.5 cm	12"	30.5 cm	7"	18 cm	4"	10 cm
AG2-1616	16"	40.5 cm	19"	48.5 cm	15" – 16"	38 – 40.5 cm	12"	30.5 cm	7"	18 cm	4"	10 cm
AG2-1618	18"	45.5 cm	21"	53.5 cm	15" – 16"	38 – 40.5 cm	12"	30.5 cm	7"	18 cm	4"	10 cm
AG2-1710	10"	25.5 cm	13"	33 cm	16" – 17"	40.5 – 43 cm	13"	33 cm	8"	20.5 cm	4"	10 cm
AG2-1713	13"	33 cm	16"	40.5 cm	16" – 17"	40.5 – 43 cm	13"	33 cm	8"	20.5 cm	4"	10 cm
AG2-1716	16"	40.5 cm	19"	48.5 cm	16" – 17"	40.5 – 43 cm	13"	33 cm	8"	20.5 cm	4"	10 cm
AG2-1718	18"	45.5 cm	21"	53.5 cm	16" – 17"	40.5 – 43 cm	13"	33 cm	8"	20.5 cm	4"	10 cm
AG2-1810	10"	25.5 cm	13"	33 cm	17" – 18"	43 – 45.5 cm	13.5"	34.5 cm	9.5"	24 cm	4"	10 cm
AG2-1813	13"	33 cm	16"	40.5 cm	17" – 18"	43 – 45.5 cm	13.5"	34.5 cm	9.5"	24 cm	4"	10 cm
AG2-1816	16"	40.5 cm	19"	48.5 cm	17" – 18"	43 – 45.5 cm	13.5"	34.5 cm	9.5"	24 cm	4"	10 cm
AG2-1818	18"	45.5 cm	21"	53.5 cm	17" – 18"	43 – 45.5 cm	13.5"	34.5 cm	9.5"	24 cm	4"	10 cm
AG2-1910	10"	25.5 cm	13"	33 cm	18" – 19"	45.5 – 48.5 cm	15"	38 cm	10"	25.5 cm	4"	10 cm
AG2-1913	13"	33 cm	16"	40.5 cm	18" – 19"	45.5 – 48.5 cm	15"	38 cm	10"	25.5 cm	4"	10 cm
AG2-1916	16"	40.5 cm	19"	48.5 cm	18" – 19"	45.5 – 48.5 cm	15"	38 cm	10"	25.5 cm	4"	10 cm
AG2-1918	18"	45.5 cm	21"	53.5 cm	18" – 19"	45.5 – 48.5 cm	15"	38 cm	10"	25.5 cm	4"	10 cm
AG2-2010	10"	25.5 cm	13"	33 cm	19" – 20"	48.5 – 51 cm	16"	40.5 cm	11.5"	29 cm	4"	10 cm
AG2-2013	13"	33 cm	16"	40.5 cm	19" – 20"	48.5 – 51 cm	16"	40.5 cm	11.5"	29 cm	4"	10 cm
AG2-2016	16"	40.5 cm	19"	48.5 cm	19" – 20"	48.5 – 51 cm	16"	40.5 cm	11.5"	29 cm	4"	10 cm
AG2-2018	18"	45.5 cm	21"	53.5 cm	19" – 20"	48.5 – 51 cm	16"	40.5 cm	11.5"	29 cm	4"	10 cm

AGILITY Aluminum Sizing Matrix




- A** – Back Shell Height
- A²** – Back Shell and Foam Height
- B** – Back Cane Width Range*
- C** – Inner Width
- D** – Flat Center Section of Shell
- E** – Length of Wings

Mid Contour

HCPCS CODES: E2615, E2620

NUMBER	A		A ²		B		C		D		E	
AG3-1410	10"	25.5 cm	13"	33 cm	13" - 14"	33 – 35.5 cm	10.5"	26.5 cm	6"	15 cm	6"	15 cm
AG3-1413	13"	33 cm	16"	40.5 cm	13" - 14"	33 – 35.5 cm	10.5"	26.5 cm	6"	15 cm	6"	15 cm
AG3-1416	16"	40.5 cm	19"	48.5 cm	13" - 14"	33 – 35.5 cm	10.5"	26.5 cm	6"	15 cm	6"	15 cm
AG3-1418	18"	45.5 cm	21"	53.5 cm	13" - 14"	33 – 35.5 cm	10.5"	26.5 cm	6"	15 cm	6"	15 cm
AG3-1420	20"	51 cm	23"	58.5 cm	13" - 14"	33 – 35.5 cm	10.5"	26.5 cm	6"	15 cm	6"	15 cm
AG3-1510	10"	25.5 cm	13"	33 cm	14" - 15"	35.5 – 38 cm	12"	30.5 cm	7"	18 cm	5"	13 cm
AG3-1513	13"	33 cm	16"	40.5 cm	14" - 15"	35.5 – 38 cm	12"	30.5 cm	7"	18 cm	5"	13 cm
AG3-1516	16"	40.5 cm	19"	48.5 cm	14" - 15"	35.5 – 38 cm	12"	30.5 cm	7"	18 cm	5"	13 cm
AG3-1518	18"	45.5 cm	21"	53.5 cm	14" - 15"	35.5 – 38 cm	12"	30.5 cm	7"	18 cm	5"	13 cm
AG3-1520	20"	51 cm	23"	58.5 cm	14" - 15"	35.5 – 38 cm	12"	30.5 cm	7"	18 cm	5"	13 cm
AG3-1610	10"	25.5 cm	13"	33 cm	15" - 16"	38 – 40.5 cm	13"	33 cm	8"	20.5 cm	6"	15 cm
AG3-1613	13"	33 cm	16"	40.5 cm	15" - 16"	38 – 40.5 cm	13"	33 cm	8"	20.5 cm	6"	15 cm
AG3-1616	16"	40.5 cm	19"	48.5 cm	15" - 16"	38 – 40.5 cm	13"	33 cm	8"	20.5 cm	6"	15 cm
AG3-1618	18"	45.5 cm	21"	53.5 cm	15" - 16"	38 – 40.5 cm	13"	33 cm	8"	20.5 cm	6"	15 cm
AG3-1620	20"	51 cm	23"	58.5 cm	15" - 16"	38 – 40.5 cm	13"	33 cm	8"	20.5 cm	6"	15 cm
AG3-1710	10"	25.5 cm	13"	33 cm	16" - 17"	40.5 – 43 cm	14"	35.5 cm	9"	23 cm	5"	13 cm
AG3-1713	13"	33 cm	16"	40.5 cm	16" - 17"	40.5 – 43 cm	14"	35.5 cm	9"	23 cm	5"	13 cm
AG3-1716	16"	40.5 cm	19"	48.5 cm	16" - 17"	40.5 – 43 cm	14"	35.5 cm	9"	23 cm	5"	13 cm
AG3-1718	18"	45.5 cm	21"	53.5 cm	16" - 17"	40.5 – 43 cm	14"	35.5 cm	9"	23 cm	5"	13 cm
AG3-1720	20"	51 cm	23"	58.5 cm	16" - 17"	40.5 – 43 cm	14"	35.5 cm	9"	23 cm	5"	13 cm
AG3-1810	10"	25.5 cm	13"	33 cm	17" - 18"	43 – 45.5 cm	15"	38 cm	10"	25.5 cm	6"	15 cm
AG3-1813	13"	33 cm	16"	40.5 cm	17" - 18"	43 – 45.5 cm	15"	38 cm	10"	25.5 cm	6"	15 cm
AG3-1816	16"	40.5 cm	19"	48.5 cm	17" - 18"	43 – 45.5 cm	15"	38 cm	10"	25.5 cm	6"	15 cm
AG3-1818	18"	45.5 cm	21"	53.5 cm	17" - 18"	43 – 45.5 cm	15"	38 cm	10"	25.5 cm	6"	15 cm
AG3-1820	20"	51 cm	23"	58.5 cm	17" - 18"	43 – 45.5 cm	15"	38 cm	10"	25.5 cm	6"	15 cm
AG3-1910	10"	25.5 cm	13"	33 cm	18" - 19"	45.5 – 48.5	16"	40.5 cm	11"	28 cm	5"	13 cm
AG3-1913	13"	33 cm	16"	40.5 cm	18" - 19"	45.5 – 48.5	16"	40.5 cm	11"	28 cm	5"	13 cm
AG3-1916	16"	40.5 cm	19"	48.5 cm	18" - 19"	45.5 – 48.5	16"	40.5 cm	11"	28 cm	5"	13 cm
AG3-1918	18"	45.5 cm	21"	53.5 cm	18" - 19"	45.5 – 48.5	16"	40.5 cm	11"	28 cm	5"	13 cm
AG3-1920	20"	51 cm	23"	58.5 cm	18" - 19"	45.5 – 48.5	16"	40.5 cm	11"	28 cm	5"	13 cm
AG3-2010	10"	25.5 cm	13"	33 cm	19" - 20"	48.5 – 51 cm	17"	43 cm	12"	30.5 cm	6"	15 cm
AG3-2013	13"	33 cm	16"	40.5 cm	19" - 20"	48.5 – 51 cm	17"	43 cm	12"	30.5 cm	6"	15 cm
AG3-2016	16"	40.5 cm	19"	48.5 cm	19" - 20"	48.5 – 51 cm	17"	43 cm	12"	30.5 cm	6"	15 cm
AG3-2018	18"	45.5 cm	21"	53.5 cm	19" - 20"	48.5 – 51 cm	17"	43 cm	12"	30.5 cm	6"	15 cm
AG3-2020	20"	51 cm	23"	58.5 cm	19" - 20"	48.5 – 51 cm	17"	43 cm	12"	30.5 cm	6"	15 cm

AGILITY Aluminum Sizing Matrix




- A** – Back Shell Height
- A²** – Back Shell and Foam Height
- B** – Back Cane Width Range*
- C** – Inner Width
- D** – Flat Center Section of Shell
- E** – Length of Wings

Max Contour


HCPCS CODES: E2615, E2620

NUMBER	A		A ²		B		C		D		E	
AG6-1416	16"	40.5 cm	19"	48.5 cm	13" - 14"	33 – 35.5 cm	11"	28 cm	6"	15 cm	8"	20.5 cm
AG6-1418	18"	45.5 cm	21"	53.5 cm	13" - 14"	33 – 35.5 cm	11"	28 cm	6"	15 cm	8"	20.5 cm
AG6-1420	20"	51 cm	23"	58.5 cm	13" - 14"	33 – 35.5 cm	11"	28 cm	6"	15 cm	8"	20.5 cm
AG6-1516	16"	40.5 cm	19"	48.5 cm	14" - 15"	35.5 – 38 cm	12.5"	32 cm	8"	20.5 cm	8"	20.5 cm
AG6-1518	18"	45.5 cm	21"	53.5 cm	14" - 15"	35.5 – 38 cm	12.5"	32 cm	8"	20.5 cm	8"	20.5 cm
AG6-1520	20"	51 cm	23"	58.5 cm	14" - 15"	35.5 – 38 cm	12.5"	32 cm	8"	20.5 cm	8"	20.5 cm
AG6-1616	16"	40.5 cm	19"	48.5 cm	15" - 16"	38 – 40.5 cm	13"	33 cm	9"	23 cm	8"	20.5 cm
AG6-1618	18"	45.5 cm	21"	53.5 cm	15" - 16"	38 – 40.5 cm	13"	33 cm	9"	23 cm	8"	20.5 cm
AG6-1620	20"	51 cm	23"	58.5 cm	15" - 16"	38 – 40.5 cm	13"	33 cm	9"	23 cm	8"	20.5 cm
AG6-1716	16"	40.5 cm	19"	48.5 cm	16" - 17"	40.5 – 43 cm	14.5"	37 cm	10"	25.5 cm	9"	23 cm
AG6-1718	18"	45.5 cm	21"	53.5 cm	16" - 17"	40.5 – 43 cm	14.5"	37 cm	10"	25.5 cm	9"	23 cm
AG6-1720	20"	51 cm	23"	58.5 cm	16" - 17"	40.5 – 43 cm	14.5"	37 cm	10"	25.5 cm	9"	23 cm
AG6-1816	16"	40.5 cm	19"	48.5 cm	17" - 18"	43 – 45.5 cm	15"	38 cm	11.5"	29 cm	9"	23 cm
AG6-1818	18"	45.5 cm	21"	53.5 cm	17" - 18"	43 – 45.5 cm	15"	38 cm	11.5"	29 cm	9"	23 cm
AG6-1820	20"	51 cm	23"	58.5 cm	17" - 18"	43 – 45.5 cm	15"	38 cm	11.5"	29 cm	9"	23 cm
AG6-1916	16"	40.5 cm	19"	48.5 cm	18" - 19"	45.5 – 48.5 cm	16.5"	42 cm	12"	30.5 cm	9"	23 cm
AG6-1918	18"	45.5 cm	21"	53.5 cm	18" - 19"	45.5 – 48.5 cm	16.5"	42 cm	12"	30.5 cm	9"	23 cm
AG6-1920	20"	51 cm	23"	58.5 cm	18" - 19"	45.5 – 48.5 cm	16.5"	42 cm	12"	30.5 cm	9"	23 cm
AG6-2016	16"	40.5 cm	19"	48.5 cm	19" - 20"	48.5 – 51 cm	17"	43 cm	12.5"	32 cm	9"	23 cm
AG6-2018	18"	45.5 cm	21"	53.5 cm	19" - 20"	48.5 – 51 cm	17"	43 cm	12.5"	32 cm	9"	23 cm
AG6-2020	20"	51 cm	23"	58.5 cm	19" - 20"	48.5 – 51 cm	17"	43 cm	12.5"	32 cm	9"	23 cm

AGILITY CARBON Sizing Matrix



- A** – Back Shell Height
- A²** – Foam Height
- B** – Outside Cane to Outside Cane
- C** – Inner Foam Width
- D** – Flat Center Section
- E¹** – Length of Wings, Shell
- E²** – Length of Wings, Foam
- F** – Shell Depth
- G** – Inner Foam Depth



Minimum Contour

NUMBER	A		A ²		B		C		D		E ¹		E ²		F		G	
AG2-1410	10"	25.5 cm	12"	30.5 cm	14"	35.5 cm	11"	28 cm	6"	15 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1413	13"	33 cm	15"	38 cm	14"	35.5 cm	11"	28 cm	6"	15 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1416	16"	40.5 cm	18"	45.5 cm	14"	35.5 cm	11"	28 cm	6"	15 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1510	10"	25.5 cm	12"	30.5 cm	15"	38 cm	12"	30.5 cm	7"	18 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1513	13"	33 cm	15"	38 cm	15"	38 cm	12"	30.5 cm	7"	18 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1516	16"	40.5 cm	18"	45.5 cm	15"	38 cm	12"	30.5 cm	7"	18 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1610	10"	25.5 cm	12"	30.5 cm	16"	40.5 cm	13"	33 cm	8"	20.5 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1613	13"	33 cm	15"	38 cm	16"	40.5 cm	13"	33 cm	8"	20.5 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1616	16"	40.5 cm	18"	45.5 cm	16"	40.5 cm	13"	33 cm	8"	20.5 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1710	10"	25.5 cm	12"	30.5 cm	17"	43 cm	14"	35.5 cm	9"	23 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1713	13"	33 cm	15"	38 cm	17"	43 cm	14"	35.5 cm	9"	23 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1716	16"	40.5 cm	18"	45.5 cm	17"	43 cm	14"	35.5 cm	9"	23 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1810	10"	25.5 cm	12"	30.5 cm	18"	45.5 cm	15"	38 cm	10"	25.5 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1813	13"	33 cm	15"	38 cm	18"	45.5 cm	15"	38 cm	10"	25.5 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm
AG2-1816	16"	40.5 cm	18"	45.5 cm	18"	45.5 cm	15"	38 cm	10"	25.5 cm	4"	10 cm	5.5"	14 cm	3"	7.5 cm	3"	7.5 cm

Mid Contour

NUMBER	A		A ²		B		C		D		E ¹		E ²		F		G	
AG3-1410	10"	25.5 cm	12"	30.5 cm	14"	35.5 cm	12"	30.5 cm	6"	15 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1413	13"	33 cm	15"	38 cm	14"	35.5 cm	12"	30.5 cm	6"	15 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1416	16"	40.5 cm	18"	45.5 cm	14"	35.5 cm	12"	30.5 cm	6"	15 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1510	10"	25.5 cm	12"	30.5 cm	15"	38 cm	13"	33 cm	7"	18 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1513	13"	33 cm	15"	38 cm	15"	38 cm	13"	33 cm	7"	18 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1516	16"	40.5 cm	18"	45.5 cm	15"	38 cm	13"	33 cm	7"	18 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1610	10"	25.5 cm	12"	30.5 cm	16"	40.5 cm	14"	35.5 cm	8"	20.5 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1613	13"	33 cm	15"	38 cm	16"	40.5 cm	14"	35.5 cm	8"	20.5 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1616	16"	40.5 cm	18"	45.5 cm	16"	40.5 cm	14"	35.5 cm	8"	20.5 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1710	10"	25.5 cm	12"	30.5 cm	17"	43 cm	15"	38 cm	9"	23 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1713	13"	33 cm	15"	38 cm	17"	43 cm	15"	38 cm	9"	23 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1716	16"	40.5 cm	18"	45.5 cm	17"	43 cm	15"	38 cm	9"	23 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1810	10"	25.5 cm	12"	30.5 cm	18"	45.5 cm	16"	40.5 cm	10"	25.5 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1813	13"	33 cm	15"	38 cm	18"	45.5 cm	16"	40.5 cm	10"	25.5 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm
AG3-1816	16"	40.5 cm	18"	45.5 cm	18"	45.5 cm	16"	40.5 cm	10"	25.5 cm	5.25"	13 cm	6.75"	17 cm	4"	10 cm	4"	10 cm



Shipping Address: 1501 S. 74th Street, Belleville, IL 62223 USA

Corporate Address: 100 N. Florida Avenue, Belleville, IL 62221 USA

(800)851-3449 | (618)277-9173 | Permobilus.com

LITG02324 Rev. 4/17 Print 4/17