



Order form

F3 Corpus



Dealer Name*:

Account #*:

City, State, Zip*:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference*:

Client Weight*:

***Required field**

Special considerations that need to be addressed here (e.g., diagnosis):

Chair Model and Wheels

Front-wheel drive F3 power base includes: Corpus seat frame, ComfortRide suspension, Enhanced Steering Performance (ESP), 4-point tie-down hardware, VoltPro charger, and R-net 120 amp electronics.

Choose Base

11332	F3 Base Corpus - F3 SPO <i>FWD F3 power base set up for one power seat function, max speed: 6 mph, weight capacity: 300 lb.</i>	\$8,768.00	K0856
108398-99-0	F3 Base Corpus - F3 MPO <i>FWD F3 power base set up for multiple power seat functions, max speed: 6 mph, weight capacity: 300 lb.</i>	\$8,768.00	K0861
11611	F3 Base Corpus - F3 SPO HD <i>FWD F3 power base set up for one power seat function, max speed: 6 mph, weight capacity: 330 lb.</i>	\$10,343.00	K0858
11612	F3 Base Corpus - F3 MPO HD <i>FWD F3 power base set up for multiple power seat functions, max speed: 6 mph, weight capacity: 330 lb.</i>	\$10,343.00	K0862
108852-99-0	Permobil Connect	No charge	

Support Wheel Options | Choose one

107330-99-0	Support Wheels for F3	No charge
11313	Omit Front Support Wheels <i>User Weight must not exceed 220lbs. Max speed will be limited to 5mph on the F3 Corpus and 6mph on the F5 Corpus. Additional restrictions may apply based on chair configuration.</i>	No charge

Chair Model and Wheels

Choose Drive Tires

108414-99-0	Flat Free Tire 3"W x 14"T Black - Standard Tread <i>Includes flat free standard tread caster wheels.</i>	No charge
108411-99-0	Pneumatic Tire 3"W x 14"T Black - Standard Tread <i>Includes flat free standard tread caster wheels.</i>	No charge
11333	Flat Free Tire 3"W x 14"T Black - Aggressive Tread <i>Includes flat free aggressive tread caster wheels and mud guards.</i>	\$514.00
108412-99-0	Pneumatic/Snow Tire 3"W x 14"T Black - Aggressive Tread <i>Not recommended for year-long use and must be ordered along with standard drive wheels. Does not include caster wheels.</i>	\$514.00
11334	Mud Guard Kit F3	\$45.00

Base Options

Choose Base Color | Choose both shroud and rim color

Standard Shroud Colors	Aqua	Sky	Violet	Flamingo	No charge
	Chili	Canary	Spring	Arctic	
	Sunset	Midnight	Gold	Bronze	
	Platinum				
Premium Shroud Colors	Carbon Fiber	Mossy-Oak	CDG Poly	Water Lily	\$223.00
	Midnight Patriot	American Liberty			
Standard Rim Colors	Aqua	Sky	Violet	Flamingo	No charge
	Chili	Canary	Spring	Arctic	
	Sunset	Midnight	Gold	Bronze	
	Platinum				
Premium Rim Colors	Carbon Fiber	Mossy-Oak	CDG Poly	Water Lily	\$55.00
	Midnight Patriot	American Liberty			

Batteries | Choose one

108999-99-0	Permobil Batteries, Grp 24 AGM GEL - Pair	\$1,095.00	E2363
109183-99-0	Permobil Batteries, Grp 34 AGM GEL - Pair	\$1,048.00	E2359

Optional Items

109170-99-0	Essentials Box with Reflectors	\$535.00	
109169-99-0	Essentials Box with Lights <i>Must also select LED Light Kit.</i>	\$872.00	
109379-99-0	LED Light Kit for F3 <i>No charge with 4 power seat function configurations.</i>	\$1,087.00	
10715	Compact Ventilator Tray - Corpus Seat <i>Custom option. 13"T, adjustable width between 8" - 13", and adjustable depth between 2" - 6.5". Additional drive and seating restrictions may apply.</i>	\$1,755.00	E1029
10787	Battery Tray for Laptop Ventilator Tray <i>Custom option. 11 3/8" W x 6 1/2" H x 1 5/8" D.</i>	\$294.00	
10788	O2 Holder for Laptop Ventilator Tray <i>Custom option. Mounts to left or right side of compact laptop ventilator tray. 7"T, adjustable width between 3.5" - 7", and adjustable depth between 3.75" - 5".</i>	\$373.00	E2208
109620-99-0	Programming Key for R-net - Permobil Version	\$575.00	

Joystick Options

Choose joystick type, handle, and mount. NOTE: For additional input devices, please refer to the Alternative Drive Systems section.

Joystick Type

11186	Permobil Joystick Module (PJSM), NC <i>2.8" color display.</i>	No charge
105246-99-0	R-net LED Remote Joystick - No Display <i>No display, no monojacks, limited to 5 profiles.</i>	No charge
107779-99-0	Color Joystick Module (CJSM) <i>High-resolution 3.5" VGA color display.</i>	\$1,294.00
105700-99-0	R-net Remote Color Joystick w/ Lights, Heavy Duty Gimbal <i>Color LCD joystick, includes 3" knob with 3/4" diameter. Not compatible with alternative joystick handles.</i>	\$1,346.00
107829-99-0	3.5mm Y-Splitter Cable for PJSM/CJSM <i>Allows 2 mechanical switches to control up to 4 separate functions through assignable buttons programming.</i>	\$39.00

Joystick Handle

10045	Cone Shaped Handle - Standard	No charge
100960-99-0	Chin Cup Handle	\$139.00 E2324
100961-99-0	T Handle	\$196.00 E2323
100962-99-0	Large Ball	\$139.00 E2323
100963-99-0	Softball	\$183.00 E2323
100964-99-0	Mushroom Handle	\$183.00 E2323
100965-99-0	Stick Handle	\$183.00 E2323
105442-99-0	1.75" Black Foam Ball	\$139.00 E2323
105846-99-0	Bodypoint 3" U-Shaped Handle w/ Flex-Shaft	\$183.00 E2323
105847-99-0	Bodypoint 4" U-Shaped Handle w/ Flex-Shaft	\$183.00 E2323
105848-99-0	Bodypoint Rubber Dome Handle	\$183.00 E2323

Joystick Mount

Fixed	Mount on Left	Mount on Right	No charge
Swing-Away	Mount on Left	Mount on Right	\$409.00 E1028
Slimline Retractable	Mount on Left	Mount on Right	\$550.00 E1028
11232	Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment and 25° of angle adjustment.</i>		\$299.00 E1028
105210-99-0	Mount Standard Joystick on Corpus Back Hinge <i>Can be chosen instead of attendant control.</i>		\$376.00
	Mount on Left	Mount on Right	
108268-99-0	Bodypoint Tri-Lock Midline Mount Assembly <i>Must select specific joystick attachment separately.</i>		\$702.00
	Mount on Left	Mount on Right	
108280-99-0	Bodypoint Tri-Lock Mount for R-net Joystick		\$292.00
108279-99-0	Bodypoint Tri-Lock Mount for Compact Joystick		\$292.00
10842	PERMOFix Mount for mo-vis/Compact Joystick <i>Permofix mount for attaching Compact and mo-vis joysticks to Tri-Lock midline mount.</i>		\$292.00

Joystick Options

Other Options

10411	Expandable Controller - R-net <i>Required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an alternative drive control.</i>	\$764.00	E2377
10513	Harness for Expandable Controller <i>Must be ordered with any expandable controller.</i>	\$678.00	E2313
109532-99-0	mo-vis Scoot Attendant Control <i>Fully proportional handlebar with configurable thumb throttles, access to wheelchair functions and status and a foldable, adjustable mounting bracket.</i>	\$4,200.00	E2331
105408-99-0	Attendant Control R-net	\$1,499.00	E2331
10582	Output Module (IOM) for R-net Control System <i>Enables user to control an external device (e.g., communication device, door opener, etc.) through the drive input (e.g., JSM, head array, etc.). May require specialty cable or wiring to be ordered separately.</i>	\$1,412.00	

Power Seating

Please Note: Power seat adjustability may be limited to prevent contact with base and/or backrest-mounted items.

Adjustable Seat Height | *If ordering power adjustable seat height, elevator travel is used to set seat heights greater than 17.5".*

11110	ActiveHeight Power Adjustable Seat Height - 12" Travel <i>Must also select power tilt when selecting this option.</i>	\$3,786.00	E2300
11117	Seat Tube Fixed	No charge	

Seat-to-Floor Height 17.5" 18.5" 19.5"

Posterior Tilt/Recline Packages | *Choose one*

11131	Corpus Seat w/ No Tilt & 180° Power Recline	\$9,555.00	E1004
11132	Corpus Seat w/ 50° Tilt	\$8,525.00	E1002
11133	Corpus Seat w/ 50° Tilt & 180° Power Recline	\$13,492.00	E1007

ActiveReach Packages | *Chest support belt is strongly recommended. Power legrest recommended for full functionality. Seating and drive capabilities may be limited in certain configurations.*

11141	5° ActiveReach	No charge	
107090-99-0	10° ActiveReach <i>ActiveHeight and power tilt required.</i>	No charge	
107091-99-0	20° ActiveReach <i>ActiveHeight, power tilt/recline and Corpus VS elevating legrests required.</i>	\$2,967.00	
11100	30° ActiveReach w/Transfer <i>ActiveHeight, power tilt/recline and Corpus VS elevating legrests required.</i>	\$3,335.00	

Seat Function Control Options | *ICS switchbox availability may be limited to certain configurations.*

10418	Seat Function Control Switch for R-net <i>Does not include operation of power seat functions through drive control.</i>	No charge	
10419	Single Seat Function Control Kit for R-net <i>Allows operation of a single power seat function through drive control.</i>	\$2,110.00	E2310
10585	Multiple Seat Function Control Kit for R-net <i>Allows operation of a two or more power seat functions through drive control.</i>	\$3,994.00	E2311
102798-99-0	Push Buttons on ICS Switchbox	No charge	
102799-99-0	Toggle Switches on ICS Switchbox	\$357.00	

Switchbox Mounting

10626	Mount ICS Switchbox onto R-net Joystick Bracket	No charge	
-------	--	-----------	--

Power Seating

104651-99-0	Mount ICS Switchbox onto Armrest Mount on Left Mount on Right	No charge
106256-99-0	Mount ICS Switchbox on Back <i>No charge when ordering Omni2 Control Module.</i> Mount on Left Mount on Right	\$172.00
105409-99-0	Do Not Install ICS Switchbox, Send as Part	No charge

Power Seating Options

104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Accommodates up to 8 momentary switches.</i>	\$1,396.00
108416-99-0	Remote Stop System <i>Consists of a 315MHz radio frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair.</i>	\$1,381.00
10750	Memory Seat Program, Corpus w/ R-net	\$144.00
106209-99-0	Legrest Movement Linked to Recline Electronically	\$55.00
10898	Enable Independent Repositioning Mode (IRM) <i>IRM allows for sequential movement of power seat functions into a saved posterior position and returns to a saved upright position.</i>	No charge

Backrest Options

Installed backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Accessible seat depth will be dependent upon backrest selection.

Corpus Backrest Kits

105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No charge
105066-99-0	Fixed Backrest Kit - No Mechanical Shear	\$44.00

Corpus Backrests | *Corpus Ergo Backs include lateral wedges and foam lumbar pads. Solid Backs do not include foam lumbar pads.*

	20" Tall <i>Power Recline will be limited to 150°</i>	23" - 28" Tall <i>Height adjustable in 1" increments</i>		
Corpus Ergo (C3GERGBCA)				
14" W	14 x 20	14 x 23-28	\$1,192.00	E2620
16" W	16 x 20	16 x 23-28	\$1,192.00	E2620
18" W	18 x 20	18 x 23-28	\$1,192.00	E2620
20" W	20 x 20	20 x 23-28	\$1,192.00	E2620
Solid Corpus Ergo				
16" W	16 x 20	16 x 23-28	No charge	
18" W	18 x 20	18 x 23-28	No charge	
20" W	20 x 20	20 x 23-28	No charge	
Set Back Height:	20"	23"	24"	25"
		26"	27"	28"

Backrest Options

Corpus Backrest Options

11335	Corpus Ergo Back - Leatherette	No charge
11336	Corpus Ergo Back - Stretch-Air	No charge
10725	Split Backrest Cushion - Corpus <i>Custom option.</i>	\$377.00
10786	Tailored Backrest Cushion - Corpus <i>Custom option. Sides of cushion are trimmed to provide additional clearance for lateral trunk supports.</i>	\$141.00
11339	Foam Lumbar Pads for Solid Back - Corpus	\$45.00
109029-99-0	BodiLink Premium 3" x 5" Lateral Supports - Pair <i>Must also select hardware.</i>	\$446.00 E0956
109030-99-0	BodiLink Premium 4" x 6" Lateral Supports - Pair <i>Must also select hardware.</i>	\$446.00 E0956
109031-99-0	Swing-Away Adj BodiLink Lateral Hardware Corpus - Pair	\$625.00 E1028

ROHO Backrests | Includes backrest shell, cushion, air inserts, UT rail for headrest and Direct Mount adapter.

	16" Tall <i>Installed height between 22" – 25"</i>	18" Tall <i>Installed height between 22" – 26"</i>	20" Tall <i>Installed height between 24" – 27"</i>		
ROHO AGILITY <i>Max Contour</i>					
14" W	14 x 16	14 x 18	14 x 20	\$874.00	E2620
16" W	16 x 16	16 x 18	16 x 20	\$874.00	E2620
18" W	18 x 16	18 x 18	18 x 20	\$874.00	E2620
20" W	20 x 16	20 x 18	20 x 20	\$874.00	E2620

Set Back Height: 22" 23" 24" 25" 26" 27"

Comfort Backrests | Includes backrest shell, cushion and Compass Power Mount.

Note: User weight may be limited based upon the size of the backrest selected.

	14" Tall <i>Installed height between 17" – 23"</i>	16" Tall <i>Installed height between 19" – 25"</i>	18" Tall <i>Installed height between 21" – 27"</i>	20" Tall <i>Installed height between 23" – 29"</i>		
Acta-Back <i>Stretch-Air</i>						
14" W		14 x 16	14 x 18		\$668.00	E2613
16" W		16 x 16	16 x 18	16 x 20	\$751.00	E2613
18" W		18 x 16	18 x 18	18 x 20	\$751.00	E2613
20" W			20 x 18	20 x 20	\$775.00	E2613
Acta-Embrace <i>Comfort-Tek</i>						
12" W	12 x 14				\$694.00	E2613
14" W		14 x 16	14 x 18		\$694.00	E2613

Comfort Backrest Options

11714	BodiLink Premium 3" x 5" Lateral Supports - Pair <i>Must also select hardware.</i>	\$446.00 E0956
11715	BodiLink Premium 4" x 6" Lateral Supports - Pair <i>Must also select hardware.</i>	\$446.00 E0956
109903-99-0	BodiLink GT2 Swing-Away Lateral Hardware - Pair <i>Must select with BodiLink lateral trunk supports on Comfort backrests.</i>	\$625.00 E1028

Aftermarket Backrest Options

105782-99-0	Tubular Backrest Frame - Corpus (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Up to several inches of seat depth may be lost when mounting this system.</i>	\$557.00
106072-99-0	Direct Backrest Frame - Corpus <i>Used to mount backrests ranging in width from 14" to 20". Some seat depth may be lost when mounting this system.</i>	\$557.00

Seat Sizing

Please select both width and depth options below.

Seat Pan Width | Choose one

Set at: 17" W 19" W 21" W 23" W No charge

Seat Pan Depth | Choose one

Kit 1, Set at: 14" 15" 16" No charge

Kit 2, Set at: 17" 18" 19" No charge

Kit 3, Set at: 20" 21" 22" No charge

Seat Cushion Options

Corpus Ergo Seat Cushions | *Corpus Ergo Air cushion features ROHO Dry Floatation technology.*

	15" D	17" D	19" D	21" D	23" D		
Corpus Ergo <i>Leatherette</i>							
17" W	17 x 15	17 x 17	17 x 19	17 x 21		\$632.00	E2605
19" W	19 x 15	19 x 17	19 x 19	19 x 21	19 x 23	\$632.00	E2605
21" W		21 x 17	21 x 19	21 x 21	21 x 23	\$632.00	E2605
23" W			23 x 19	23 x 21	23 x 23	\$632.00	E2605
Corpus Ergo <i>Stretch-Air</i>							
17" W		17 x 17	17 x 19	17 x 21		\$632.00	E2605
19" W		19 x 17	19 x 19	19 x 21	19 x 23	\$632.00	E2605
21" W		21 x 17	21 x 19	21 x 21	21 x 23	\$632.00	E2605
23" W			23 x 19	23 x 21	23 x 23	\$632.00	E2605
Corpus Ergo Air <i>C3GRHOSCA</i>							
17" W		17 x 17	17 x 19	17 x 21	17 x 23	\$667.00	E2622
19" W		19 x 17	19 x 19	19 x 21		\$667.00	E2622
21" W			21 x 19	21 x 21		\$667.00	E2622

ROHO Neoprene Seat Cushions

	8 cells 14.75" D	9 cells 16.5" D	10 cells 18.25" D	11 cells 20" D	12 cells 21.75" D	13 cells 23.5" D		
Quadtro Select <i>High Profile</i>								
9 cells / 16.5" W	9 x 8	9 x 9	9 x 10	9 x 11	9 x 12		\$532.00	E2624
10 cells / 18.25" W		10 x 9	10 x 10	10 x 11	10 x 12	10 x 13	\$532.00	E2624
11 cells / 20" W		11 x 9	11 x 10	11 x 11	11 x 12	11 x 13	\$532.00	E2624
12 cells / 21.75" W			12 x 10	12 x 11	12 x 12	12 x 13	\$532.00	E2624
13 cells / 23.5" W		13 x 9	13 x 10	13 x 11	13 x 12	13 x 13	\$592.00	E2625
Single Compartment <i>High Profile, Sensor-Ready (* = non-SR)</i>								
9 cells / 16.5" W	9 x 8	9 x 9	9 x 10	9 x 11	9 x 12*		\$532.00	E2622
10 cells / 18.25" W		10 x 9	10 x 10	10 x 11	10 x 12	10 x 13*	\$532.00	E2622
11 cells / 20" W		11 x 9	11 x 10	11 x 11	11 x 12	11 x 13*	\$532.00	E2622
12 cells / 21.75" W			12 x 10	12 x 11	12 x 12	12 x 13*	\$592.00	E2622
13 cells / 23.5" W			13 x 10	13 x 11			\$592.00	E2623

Seat Cushion Options

ROHO Hybrid Seat Cushions

Hybrid Select	14" D	15" D	16" D	17" D	18" D	20" D		
14" W	14 x 14		14 x 16				\$585.00	E2624
15" W		15 x 15	15 x 16	15 x 17			\$585.00	E2624
16" W		16 x 15	16 x 16	16 x 17	16 x 18	16 x 20	\$585.00	E2624
17" W			17 x 16	17 x 17	17 x 18	17 x 20	\$585.00	E2624
18" W			18 x 16	18 x 17	18 x 18	18 x 20	\$585.00	E2624
20" W			20 x 16	20 x 17	20 x 18	20 x 20	\$585.00	E2624

Comfort Seat Cushions | *User weight may be limited based upon the size of the cushion selected.*

	14" D	16" D	18" D	20" D	22" D		
Embrace Anti-Thrust <i>Comfort-Tek</i>							
12" W	12 x 14					\$465.00	E2607
14" W	14 x 14	14 x 16				\$465.00	E2607
M2 Anti-Thrust <i>Comfort-Tek</i>							
14" W	14 x 14	14 x 16				\$474.00	E2607
16" W		16 x 16	16 x 18	16 x 20		\$474.00	E2607
18" W		18 x 16	18 x 18	18 x 20	18 x 22	\$474.00	E2607
20" W		20 x 16	20 x 18	20 x 20	20 x 22	\$474.00	E2607
22" W			22 x 18	22 x 20	22 x 22	\$508.00	E2608

Cushion Options

107505-99-0	ROHO Smart Check	\$185.00
11113	Incontinence Cover for Ergo Cushion <i>Sent as a part.</i>	\$203.00
11173	Incontinence Cover for ROHO Cushion <i>Sent as a part.</i>	\$82.00
105118-99-0	Omit Seat Cushion	No charge

Armrest Options

Corpus armrests are infinitely adjustable from 8" - 15", measured from the seat pan to the top of the armrest pad.

Armrest Set-Up | *Width must be wider than backrest to ensure proper fit.*

106543-99-0	Narrow Armbar Kit (13"W - 19"W) <i>For 13" width, armrests will be mounted to inside of bracket. Set Width Between Armrests to:</i> 13" 15" 17" 19"	No charge
106544-99-0	Wide Armbar Kit (19"W - 23"W) <i>Set Width Between Armrests to:</i> 19" 21" 23"	No charge
105069-99-0	Height Adjustable Armrests - NC <i>No charge when ordering power tilt and/or recline. Set Armrest Height to:</i> 8" 9" 10" 11" 12" 13" 14" 15"	No charge
105081-99-0	Height Adjustable Armrests <i>Set Armrest Height to:</i> 8" 9" 10" 11" 12" 13" 14" 15"	\$391.00
105071-99-0	Medial Swivel <i>Medial swivel cannot be enabled when selecting fixed option. Must select fixed option when ordering arm trough adapters.</i> Left Fixed Left Enabled Right Fixed Right Enabled	No charge

Armrest Options

Choose Armrest Pads | *Must select left and right. Priced per side.*

CL Pads	Corpus Leatherette						No charge
	13" Left	16" Left	18" Left	13" Right	16" Right	18" Right	
CF Pads	Corpus Black Fabric						No charge
	13" Left	16" Left	18" Left	13" Right	16" Right	18" Right	
CG Pads	Corpus Gel						\$202.00
	13" Left	16" Left	18" Left	13" Right	16" Right	18" Right	
MFS	Moldable Forearm Support						\$282.00
	Left	Right					
MCHS	Contoured Hand Support - Medium						\$140.00
	<i>Used with Moldable Forearm Support.</i>						
	Left Medium	Right Medium					
MLHS	Contoured Hand Support - Large						\$147.00
	<i>Used with Moldable Forearm Support.</i>						
	Left Large	Right Large					
ATA	Arm Trough Adapters						\$377.00
	1" Pair	7/8" Pair	Install	Send as Parts			

Optional Items | *Priced per side/support.*

Armrest Pouch	Left	Right	\$267.00
Armrest Supports	Left Outer	Left Inner	Right Outer
	Right Inner		\$385.00
Elbow Supports (Must also select hardware)	Left Elbow Support	Right Elbow Support	\$191.00
Elbow Support Hardware	Left	Right	\$280.00

Legrest Options

Legrest Type | Complete column for legrest type chosen

	Corpus <i>Includes mechanical articulation Footplate range is 6" - 17"</i>	Corpus VS <i>Only available with > or = 20° ActiveReach; Includes power articulation Footplate range is 9" - 17"</i>		
Type	Power	Power	\$4,102.00	E1012
	Manual		No charge	
Footplate Style	2pc Wide - 9"D x 7.5"W	2pc Wide - 9"D x 7.5"W	No charge	
	2pc Narrow - 9"D x 6"W	2pc Narrow - 9"D x 6"W	No charge	
	1pc Narrow - 9.75"D x 12.5"W	1pc Wide - 9"D x 16.5"W	\$294.00	
Footplate Height <i>Height measured from seat pan with legrest at 90°</i>	Set footplate height at:	Set footplate height at:		
Calf Supports				
B Pads	3.5"W x 4.5"T	3.5"W x 4.5"T	No charge	
J Pads	5"W x 6.5"T	5"W x 6.5"T	No charge	
A Pads	5.5"W x 6"T	5.5"W x 6"T	No charge	
H Pads (curved)	7"W x 7.5"T	7"W x 7.5"T	\$419.00	
D Pads	8"W x 8"T	8"W x 8"T	\$419.00	
Knee Support <i>Both items required for Corpus VS Legrest</i>				
		Single Post Knee Support/Abductor	\$292.00	E0957
		Adj/Removable Knee Support Hardware	\$306.00	E1028

Footplate Options

106530-99-0	Stealth Gel Pads - Narrow Footplates - Pair	\$470.00
106225-99-0	Stealth Gel Pads - Wide Footplates - Pair	\$481.00

Accessories

Belt Options | Must select one of the following

104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No charge
104520-99-0	Bodypoint Padded Hip Belt - Medium	\$167.00
104514-99-0	Bodypoint Non-Padded Hip Belt - Large	\$142.00
104516-99-0	Bodypoint Padded Hip Belt - Large	\$223.00
100519-99-0	Positioning Belt Retractable	\$332.00
100600-99-0	Positioning Belt Retractable Long	\$345.00
108267-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Medium	\$224.00
108266-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Large	\$236.00

Head Support Options

107470-99-0	Permobil 7"W Headrest <i>Must also select headrest mounting hardware.</i>	\$434.00	E0955
107469-99-0	Permobil 10"W Headrest <i>Must also select headrest mounting hardware.</i>	\$434.00	E0955
107468-99-0	Permobil 14"W Headrest <i>Must also select headrest mounting hardware.</i>	\$434.00	E0955
107480-99-0	Adj/Removable Headrest Link Hardware - Standard	\$299.00	E1028
109178-99-0	Extra Link for Link Headrest Hardware <i>Sent as a part.</i>	\$55.00	
108721-99-0	Adj/Removable BodiLink Head Support Hardware - MAX 8 <i>Required with more than 2 swing-away headrest accessories or chin boom.</i>	\$365.00	E1028
109624-99-0	Extra 4" Link for BodiLink MAX 8 Hardware <i>Sent as a part.</i>	\$80.00	

Accessories

11305	Swing-Away Egg Switch Assembly Mount on Upper Left Mount on Upper Right	\$640.00	E1028
11306	Swing-Away Microlight Assembly Mount on Upper Left Mount on Upper Right	\$730.00	E1028
11307	Swing-Away Facial Pad Assembly 4"L x 2"T Mount on Upper Left Mount on Upper Right	\$519.00	E1028
11308	Swing-Away Facial Pad Assembly 6"L x 2"T Mount on Upper Left Mount on Upper Right	\$519.00	E1028
11276	Provide Headrest Receiver for AGILITY Backrest <i>Select for AGILITY backs and third party backs ordered off chair.</i>	No charge	
106331-99-0	Permobil 10"W Headrest, Leatherette <i>Only available with Corpus backrests. Must also select hardware.</i>	\$397.00	E0955
10552	Adj/Removable Headrest Hardware	\$306.00	E1028
107499-99-0	Universal Headrest Adapter	\$302.00	

Leg Support Options | Lateral pelvic thigh supports (LPTS)

103070-99-0	UT Amp. Support Pad (H), Mesh Cover - 1 Each <i>7"W x 7"T curved pad. Must also select hardware.</i>	\$181.00	E1020
103071-99-0	UT Amp. Support Pad (D), Mesh Cover - 1 Each <i>8"W x 8"T flat pad. Must also select hardware.</i>	\$181.00	E1020
10553	Adj Removable - Amp Support Hardware - 1 Each Mount on Left Mount on Right	\$306.00	
109040-99-0	BodiLink Premium LPTS Pad (B), Stretch-Air - Pair <i>4.5"W x 3.5"T pad. Must also select hardware.</i>	\$359.00	E0953
109043-99-0	BodiLink Premium LPTS Pad (A), Stretch-Air - Pair <i>6"W x 5.5"T pad. Must also select hardware.</i>	\$359.00	E0953
109041-99-0	BodiLink Premium LPTS Pad (C), Stretch-Air - Pair <i>8"W x 3.5"T pad. Must also select hardware.</i>	\$359.00	E0953
109042-99-0	BodiLink Premium LPTS Pad (G), Stretch-Air - Pair <i>12"W x 3.5"T pad. Must also select hardware.</i>	\$359.00	E0953
109044-99-0	Adj/Removable BodiLink LPTS Hardware Med Corpus - Pair <i>Height adjustable.</i>	\$689.00	E1028
109194-99-0	Adj/Removable Permobil Lateral Thigh/Hip Hardware - Pair	\$602.00	E1028
106762-99-0	Abductor Kit Corpus - Leatherette <i>Not compatible with Corpus VS Legrest. Must also select hardware.</i>	\$423.00	E0957
10555	Adj Removable - Abductor Hardware - 1 Each	\$306.00	E1028

Chest Support Options

105840-99-0	Bodypoint Monoflex Chest Support Belt Kit - Medium <i>18.25"L x 4"W Pad Size.</i>	\$204.00	E0960
105839-99-0	Bodypoint Monoflex Chest Support Belt Kit - Large <i>21.5"L x 4.5"W Pad Size.</i>	\$204.00	E0960
105850-99-0	Bodypoint Monoflex Chest Support Belt Kit - X-Large <i>24.5"L x 5"W Pad Size.</i>	\$204.00	E0960

Other Options

105260-99-0	Transfer Handles	\$550.00	
105292-99-0	Push Handles, Corpus Seat	\$489.00	
106210-99-0	RAM Self-Leveling Cup Holder	\$93.00	
106806-99-0	RAM X-Grip Phone Holder	\$93.00	
107100-99-0	RAM X-Grip 10" Tablet Holder w/ Swing-Away Mount	\$467.00	

Accessories

107189-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. No charge with 4 power seat function configurations.</i>	\$231.00	
	Mount on Left Armrest Mount on Left Seat Rail Mount on Right Armrest Mount on Right Seat Rail		
105805-99-0	Medical Necessities Bag Hooks	\$111.00	
102059-99-0	Medical Necessities Bag	\$321.00	
105262-99-0	Joystick Cutout Tray - 16" W <i>Must also select hardware.</i>	\$438.00	E0950
105263-99-0	Joystick Cutout Tray - 18" W <i>Must also select hardware.</i>	\$438.00	E0950
100578-99-0	Small Angle Adjustable Tray <i>Must also select hardware.</i>	\$339.00	E0950
	Mount on Left Mount on Right		
10556	Adj Removable - Tray Hardware - 1 Each	\$306.00	
108230-99-0	Permobil Adjustable Tray - 13"W <i>Must also select hardware.</i>	\$284.00	E0950
108228-99-0	Permobil Adjustable Tray - 15"W <i>Must also select hardware.</i>	\$284.00	E0950
108227-99-0	Permobil Adjustable Tray - 17"W <i>Must also select hardware.</i>	\$284.00	E0950
108229-99-0	Permobil Adjustable Tray - 19"W <i>Must also select hardware.</i>	\$284.00	E0950
Permobil Adjustable Tray Hardware	Left Side Adjustable Right Side Adjustable	\$208.00	
105224-99-0	Oxygen Holder, Corpus	\$301.00	E2208
105223-99-0	Crutch Holder, Corpus	\$301.00	E2207

Alternative Drive Systems

Omni2 Control Module required for all alternative drive systems.

Alt Proportional Joysticks

109525-99-0	mo-vis Micro Joystick for Omni <i>Operating Force = 8 g (0.28 oz). Neutral to Full Deflection = 0.13" (3.3 mm).</i>	\$3,200.00	E2312
109527-99-0	mo-vis Multi Joystick for Omni <i>Operating Force = 50 g (1.8 oz). Neutral to Full Deflection = 0.31" (8 mm).</i>	\$3,200.00	E2312
103894-99-0	Compact Joystick Lite for Omni - (9-pin connector) <i>Operating Force = 50 g (1.8 oz). Neutral to Full Deflection = 0.7" (17mm).</i>	\$2,003.00	E2373
103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 220 g (7.8 oz). Neutral to Full Deflection = 1.1" (28mm).</i>	\$1,771.00	E2373

Alt Proportional Joystick Mounts | *Please see Alt Drive Accessories section for additional mounting options.*

107857-99-0	Swing-Away Chin Boom for Micro/Multi Joystick - Left	\$465.00	E1028
107947-99-0	Swing-Away Chin Boom for Micro/Multi Joystick - Right	\$465.00	E1028
107946-99-0	Swing-Away Chin Boom for Compact Lite Joystick - Left	\$465.00	E1028
107948-99-0	Swing-Away Chin Boom for Compact Lite Joystick - Right	\$465.00	E1028
104690-99-0	Alt Joystick Swing-Away Mount Kit - Armrest	\$310.00	E1028
	Mount on Left Mount on Right		
107135-99-0	Alt Joystick / Omni Combo Swing-Away Mount Kit - Armrest	\$586.00	E1028
	Mount on Left Mount on Right		

Alternative Drive Systems

104143-99-0	Alt Joystick Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	\$966.00	E1028
10752	Chin Control Harness - Large <i>For neck width up to 7½".</i>	\$687.00	E1028
10753	Chin Control Harness - Small <i>For neck width up to 5½".</i>	\$687.00	E1028

Alt Proportional Joystick Accessories

104904-99-0	Dual Switches for Compact Joystick - Fixed Mount Mount Left of JS Mount Right of JS Mount front of JS Mount back of JS	\$288.00	
104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mount	\$518.00	

Total Control Head Array

109558-99-0	Permobil Total Control Head Array <i>Features an occipital pad sensor and two lateral arms configured with proximity sensors. Must select swing away hardware and mounting hardware below to complete assembly.</i>	\$4,819.00	E2330
11575	Swing-Away Hardware for Side Sensor - Left	\$263.00	E1028
11576	Swing-Away Hardware for Side Sensor - Right	\$263.00	E1028
11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the notes section.</i>	No charge	
109563-99-0	Adj/Removable BodiLink HA Hardware - MAX 8	\$365.00	E1028
109564-99-0	Adj/Removable BodiLink HA Hardware - MAX 8 Flip Down <i>Flip-down joint that allows head array to swing back 120°.</i>	\$501.00	E1028
10976	(Optional) Add Left Lateral Arm w Proximity Sensor	\$688.00	E1028
10977	(Optional) Add Left Lateral Arm w Egg Switch	\$563.00	E1028
10978	(Optional) Add Left Lateral Arm w Pad (No Sensor or Switch)	\$459.00	E1028
10979	(Optional) Add Right Lateral Arm w Proximity Sensor	\$688.00	E1028
10980	(Optional) Add Right Lateral Arm w Egg Switch	\$563.00	E1028
10981	(Optional) Add Right Lateral Arm w Pad (No Sensor or Switch)	\$459.00	E1028

ASL FUSION Head Array

109582-99-0	ASL FUSION Proportional and Digital Head Array <i>Includes occipital (ASL 935) and lateral (ASL 933) pads, and pulse mechanical switch (ASL 320BK) for reset/mode changes. Must select hardware listed below to complete assembly.</i>	\$6,715.00	E2328
11580	Swing-Away Hardware for Lateral Sensor - Left	\$342.00	E1028
11581	Swing-Away Hardware for Lateral Sensor - Right	\$342.00	E1028
11582	Swing-Away Hardware for Reset/Mode Switch - Right	\$342.00	E1028
11583	Adj/Removable Headrest Mounting Hardware <i>Includes ASL 606 and BodiLink receiver.</i>	\$342.00	E1028
109583-99-0	ASL FUSION Display, Programmer and Attendant Control w/ Mount <i>Includes ASL 165 and quick release mount.</i>	\$730.00	K0108

Sip & Puff Systems

10573	Sip & Puff System for Omni	\$1,937.00	E2325, E2326
10574	Sip & Puff System with Fiber Optic Opti-Stop	\$2,577.00	E2325, E2326
10762	Adjustable Mounting Hardware - Sip & Puff <i>Required with any sip & puff system.</i>	\$310.00	

Alternative Drive Systems

Proportional Foot Control

10767	Proportional Foot Control <i>Not available with ActiveReach packages.</i> Mount on Left Mount on Right	\$3,690.00	E2328
10766	Adjustable Mounting Hardware - Foot Control <i>Required with proportional foot control.</i>	\$310.00	

Alt Drive System Setup

Omni2 Options | *Required for all alternative drive systems.*

108500-99-0	Omni2 Control Module for R-net - Corpus <i>3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology. Must select mounting hardware below. Specify alternative input device being used:</i>	\$2,061.00	
103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet Mount on Left Mount on Right	\$310.00	

Alt Drive Accessories

PERMOfix Mounting Solutions

106655-99-0	PERMOfix UniTrack Fixed Mount Kit	\$312.00	
106656-99-0	PERMOfix UniTrack Removable Mount Kit	\$549.00	
106657-99-0	PERMOfix UniTrack Swing-Away Mount Kit	\$786.00	
105907-99-0	Proximity Sensor + Swivel Assembly <i>Swivel assembly allows sensor to be mounted with PERMOfix. Since Proximity Sensors require power, CANNOT be used as On/Off switches.</i>	\$397.00	E1028
105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix.</i>	\$268.00	
105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix.</i>	\$162.00	
106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Microlight, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix.</i>	\$142.00	

Switches and Cables

101769-99-0	Microlight Switch Qty:	\$217.00	
101767-99-0	Buddy Button Qty:	\$183.00	
106293-99-0	Egg Switch Black Qty:	\$183.00	
101771-99-0	Mini Cup Switch Qty:	\$183.00	
107500-99-0	Proximity Switch Qty:	\$723.00	
10772	Dome Switch, Backrest Mount w/ Gooseneck Qty: <i>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights.</i> Mount on Left Mount on Right	\$573.00	
10846	Dome Switch, Backrest Swing-Away Mount w/ Gooseneck Qty: <i>Swing-away, mounts at top of backrest shell. Custom charge may be applied if ordered with short backrest heights.</i> Mount on Left Mount on Right	\$847.00	
105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	\$36.00	
105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	\$39.00	
105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	\$47.00	
105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	\$51.00	
105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	\$39.00	

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- The client's name or code and the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.