

Adjusted by:

Tel:

Save the form to the desktop and open it with Adobe Reader.

E-mail:

Date:

The price changes depending on the choices made in the form £:

Standard configuration includes; chassis, seat, calf strap, backrest lock, standard colour option, backrest, carbon rear axle, castor and fork, titanium foot rest or foldable footplates, push to lock S3 brake, Panthera crosslaced steel spoke rear wheel, titanium pushrims and Schwalbe right Run tires.

Choice required: Max. user weights : Seat width **27-42 cm 100 kg**, Seat width **45-50 cm 150 kg**, (S3 Large) Seat width **39-50 cm 150 kg**

Optional **Upcharge /Add-price** **Marked yellow** Click on "**blue underlined text**" for more info or image

Original design
1) Chassis

Alternative options

Other options

- [S3](#) **Seat width 33-50** cm, Seat depth 40 cm
 - [S3 Short](#) **Seat width 30-45** cm, Seat depth 27-33 cm (seat width 30-33), 35 cm (seat width 36-45)
 - [S3 Short Low](#) **Seat width 33-39** cm, Seat depth 35 cm
 - [S3 Large](#) **Seat width 39-50** cm, Seat depth 45 cm
 - [S3 Long](#) **Seat width 33-50** cm, Seat depth 45 cm
 - [S3 0°](#) **Seat width 33-50** cm, Seat depth 35, 37.5, 40, 42.5 cm
 - [S3 Short Abd](#) **Seat width 27-33** cm, Seat depth 27-33 cm
- Clear

2) Chassis colour

- [Black](#) - Anodic Black
 - [Gray](#) - Metallic gray
- Clear
- Optional colour** - A delivery time of 20 days for optional colour. Fill in art.no. in the square

3) Seat width

Check allowed seat width at "1) Chassis"

Clear	27 cm	30 cm	33 cm	36 cm	39 cm	42 cm	45 cm	50 cm
-------	-------	-------	-------	-------	-------	-------	-------	-------

4) Seat depth

Clear	35 cm	37.5 cm	40 cm	42.5 cm	45 cm	Option only needed for S3 0°
-------	-------	---------	-------	---------	-------	------------------------------

5) Back rest

Height (HH) [INFO](#)

Standard 523BBHH VER2	20 cm	35 cm	Tapered 525BBHH VER2	INFO	Tapered 6 cm
	25 cm	40 cm	Curved 524BBHH VER2		Flared 3 cm
Clear	30 cm	45 cm	Frame from solid backrest		Curved/Tapered
Clear	Clear				

6) Rear axle

Carbon std 33023BB	Carbon +15 mm 33023BC*	* seat width BB+1	Other rear axle or seat height
Clear	Use this rear axle when choosing arm support at Accessories and/or increased seat height		Increase the seat height

7) Castor and fork

120 mm	90 mm
150 mm	150 mm
Clear	

8) Footrest

V-shape for S (titanium) 36001BB	V-shape extended, reinforced (titanium) * 36111BB	Carbon footplate 65040BB
Small heel bar for S series 3600134 (width 33)	* std for S3 Large	
Clear	Foldable footplates; art.no.	
	Heel strap for footplates	
	Plastic plate for footrest 36010BB	

9) Brake

Brake S3 4630000	One hand brake, right	art.no.	Assistant brake *
Clear	One hand brake, left	Art.no.	
			* Not for short / low chassis

All prices for rearwheels are for pair

10) Rear wheel - without pushrim

[Rear wheel sizes](#)

[24" std](#) [art.no.](#)

Clear

22" [standard](#)

24" [Spinergy X](#) black spokes

25" [Spinergy X](#) black spokes

26" [Spinergy X](#) black spokes

24" [Spinergy X](#) white spokes

25" [Spinergy X](#) white spokes

26" [Spinergy X](#) white spokes

[One arm drive](#)

[Rear wheel 24x2" MTB](#)
Titanium pushrims

Tyre - choose tyre

[Schwalbe Right Run](#) 22", 24", 25", 26"

Clear

[Rugged 24"](#)

[Schwalbe Marathon Plus](#), 24", 25", 26"

[Schwalbe One](#) 24", 25"

All prices for tyres are for pair

[PU 24" \(massive\)](#)

reTyre does not fit PU tyres

[reTyre skins 24"](#)

[reTyre skins 25"](#)

[reTyre skins 26"](#)

Clear

11) Pushrim

[Aluminium](#)

1902430 (x2 / wheelchair)

Clear

[Paragrip](#)

19051DDBT x 2

[Tetragrip](#)

19052DDBT x 2

All prices for pushrims are for pair

[Maxgrip 22", 24"](#)

1903060 x 2 (60=24", 80=22"
(contains latex)

[Titanium](#)

190DD11 (x2 / wheelchair)

Accessories

Side guard [art.no.](#)

[Side guard S3](#)

[video](#)

Clear

Arm support [art.no.](#)

[Arm rest S3](#) *

[video](#)

* 26" std for 24" wheels for Short Low
26" std for Assistant brake

Side guard Carbon [art.no.](#)

[Without fender](#)

[Small fender](#)

[Large fender](#)

Push handles [Video](#)

[Height adjustable S3](#)

4252060 VER2

[Height adjustable S3 +10 cm](#)

4252063 VER2

Max. userweight 100 kg

Clear

Push handle foldable

[Foldable S3](#)

4245000

[Video](#)

Push bar

[Push bar S3](#)

42400BB VER2

Max. userweight 100 kg

Anti-Tip [Video](#)

[Anti-tip](#) std seat height

4353000

[Anti-tip](#) S3 Short Low

4355000

[Anti-tip](#) S3 Large

4353080

Clear

[Other seat height](#)

Art.no.

Other

[Tube protection, front](#)

9900251

[Tube protection, rear](#)

9900251

[Crutch holder S3](#)

9900301 (x1/ crutch)

[Backrest wedge](#)

4710012

[Freewheel](#)

Cushion [Wedge cushion](#)

[2,5 cm](#)

572BB40

[Wedge cushion](#)

57710BB

[5 cm](#)

575BB40

Clear

Spoke protectors

[24" logo](#)

4500000

[22" logo](#)

4500100

[24" no logo](#)

4500011

[22" no logo](#)

4500111

Clear

Spoke protectors with theme

[22" Splash](#)

4500322

[24" Splash](#)

4500324

[22" Car rim](#)

4500122

[24" Car rim](#)

4500124

[22" Wild animals](#)

4500422

[24" Wild animals](#)

4500424

[22" Smiley](#)

4500222

[24" Smiley](#)

4500224

Clear

Chose art.no. for Freewheel according to the footrest height of the wheelchair
Read more about Freewheel on the [homepage](#).

SPLASH CAR RIM WILD ANIMALS SMILEY

Set-up adjustments

Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard = 20 mm

A= _____mm

Footrest height - Foot plates

Note measure (A) and (B) for desired footrest

A= _____mm B= _____mm

Rear axle adjustment

Note measure (A) for desired balance
Standard = 60 mm. at raised rear axle = 51 mm
Standard = 40 mm for S3 Short, Short Low
Standard = 31 mm Short ABD (6 degree camber)

A= _____mm

Extra