

Room to grow

Medium chair at its smallest dimensions

- 1 Backrest height **15½"**
- 2 Distance between laterals **6½"**
- 3 Seat depth **11"**
- 4 Seat width **8½"**
- 5 Armrest height above seat **7"**
- 6 Seat height above footboard **7¾"**




Medium chair at its largest dimensions

- 1 Backrest height **19"**
- 2 Distance between laterals **12"**
- 3 Seat depth **16"**
- 4 Seat width **14"**
- 5 Armrest height above seat **10½"**
- 6 Seat height above footboard **15¼"**




Medium standard base
R840 Rifton Activity Chair

User dimensions (inches)	
Height	40 – 56
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.	
Item dimensions (inches)	R840 medium Standard base
Frame width	short legs: 23¼ long legs: 24½
	short legs w/ casters: 24½ long legs w/ casters: 25½
Seat height above floor	short legs: 12½ – 16½ long legs: 19 – 23
	short legs w/ casters: 16 – 20 long legs w/ casters: 19 – 23
Seat angle (tilt-in-space) – degrees	15° forward, 15° back
Backrest angle – degrees	5° forward, 20° back
Footboard knee angle – degrees	45° – 110°
Seat height above footboard	7¾ – 15¼
Seat width with hip guides without hip guides	8½ – 11½ 14
Seat depth	11 – 16
Armrest height above seat	7 – 10½
Distance between laterals	6½ – 12
Backrest height	15½ – 19
Headrest height above seat	17½ – 24½
Max. working load (lb)	150



Medium hi/lo base
R850 Rifton Activity Chair

User dimensions (inches)	
Height	40 – 56
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.	
Item dimensions (inches)	R850 medium Hi/lo base
Frame width	27½
Seat height above floor	12 – 23
Seat angle (tilt-in-space) – degrees	15° forward, 25° back
Backrest angle – degrees	5° forward, 20° back
Footboard knee angle – degrees	45° – 110°
Seat height above footboard	7¾ – 15¼
Seat width with hip guides without hip guides	8½ – 11½ 14
Seat depth	11 – 16
Armrest height above seat	7 – 10½
Distance between laterals	6½ – 12
Backrest height	15½ – 19
Headrest height above seat	17½ – 24½
Max. working load (lb)	150

Retail order form - January 2023

Use dimension chart to select appropriate size. All Rifton products are made to order and are non-returnable and non-exchangeable.



* Required		Standard R840 medium	Hi/lo R850 medium
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	Adjustable	R841 \$956	R851 \$1,872
	Adjustable w/ back spring	R842 \$1,100	R852 \$2,045
* Base	Without spring	R843 \$557	R853 \$4,025
	With spring	R844 \$669	N/A
* Arm supports (pair)	Armrests	R868 \$391	R868 \$391
	Forearm prompts	R869 \$977	R869 \$977
* Adjustable legs	Short	R893 \$88	N/A
	Long	R894 \$131	N/A
	Short with casters	R895 \$472	N/A
	Long with casters	R896 \$514	N/A
* Pads	Tan	R835 \$391	R835 \$391
	Blue	R836 \$391	R836 \$391
	Red	R837 \$391	R837 \$391
	Pink	R974 \$391	R974 \$391
	Purple	R975 \$391	R975 \$391
	Green	R976 \$391	R976 \$391
Headrest	Contoured	R897 \$382	R897 \$382
	Flat	R898 \$352	R898 \$352
	Adjustable winged	R881 \$472	R881 \$472
Laterals (pair)	Small	R828 \$308	R828 \$308
	Small with chest strap	R829 \$391	R829 \$391
	Large	R838 \$340	R838 \$340
Push handles (recommended for all chairs with wheels)		R878 \$251	Included
Tray (requires armrests)		R668 \$433	R668 \$433
Handhold		R649 \$39	R649 \$39
Chest strap (wide)		R863 \$131	R863 \$131
Butterfly harness	Regular	R864 \$194	R864 \$194
	Slim cut	R882 \$194	R882 \$194
Thigh belt		R865 \$131	R865 \$131
Pelvic harness		R866 \$251	R866 \$251
Hip guides (pair)		R867 \$331	R867 \$331
Abductor	Small	R884 \$194	R884 \$194
	Large	R891 \$228	R891 \$228
Adductors (pair)		R892 \$308	R892 \$308
Leg prompt		R858 \$352	R858 \$352
Footboard (recommended for chairs with long legs)		R848 \$699	Included
Ankle straps (pair)		R899 \$131	R899 \$131
Sandals (pair)		R954 \$151	R954 \$151
Wedges (pair, require sandals)		R698 \$110	R698 \$110
Lumbar and seat support kit		R887 \$131	R887 \$131
Backrest filler pad		R849 \$194	R849 \$194
Whitmyer headrest adapter		R804 \$31	R804 \$31
Total		R840 \$	R850 \$