



Permobil Paediatric Equipment

permobil

Explorer Mini

The Explorer Mini is designed to introduce self-initiated movement for children with mobility impairments aged 12-36 months.



Promotes stable, upright posture

Explorer Mini	
Max speed	2.4 km/h
Base width	19"
Base length	25"
Weight limit	15.9 kg
Driving range	5.6 km

Height adjustable

Multiple weight-bearing surfaces

Snooza Chair

The all-foam Snooza Chair with its planar angles uses gravity to keep the occupant safe, providing good postural stability while they are out of their wheelchair.

Snooza Chair	
Size options	4
Overall width	18" - 23.5"
Overall length	24" - 36"
Seat width (hip)	8" - 12"
Seat depth	10" - 18"



Convenient carry handles

A range of accessory options



Multiple colour choices

K300 PS Jr

A versatile paediatric power chair, the expandable seat pan allows for a custom fit no matter how fast a child grows.

Easy to grow with the child

Power seat elevator and power tilt

Flip back armrests for transfers

K300 PS Jr	
Max speed	8 km/h, 10 km/h
Base width	24"
Base length	40"
Weight limit	75 kg
Driving range	25 km

M300 PS Jr

The M300 PS Jr is a mid-wheel drive variant, boasting similar adjustability to the K300 PS Jr with a tighter turning radius.

Virtually unlimited adjustability

Tight turning radius for indoor use

M300 PS Jr	
Max speed	10 km/h
Base width	24"
Base length	36.5"
Weight limit	75 kg
Driving range	20 km

Koala

Designed to grow with a child for years to come, the Koala is the ideal paediatric wheelchair for toddlers. Featuring a rugged front-wheel drive system.

Koala	
Max speed	6.5 km/h
Base width	23"
Base length	32"
Weight limit	35 kg
Driving range	24 km



Optional remote stop

Integrated elevation plus tilt feature

Ideal for toddlers

Optional footplate bumper

K450 MX

The all-wheel suspension system of the K450 MX provides a smooth ride while its unique seat-to-floor function is perfect for exploring and socialising.

K450 MX	
Max speed	8 km/h
Base width	25"
Base length	35"
Weight limit	60 kg
Driving range	27 km



Growable

All-wheel suspension

Powered seat-to-floor function



Height-adjustable and easily removable push handle

Fully adjustable and ready for growth

TiLite Twist

Designed from the ground up for kids, the lightweight Twist encourages exploration, independence, and playfulness.

TiLite Twist

Seat width	8" - 15"
Seat depth	8" - 16"
Weight limit	75 kg
Transport weight	5.4 kg



Maintains access to wheels throughout growth

Carbon fibre frame panels

Unique footrest adjustment for ultimate fit

TiLite Pilot

Ultra lightweight with built-in forward growth adjustability. Allows for 2" of growth in seat width and 3" in seat depth.

TiLite Pilot

Seat width	8" - 14"
Seat depth	9" - 15"
Weight limit	54 kg
Transport weight	4.9 kg

F5 Corpus VS

Suitable for older kids and adults, the F5 VS empowers the user with face-to-face interactions and improved health outcomes.

F5 Corpus VS

Max speed	10 km/h
Base width	25.5"
Base length	43"
Weight limit	136 kg
Driving range	32 km

Ergonomic Corpus®
Seating System

Fully programmable
standing sequences



Our personal approach makes all the difference.

Here at Permobil, we take a personal approach to mobility and seating. This means that every solution we create is inspired by the real lives of people with disabilities and their families. Each of our employees strive every day to solve the challenges of mobility, positioning, comfort and safety. These are challenges that our users face every day.

It also means that, no matter what our users are striving for, we are there to support them using the latest technologies and clinical evidence.

We believe that, “Every person with a disability has the right to have his or her needs compensated as far as possible by aids with the same technical standards as those we all use in our everyday lives” - Dr. Per Uddén, Permobil Founder.

For more information, visit permobil.com or contact us using the details below.

permobil

Permobil Australia | 1300 845 483 | permobil.com | info.au@permobil.com

Follow us!



PermobilAus



@permobilaustralia

Paediatric equipment AU Rev C